Betriebsanleitung  
Operating instructions  
Notice d'instructions

Hei-VAC Control

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Seite 4 – 39

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Translation of the original instructions

Page 44 – 79

Traduction de la notice originale

Page 84 – 119

Zertifikate/Certifications

# Contents

## Introduction

About this document.....	44
Typographic conventions.....	44
Copyright protection.....	44

## General notes

Directives applied, product certification.....	45
California Residents.....	45
Copyright notice.....	45
Residual risk.....	45
Intended use.....	45
Reasonably foreseeable misuse.....	45
Compliant use.....	46
Transportation.....	46
Storage.....	46
Acclimatization.....	46
Permissible ambient conditions.....	46

## Safety

General safety instructions.....	46
Electrical safety.....	46
Operational safety.....	47
Occupational safety.....	47
Personal protective equipment (PPE).....	47
Environmental protection.....	47
Biohazard.....	48
Other regulations.....	48

## Device description

Mechanical design.....	49
Control elements.....	49
Connections.....	49
Principle of operation.....	50
Display and signaling.....	50
Display icons.....	51
Signal sounds.....	53
Key combinations.....	53

## Commissioning

Set up the device.....	54
Table top version .....	54
Mounting on the Hei-VAP rotary evaporator .....	54
Power supply .....	56
Prepare wall power supply plug .....	56
Connect to power supply.....	56
Vacuum connection.....	57

## Operation

Initial start-up .....	58
Switching the device on/off .....	58
Operation .....	59
Menu navigation .....	59
Function Starting/stopping a program .....	59
Operating modes.....	60
Pump down .....	60
Vac control .....	61
Change operation mode.....	62
Programs .....	62
Configuration .....	63
Readjustment .....	63
RS-232 .....	63
Sensors.....	64
Display .....	64
Autostart.....	65
Factory settings.....	65
Differential pressure measurement .....	65
Graphic (curve) .....	66
Venting .....	66

## Troubleshooting

Troubleshooting.....	67
Error display .....	67
Error messages .....	68
Error of external components .....	71
Reset .....	71

## Appendix

Technical specifications .....	72
Wetted materials.....	73
Interface commands .....	73
Pin assignment (RS232).....	73
Scope of delivery.....	74
Accessories .....	74
Device service.....	75
General cleaning instructions.....	75
Clean venting valve (air admittance valve) .....	75
Clean controller pressure sensor.....	75
Sensor readjustment.....	75
Adjustment at atmospheric pressure .....	76
Adjustment under vacuum .....	76
Adjustment at a reference pressure .....	77
Repairs .....	77
Maintenance.....	77
Disposal.....	78
Warranty statement.....	78
Contact details .....	78
Certificate of decontamination.....	79

## About this document

This operating instructions manual describes the features and operation of vacuum controllers, type Hei-VAC Control.

The operating instructions are an integral part of the described device!

## Typographic conventions

Standardized symbols, signal words and highlighting elements are used in this document to warn of hazards and to identify important information or rather special text contents.

Symbol	Signal word / explanation
	<p>Warning symbols in combination with a signal word indicate dangers:</p> <p><b>DANGER</b> Indicates an immediate dangerous situation. Failure to respect the indications will result in death or serious injury.</p> <p><b>WARNING</b> Indicates a potential danger. Failure to respect the indications will result in serious injuries.</p> <p><b>CAUTION</b> Indicates a potential hazard which, if not avoided, damage to property and minor to moderate injuries can occur.</p>
	<p>Mandatory signs are used to indicate important and useful information on handling a product.</p> <p>This information is used to ensure operational safety and to maintain the value of the product.</p>
[GUI]	<p><b>Parameter</b> designations, <b>display texts</b>, and <b>device labels</b> are typographically highlighted in the text and in tables to facilitate their assignment on the device.</p>
→	<p>The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.</p>

## Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Instruments GmbH & Co. KG. Any violation is subject to compensation for damage.

## Directives applied, product certification

### CE marking

The device meets all requirements of the following directives:

- Low voltage Directive 2014/35/EU
- EMC Directive 2014/30/EU

## California Residents

Important information for California residents regarding Prop 65. Please visit [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov) for more information.

## Copyright notice

The software implemented in this product is protected by copyright laws. The rights holder is Heidolph Instruments GmbH & Co. KG, Walpersdorfer Straße 12, 91126 Schwabach/Germany. Any open source software components in this software are excluded from our copyright. Further information is available in the service area on our website [www.heidolph-instruments.com](http://www.heidolph-instruments.com).

Violations of copyright (e.g. unauthorized use or modification of the software) may result in civil claims (e.g. omission, compensation for damages) and/or entail criminal penalties within the meaning of the legislation of the destination country.

## Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

## Intended use

Hei-VAC Control is a laboratory instrument for measuring and/or controlling a vacuum in installations defined for this purpose.

Any other use of this device is not considered as intended!

Due to its design, the device in its delivery condition may only be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

## Reasonably foreseeable misuse

Additional measures may be necessary, and/or specific directives and safety regulations may have to be observed for use under conditions or for purposes deviating from the intended use (see, for example, section „Other regulations“ on page 48). Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

## Compliant use

The user is generally responsible for evaluating the conformity of his application and, if necessary, for taking additional measures.

## Transportation

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device. Keep the original packaging in a dry and protected place for later use.

## Storage

Always store the device in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

## Acclimatization

After each transportation and after storage under critical climatic conditions (e.g. high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

## Permissible ambient conditions

The device is designed for indoor use only. The device is **NOT** suitable for outdoor use! The device is **NOT** suitable for use in potentially explosive areas!

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

## General safety instructions

- Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.
- Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself, the power supply cord and, where applicable, on connected devices or the supply connections.
- If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.
- Only use the device in accordance with the regulations on intended use („Intended use“ on page 45).

## Electrical safety

- Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.
- Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).
- Only use the short circuit protected wide range power supply provided with the device.
- The device must only be supplied with power from a properly grounded mains socket-outlet.

- Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Instruments.
- Always disconnect the device from the power supply system before carrying out any maintenance, cleaning or repair work.

### **Operational safety**

- Do not make any unauthorized changes or modifications to the device!
- Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!
- Rectify malfunctions or faults on the device immediately.
- Switch off and disconnect the device from the power supply, preventing reconnection, if it is not possible to eliminate the malfunction or rectify the fault immediately.
- Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

### **Occupational safety**

- Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.
- Depending on the process, an explosive mixture can be formed in the installations or other dangerous situations could occur. Thus, never control critical processes unattended!

### **Personal protective equipment (PPE)**

The operating company must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The corresponding instruction of the personnel is solely within the operating company's responsibility.

### **Environmental protection**

When processing environmentally hazardous substances, take appropriate measures to avoid hazards to the environment.

The evaluation of corresponding measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

## Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking the device with a biohazard warning symbol.

The evaluation of corresponding measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

## Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Instruments.

The operator is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

## Mechanical design

### Control elements

#### 1 LC display

[VENT] button

- Keystroke < 2 sec = momentarily venting, control continues.
- 2 ▪ Keystroke > 2 sec = venting to atmospheric pressure (max. 1050 mbar), control stops.
- Keystroke while venting = venting stops.

[MODE] button

- 3 ▪ During operation: To switch from [Pump down]/[Auto mode] to [Vac control].
- In Standby mode: open [Functions] selection.

[START/STOP] button


- 4 ▪ Activate or deactivate controller function.
- Start program

[Standby] button

- 5 ▪ As soon as operating voltage is applied on the connected pump, the device can be activated or deactivated using this button

- 6 Rotary switch to select and activate the device functions and to set the operating parameters.

### Connections

#### 1 Port for Heidolph vacuum valve

#### 2 Port for Heidolph vacuum valve

#### 3 Connection socket power supply

#### 4 serial interface RS 232 C (Sub-D)

#### 5 Hose nozzle, vacuum connection

#### 6 Venting, connection for external venting, e. g., inert gas

#### 7 Rating plate

#### 8 Spring clip as fixation for built-in version or foot for table top version

## Principle of operation

The controller was designed for applications requiring a regulated vacuum. The device work with a two-point control mode to switch an in-line isolation valve.

The controller can be freely programmed: ten program memory locations are available for processes, whereas every program consists of up to ten program steps (time and pressure) with control functions such as venting, pump down and ramp function.

The controller measures the relative pressure with regard to a reference sensor (VSK 3000).

The vacuum process is managed via connected vacuum pumps as well as in-line isolation- and/or air admittance valves: process vacuum, cooling water and venting on demand.

Valves and/or vacuum pumps are necessary to operate the controller. Without those components the controller can only be used as vacuum measurement device.

A venting valve and a pressure sensor (= capacitive ceramic diaphragm vacuum sensor) are installed in the controller.


The pressure sensor is chemically highly resistant with a high measurement accuracy, regardless of the gas type.

When switched on, the controller verifies automatically the actual system configuration (connected components).

Operating elements are the rotary switch control, buttons on the control panel and full text menus on the display.

## Display and signaling

After booting, the actual pressure and the last activated operation mode appear on the controller display:

Status bar (title bar)

Operating mode pump down, vacuum controller, program

1 Process time: hh:mm:ss (only during running process)

2 Bar graphic: graphical display of actual pressure




3 Numerical value: actual pressure = digital pressure display

VAC: Vacuum

4 mbar: Pressure unit according to pre-setting (mbar, Torr, hPa)

ATM: Atmospheric pressure

## Display icons

Symbol	Meaning
	Vacuum control is running (animation)
<b>00:00:00</b>	Process time; runtime vacuum control (hh:mm:ss)
	Pump is running; in combination with percentage sign = motor speed (only for VARIO systems)
	Venting valve is active, open (VENT); Flashing cycle: continuous venting switched on.
	Coolant valve switched on, open
	In-line valve switched on, open
	Emission condenser (Peltronic) connected
	Level sensor activated (only when level sensor is connected)
	Pump down – continuous pumping
	Pump down: lower pressure limit reached VACUU-LAN: pump down to set pressure Vac control: for 2-point control – pump down to set pressure

## Device description

VACUU-LAN: pressure increase to switch on pressure  
Vac control: preset maximum exceeded

VARIO control: pump down to set point  
Auto mode: Searching for boiling point and adjusting the process pressure while boiling pressures are changing.

VARIO control: Set pressure reached  
Auto mode: Boiling point reached and adjusting the process pressure.

2-point control: pressure in hysteresis, pump switched on

2-point control: pressure in hysteresis, pump switched off

Turbo mode switched on (for VARIO@ pump in combination with turbo-molecular pump).

Clock

- Program completed (The clock icon keeps flashing until the Start/Stop button has been pressed to acknowledge the end of program.)
- VACUU-LAN mode: delay time elapses

Lock – operation locked

HI mode for pump down = optimum speed for the respective pressure.

Percentage value for pump down motor speed.

Set value for Vac control.

Flashing: Warning!

## Signal sounds

Setting Sound On in the controller configuration is required to hear the audio signals.

<p>1x  </p>	<p>Short beep for each keystroke.</p>
<p>2x  </p>	<p>Audio warning for error indication. In short intervals a number of warning beeps are to be heard. This audio warning is active until error clearance or reset. See also section „Troubleshooting“ on page 67.</p>

## Key combinations

Different special menus and special functions can only be accessed through key combinations:

### Controller in Standby mode

- Press and hold the rotary switch, press [Standby]:
  - Call up menu [Language selection]
  - Call up menu [Pressure unit]
- Press and hold [VENT], press [Standby]:
  - Call up menu [Function]

### Menu [Function] opened

- Press and hold [MODE], press rotary switch:
  - Activate Vacuubus configuration (frame marker).
  - Transmit parameters from selection [Vario init], confirm

### Mode [Vac control]

- Keep pressed [Rotary switch] and turn it
  - Adjust set vacuum quickly

### Mode [Pump down]

- Keep pressed [Rotary switch] and turn it
  - Adapt motor speed (only in combination with VARIO®)

### Set up the device

#### CAUTION

- Measuring error due to an obstructed vacuum line. Prevent overpressure > 1060 mbar (> 795 Torr) in the piping system.
- Condensate can falsify the measurement. Position the vacuum hose in such a way that condensate cannot flow towards the controller and its vacuum sensor. No liquid should accumulate inside the vacuum hose. Install vacuum hoses in such a way that condensate cannot flow into the controller.
- Particles, liquids or dust may not enter the controller. Install suitable separators and/or filters at the intake of the vacuum system. Appropriate filters are for example chemically resistant, and resistant to clogging.
- Installation and operation in areas where explosive atmospheres can occur is not allowed.
- The use of venting valves is only permitted if thereby explosive mixtures normally do not occur in the interior of the controller or, if they do occur, are likely to do so only infrequently and for a short period. If necessary vent with inert gas.

#### Table top version

The table top controller can be installed and connected directly on top of the work bench, (e.g. on the laboratory table).

The table top version is supplied with a hose nozzle. It should be positioned in a way that the connected vacuum hose cannot kink.

#### Mounting on the Hei-VAP rotary evaporator

The Hei-VAC Control can be directly mounted on a Hei-VAP rotary evaporator.

- Screw the pre-assembled vacuum valve on the rear panel of the evaporator (two screws):


## Commissioning

- Mount the support rod on the evaporator and fix the controller on the upper part of the support rod. Be sure to check and feel snapping the support rod into the recess on the back of the controller.
- Connect the wall-plug power supply to the controller [1].
- Connect the cable of the vacuum valve to the controller [2].

- Tube connections:
  - Connect the upper connection of the vacuum valve to the upper connection of the condenser [3].
  - Connect the lower connection of the vacuum valve to the vacuum source [4].
  - Connect the measuring and venting tube of the controller to the second upper connection of the condenser [5].

## Power supply

Wall power supply plug: short-circuit-protected multi-voltage power supply with integrated overload

protection and country-specific adapters:

### Prepare wall power supply plug

Take the wall power supply kit and the adapters out of the packaging. Select the mains plug that fits to your mains socket. Connect the mains plug to the metal contacts of the wall power supply plug and slide the mains plug until it locks.

#### Remove mains plug

Press the locking knob on top of the wall power supply plug and remove the mains plug from the power supply.

#### CAUTION

Please install the power supply line of the controller in such a way, that no damage is caused to the cable due to sharp edges, chemicals or hot surfaces.

---

### Connect to power supply

Connect the wall power supply to the controller. Plug the wall power supply into the mains socket. As soon as the green LED on the wall power supply lights up, the device is ready for use.

## Vacuum connection

### CAUTION

Always use vacuum hoses that are designed for the specific vacuum area. Make sure that the vacuum hoses are correctly dimensioned:

- Measure the hose length taking into account the maximum shrinkage.
- Keep all hose lengths as long as necessary and as short as possible.

Before switching on the device, make sure that all hose connections are safely and gas-tight fixed.

Vacuum hoses may not be kinked!

Check all hose connections for leaktightness and tight fit at regular intervals.

Clean all hoses and hose connections at regular intervals or immediately when there are obvious contaminations to avoid damages to components in contact with media due to penetrating liquids/fluids.

Please note the following limits:

- Maximum admissible pressure at vacuum sensor: 1,5 bar (absolute).
- Measuring limit of the controller: approx. 1060 mbar.

Please note that filters can affect the measurement and control functions of the controller.


Depending on design and installation, the controller provides several options for integration into the vacuum system:

Example 1: vacuum hose made of caoutchouc, directly plugged on the hose nozzle:


Example 2: vacuum hose made of PTFE plugged on hose nipple, fixed with union nut:

## Initial start-up

During the initial start-up, the user language and the required display unit have to be selected:

- Make sure that the controller is supplied with the necessary operating voltage (main device on).
- Press and hold the rotary switch and press the Standby button on the controller.
  - The language selection opens.
- Select the required language with the rotary switch and confirm the entry (press rotary switch one time).
  - The display unit selection opens.
- Select the required display unit with the rotary switch and confirm the entry (press rotary switch again one time).

- The initial screen with company logo and firmware version are shown for approximately two seconds. Afterwards, the pressure display appears and the device is ready for use.

## Switching the device on/off

When switching on the connected pump, the device is also switched on or off. The initial screen with company logo and firmware version are shown for approximately two seconds. Afterwards, the pressure display appears.

## Operation

The handling of the controller is menu-driven via the control buttons and the rotary switch of the device.

### Menu navigation

Press the [MODE] button to open the [FUNCTION] selection.

Turn the rotary switch to navigate between the menu items.

- The selected menu item turns black (in the example: menu item [back]).

Afterwards, press the rotary switch to activate the marked mode or to open the configuration menu or to turn back to the higher level menu (select [back]).

Existing submenus are highlighted with three points. Use the rotary switch to mark the required entry and press the rotary switch to open the submenu.

Values that have to be adjusted appear on the right display side:

Each time you press the rotary switch, the marking jumps to the next selection point.

- By turning the rotary switch clockwise or counterclockwise, you can adjust editable values.
- Confirm each entry by pressing the rotary switch.

Function
<b>Pump down</b>
Vac control
Auto mode
Program
VACUULAN
Configuration
----- Back -----

RS-232...
Sensors...
Display...

Display
Brightness 100 %
Contrast 25 %
Sound On

All editable values can be adjusted during a running process as well as when control is stopped.

Modifications are immediately applied and will not be reset when the device is switched off!

In the submenu [Sensors], the display returns to the previous menu only after the selection of a sensor.

When inactive, the display returns back to the pressure display after approx. 20 seconds.

### Function start/stop program

To start/stop the selected function (mode) or the loaded program, use the [Start/Stop] button.

The controller works in delivery status with the default settings of the factory setting.

## Operating modes

The controller offers different modes. To select the required mode, press the [MODE] button one time.

A selected operation mode can be adapted and optimized for the process through the corresponding operation menu. The settings in an operation menu include mainly: motor speed, set vacuum or time presettings.

### Call up menu of a mode and adjust values

- Press the [MODE] button and select the required mode with the rotary switch.
- Press the rotary switch to open the corresponding submenu.
- Turn the rotary switch to mark a menu item.
- Press the rotary switch again to adjust the value in the Edit mode.
- Confirm the input value by pressing the rotary switch again.

## Pump down

Continuous pump down with pressure and time presettings. VARIO®: pump down with adjustable motor speed (pumping

speed) and continuous speed control.

Pump down <span style="float: right;">978mbar</span>	
Speed	100 %
Minimum	Off
Duration	Off
----- Graphic -----	
----- Back -----	

Parameters	Meaning
Minimum (mbar, Torr, hPa)	Vacuum set point that shall be reached by pumping down. Once reached, the controller switches off the vacuum pump or closes the in-line valve. Adjustment Range: Off; 1-1060
Duration (min)	Presetting process runtime from Start. Adjustment Range: Off; 1-1440

If Minimum and Duration are set to OFF, pump down has to be stopped by pressing the START/STOP button!


### Application example cabinet dryer

Set Minimum to a vacuum value below boiling pressure. The controller will switch the vacuum pump off, once the liquid has completely evaporated.

## Vac control

In this mode, the device controls to a preset vacuum value.

Parameters	Meaning
Set vacuum (mbar)	Setting for lower vacuum level for 2-point control or precisely for VARIO pump. Adjustment range: Turbo <sup>1)</sup> ; 1–1060
Hysteresis (mbar, Torr, hPa)	Only for VMS+vacuum pump or with in-line valve: Control range for 2-point control. Adjustment range: Auto <sup>2)</sup> ; 1–300
Maximum (mbar, Torr, hPa)	Setting for upper vacuum level. Once reached, the vacuum control switches off. Adjustment Range: Off; 1–1060
Duration (min)	Presetting process runtime from Start. Adjustment Range: Off; 1–1440

<sup>1)</sup> Turbo mode: auto-adapting vacuum control for best ultimate vacuum, best initial vacuum for operation with a turbomolecular pump.

<sup>2)</sup> Hysteresis value Auto

Set vacuum (mbar)	5	10	50	80	100	200	500	700	900	1000
Hysteresis (mbar)	2	2	5	8	9	17	40	55	71	78

### Adjust set vacuum during operation

In mode [Vac control], the set vacuum can be adapted during running operation:

- Adapt set vacuum - quick tuning
  - Keep pressed the rotary switch and turn it to the right: Increase set vacuum (venting).
  - Keep pressed the rotary switch and turn it to the left: Reduce set vacuum (pump activated).
  - Controller controls to the new set vacuum which is displayed while releasing the rotary switch.
- Adapt set vacuum - fine tuning
  - 1 detent = 1 pressure value (mbar, Torr, hPa)
  - Controller controls to new set vacuum.

#### Application example filtration

Set the set vacuum higher than the boiling pressure of the liquid and set Maximum value even a little bit higher. If the filter runs dry or if the filter is fractured, the pressure will increase and the control will be stopped automatically.

## Change operation mode

During operation, you can switch between the modes [Pump down] and [Vac control].


However, switching the mode during operation via the Mode button works only temporarily: after stopping, the controller switches back to its primary mode.

To switch to a mode permanently, the running process must be stopped and restarted after selection of the new mode!

- Press the [MODE] button to switch directly to the [Vac control] mode from the [Pump down] mode:
  - Semiautomatic distillation. Recommended for applications for which the process vacuum is still to be determined. Firstly, the vacuum pump is pumping down quickly in mode Pump down. As soon as the required process vacuum has been reached, e. g. boiling vacuum, this vacuum can be maintained by switching to Vac control. The actual pressure is adapted as the required set vacuum in the controller.
- Press the [MODE] button to switch directly to the [Pump down] mode from the [Vac control] mode:
  - With a connected VARIO® pump, a controller working in Auto mode will detect and track the boiling point automatically. The vacuum will be adapted continuously to the process. If a particular process vacuum is required, the mode can be switched back again to Vac control.

## Programs

Ten programs with vacuum and time presetting can be stored and edited in the device memory. Before selecting and activating a program, the running process has to be stopped!

Parameters	Meaning
Edit	Edit default values for process sequences of the actual program or edit an existing program.
Open	Load program out of one of the 10 memory slots.
Store:	Store the program under the selected number (memory capacity for up to 10 programs)


For further descriptions of the program functions, please contact Heidolph Instruments' technical service.

## Configuration

In the [Configuration] menu, different device parameters can be adjusted. Furthermore, the vacuum sensor can be readjusted and the device can be reset to the factory settings.


## Adjustment

This menu is used to adjust the vacuum sensor.

## RS-232

This menu is used to adjust different interface parameters.


Parameters	Selection	Meaning
Baud	▪ 19200	Default setting for transmission speed. The baud rate of data transfer of transmitter and receiver must correspond.
	▪ 9600	
	▪ 4800	
	▪ 2400	
Parity	▪ 8-N-1	Default setting for parity check, a method for error detection.
	▪ 7-O-1	
	▪ 7-E-1	
Handshake	▪ RTS-CTS	Preset for continuous data transmission without loss – data flow control.
	▪ Xon-Xoff	
	▪ None	
Remote	▪ Off	Control commands not enabled, only queries possible.
	▪ On	Connection for communication via RS 232 interface enabled.

When selecting [Remote On] the controller itself is only operable via an external device. All control elements of the controller are locked, except the [Standby] button.

VACUU-CONTROL® detects the [Remote On/Off] state automatically and retains that setting.

## Sensors


In the menu [Sensors], all connected sensors are listed. The internal sensor is generally displayed as [Sensor]. The entries of external sensors are listed as [Sensor type] with the corresponding sensor address.

The display switches automatically to the previous menu when selecting a sensor with the rotary switch.

## Display


In this menu, different device parameters can be adjusted and the menu language can be selected:

Selection	Meaning
Brightness	Adjust intensity of the display backlight brightness in a range of 0–100 %.
Contrast	Adjust contrast of the display in a range of 0–100 %.
Warning sound	Off – Switch off keystroke sound and warning sound. On – Switch on keystroke sound and warning sound.
Unit	Preset pressure unit for user interface (mbar/Torr/hPa).
Languages	Set user language for the user interface.

## Autostart

In this menu, the behavior of the device upon return of the power supply after a power interruption is defined:

Before activating the function ensure that no hazardous situations may occur due to the automatic restart of the process. Check in advance whether the Autostart feature can be used safely with the intended application.

### Operating behavior when selecting [Autostart on]:


In case of mains return after a power interruption, when the function [Autostart] is activated, the device restarts with the last set operating parameters without the need for intervention on the part of the operator.

### Operating behavior when selecting [Autostart Off]:

In case of mains return, when the function [Autostart] is deactivated, the device remains in [Stop] state and has to be restarted by the user with the [Start/Stop] button.

## Factory settings

Use this menu to reset the device to the factory settings. Please note that, doing so, all user data (programs) will be irretrievably deleted:

Selection	Meaning
Load	Load default settings.
Cancel	Leave menu item without default setting.


## Differential pressure measurement

For detailed descriptions about the Function menu with address assignment or about differential pressure measurement, please contact Heidolph Instruments' technical service.

## Graphic (curve)

In addition to the basic display with bar graphic, a real-time graphic with actual pressure/time curve can be accessed under [Graphic].

- During operation, press the rotary switch shortly.
- Mark the point [Graphic] and press the rotary switch to confirm the selection:

### Title bar

- 1 Mode, vacuum set point; for Vac control or Auto mode, actual vacuum value (actual pressure), elapsed process time

### Axis – pressure

- 2 Unit according to pre-settings (mbar, Torr, hPa)

### Axis – time

- 3 Continuous, automatically scaling time (minute, hour)

- Press the rotary switch again to exit the graphic (curve).

## Venting

The system can be vented at any time via the [VENT] button.

### DANGER

Depending on the process, venting the system can cause the formation of explosive mixtures.

Do not vent processes with ambient air when the formation of explosive mixtures can be caused.

If necessary vent with inert gas (max. 1.2 bar, absolute).

Certain processes may cause overpressure!

- Press the button one time to vent the system shortly.
  - Venting impulse, venting valve opens momentarily → short-term pressure increase.
- Keep pressed the button for more than two seconds to vent the system until reaching the atmospheric pressure (max. 1060 mbar).
  - Icon for venting valve is flashing, venting valve opens → continuous pressure increase until atmospheric pressure → venting valve closes, controller stops.
  - Continuous venting stops when pressing [VENT] button again.

## Troubleshooting

### Error display

The major symbol for fault indication is the warning triangle. Additionally displayed icons and sounds refer to the cause of fault:

Flashing: Warning!

Where applicable with

- flashing component icon
- warning sound (only when switched on)
- flashing backlight

In combination with number = Vacuubus address of the component which is defective

Display icons	Error message	Signal sound
↓ <b>1087.1</b>	Limit pressure reached	1x
	Process time elapsed	1x
	Venting valve	2x
	In-line suction valve	3x
	Coolant valve	4x
	External sensor removed or defective Internal sensor defective	5x 7x
	Vario pump	6x
	VACUULAN process pressure not reached within 99 hours.	8x

	Digital I/O module: <ul style="list-style-type: none"> <li>Fault indicator triggered</li> <li>Fault special configurations</li> </ul>	9x
	Level sensor triggered; flask full	10x
	Emission condenser Peltronic (too hot)	11x
	Analog I/O module:	12x

A defective I/O module, which is configured as a remote module, does not trigger a warning alert. Control stops. Alert display by the flashing warning triangle.

## Error messages

If a fault cannot be rectified with the described suggestions, please contact an authorized sales representative or our technical service (see section „Contact details“ on page 78).

Message	Possible cause	Remedy
Sensitive process not controllable	Motor speed too high Pumping speed too high	Reduce motor speed
Frequent error messages of connected components	Several controllers are connected. Several VACUU-BUS components of the same type are using the same address.	Use only one controller for one VACUU-BUS system. Only in Function/Vacuubus menu: Configure addresses.
VENT button does not work. Internal air admittance valve cannot be triggered.	Venting function deactivated. Setting of int.Air V in menu Function is switched Off or Auto. External air admittance valve is connected. and/or External vacuum sensor is connected.	Check why Venting is deactivated. Check if Venting by internal air admittance valve can be used without risk. Venting safe? Enable the function in menu Function. Adjustment: Auto or On.
Internal air admittance valve does not switch.	Air admittance valve soiled.	Clean venting valve (air admittance valve)
Function or menu item cannot be used	Function or menu item possibly only usable with short-cut (key combination).	Press the correct key combination; function of buttons and key combinations
Vario pump icon flashes	VARIO pump and VMS are both connected at the same time. VARIO pump defective. VMS defective or cable is not connected. Cable break of the connection cable.	Remove VMS from VARIO pump and restart controller. Check VARIO pump for defective parts. Check VMS for defective parts. Check cable connections. Replace defective parts.

Air admittance valve icon flashes	<p>External air admittance valve is removed.</p> <p>Plug disconnected</p> <p>External air admittance valve is defective.</p>	<p>Check the connection of the external air admittance valve.</p> <p>Check plug connection.</p> <p>Replace defective parts.</p> <p>Use internal air admittance valve.</p> <p>Reconfiguration without air admittance valve.</p>
In-line suction valve icon flashes	<p>In-line suction valve is removed.</p> <p>Plug disconnected.</p> <p>In-line suction valve is defective.</p>	<p>Check the connection of the in-line suction valve.</p> <p>Check plug connection.</p> <p>Replace defective parts.</p> <p>Reconfiguration without in-line suction valve.</p> <p>Switch-off the controller; On/Off button.</p> <p>Remove in-line suction valve and switch on controller again.</p>
Coolant valve icon flashes	<p>Coolant valve is removed.</p> <p>Coolant valve is defective.</p>	<p>Check the connection of the coolant valve.</p> <p>Replace defective parts.</p> <p>Reconfiguration without coolant valve.</p>
Level sensor icon flashes	<p>Level sensor is triggered.</p> <p>Level sensor is removed.</p> <p>Level sensor was triggered with empty flask.</p> <p>Cable break.</p> <p>Level sensor defective.</p>	<p>Empty flask/catch pot.</p> <p>Check position of level sensor.</p> <p>Adjust level sensor or delete the sensor from controller (by loading default).</p> <p>Check plug connection.</p> <p>Replace defective parts.</p>
Peltronic icon flashes	<p>Peltronic emission condenser too hot.</p> <p>Plug disconnected.</p>	<p>Let the Peltronic emission condenser cool down.</p> <p>Check plug connection.</p>
Title bar without text	<p>No controllable device connected (In-line suction valve, VMS, VARIO pump).</p>	<p>Check device connections and cable.</p> <p>Replace defective parts.</p> <p>Connect a controllable device to the controller.</p> <p>Use the controller as measuring gauge.</p>
No key reaction – only On/Off, PC icon displayed	<p>Remote switched On.</p> <p>Controller only controllable via external device.</p>	<p>Switch off Remote (switch off, when switching on, press rotary switch shortly, in the Configuration menu, select RS232 and switch Remote to Off).</p> <p>Control controller via end device.</p>
No reaction to key actuation	<p>Controller defective</p>	<p>Contact service and return device for repair.</p>
No display	<p>Controller switched off/defective.</p> <p>Power supply disconnected.</p> <p>Power supply not correctly connected.</p> <p>Mains voltage failure.</p> <p>Cable break.</p>	<p>Switch on controller; On/Off button.</p> <p>Check plug connection and wall power supply for correct connection.</p> <p>Replace defective parts.</p> <p>Contact service and return device for repair.</p>

Blank display	<p>Too many devices connected, e. g., valves.</p> <p>Short circuit of a connected device.</p> <p>Short circuit at RS232 interface.</p> <p>Controller defective.</p>	<p>Power input of all connected devices may not exceed the maximum power consumption: controller with wall power supply max. 30 W, controller + VARIO max. 25 W.</p> <p>Replace defective parts.</p> <p>Check RS232 plug connection.</p> <p>Contact <a href="#">service</a> and return device for repair.</p>
Incorrect pressure display	<p>Humidity inside the vacuum sensor.</p> <p>Vacuum sensor soiled.</p> <p>Vacuum sensor not adjusted.</p> <p>Vacuum sensor not correctly adjusted.</p>	<p>Identify and remove source of humidity.</p> <p>Dry the vacuum sensor, e. g., by pumping down.</p> <p>Clean the vacuum sensor.</p> <p>(Re)adjust the vacuum sensor.</p>
Digital pressure gauge flashes	<p>Pressure display flashing with 0.0: Vacuum adjustment not correctly carried out.</p> <p>Pressure display flashing with 1060: Overpressure! Pressure &gt; 1060 mbar.</p>	<p>(Re)adjust internal or external vacuum sensor.</p> <p>Risk of bursting: Discharge the system immediately by venting.</p>
No digital pressure display.	<p>External vacuum sensor is defective.</p> <p>External vacuum sensor is removed.</p> <p>Internal vacuum sensor defective.</p>	<p>Replace defective parts.</p> <p>Reconnect external vacuum sensor.</p> <p>Contact <a href="#">service and return device for repair</a>.</p>
Sensors submenu is permanently displayed	<p>No sensor selected in the submenu Sensors.</p>	<p>Select the required sensor by turning the rotary switch.</p>
After loading defaults Language selection appears	<p>Special factory settings have been loaded.</p>	<p>Set language and pressure.</p> <p><b>IMPORTANT:</b> Check if the loaded default settings are suitable for your vacuum apparatus.</p>
Error I/O module	<p>Plug disconnected.</p> <p>An error occurred in the system, the I/O module passed the error alert to the controller.</p>	<p>Check plug connection.</p> <p>Remedy external fault.</p>
VSP sensor displays wrong values	<p>VSP sensor configured as VSK.</p>	<p>Use menu Function/Vacuubus to reconfigure the sensor correctly as VSP.</p>
Controller in operation, pressure display flashes	<p>VSK sensors are measuring negative difference pressure.</p>	<p>Select another vacuum sensor in menu Sensors.</p>
Process time elapsed	<p>All program steps of a program are completed.</p> <p>Program end reached.</p>	<p>Acknowledge (confirm) indication by pressing the Start/Stop button.</p>
Clock icon flashes	<p>Process time elapsed.</p>	<p>Acknowledge end of process by pressing the Start/Stop button.</p>
Pump down stops, flashing arrow down icon	<p>Pressure below preset minimum value.</p>	<p>Acknowledge message by pressing the Start/Stop button.</p> <p>If possible readjust default value (min.).</p>
Vac control stops, flashing arrow up icon	<p>Preset maximum value exceeded.</p>	<p>Acknowledge message by pressing the Start/Stop button.</p> <p>If possible readjust default value (max.).</p>
Program -	<p>Program not yet stored.</p>	<p>Store program under a free program number.</p>

### Error of external components

Error messages for defective external components such as in-line suction valve, vacuum sensor, etc. cannot be reset. Replace defective accessories or send defective accessories for repair to your authorized dealer or our service.

### Reset

#### Auto reset

The following error indications will be reset automatically with remedy:

- Overpressure
- Process time elapsed
- Limit pressure reached
- Error air admittance valve
- Error Peltronic

#### Active reset

Several error indications need to be reset with remedy. Depending on the fault severity different actions are required.

→ Press the Start/Stop button to reset the following error indications:

- In-line suction valve error
- Coolant valve error
- External vacuum sensor removed
- I/O module activated Error indication
- Level sensor triggered

→ Load the factory settings to reset the following

error indications:

- I/O module Remote module, missing set value presetting
- Level sensor removed and disconnected

Save stored programs before loading the default settings.

## Technical specifications

<b>General device data</b>	
Model	<b>Vacuum control Hei-VAC Control</b>
Device dimensions (W × H × D)	123 × 124 × 83 mm (5 × 5 × 3.5 in.)
Dimensions of version with foot (W × H × D)	144 × 124 × 115 mm (6 × 5 × 4.5 in.)
Device weight	440 g (0.97 lb)
Weight (with foot)	570 g (1.3 lb)
Display	LC display with brightness control
Pressure sensor	Integrated (ceramic diaphragm/aluminium oxide, capacitive, independent from gas type, absolute pressure)
Protection class (front side)	IP42 (IP42)
<b>Electrical data</b>	
Rated voltage	24 VDC (±10 %)
Power (max.)	3.4 W
Admissible electricity of the connected valves	4 A
Interfaces	RS232, SUB-D, 9-pin
<b>Permissible ambient conditions</b>	
Operating temperature	10 °C – 40 °C (50 °F – 104 °F)
Storage/transport temperature	-10 °C – 60 °C (14 °F – 140 °F)
Rel. Humidity	30 – 85 %, without condensation
Altitude, max.	3000 m above sea level
<b>Vacuum data</b>	
Measuring range, absolute	1080 – 0.1 mbar (810 – 0.1 Torr.)
Control range	1060 – 0.1 mbar (795 – 0.1 Torr.)
Resolution	0.1 mbar (0.1 Torr.)
<b>Vacuum sensor VSK 3000</b>	
Max. admissible pressure, absolute	1.5 bar (1125 Torr.)
Max. admissible media temperature (gas)	Temporarily: 80 °C (176 °F) Continuous operation: 40 °C (104 °F)
Measurement uncertainty	< ±1 mbar (< ±0.75 Torr.)
Temperature coefficient	< ±0.07 mbar / K (< ±0.05 Torr. / K)
<b>Venting</b>	
Max. admissible pressure, absolute	1.2 bar (900 Torr.)
Connections	Hose nozzle for tube, inner diameter $\varnothing = 4 - 5$ mm

## Wetted materials


Component	Material
Vacuum connection, hose nozzle	PP
Sensor	Aluminum oxide ceramic
Sensor housing	PPS/Glass fiber
Sensor seal	Chemically resistant fluoroelastomer
Venting valve seal	FFKM

## Interface commands

If necessary, please contact Heidolph Instruments' technical service.

## Pin assignment (RS232)

Pin	Designation	Function
1	DCD	-
2	RxD	Received data
3	TxD	Transmission data
4	DTR	+10 V
5	GND	Mass
6	DSR	-
7	RTS	Transmission request
8	CTS	Ready to send
9	RI	+5 V (bluetooth, remote control)

## Scope of delivery

Item	Quantity	Product no.
Controller Hei-Vac control	1	11-001-591-28
Wall power supply plug 30 W, 24 V; with all adapters	1	11-300-008-12
Vacuum valve, pre-assembled	1	569-00080-00
Foot CVC 3000	1	11-300-008-11
Operating instructions	1	01-005-005-87
Warranty registration	1	01-006-002-78

## Accessories

Detailed information on the available accessories for your device variant can be found on our website at [www.heidolph-instruments.com](http://www.heidolph-instruments.com).

If required, contact an authorized dealer or our technical service, see section „Contact details“ on page 78.

## Device service

### DANGER

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

When cleaning, avoid the penetration of liquids.

## General cleaning instructions

Wipe all surfaces of the device with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

### CAUTION

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

## Clean venting valve (air admittance valve)

- Apply slight overpressure of dry air or inert gas to the vacuum port.
- Press the [VENT] button several times until gas escapes through the venting port.
- Repeat this procedure until you hear the clicking of the valve and a gas stream is noticeable at the venting port.

## Clean controller pressure sensor

- Fill a small amount of solvent via the vacuum port in the controller, e. g., white gas.
- Let the solvent react for a few minutes.
- Drain the solvent.
  - Dissolved substances or discolorations in the solvent are possible.
- Repeat this procedure until no more pollutants are in the solvent.
- Leave the controller until its interior has dried.
- Readjust the pressure sensor.

## Sensor readjustment

For readjustment the reference pressures need to be known with certainty. In the pressure range 20 – 700 mbar (15 – 525 Torr) no adjustment is possible.

Check the accuracy of the pressure sensor in case of irregularities in the pressure display. Readjust the sensor in two steps: at atmospheric pressure and under vacuum.

Do not adjust at atmospheric pressure, if the pressure at the location of the device is not exactly known (pay attention to height above sea level)!

Any kind of pollution of the vacuum system, e. g., oil, substances, or humidity could falsify the adjustment. Clean polluted sensors before readjustment!

## Adjustment at atmospheric pressure

An adjustment at atmospheric pressure is only possible if the pressure is higher than  $> 700$  mbar.

- Vent the controller or the external pressure sensor VSK 3000.
- Make sure that the pressure sensor (internal or external) is really at atmospheric pressure.
- Determine the exact atmospheric pressure of your location, e. g. by barometer, inquiry at the meteorological office or the airport.
- Call up the menu Configuration.
- Turn the rotary switch and place the bar marking on Adjustment.
- Press the rotary switch.
  - Marking jumps to numeric value.
- Adjust the exactly determined local atmospheric pressure by turning the rotary switch.
- Press the rotary switch.
  - Pressure sensor adjusted to atmospheric pressure.

## Adjustment under vacuum

An adjustment at atmospheric pressure is only possible if the pressure is lower than  $< 20$  mbar.

- Evacuate the controller or the external pressure sensor VSK 3000 to a pressure  $< 0,1$  mbar.

Adjustment under vacuum with an actual pressure higher than  $< 0,1$  mbar reduces the accuracy of the measurement. If the pressure during adjustment is significantly higher than  $> 0,1$  mbar, the device is not optimally adjusted and has to be adjusted to a reference pressure.

- Call up the menu Configuration.
- Turn the rotary switch and place the bar marking on Adjustment.
- Press the rotary switch.
- Marking jumps to numeric value.
- Adjust the pressure value to 0 by turning the rotary switch.
- Press the rotary switch.
  - Pressure sensor adjusted to vacuum.

The readjustment of a VSP 3000 can only be carried out in warmed-up state. Adjustment is not possible during the warm-up time.

Use a high vacuum pump for the adjustment of a VSP sensor.

After connection to power supply and after the pressure has reached  $< 10-3$  mbar, wait 20 minutes before adjusting the VSP sensor.

Carry out the adjustment in the same order as described above for VSK.

## Adjustment at a reference pressure

Instead of adjustment under vacuum to a pressure  $< 0,1$  mbar, the adjustment to a precisely known reference pressure within the range of 0 – 20 mbar is possible.

- Evacuate the controller or the external pressure sensor VSK 3000 to a pressure between 0 – 20 mbar.
- Call up the menu Configuration.
- Turn the rotary switch and place the bar marking on Adjustment.
- Press the rotary switch.
- Adjust the pressure value to the actual reference pressure by turning the rotary switch.
- Press the rotary switch.
  - Pressure sensor adjusted to reference vacuum.

The measurement uncertainty of the reference pressure will directly affect the measurement uncertainty of the controller.

If the ultimate vacuum of a diaphragm pump is adjusted and if the pressure is not determined with a precise vacuum gauge, a measurement error may eventually occur; especially if the diaphragm pump does not reach the ultimate vacuum, e.g., due to condensate, failure or soiling of the valves or leaks.

## Repairs

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim.

The owner is generally liable for damage caused by unauthorized repairs.

If repairs are required, contact an authorized dealer or our technical service, see section „Contact details“ on page 78.

Include the completed certificate of decontamination with every device return, see „Certificate of decontamination“ on page 79.

## Maintenance

There are no user-serviceable components in the unit housing. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example), contact our technical service, see section „Contact details“ on page 78.

## Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
- Properly remove and dispose of residues of substances that are hazardous to health, the environment and biohazardous!

## Warranty statement

Heidolph Instruments GmbH & Co. KG provides a three-year warranty for material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage [www.heidolph.com](http://www.heidolph.com).

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

## Contact details

### Heidolph Instruments GmbH & Co. KG

Technical Service  
Walpersdorfer Str. 12  
D-91126 Schwabach/Germany

Email: [service@heidolph.de](mailto:service@heidolph.de)

### Representations

To find your local Heidolph distributor please visit [www.heidolph.com](http://www.heidolph.com)

## Certificate of decontamination

Enclose the certificate of decontamination, duly completed, with your device return.  
Submissions without a certificate of decontamination cannot be processed!

# CERTIFICATE OF DECONTAMINATION

IN CASE OF RETURNS

Please fill in the required fields.

**Note: The sender must package the goods properly and appropriately for transport.**

Heidolph Instruments GmbH & Co. KG  
Walpersdorfer Straße 12  
91126 Schwabach, Germany  
Phone: +49 (0) 9122 9920-380  
**Fax: +49 (0) 9122 9920-19**  
Email: service@heidolph.de

### SENDER

Name \_\_\_\_\_ First name \_\_\_\_\_  
Company/institution \_\_\_\_\_ Department \_\_\_\_\_  
\_\_\_\_\_ Workgroup \_\_\_\_\_  
Address \_\_\_\_\_  
ZC/City \_\_\_\_\_  
Country \_\_\_\_\_ Phone \_\_\_\_\_  
Email \_\_\_\_\_

### DEVICE DETAILS

Article number \_\_\_\_\_ Serial no. \_\_\_\_\_  
Ticket number \_\_\_\_\_  
Reason for sending in \_\_\_\_\_  
\_\_\_\_\_

**Has the device been cleaned, decontaminated/disinfected?** Yes No (Please mark as applicable)  
If yes, which measures were carried out?  
\_\_\_\_\_

**Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous?** Yes No (Please mark as applicable)  
If yes, with which substances did the device come into contact?  
\_\_\_\_\_

### LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product (from serial number on 200067811) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

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## Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,  
Heidolph Instruments GmbH & Co. KG  
Walpersdorfer Straße 12  
91126 Schwabach / Germany

Product: Laboratory evaporator  
Model:

Hei-VAP Industrial xxx

518-xxxxx-xx-x

**Description:**

Rotary Evaporators with heating bath and human machine interface

The object of the declaration described above is in conformity with the relevant  
UK Statutory Instruments (and their amendments):

*2008 No. 1597*

*The Supply of Machinery (Safety) Regulations*

*2008 2016 No. 1091*

*The Electromagnetic Compatibility Regulations 2016*

*2012 No. 3032*

*The Restriction of the Use of Certain Hazardous Substances in  
Electrical and Electronic Equipment Regulations 2012*

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):

ProductIP ( UK ) Ltd.  
8, Northumberland Av.  
London WC2N 5BY

Signed for and on behalf of Heidolph Instruments GmbH & Co. KG  
Walpersdorfer Straße 12, 91126 Schwabach / Germany

Schwabach, 01.09.2021

Wolfgang Jaenicke  
Managing Director

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices by Heidolph Instruments GmbH + Co. KG fulfill the requirements of the EU directive 2011 / 65 / EU (RoHS) and the amended of directive 2015 / 863.

All devices are compatible with the requirement of the RoHS for the following materials:

Max. 0,01% of the weight in homogeneous material for cadmium and max. 0,1 % of the weight in homogeneous material for lead, mercury, hexavalent chromium, polybrominated biphenyl, polybrominated diphenyl ether, Di (2-ethylhexyl) phthalate, butyl benzyl phthalate, dibutyl phthalate, diisobutyl phthalate.

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 22.07.2019

Stefan Peters  
Vice President Marketing, Innovation & Technology

Marcell Sarré  
Vice President Quality Management & Technical Service

## China RoHS DECLARATION OF CONFORMITY

Heidolph Instruments GmbH & Co.KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

### Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Instruments GmbH & Co.KG products is 25 years.

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This table is created according to SJ/T 11364-2014

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**Note:** Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
- O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
- X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Instruments GmbH & Co.KG may enter into further devices or can be used together with other appliances.

With these products and appliances in particular, Heidolph Instruments GmbH & Co.KG will not take responsibility for the EFUP of those products and appliances.

Schwabach, 06.08.2021

Wolfgang Jaenicke  
Chief Executive Officer CEO

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Instruments GmbH & Co. KG PRODUCTS
© Heidolph Instruments GmbH & Co. KG