
Hei-PROCESS

Process Solutions

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Quick Start Guide
Guide d'initiation rapide

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Seite 4 – 12

Quick start guide

Page 16 – 24

Guide d'initiation rapide

Page 28 – 36

Inhalt

Einleitung

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

The screenshot shows the login interface for the OSPIN Modular Bioprocessing system. At the top, there is a red circular logo with the text 'OSPIN' and 'Modular Bioprocessing' to its right. Below the logo, there are two input fields: 'username or email' and 'password'. The 'username or email' field contains a person icon, and the 'password' field contains a lock icon. Below these fields is a red button labeled 'Login'. At the bottom, there are two links: 'No account yet? Sign up!' and 'Forgot password'.

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The screenshot shows the sign-up interface for the OSPIN Modular Bioprocessing system. At the top, there is a red circular logo with the text 'OSPIN' and 'Modular Bioprocessing' to its right. Below the logo, there are five input fields: 'username', 'email', 'password', 'confirm password', and 'phone number (optional)'. The 'username' field contains a person icon, the 'email' field contains an envelope icon, and the 'password' and 'confirm password' fields contain lock icons. Below these fields is a red button labeled 'Sign Up'. At the bottom, there is a link: 'Back to Login'.

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The screenshot shows the OSPIN user interface. On the left is a sidebar with navigation items: Devices, Profile, Alerts, and Messages. The main content area is divided into three sections: 'User Information' (Username, Name, Organization, Email, Phone, Change Password), 'Licences' (No licences found, Claim Licence button), and 'Delete Account' (Warning message, Delete Account button). A 'Claim Licence' dialog box is open on the right, with a 'Licence Key' input field and 'Close' and 'Confirm' buttons.

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The screenshot shows the 'Request Password Reset' form. At the top is the OSPIN logo and the text 'Modular Bioprocessing'. Below is a text input field labeled 'username or email' with a question mark icon. A red button labeled 'Request Password Reset' is positioned below the input field. Underneath the button are the links 'Back to Login' and 'I already received a code'.

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The screenshot shows a blue circular icon with a white exclamation mark, indicating an error or warning message.

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The screenshot displays the OSPIN Showroom Gateway interface. On the left, there is a navigation menu with options: Devices, Profile, Alerts, and Messages. The main content area features a search bar and a list of individual devices and setups.

Individual Devices:

- Hei-Chill:** ready for process
- Hei-Connect:** Mix & Heat #1 started by flooper on 4/3 Bar 11:27, 0h, 14m
- Hei-FLOW Ultimate 120:** ready for process
- Hei-TORQUE Ultimate 200:** Stirring #5 started by flooper on 4/3 Bar 11:56, 0h, 15m
- Hei-VAP:** ready for process

Setups:

- Automated Dosing Controller:** ready for process
- Basic Chemical Reactor:** blocked by Hei-TORQUE Ultimate 200
- Hei-Connect + Hei-FLOW Ult...:** blocked by Hei-Connect

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The screenshot displays the Felix' Internal Gateway interface. On the left, there is a navigation menu with options: Devices, Profile, Alerts, and Messages. The main content area features a 'My Processes' section with a 'Create Process' button and a grid of process cards.

My Processes:

- Clone_Test_3:** Finished Feb 02 2023 14:05, Duration: 2m 0s, finished
- Process template 1:** Created Jan 05 2023 09:33, executable
- Longterm_Test:** Finished Nov 02 2022 08:05, Duration: 15h 0m 0s, finished
- pH Control & Mixing Bolus test:** Finished Jun 16 2022 16:18, Duration: 19m 14s, finished
- Hei-TORQUE Ultimate 200 Clone_Testing:** Finished Mar 01 2023 11:35, Duration: 2m 0s, finished
- Hei-TORQUE Ultimate 200 Testing:** Finished Mar 01 2023 11:32, Duration: 2m 0s, finished
- Basic Chemical Reactor Clone_Clone_MyProcess...:** Finished Feb 20 2023 13:02, Duration: 26m 52s, finished
- Basic Chemical Reactor Clone_MyProcess-Feb 27...:** Duration: 24s, finished
- Basic Chemical Reactor MyProcess-Feb 27 2023 1...:** Finished Feb 27 2023 16:31, Duration: 18s, finished
- Hei-Connect + Hei-FLOW Ultimate 600 Test_1:** Created Feb 27 2023 14:18, executable
- pH Control & Mixing Clone_pHTest19_300ml dil...:** Created Feb 24 2023 12:17, executable
- Basic Chemical Reactor MyProcess-Feb 23 2023 1...:** Created Feb 23 2023 15:10, executable
- Hei-TORQUE Ultimate 200 MyProcess-Feb 17 2023 1...:** Finished Feb 17 2023 12:17, Duration: 2m 19s, finished
- Hei-TORQUE Ultimate 200 MyProcess-Feb 17 2023 1...:** Finished Feb 17 2023 12:10, Duration: 4m 0s, finished
- Hei-TORQUE Ultimate 200 MyProcess-Feb 17 2023 1...:** Finished Feb 17 2023 12:02, Duration: 5m 0s, finished
- Hei-TORQUE Ultimate 200 MyProcess-Feb 17 2023 1...:** Finished Feb 17 2023 11:52, Duration: 1m 34s, finished

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Contents

Introduction

About this document	16
Typographic conventions	16
California Residents	16
Copyright protection	16

General notes

Notes on the product	17
Intended use	17
Compliant use	17
Safety information	17
Other regulations	17

Basics

General information	18
Logging into the system	18
OSPIN App	19

Basic functions

Home page	21
Sidebar	21
Devices view	21

Appendix

Contact details for Heidolph Instruments GmbH	24
Contact details for OSPIN GmbH	24

About this document

This quick start guide describes the requirements and procedure for using the Hei-PROCESS application via an OSPIN gateway and the OSPIN software.

Typographic conventions

Standardized symbols, signal words and highlighting elements are used in this document to warn of hazards and to identify important information or rather special text contents.

Symbol	Signal word / explanation
	Mandatory signs are used to indicate important and useful information on handling a product. This information is used to ensure operational safety and to maintain the value of the product.
[GUI]	User interface labels are highlighted typographically in the text and in tables to make them easier to assign.
→	The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.

California Residents

Important information for California residents regarding Prop 65. Please visit www.P65Warnings.ca.gov for more information.

Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Instruments GmbH & Co. KG. Any violation is subject to compensation for damage.

Notes on the product

The product described has been designed to the state-of-the-art standards and recognized safety rules current at the time of development. The product described does not create any direct hazards.

Intended use

The product described has been designed solely for the implementation of process solutions for laboratory equipment. Any other use is deemed to be not as intended!

The product described may only be used by authorized and competent personnel. Training and qualification of the operating personnel as well as ensuring they act responsibly are the sole responsibility of the operator!

Compliant use

The user is generally responsible for evaluating the conformity of his application and, if necessary, for taking additional measures.

Safety information

Before using the product described, familiarize yourself with all safety regulations and occupational guidelines applicable in the place of use, and all safety regulations and guidelines for the activated devices, and follow these at all times (operating instructions, guidelines, directives, legal provisions, etc.).

If there is missing or misleading information on the product or on occupational safety, contact the responsible safety specialist or our technical service.

Contact our technical service in case of non-remediable and/or recurring faults (for contact details see section „Contact details for Heidolph Instruments GmbH“ on page 24).

The owner/operating company is solely responsible for all damage resulting from unauthorized changes to the products described or disregard of the safety instructions, hazard warnings or manufacturer's instructions!

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

General information

Hei-PROCESS is a user-friendly, cloud-based application that connects laboratory equipment, enabling control and optimization of laboratory processes.

The complexity and functional scope of the application depends on the individually selected solution (Hei-PROCESS Core or Hei-PROCESS Expert).

The laboratory equipment is connected via a so-called OSPIN gateway.

Processes are created, monitored and controlled through any device with internet browser and internet connection.

This document describes the procedure for logging into the system and the basic system functions.

All software functions – where necessary – are described in detail in the extensive online help of OSPIN in the browser window. In addition, a support function is available to the user, not only for the OSPIN but also for the Heidolph system components.

Detailed information on the hardware link is given in the operating instructions for the OSPIN gateway and the operating instructions for the Heidolph devices.

Logging into the system

After the hardware installation, use the login data provided by OSPIN to log into the system. These data are sent to the user after acquiring a license from Heidolph Instruments:


- Make sure that all the required and supported devices are connected to the assigned OSPIN gateway via the respective RS232 port.
- Make sure that all OSPIN gateways involved are connected to the internet and all connected laboratory equipment is switched on.

OSPIN App

The control of the OSPIN gateway and the connected devices or setups is always realized via the OSPIN app. Start the internet browser on your device and open the URL <https://ospin-app.com>.

Registered user

- As registered user, enter your assigned user name or the email address used to sign up, and your password and click [**Login**] to log into the system:

OSPIN Modular Bioprocessing

username or email ⓘ


password

Login

No account yet? Sign up! Forgot password

New user

- Click [**No account yet? Sign up!**] to create a new user account. You can create an unlimited number of accounts without any extra costs!
- In the following screen, store all required information:

OSPIN Modular Bioprocessing

username

email

password

confirm password

phone number (optional)


Sign Up

Back to Login

Claim license key

To associate the OSPIN Gateway with your user account, you need to claim the provided license key: This will result in making you the device owner and allow you to manage the access for other users to it.


→ License claiming is possible when navigating to the **[Profile]** section via the side bar:

- Click on **[Claim Licence]** in the section **[Licences]**.
- Enter your license key in the dialog window **[Claim Licence]** and confirm by clicking on **[Confirm]**.
- If you use a Hei-PROCESS Expert setup you additionally need to register the device ID which can be found on the label of the OSPIN Gateway.

Reset password

→ Click **[Forgot password]** to reset your password if necessary:

→ Follow the instructions on the screen.

In case of an incorrect input you will be prompted to enter your user information again. If, despite correct data input, you are unable to log into the system, please contact your system administrator.

Home page

After logging in successfully, the application Home page appears in the browser.

Sidebar

The Sidebar on the left of the screen provides the following options:

	
Devices 	→ Devices: Open [Devices view], see figure of Main window above.
Profile 	→ Profile: open [Profile data] view. All user-specific data are displayed in this view. Also in this same view, you can claim your license.
Alerts 	→ Alerts: Open [Alerts] view. All alert messages received are listed in this view.
Messages 	→ Messages: Open [Messages] view. All messages received are listed in this view.
Help 	→ Help: Open the online help. The online help provides extensive help on use of the OSPIN software. The online help is updated continuously, so that the most up-to-date information is always available to the user. Use the following link to open the online help: https://ospin.refined.site/
Changelog 	→ Changelog: Open the change log. Use the change log to trace all updates to the OSPIN software.
Logout 	→ Logout: Log out the user

Devices view

The devices view [**Devices**] shows the registered devices with their respective availability status as so-called device widgets. In each device widget you can find the COM port through which the device – or multiple devices in case of a setup – is connected to the gateway.

Functions

The screenshot displays the OSPIN Showroom Gateway interface. At the top, there is a search bar and a title 'OSPIN Showroom Gateway' with the location 'Heidelberg Instruments'. Below this, there are filter options: 'Show disconnected devices & setups' and 'Aggregate devices & setups by type'. The main content is divided into two sections: 'Individual Devices' and 'Setups'.

Individual Devices:

- Hei-Chill:** ready for process
- Hei-Connect:** Mix & Heat #1 started by flooper on 4/18/23 11:27, 0h, 14m
- Hei-FLOW Ultimate 120:** ready for process
- Hei-TORQUE Ultimate 200:** Stirring #5 started by flooper on 4/18/23 11:56, 0h, 15m
- Hei-VAP:** ready for process

Setups:

- Automated Dosing Controller:** ready for process
- Basic Chemical Reactor:** blocked by Hei-TORQUE Ultimate 200
- Hei-Connect + Hei-FLOW Ult...:** blocked by Hei-Connect

On the left side, there is a navigation menu with 'Devices', 'Profile', 'Alerts', 'Messages', 'Help', 'Changelog', and 'Logout'.

→ If you click on a device widget an overview of all processes associated with this device or setup will be shown.

The screenshot displays the Felix' Internal Gateway interface. At the top, there is a title 'Felix' Internal Gateway' and navigation tabs for 'My Processes', 'Logs', 'Setup', and 'Access'. Below this, there is a '+ Create Process' button and a pagination indicator 'showing 1 - 40 out of 166'.

The main content is a grid of process cards, each showing a device icon, a process name, a status (finished or executable), and a duration. The cards are organized by device type:

- Hei-VAP:** Clone_Test_3 (finished), Longterm_Test (finished), Bolus test (finished)
- Hei-TORQUE Ultimate 200:** Clone_Testing (finished), Testing (finished)
- Basic Chemical Reactor:** Clone_Clone_MyProcess... (finished), Clone_MyProcess-Feb 27... (finished)
- Hei-Connect + Hei-FLOW Ultimate 600:** Test_1 (executable)
- pH Control & Mixing:** Clone_pHTest19_300ml dil... (executable)
- Basic Chemical Reactor:** MyProcess-Feb 23 2023 1... (executable)
- Hei-TORQUE Ultimate 200:** MyProcess-Feb 27 2023 1... (finished), MyProcess-Feb 17 2023 1... (finished), MyProcess-Feb 17 2023 1... (finished), MyProcess-Feb 17 2023 1... (finished)

On the left side, there is a navigation menu with 'Devices', 'Profile', 'Alerts', 'Messages', 'Help', 'Changelog', and 'Logout'.

The following options are available:

- **My Processes**: Open list of the defined processes
- **Log**: Open list of the recorded events
- **Setup**: Open overview of the defined process and device parameters
- **Access**: Open overview of the access rights
- Button [**Create Process**]: create new device process, follow the instructions on the screen. Each new process is displayed as an additional function button. Click on these function buttons to open the screen for editing the process.

Contact details for Heidolph Instruments GmbH

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

The login form features the OSPIN logo and the text 'Modular Bioprocessing' at the top. It contains two input fields: 'username or email' with a user icon and a help icon, and 'password' with a lock icon. A red 'Login' button is positioned below the fields. At the bottom, there are two links: 'No account yet? Sign up!' and 'Forgot password'.

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The sign-up form features the OSPIN logo and the text 'Modular Bioprocessing' at the top. It contains five input fields: 'username' with a user icon, 'email' with an envelope icon, 'password' with a lock icon, 'confirm password' with a lock icon, and 'phone number (optional)' with a phone icon. A red 'Sign Up' button is positioned below the fields. At the bottom, there is a link: 'Back to Login'.

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The screenshot shows the 'Request Password Reset' form. At the top is the OSPIN logo and 'Modular Bioprocessing'. Below is a text input field labeled 'username or email' with a help icon. A 'Request Password Reset' button is below the input field. Underneath are links for 'Back to Login' and 'I already received a code'.

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Fonctions

The screenshot displays the OSPIN Showroom Gateway interface. At the top, there is a search bar and a title 'OSPIN Showroom Gateway' with the subtitle 'Heidolph Gateway' and 'Heidolph Instruments'. Below this, there are two radio buttons: 'Show disconnected devices & setups' (selected) and 'Aggregate devices & setups by type'. The main content is divided into two sections: 'Individual Devices' and 'Setups'.

Individual Devices:

- Hei-Chill:** ready for process
- Hei-Connect:** Mix & Heat #1 started by flowser on 4/18/23 11:27, 0h,14m
- Hei-FLOW Ultimate 120:** ready for process
- Hei-TORQUE Ultimate 200:** Stirring #5 started by flowser on 4/18/23 11:56, 0h,15m
- Hei-VAP:** ready for process

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The screenshot displays the Felix' Internal Gateway interface. At the top, there is a title 'Felix' Internal Gateway' and a navigation menu with 'My Processes', 'Logs', 'Setup', and 'Access'. Below this, there is a '+ Create Process' button and a pagination bar showing '1 - 40 out of 166'.

The main content is a grid of process cards, each representing a different process or setup. Each card includes a device icon, a title, a status (finished or executable), a duration, and a 'Created' date. The cards are arranged in a grid with columns and rows.

Process Cards:

- Hei-VAP Clone_Test_3:** Finished Feb 02 2023 14:05, Duration: 2m 0s, finished
- Hei-VAP Process template 1:** Created Jan 05 2023 09:33, executable
- Hei-VAP Longterm_Test:** Finished Nov 02 2022 08:05, Duration: 15h 0m 0s, finished
- pH Control & Mixing Bolus test:** Finished Jun 16 2022 16:18, Duration: 19m 14s, finished
- Hei-TORQUE Ultimate 200 Clone_Testing:** Finished Mar 01 2023 11:35, Duration: 2m 0s, finished
- Hei-TORQUE Ultimate 200 Testing:** Finished Mar 01 2023 11:32, Duration: 2m 0s, finished
- Basic Chemical Reactor Clone_Clone_MyProcess-...:** Finished Feb 20 2023 13:02, Duration: 26m 52s, finished
- Basic Chemical Reactor Clone_MyProcess-Feb 27...:** Duration: 24s, finished
- Basic Chemical Reactor MyProcess-Feb 27 2023 1...:** Finished Feb 27 2023 16:31, Duration: 18s, finished
- Hei-Connect + Hei-FLOW Ultimate 600 Test_1:** Created Feb 27 2023 14:18, executable
- pH Control & Mixing Clone_pHTest19_300ml dil...:** Created Feb 24 2023 12:17, executable
- Basic Chemical Reactor MyProcess-Feb 23 2023 1...:** Created Feb 23 2023 15:10, executable
- Hei-TORQUE Ultimate 200 MyProcess-Feb 17 2023 1...:** Finished Feb 17 2023 12:17, Duration: 2m 19s, finished
- Hei-TORQUE Ultimate 200 MyProcess-Feb 17 2023 1...:** Finished Feb 17 2023 12:10, Duration: 4m 0s, finished
- Hei-TORQUE Ultimate 200 MyProcess-Feb 17 2023 1...:** Finished Feb 17 2023 12:02, Duration: 5m 0s, finished
- Hei-TORQUE Ultimate 200 MyProcess-Feb 17 2023 1...:** Finished Feb 17 2023 11:52, Duration: 1m 34s, finished

On the left side, there is a navigation menu with 'Devices', 'Profile', 'Alerts', 'Messages', 'Help', 'Changelog', and 'Logout'.

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Coordonnées Heidolph Instruments GmbH
Coordonnées OSPIN GmbH
© Heidolph Instruments GmbH & Co. KG