# Hei-CHILL

Chiller

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Operating instructions

Hei-CHILL 250  
Hei-CHILL 350  
Hei-CHILL 600  
Hei-CHILL 1200

Translation of the original instructions  
Page 4 – 45

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<b>About this Document</b>	<b>4</b>
Version and variants .....	4
About this manual .....	4
Symbols and keywords.....	5
<b>Safety Instructions</b>	<b>6</b>
General safety instructions.....	6
EU Declaration of Conformity.....	6
Intended use.....	6
Installation .....	7
Qualifications of personnel .....	7
Operating company's obligations.....	8
Installation site .....	8
Modifications to the device.....	8
Safety of personnel.....	8
Safety during use .....	9
Disposal .....	9
<b>Device Description</b>	<b>10</b>
Device overview .....	10
Device overview Hei-CHILL 250 / 350 / 600 / 1200 .....	10
Hei-CHILL control panel.....	11
<b>Start Up</b>	<b>12</b>
Setting up the device .....	12
Connecting/disconnecting the power cord .....	13
Switching the device on and off .....	13
Filling with thermal fluid .....	14
Filling for the first time.....	15
Menu navigation.....	17
Entering the temperature setting .....	18
Setting the temperature limits.....	18
Setting the temperature setpoint .....	20
<b>Operation</b>	<b>21</b>
Overview of menu structure.....	21
Temperature control operation.....	22
Interfaces.....	23
RS 232 interface.....	23
Interface commands .....	24
Cable and interface test for RS232 .....	25
Alarm output.....	26
Autostart.....	27
Restoring the factory settings.....	28

<b>Cleaning and Maintenance</b>	<b>29</b>
Cleaning.....	29
Cleaning the device.....	29
Cleaning the air filter.....	29
Maintenance .....	30
Maintenance .....	30
Maintenance intervals.....	31
Checking the thermal fluid .....	31
<b>Troubleshooting</b>	<b>32</b>
Malfunctions and troubleshooting.....	32
<b>Assembly</b>	<b>33</b>
Electrical connection .....	33
Peripheral connections .....	34
Overview of connections.....	34
Connecting the alarm output.....	34
Temperature control hoses.....	35
Connecting the chiller to the external consumer.....	36
Calibrating the temperature sensor .....	37
<b>Disassembly and Storage</b>	<b>38</b>
Dismantling, transportation and storage .....	38
Emptying the device .....	38
Transportation and storage .....	38
<b>Accessories and Spare Parts</b>	<b>39</b>
Scope of delivery.....	39
Accessories.....	39
<b>Attachments</b>	<b>40</b>
Technical data .....	40
Thermal fluids .....	41
Cooling unit .....	42
<b>Service</b>	<b>43</b>
Contact / Technical Service .....	43
Warranty.....	44
Certificate of decontamination.....	45

# About this Document

## Copyright

This operating manual is protected by copyright and is exclusively intended for the internal use of the purchaser.

Handing over this operating manual to third parties, reproducing it in any type or form – in part or in full – and utilizing and/or communicating its contents is prohibited without prior written approval from the manufacturer.

Any violations will result in an obligation to provide compensation for damages. Other claims remain unaffected.

## ➤ Version and variants

### Version

This manual describes the function, operation and maintenance of the circulation chillers High-CHILL 250, High-CHILL 350, High-CHILL 600 and High-CHILL 1200.

Version	Alteration date
1.0	11/2022

### Variants

The device is available in different variants. Certain features or functions are only available in specific product variants. These variants are described in this operating manual.

## ➤ About this manual

This operating manual is an integral part of the device described here.

- ➔ Carefully read through this operating manual and observe all of the safety instructions and warnings notices.
- ➔ Ensure that every user has carefully read the operating manual before using the device for the first time.
- ➔ Store this operating manual so that it is accessible to all users at all times.
- ➔ Pass the operating manual on to any subsequent owner.

The current version of this operating manual can also be found at:  
[https://heidolph-instruments.com/en/service/downloads/operation-manuals#Rotary Evaporators/Chillers \(High-CHILL 250-1200\)](https://heidolph-instruments.com/en/service/downloads/operation-manuals#Rotary%20Evaporators/Chillers%20(High-CHILL%20250-1200))

In addition, also observe the local regulations at the installation site.

## ➤ Symbols and keywords

Standardized symbols and signal words are used in this operating manual to warn against any dangers and provide important instructions. These instructions must be strictly observed to avoid accidents and damage.

The following terms and basic symbols are used:

Symbol	Additional signal words / explanation
<b>Warning signs</b> 	The yellow triangle indicates hazardous situations. It is used in combination with the following signal words: <b>DANGER:</b> Indicates a hazardous situation which, if not avoided, will result in serious injury or death. <b>WARNING:</b> Indicates a hazardous situation which, if not avoided, may result in serious injury or death. <b>CAUTION:</b> Indicates a hazardous situation which, if not avoided, may result in property damage and minor or moderate injury.
<b>Prohibitory signs</b> 	Prohibited: The red circle indicates a situation that should be avoided under all circumstances and which, if not avoided, may result in serious injury or death.
<b>Mandatory signs</b> 	Must be observed: The blue circle indicates important information. Please observe this information to avoid any property damage.

### Other symbols used:

Symbol	Description
→	Handling instruction, action required
✓	Result of action
▪	List of information
▪	
▪	
a.	List of variants
b.	
c.	

## ➤ General safety instructions

### EU Declaration of Conformity

This device complies with the following EC-Directives\*:

- 2006/42/EC Machinery Directive
- 2014/30/EU Electromagnetic Compatibility Directive

\*Also see the attachment "EU Declaration of Conformity".

The device has been manufactured according to state-of-the-art technology and in compliance with recognized safety regulations. However, risks may still arise during installation, operation and maintenance.

➔ Please ensure that the operating manual is available at all times.

#### The device may only be used under the following conditions:

- ➔ Only operate the device if it is in a technically perfect condition.
- ➔ Only operate the device if you are properly aware of the risks and required safety measures.
- ➔ Only operate the device in accordance with the instructions given in this operating manual.
- ➔ If any information is ambiguous or missing, ask your superior or contact the manufacturer.
- ➔ Do not operate anything on the device without authorization.
- ➔ Only use the device in accordance with its intended use.

### Intended use

The device is intended for use by trained and authorized personnel only.

- |  |   |
|--|---|
| The device is suitable for the following uses: | <ul style="list-style-type: none"> <li>• Cooling of rotary evaporators: Hei-VAP Core, Hei-VAP Expert, Hei-VAP Ultimate</li> <li>• Tempering and delivering non-combustible heat transfer liquids in a closed circuit</li> <li>• Economical cooling as alternative to tap water systems</li> </ul> |
|--|---|

- |   |  |
|---|--|
| The device may be operated in research laboratories, other laboratories and production facilities in the following sectors: | <ul style="list-style-type: none"> <li>• Chemistry</li> <li>• Pharmacy</li> <li>• Biology</li> <li>• Environmental analytics</li> <li>• Basic research</li> <li>• Similar research laboratories</li> </ul> |
|---|--|

### Reasonable foreseeable misuse

Any use which deviates from the device's intended use is considered to be improper. The manufacturer is not liable for damage that occurs as a result. The risk is borne by the operator alone.

Amongst other things, the following uses are considered to be foreseeable misuse:

- Operation of the device without heat transfer liquid
- Incorrect connection of tubes
- Placement of the device on a tabletop surface, only permitted for Hei-CHILL 250 and Hei-CHILL 350
- Setting an incorrect pump pressure
- Operation of the device with a glass reactor without overpressure protection
- Chilling foodstuffs

# Installation

### Electrical safety

- The device may only be connected up if the available power supply voltage corresponds to the information stated on the rating plate for the device.
- The power supply connection must be easily accessible at all times.
- An RCD circuit breaker (residual current circuit breaker) must be fitted to the electrical system in the building so that the device is disconnected from the power circuit in the event of a fault.
- Repairs may only be performed by electricians authorised by Heidolph Instruments.
- Never operate the device with a damaged power cord.
- Always turn the device OFF and disconnect the power cord before carrying out any maintenance or repair work.

### Thermal fluids

- The device is exclusively designed for nonflammable thermal fluids in Class I according to DIN 12876-1.
- Each thermal fluid covers a specific temperature range. This temperature range must match the temperature range of your application. If these temperatures are exceeded or undercut, serious hazards to health or environment can be caused (also see chapter "Attachment, Technical data", section "Thermal fluids").
- We recommend to use the Heidolph thermal fluid. This has been tested and approved by Heidolph Instruments.
- The safety data sheet must be consulted for the intended use of the device.

### Tubes

- Besides material durability please note, that permissible temperature range and maximum pressure must match with your application when selecting tubes for connecting your periphery devices.
- We recommend to use the Heidolph tubes. These have been tested and approved by Heidolph Instruments.

# Qualifications of personnel

- The device may only be operated by trained personnel.
- The device may only be operated by persons who have been instructed and supervised in its proper use by trained specialist personnel.
- The device may only be operated by specialist personnel who are above the legal minimum age.
- Other persons may only work on the device under the constant supervision of experienced and trained specialist personnel.
- This operating manual must be read and understood by all persons working with the device.
- Personnel must receive safety training that ensures responsible and safe working practices.

# ➤ Operating company's obligations

## Installation site

- The device must be installed in a suitable location.
- The device must be installed on a firm and stable surface.
- Ensure that the device and all of its components and connections are easily accessible at all times.
- The front and side vents of the device must be a minimum of 40 cm (16 inch) away from walls or vertical surfaces so air flow is not restricted.
- It is not permitted to place or store any objects such as accessories, tools or chemicals within this clearance area during operation.
- All screw connections must be securely tightened.
- It is not permitted to operate the device near to highly flammable or explosive substances.
- Only operate the device (depending on the substances being used) in conjunction with an extractor hood (see DIN EN 14175 and DIN 12924).
- The device is designed for indoor use only and under the following ambient conditions:

Ambient temperature	5 - 31 °C at 80 % relative humidity 32 - 40 °C decreasing linearly to a maximum 50 % relative humidity
Installation altitude	0 - 2,000 m above sea level
Contamination level	2
Overvoltage category	II
Permissible supply deviations	± 10 %

- If the device is operated in corrosive atmospheres, the service life of the device will decrease based on the concentration, duration and frequency of the exposure to the corrosive atmosphere e.g. concentrated hydrochloric acid (HCl).

## Modifications to the device

- It is not permitted to make any unauthorised modifications or changes to the device.
- Do not attach or install any parts that have not been approved by the manufacturer.
- Unauthorized modifications or changes will void the EC Declaration of Conformity for the device and operation of the device will no longer be permitted.
- The manufacturer is not liable for any damage, dangers or injuries that result from unauthorized modifications and changes or due to the non-observance of the instructions in this manual.

## Safety of personnel

- ➔ Ensure that the device is only operated by qualified specialist personnel and trained employees.
- ➔ Observe the following instructions to avoid any personal injuries and property damage:
  - Laboratory regulations
  - Accident prevention regulations
  - Hazardous Substance Act
  - Other generally accepted occupational health and safety regulations
  - Local regulations

### ➤ Safety during use

- Wear the appropriate clothing when working on the device (protective clothing, protective glasses and, if necessary, safety gloves).
- Do not use the device in potentially explosive areas. The device is not protected against explosion. There is no explosion or ATEX protection available.
- Do not carry out work with naked flames in the vicinity of the device (risk of explosion).
- Do not operate or assemble devices in the vicinity which are emission or radiation sources (electromagnetic waves) for the frequency range ( $3 \cdot 10^{11}$  Hz to  $3 \cdot 10^{15}$  Hz).
- Do not operate or assemble appliances in the vicinity of the device which constitute emission or radiation sources for ionizing radiation or in the ultrasonic range.
- Do not operate the device where adiabatic compression or shock waves might occur (shock wave combustion).
- Attach all tube connections securely to ensure safe operation.
- Connect all cables and tubes without kinks and locate them outside the operating and danger zone.
- Do not spill liquids over the device or any parts of it.
- Remove any accidentally spilled liquids immediately.
- Eliminate errors immediately.
- Avoid putting pressure on the display when you are not operating the device.
- Do not use abrasive material to clean the surface. Only wipe with damp cloth.
- Always switch the device OFF after use.

### ➤ Disposal

#### Refrigerant

- The device may only be disposed of by a qualified refrigeration technician.
- Dispose of the refrigerant in accordance with EC Directive 215/2067 in combination with 517/2014/EC. (Type and quantity of the refrigerant can be seen on the rating plate; also see chapter „Attachement, Technical data“, section „Refrigerant unit“.)

#### Device and packaging

- Check the device components for hazardous substances and solvents.
- Clean all components before disposal.
- Dispose of the device in accordance with the relevant national regulations.
- Dispose of the packaging material in accordance with the appropriate national regulations.


## ➤ Device overview

### Device overview Hei-CHILL 250 / 350 / 600 / 1200

# Hei-CHILL control panel

### ➤ Setting up the device

**Danger:**  
**Risk of poisoning!**

Toxic vapours may be emitted depending on the thermal fluid being used and the type of operation.

Ensure that these vapours are extracted adequately.

**Warning:**  
**Risk of crushing on the Hei-CHILL 600 and Hei-CHILL 1200**

The chillers are set up on the floor and fitted with castors to make them easier to manoeuvre. As a result, the device could roll away and possibly tip over.

Set up the device on a level, non-slip surface with sufficient load-bearing capacity!

Secure the device against rolling away by applying the locking brake!

**Caution:**  
**Possible damage to the cooling system**

The liquid from the compressor may leak into the circulation system for the thermal fluid if the device is not set up in an upright position.

The device may not be tilted. Only set up, move and operate the device in an upright position.

**Caution:**  
**Damage to the device**

The device may overheat in the event of insufficient airflow.

Maintain a minimum distance of 40 cm between the ventilation grille and walls or other vertical surfaces.

- The chiller can be set up lower than the apparatus to be cooled. In this case the pump must be switched on and all tubes must be checked for tightness. Otherwise thermal liquid may overflow or discharge.
- As long as the process is operated as a closed circuit, there is no danger of the device overflowing when topping up the thermal fluid.
- You can position multiple devices next to one another.

- ➔ Set up the device on a solid, smooth, level and temperature-resistant surface.
- ➔ Clean the installation surface and the feet or castors on the device at regular intervals.
- ➔ Apply the locking brake on those devices fitted with castors to prevent the device rolling away.
- ➔ Do not place any heavy objects on the device.

#### **Hei-CHILL 250 and Hei-CHILL 350**

The devices can be set up on an even surface.

- ➔ Carry the device by gripping underneath the device.


## Start-Up

### Hei-CHILL 600 and Hei-CHILL 1200

The devices may only be set up on the floor.


- Release the locking brake on the castors by pressing the lever downwards to [Off].
- Apply the locking brake to the castors by pressing the lever upwards to [On].

## ▶ Connecting/disconnecting the power cord


### Connecting the power cord

The power cord has a three-pole connector.

The device socket is located on the rear of the device.

- ✓ The device is switched off.
- Insert the power cord connector into the device socket.
- ✓ The connector has been securely inserted when it has been pushed in up to the stop.
- Now insert the power cord plug into the power socket.

### Disconnecting the power cord

- First, remove the plug from the power socket.
- Then disconnect the power cord connector from the device.


## ▶ Switching the device on and off

### Switching on the device

- Press the toggle switch upwards to "I".
- ✓ The device is switched on.

### Switching off the device

- Press the toggle switch downwards to "O".
- ✓ The device is switched off.

### ➤ Filling with thermal fluid

**Danger:**  
**Risk of fire**

If the device is filled with a thermal fluid with a very low flash point, it could catch fire under certain circumstances.

Only use thermal fluids with a flash point that is at least 20 K above the temperature range for the application. Do not use any ethanol as an additive because it is a highly flammable substance.

**Warning:**  
**Risk of electric shock**

If the device is overfilled, thermal fluid may overflow into the device and cause the device to short circuit.

Ensure that the device is not overfilled. Carefully monitor the fill level indicator and the thermal expansion of the thermal fluid.

**Warning:**  
**Risk of electric shock**

If the device is filled in a careless manner, thermal fluid may splash into the device and cause the device to short circuit.

Use a funnel for the filling process. Avoid splashing the thermal fluid.

**Warning:**  
**Risk of poisoning! Risk of injury!**

When handling chemically prepared thermal fluids, any splashes could cause injuries to areas of bare skin.

Always wear safety goggles, protective gloves and protective clothing when filling the chiller.

**Caution:**  
**Possible damage to the device**

If you operate the device without thermal fluid, the cooling circuit pump could overheat and thus damage the device.

Only operate the device if it has been sufficiently filled with thermal fluid.

The thermal fluid from Heidolph is suitable for all standard applications. It covers a temperature range from minus 30 °C to plus 90 °C. (See chapter "Spare parts and accessories, Accessories")

## Filling for the first time


Depending on the version of the device, the reservoir for the thermal fluid can hold between 2-14 litres (see chapter "Attachments, Technical data"). The reservoir on the device is filled with thermal fluid from above.

Filling of the device is supported by a software program. If the fill level is too low:

- the pump and cooling unit will not start
- FILL will appear on the display and the fill level indicator will light up after the device is switched on.


- Connect up your external consumer (see chapter "Assembly, Peripheral connections").
- Close the drain screw by turning it to the right as far as it will go.

Drain screw

- Attach a suitable hose with an inner diameter of  $\varnothing$  8 mm to the overflow pipe on the device.
- Place the hose into a suitable container to capture any overflowing thermal fluid.

Overflow pipe


If air penetrates into the temperature control circuit while filling the device – e.g. due to an open or damaged bleed valve – the reservoir and also the connected consumer could empty.

Ensure that you use a sufficiently large container for the overflow pipe.

- Switch the device ON.
- ✓ A signal tone sounds.
- ✓ The software version will be shown on the display, followed by "FILL".
- Carefully remove the cover for the filling pipe by pulling it upwards (no screw connection).
- Carefully fill the reservoir with thermal fluid via the filling pipe with the aid of a funnel.

Filling pipe for reservoir

Sealing cap for reservoir

## Start-Up

- The maximum fill level for the device is shown by the upper arrow.
- The minimum fill level for the device is shown by the lower arrow.

The fill level indicator should not rise above the maximum filling level.


→ Observe the fill level indicator and fill the device up to the maximum fill level.

→ Press the confirmation button when the maximum fill level has been reached or sufficient fluid has been added.

- ✓ The pump starts.
  - ✓ The current temperature is shown on the display.
  - ✓ The fill level will now fall because the external consumer is being filled.
- Carefully refill the thermal fluid.
- ✓ If the fill level falls too far, the device will automatically switch to FILL mode and the pump and cooling unit will be switched off.
- Continue to fill the device until the fill level indicator no longer falls.

✓ The reservoir has been filled sufficiently when the fill level settles between the minimum and maximum marks on the fill level indicator.

→ Close the reservoir by carefully pushing the cover into the filling pipe.

### ➤ Menu navigation

The menu is divided into a main menu and 2 submenus.

Use the confirmation button to access the next menu and also to save changed values where relevant.


Use the arrow buttons to move up and down through the parameters in the menu and to change the values set for the parameters in the submenu.


#### Main menu

The current temperature of the device is shown here.

#### Submenus 1 and 2

You can select the different parameters and change the values where relevant here.

- Switch the device ON.
- Press the confirmation button.
- ✓ You will call up submenu 1.
- Press the up  or down  arrow buttons.
- ✓ The different parameters will be displayed in sequence (see chapter "Operation, Overview of menu structure").
- Press the confirmation button to select a parameter.
- ✓ You will call up submenu 2.
- ✓ The set value for the selected parameter will be displayed.
- Press the up  or down  arrow buttons.
- ✓ The set value will be changed upwards or downwards.
- Press the confirmation button to confirm the value.
- ✓ The set value will be changed upwards or downwards.

## ➤ Entering the temperature setting

Before placing the chiller into operation, the desired temperature setting needs to be entered depending on the thermal fluid being used.

**Warning:**  
**Risk of poisoning! Risk of injury!**

Each individual thermal fluid covers a certain temperature range. If the temperature exceeds or falls below this temperature range, it could damage the device and cause the thermal fluid to escape. There is a risk of slipping and poisoning in the event of escaping thermal fluids.

Check the permissible temperature range for the thermal fluid being used and carefully set the temperature limits on the device.

The desired cooling temperature is maintained by defining the following two variables:

- **Temperature limit:**  
This value defines the upper and lower maximum values that must not be exceeded. If the temperature exceeds or falls below these values, a warning signal will sound.
- **Temperature setpoint:**  
This value defines the cooling temperature for the thermal fluid. Until this setpoint is achieved, the current temperature will fluctuate above and below this setpoint and thus also between the upper and lower temperature limits.

## Setting the temperature limits

Both values for the temperature limits are dependent on the thermal fluid being used.

**Possible setting range:**

Lower temperature limit	"Lo" -12° C to 41° C
Upper temperature limit	"Hi" -11° C to 42° C

**Sensible temperature limits are:**

Heidolph thermal fluid Temperature limits: (also see chapter "Spare parts and accessories, Accessories")	"Lo" min. -12 °C "Hi" max. 42 °C
--	-------------------------------------


## Start-Up

### Setting the upper temperature limit

- ✓ The device is switched on.
- Press the confirmation button.
- Press the down arrow button  $\blacktriangledown$  once.
- ✓ "Hi" appears on the display.
- Press the confirmation button.
- ✓ The display flashes.
- Press the up  $\blacktriangle$  or down  $\blacktriangledown$  arrow buttons to set the desired value.
- Confirm the desired value with the confirmation button.
- ✓ The set value is saved.
- ✓ You return to the main menu.

Upper temperature limit

If you do not press any buttons for more than 4 seconds, the changed values will be saved. You return automatically to the main menu.

### Setting the lower temperature limit

- ✓ The device is switched on.
- Press the confirmation button.
- Press the down arrow button  $\blacktriangledown$  twice.
- ✓ "Lo" appears on the display.
- Press the confirmation button.
- ✓ The display flashes.
- Press the up  $\blacktriangle$  or down  $\blacktriangledown$  arrow buttons to set the desired value.
- Confirm the desired value with the confirmation button.
- ✓ The set value is saved.
- ✓ You return to the main menu.

Lower temperature limit

### Setting the temperature setpoint

The temperature range for the setpoint is automatically limited by the temperature limits. It can be set to a maximum of 2 °C below the upper temperature limit and 2 °C above the minimum temperature limit.

In the case of a temperature setpoint or current temperature < 5 °C, the yellow LED on the device will light up. It warns you about the possible incorrect use of the thermal fluid and resulting damage to the device.

- Press the confirmation button.
- Use the arrow button to select the submenu for entering the setpoint "Set".
- Set the desired temperature setpoint.
- Save the temperature setpoint with the confirmation button.

Entered setpoint

If the temperature setpoint you have entered is outside of the defined temperature limits, a warning signal will sound. The value will not be saved.

## ➤ Overview of menu structure

Current temperature

11.0

Enter setpoint

SET

Setpoint

10.0

Set upper limit

H.1

Upper limit

40.0

Set lower limit

Lo

Lower limit

-10.0

Activate / deactivate  
autostart

A.Std

Autostart activated

Auto

Autostart deactivated

Off

Configure alarm output

AL.Ou.

Alarm for: Error and alarm

E.A.

Alarm for:  
Error, alarm and warning

EALW

Adjust temperature sensor

CAL

Temperature sensor  
adjusted

done

Restore factory settings

DEF

Factory setting restored

done

End of menu

End

## ➤ Temperature control operation

### **Caution:** **Possible damage to the device**

If the ventilation openings are not kept free, the motors for the pump and compressor could overheat.

Never cover the ventilation grille. Ensure that there is always a distance of 40 cm between the device and wall or other devices.

- ✓ All peripheral devices are securely connected and the refrigerant can circulate freely.
- ✓ The fill level indicator for the thermal fluid lies between the minimum and maximum levels.
- Switch the device ON.
- ✓ A signal tone sounds.
- ✓ The software version will be displayed, followed by the current temperature.
- ✓ If the fill level is sufficient, "StbY" for standby operation will flash for around 4 seconds on the display.
- Press the confirmation button "Select".
- ✓ The pump for the thermal fluid circuit will start.
- ✓ Depending on the ambient temperature and the temperature setpoint, the cooling unit will start up after a short time.
- ✓ As soon as the cooling unit is active, the blue cooling LED will light up.
- Monitor the fill level indicator and refill the device with thermal fluid if necessary.

### **Caution:** **Possible damage to the device**

If you set the chiller to cool down to a minus temperature, it is possible that the thermal fluid will freeze.

If you want to cool the device down to minus temperatures, ensure that you are using an appropriate thermal fluid.(also see chapter "Spare parts and accessories, Accessories")

If there is a large difference in temperature between the chiller and the external consumer, the temperature drop between the two devices could be so high that it will not be possible to maintain the cooling output.

In this case, increase or decrease the setpoint for the cooling temperature accordingly or ensure that there is a lower volume of items to be cooled in your external consumer.

## › Interfaces

**Warning:**  
**Risk of electric shock!**

If the voltages received by the interface inputs and outputs are too high and not sufficiently insulated, metal components such as e.g. the housing could become live in the event of a fault.

Safely isolate small voltage inputs and outputs above 25 V AC or 60 V DC according to DIN EN 61140 and use double or reinforced insulation according to DIN EN 60730-1 or DIN 60950-1.

Exclusively use shielded connecting cables. Connect the shielding to the connector housing.

**Caution:**  
**Possible damage to the device**

If the connections on the device are not covered, liquids could reach the live components. Liquids penetrating into the device could result in faults at the interfaces and cause a short circuit.

The interface connections on the rear of the device must always be kept covered.


The display can be used to set which signal is issued via the interface in each particular fault situation.

## RS 232 interface

The Heidolph rotary evaporators Hei-VAP Expert and Hei-VAP Ultimate can communicate with the chiller and control the cooling process via the RS 232 interface (see operating manual for the rotary evaporator).

The RS 232 interface and the alarm output are located on the left rear of the device above the power connection.

## Interface commands

*Write commands	Response from chiller	Meaning
OUT_SP_00_XX.X	OK	Temperature setpoint thermal fluid
OUT_SP_04_XX.X	OK	Flow line temperature: Limit upper value
OUT_SP_05_XX.X	OK	Flow line temperature: Limit lower value
START	OK	Switch device to standby
STOP	OK	Switch device to standby (pump, cooling unit off)

*Read command	Response from chiller	Meaning
IN_PV_00	XX.X	Request bath temperature (flow line temperature)
IN_SP_00	XX.X	Request temperature setpoint
IN_SP_04	XX.X	Request flow line temperature limit "Hi"
IN_SP_05	XX.X	Request flow line temperature limit "Lo"
TYPE	X	Request device type:
	1	Hei-CHILL 250
	2	Hei-CHILL 350
	3	Hei-CHILL 600
	4	Hei-CHILL 1200
VERSION	XX.X	Request software version number
STATUS	X	Request device status:
	0	OK
	-1	Fault
STAT	XXXXXXXX	Request fault diagnosis:
	0000000	No fault
	I000000	Fault (hardware)
	0I00000	Alarm: Pump is overheated
	00I0000	Alarm: PCB (board) too hot
	000I000	Alarm: Fill level for thermal fluid too low
	0000I00	Warning: Fill level for thermal fluid too low
	00000I0	Warning: Temperature above setpoint
	000000I	Warning: Temperature below setpoint
SERIAL	SXXXXXXXXX	Serial number analogue to rating plate

\*Interface parameters RS 232: 9600 Baud, No Parity, 8 Bit, 1 Stopbit

## Operation

Error message	Description
ERR_2	Incorrect entry (e.g. data buffer overflow)
ERR_3	Incorrect command
ERR_5	Syntax error in value
ERR_6	Impermissible value
ERR_32	The upper temperature limit is lower or the same as the lower temperature limit.

- The RS 232 interface can be operated with or without a hardware handshake (RTS/ CTS).
- Do not send commands in a package, pause of at least 0.1 seconds between commands
- End every command with CR (Carriage Return, Hex=0D) or LF (Linefeed, Hex=0A)
- The response from the chiller always ends with CRLF
- Date format XX.X or -XX.X (X stands for a number)
- As well as "\_", " " (blank space) is also permissible

## Cable and interface test for RS232

Signal	Computer				Thermostat		Signal
	9-Pole Sub-D socket		25-Pole Sub-D socket		9-Pole Sub-D socket		
	with a hardware handshake	without a hardware handshake	with a hardware handshake	without a hardware handshake	with a hardware handshake	without a hardware handshake	
RxD	2	2	3	3	2	2	TxD
TxD	3	3	2	2	3	3	RxD
DTR	4	-	20	-	4	-	DSR
Signal Ground	5	5	7	7	5	5	Signal Ground
DSR	6	-	6	-	6	-	DTR
RTS	7	-	4	-	7	-	CTS
CTS	8	-	5	-	8	-	RTS

### With a hardware handshake:

- If connecting a thermostat to the PC, use a straight-through connection (1:1) and not a null modem cable.
- The RS232 interface can be directly connected to the PC with a straight-through (1:1) cable.

### Without a hardware handshake:

- Set the corresponding operating mode on the PC.
- Use shielded connecting cables.
- Connect the shielding to the connector housing.
- Galvanically isolate the cables from the rest of the electronics.
- Do not connect up unused pins.

## Alarm output

The device differentiates between three different fault messages.

- **Alarm:**  
An alarm is triggered if a fault is relevant to safety. The corresponding components on the device will be switched off. The device switches to standby operation.
- **Error:**  
An error indicates a fault that was caused by a hardware error. The device switches to standby operation.
- **Warning:**  
Warnings indicate faults that are not relevant to safety. The device continues to operate.

All fault messages are issued as an error on the display (see chapter "Troubleshooting, Malfunctions and troubleshooting").

In addition, the faults can be output to a connected device via the alarm output (also see chapter "Assembly, Peripheral connections", section "Alarm output").

The device is configured in the factory to always output safety-relevant faults to the alarm output. You can also define whether the fault should be output in the event of the following additional faults.

### Defining fault messages for the alarm output

- ✓ The device is switched on.
- Call up the alarm output submenu "AL.Ou."
- Call up the submenu with the confirmation button.
- Select one of the two following options in the submenu:
- ✓ The alarm output is activated for:
  - an error
  - an alarm
- ✓ The alarm output is activated for:
  - an error
  - an alarm
  - a warning

Alarm output

Errors and alarms option

Option with additional warnings

- Save the setting with the confirmation button.
- You return to the main menu.

The setting needs to be saved using the confirmation button within 4 seconds. Otherwise, the setting will not be saved.

## Autostart

The device switches off in the event of a power failure.

The device will be switched on again as soon as the device is connected to the power grid again. However, the pump for the thermal fluid circuit needs to be restarted.

If the device is connected to one of the two rotary evaporators Hei-VAP Expert and Hei-VAP Ultimate, these two rotary evaporators will control how the device is restarted.

You can define in the "Autostart..." submenu whether the device should automatically start in temperature control operation after a power failure.

**Warning:**  
**Risk of explosion! Risk of frostbite! Risk of poisoning!**  
**Risk of injury! Risk of slipping!**

If the device automatically resumes full operation after a power failure, loose screw connections or defective glass containers could result in the uncontrolled escape of liquids. The liquids that may escape as a result pose a risk of frostbite, poisoning and slippage. Burst glassware also poses a significant risk of injury and can result in exothermic reactions.

If you have activated the autostart function, do not leave the chiller and connected peripheral devices unsupervised under any circumstances.


Check that the hose connections are secure at regular intervals. Only work with flawless glass containers.

### Activating / deactivating the autostart

- ✓ The device is switched on.
- Call up the autostart submenu "A.Std".
- Call up the submenu with the confirmation button.
- Select one of the two following options in the submenu:
- ✓ Autostart is activated.
- ✓ The device will resume full operation after a power failure.

Autostart

Autostart activated

Autostart deactivated

- ✓ Autostart is deactivated.
- ✓ The device will be switched to standby operation after a power failure.
- ✓ Temperature control operation will need to be restarted.

## ➤ Restoring the factory settings

All settings that you have changed for the parameters in the submenu can be restored to the original factory settings under the factory settings menu item.

- ✓ The device is switched on.
- Call up the factory settings submenu "dEF".
- Briefly press the confirmation button once.
  
- Press the confirmation button again and hold it down for around 3 seconds.
- ✓ "donE" will be shown on the display.

Factory settings

- ✓ You return to the main menu.
- ✓ The factory settings have been restored:
  - The temperature limits will have been reset to 45 °C and 5 °C.
  - The faults that are output via the alarm output will be reset to alarms and errors.

## ➤ Cleaning

### Cleaning the device

The housing and surfaces of the device can be cleaned using a damp cloth and, if required, a mild soap solution.

**Warning!**  
**Risk of electric shock!**

Liquids penetrating into the device could cause the device to become live if it is still connected to the electricity supply.

Before cleaning the device:

- Switch off the device.
- Disconnect the power cord plug.

**Caution:**  
**Damage to the device**

Liquids penetrating into the device could damage the electronics. Only wipe the surfaces with a damp cloth and not a wet cloth.

**Caution:**  
**Damage to surfaces**

The surfaces could be damaged by improper cleaning.

Do not use the following under any circumstances:

- Chlorine bleach or cleaning agents containing chlorine
- Substances containing solvents (e.g. acetone)
- Ammonia
- Abrasive cleaning agents such as cleaning wool, scouring agents or cleaning agents with metallic particles

### Cleaning the air filter

As well as the mechanical fixtures, the ventilation grille in front of the air filter is also held in place with magnets.

- ✓ The device is switched off and disconnected from the power grid.
- ➔ Reach under the ventilation grille using both hands.
- ➔ Lightly push the ventilation grille upwards and then pull it out of the device at the bottom.
- ➔ Brush or vacuum the air filter and, if necessary, the ventilation grille.
- ➔ Push the ventilation grille upwards so that the holes are above the two retaining pins.
- ➔ Carefully push the ventilation grille onto the device again at the bottom.

## ➤ Maintenance

### Maintenance

---

**Danger!**  
**Risk of poisoning! Risk for the environment!**

Improper repairs to the internal cooling circuit may result in leaks and the refrigerant could escape in an uncontrolled manner. This could result in employees being poisoned and also serious environmental damage.

Repair and maintenance work on the internal cooling circuit may only be carried out by specialist personnel authorised by Heidolph Instruments.

---

**Warning!**  
**Risk of electric shock! Risk of injury!**

You could come into contact with live or moving parts when carrying out repairs. This could result in electric shocks, collisions, cutting or crushing

The device must be disconnected from the power grid before all maintenance work.

Repairs may only be performed by specialist personnel authorised by Heidolph Instruments.

---

**Warning!**  
**Risk of burns! Risk of frostbite**

Touching hot or cold device components, accessories and thermal fluid could result in burns and also frostbite.

Always wear protective gloves and protective clothing when carrying out maintenance work.

Allow the device components, accessories and thermal fluid to reach room temperature before touching them.

---

If the device has come into contact with hazardous materials, the device needs to be decontaminated before completing any maintenance work.

---

### Maintenance intervals

The maintenance intervals stated in the following table must be observed. The following maintenance tasks are obligatory before all prolonged unsupervised operation.

<b>Interval</b>	<b>Maintenance work</b>
Daily	Check the drain screw by inspecting it from the outside
Monthly	<ul style="list-style-type: none"><li>• Check the external hoses for material fatigue</li><li>• Clean the air filter</li></ul>
Six monthly	Check the thermal fluid

Any repairs that may be required may only be performed by a specialist authorized by Heidolph Instruments.

In this case, contact your Heidolph Instruments dealer or a representative of Heidolph Instruments.

### Checking the thermal fluid

Contaminated or degenerated thermal fluid must be replaced. Continuing to use the thermal fluid in this case is only permitted after obtaining corresponding test results. The thermal fluid must be checked in accordance with DIN 51529.

## Troubleshooting

### Malfunctions and troubleshooting

**Alarms** are relevant to safety. The device will be switched to standby operation.

- The red warning LED flashes.
- A short acoustic signal will sound rhythmically around 5 minutes after the start of the fault.

Alarm	Possible cause	Remedy
	Fill level for thermal fluid below minimum level	<ul style="list-style-type: none"> <li>• Switch off the device</li> <li>• Refill thermal fluid</li> <li>• If necessary, check leak tightness of connections</li> </ul>
	Pump is blocked.	<ul style="list-style-type: none"> <li>• Switch off the device</li> <li>• Reduce the viscosity of the thermal fluid</li> <li>• Remove any possible contamination of the thermal fluid</li> </ul>
	The temperature of the electronics is higher than 75 °C.	<ul style="list-style-type: none"> <li>• Switch off the device and allow it to cool down</li> <li>• Check the air intake and outlet for free circulation of air</li> <li>• If necessary, reduce the application requirements to match the cooling output</li> </ul>

**Warnings** are not relevant to safety. The device continues to operate. A prolonged acoustic signal will sound that is then repeated periodically.

Warnings	Possible cause	Remedy
	Fill level for thermal fluid below minimum level	<ul style="list-style-type: none"> <li>• Switch off the device</li> <li>• Refill thermal fluid</li> <li>• If necessary, check leak tightness of connections</li> </ul>
	Upper temperature limit has been exceeded	<ul style="list-style-type: none"> <li>• Switch off the device</li> <li>• Check if the application is generating heat</li> <li>• Check the setpoint</li> </ul>
	Temperature has fallen below lower limit	<ul style="list-style-type: none"> <li>• Switch off the device</li> <li>• Check cooling application</li> <li>• Check the setpoint</li> </ul>

→ After resolving the cause of the fault, restart the device and press the confirmation button.

✓ The warning or alarm has been resolved.

**Error:** If the fault occurs again after restarting the device or you receive the error message "EXXX", note down the message shown on the display and contact Heidolph Instruments or your local Heidolph Instruments dealer.

## ➤ Electrical connection

Repairs to the device must **only be carried out by a qualified electrician approved by Heidolph Instruments**. Improper repairs could result in serious dangers. Please contact your specialist dealer if the device needs to be repaired.

### Warning: Risk of electric shock

The device must be connected to an earthed power socket.

The rated voltage for the device and the supply voltage must be the same. The rated voltage for the device can be found on the rating plate on the rear of the device.

The device is earthed via the supplied power cord when delivered. If the original plug is replaced, it is essential that the new plug is connected to a protective earth!

### Colour coding for the electrical connection:

Europe		North America	
<b>GREEN/ YELLOW</b>	PE: Protective earth <b>(earth)</b>	<b>GREEN</b>	PE: Protective earth <b>(earth)</b>
<b>BLUE</b>	N: Neutral conductor	<b>WHITE</b>	N: Neutral conductor
<b>BROWN</b>	P: Phase	<b>BLACK</b>	P: Phase

The device is connected to the power supply with the supplied power cord. The device socket is located on the rear of the device.


If you want to use the device in a country with a different type of plug:

- The supplied plug may only be replaced by a qualified electrician.
- If you want to use an adapter, it must be earthed and approved in accordance with local regulations.
- The cable, adapter and plug must correspond at least to the power stated on the rating plate.

## Peripheral connections

### Overview of connections

### Connecting the alarm output

The alarm output can be used to output faults to an external device via a potential-free contact.

#### Warning: Risk of electric shock!

If the voltages received by the interface inputs and outputs are too high and not sufficiently insulated, metal components such as e.g. the housing could become live in the event of a fault.

Safely isolate small voltage inputs and outputs above 25 V AC or 60 V DC according to DIN EN 61140 and use double or reinforced insulation according to DIN EN 60730-1 or DIN 60950-1.

Exclusively use shielded connecting cables. Connect the shielding to the connector housing. Cover unused plug connections with protective caps.


## Assembly

A maximum 30 V DC at 0.2A may be connected to the alarm output.

### Sleep state

If the device is switched off or there is a fault after switching it on, the alarm output is in a sleep state:

- Pin 1 and 2 are open.
- Pin 3 and 2 are closed.

### Good state

If the device is switched on and operating without any faults, the alarm output is in a good state:

- Pin 1 and 2 are closed.
- Pin 3 and 2 are open.

## Temperature control hoses

### Warning: Risk of frostbite!

If you use unsuitable hoses at very low cooling temperatures, the thermal fluid may escape and you could experience frostbite if you come into contact with it.

Only use hoses with a temperature resistance that corresponds to the operating temperature range of the device.

### Warning: Risk of frostbite!


If you come into contact with cold hoses, you could experience frostbite if you are working with very low cooling temperatures.

Only use insulated hoses at operating temperatures  $< 0\text{ }^{\circ}\text{C}$ .

All hose connections are located on the top rear of the device.

- Pump flow line
- Pump return
- Overflow pipe

All hose connections are fitted with a hose connector suitable for hoses with an internal diameter of  $\varnothing 8\text{ mm}$ .

Carefully select hoses and connections for the chiller depending on the fluid, temperature and pressure of your application. (also see chapter "Scope of delivery and accessories, Accessories")

The user is responsible for any damage that results from improper connections and the selection of unsuitable hoses.

## Connecting the chiller to the external consumer

### Warning:

#### Risk of electric shock! Risk of frostbite!

If the hoses between the chiller and the external consumer are not leak-tight and connected without kinks, any escaping thermal fluid could result in frostbite and electric shock in the event of contact with live components.

- Lay all hoses without kinks and protected against kinks.
- Check the leak tightness of all hose connections and hoses.
- Secure all hoses against slippage.
- Maintain a clearance distance between hoses and power connections or live components.
- Check the leak tightness of your consumer. It is not permitted to connect any open consumers.

### Warning:

#### Risk of explosion! Risk of frostbite! Risk of poisoning!

#### Risk of injury! Risk of slipping!

If the hoses are improperly connected or glassware is exposed to high pressure, the hoses and glassware could burst due to overpressure.

The liquids that may escape as a result pose a risk of frostbite, poisoning and slippage. Burst glassware also poses a significant risk of injury and can result in exothermic reactions.

The maximum pressure at the outlet from the chiller is stated in the Attachments (see chapter "Attachments, Technical data").

If necessary, reduce the pressure in accordance with the connected peripheral devices by fitting a pressure regulator upstream at the outlet from the chiller.

Ensure that all hose connections are leak-tight and secured against kinks and slippage to avoid overpressure.

If air penetrates into the external consumer and it is higher than the chiller, it could cause the external circuit to empty into the chiller when the pump is switched off.

This could cause the chiller to overflow.

- Shorten the temperature control hoses as much as possible while ensuring that they can be laid without any kinks and without any tension.
- Secure the temperature control hoses against slippage with the aid of hose clamps.

### Calibrating the temperature sensor

The internal thermometer in the device is precisely calibrated upon delivery.

If you detect a possible temperature deviation when checking the thermal fluid using a reference thermometer, the measurement values on the temperature sensor can be adjusted.

For this purpose, the reference thermometer needs to be fitted to the flow line on the device.

The factory calibration will be overwritten when an adjustment is made and cannot be restored! Only use a reference thermometer with a high level of precision for the calibration process.

- ✓ The chiller is switched off and disconnected from the power grid.
- ✓ Allow the device and thermal fluid to warm up or cool down to room temperature.
- Attach a T-piece suitable for the interior diameter of the hose ( $\varnothing$  8 mm) between the connection for the pump flow line and your connection hose.
- Slide your reference thermometer into the still open end of the T-piece and fix the thermometer in place.
- ✓ The calibrated reference thermometer has been fitted in the flow line in such a way that it will ensure a reliable measurement process.

- Switch the device on and start temperature control operation.
- Call up the calibrate submenu "CAL".
- Enter the temperature read from the reference thermometer.

Adjust temperature sensor

- Then press and hold down the confirmation button for around 3 seconds.
- ✓ "donE" will be shown on the display.
- ✓ The new value has been adopted.
- Switch off the chiller and disconnect it from the power grid.

Temperature sensor adjusted

- Remove the reference thermometer and the T-piece.
- Immediately reconnect the open end of the hose to the flow line.
- ✓ The cooling circuit has been closed to ensure reliable operation.
- ✓ The device is ready to use.

## ➤ Dismantling, transportation and storage

### Emptying the device

#### **Warning!** **Frostbite**

You could suffer frostbite if you come into contact with the thermal fluid while emptying it.

Switch off the device at an early stage before emptying it so that the thermal fluid can warm up to room temperature.

Always wear protective gloves and protective clothing when emptying the device.

- Switch the device OFF.
- Allow the device and the thermal fluid to cool down or warm up to room temperature.
- Position a container with a sufficient capacity directly underneath the drain screw.

The thermal fluid will immediately flow out of the device when the drain screw is opened.

- Turn the drain screw in an anticlockwise direction to open it.
- Observe the regulations for the disposal of used thermal fluid (also see chapter "Safety, Disposal").

### Emptying the refrigerant

#### **Danger!:** **Risk of injury! Risk of poisoning! Risk of serious environmental damage!**

The refrigerant unit is under pressure (19 - 28 bar) and contains substances harmful to the environment. Improper emptying of the refrigerant could result in a serious risk of poisoning and also risks to the environment.

The refrigerant in the internal cooling circuit may only be emptied by a qualified refrigeration technician.

### Transportation and storage

- Only store and transport the device and its components after emptying, cleaning and dismantling it down to the basic device.
- When lifting the device, wear protective gloves and grip underneath the device.
- Store and transport the device and its components in the original packaging or in another suitable container to prevent any damage.
- Carefully seal the packaging against any unauthorized or unintentional opening.
- Store the device in a dry and frost-free location.

Improper storage and transport may result in damage to the device and the mechanical components in the device.

Avoid jolts and vibrations when transporting the device.

## ➤ Scope of delivery

Component	Quantity	Product number
Hei-CHILL 250*	1	591-00600-00
Hei-CHILL 350*	1	591-00610-00
Hei-CHILL 600*	1	591-00620-00
Hei-CHILL 1200*	1	591-00630-00
Operating manual	1	01-005-005-83
Warranty registration	1	01-006-002-78
Power cord	1	Country-specific

\* The product number is valid for 230 V devices in the european community. Please contact Heidolph Instruments for further product numbers on country specific varieties.

## ➤ Accessories

Component	Temperatur range °C	Quantity	Product number
Tube with fabric reinforcement, by the meter, not insulated, maximum pump pessure < 16 bar	- 20 to 60	-	591-38000-00
Thermal fluid (10 L)*	- 30 to 90	1	569-00611-00
Thermal fluid (20 L)*	- 30 to 90	1	569-00610-00

\* The safety data sheet for the thermal fluid can be found at: <https://heidolph-instruments.com/en/service/downloads/operation-manuals#Rotary-Evaporators/Chillers>

Additional accessories can be found in our general catalogue or at: [www.heidolph-instruments.com](http://www.heidolph-instruments.com)

## ➤ Technical data

Standard supply voltage	1 x 230 V (50 Hz) or 1 x 115 V (60 Hz)			
Overvoltage category	II			
Protection class (IEC 61140)	1 ⊕			
Protection class (IEC 60529)	IP 32			
Sound pressure level (dB(A)) (in accordance with IEC 61010)	< 70			
Alarm output connection	max. 30 V DC, 0.2 A			
Contamination level	2			
Distance to surroundings (front and rear) (cm)	40			
Classification according to DIN 12876-1 (class/designation)	I/NFL			
Display	7-segment, LED			
Refrigerant	R-134a			
GWP <sub>(100a)</sub> *	1430			
Operating temperature range (°C)	-10 to 40			
Temperature stability (K)	±0.5			
Storage temperature range (°C)	5 - 40			
Transportation temperature range (°C)	-20 to 60			
<b>Model</b>	<b>Hei-CHILL 250</b>	<b>Hei-CHILL 350</b>	<b>Hei-CHILL 600</b>	<b>Hei-CHILL 1200</b>
Cooling	Air-cooled	Air-cooled	Air-cooled	Air-cooled
Cooling output (W)	250	350	600	1200
Rated power intake 230 V / 115 V (W)	230 / 230	500 / 500	700 / 750	1,150 / 1,100
Thermal fluid fill volume min - max (L)	2 - 4	4 - 7	4 - 8	7 - 14
Max. flow rate (L, H <sub>2</sub> O 20°C)	9	9	20	20
Max. flow pressure (bar, H <sub>2</sub> O 20°C)	0.35	0.35	1.3	1.3
Dimensions (B x T x H) (mm)	200 x 350 x 465	240 x 400 x 500	350 x 480 x 595	450 x 550 x 650
Weight (kg)	28	36	52	64

\* Time horizon 100 years – according to IPCC IV (also see following section "Cooling unit").

## Thermal fluids

Each individual thermal fluid covers a recommended temperature range and must be suitable for the temperature range for your application.

→ Never use contaminated or degenerated thermal fluid.

The thermal fluid becomes more viscous at the lower limit of the temperature range and influences the temperature stability, pump output and cooling output. Especially in water, ice may form that could destroy the device.

Vapours and odours may develop at the upper limit of the range.

Therefore, only fully utilise the temperature range where required.

### Permissible recommended thermal fluids

Heidolph designation	Chemical designation	Temperature range in °C	Viscosity (kin) in mm <sup>2</sup> /s (at 20 °C)	Viscosity (kin) in mm <sup>2</sup> /s at temperature
Thermal fluid	Monoethylene glycol / water	-30 to 90	4	50 to -25 °C

\* Safety data sheets for the thermal fluid can be requested at any time if required.

#### Caution:

#### Damage to the device, damage to the environment

If water is used as the thermal fluid, incorrect additives could damage the device and also the environment if the additives escape from the device.

Do not use the following thermal fluids or additives under any circumstances:

- Distilled, deionised, fully demineralised water or sea water (due to their corrosive properties)
- Water containing iron (due to rusting)
- Untreated river water (due to the formation of algae)
- Inorganic salts (due to flocculation)
- Chlorine or ammonia as additives (e.g. in cleaning agents and disinfectants)

The following water is permissible as a thermal fluid:

- Water free of impurities
- Water with a pH value between 6.0 and 8.5
- Ultra pure water and distilled water after the addition of 0.1 g sodium carbonate (Na<sub>2</sub>CO<sub>3</sub>) per litre of water

The proportion of alkaline earth ions in the water must lie between 0.71 mmol/L and 1.42 mmol/L (corresponds to 4.0 and 8.0 °dH). Hard water will cause limescale to develop in the device.

In the case of the Heidolph thermal fluid, the water content will fall after prolonged operation at high temperatures and the mixture will become flammable (flash point 128 °C).

Use a hydrometer to check the mixing ratio.

## Cooling unit

**Danger!:**  
**Risk of injury! Risk of poisoning! Risk of serious environmental damage!**

The refrigerant unit is under pressure (19 - 28 bar) and contains substances harmful to the environment. Improper emptying of the refrigerant could result in a serious risk of poisoning and also risks to the environment.

The refrigerant in the internal cooling circuit may only be emptied by a qualified refrigeration technician.

The compressor in the cooling unit is hermetically sealed and fitted with overload protection that responds to the temperature and power consumption of the compressor.

### Cooling output in kW at ambient temperatures in °C

Cooling output kW	Hei-CHILL 250	Hei-CHILL 350	Hei-CHILL 600	Hei-CHILL 1200
at 20 °C	0.25	0.35	0.6	1.2
at 10 °C	0.2	0.27	0.5	1.05
at 0 °C	0.15	0.20	0.36	0.75
at -10 °C	0.09	0.12	0.15	0.4

\* The cooling output is measured with the thermal fluid maintained at a certain temperature. The ambient temperature for the measurements is 20 °C.

### Refrigerant and filling level

The cooling unit contains refrigerant. Observe the special regulations for its disposal (see here chapter "Safety, Disposal").

The refrigerant R-134a contains fluorinated greenhouse gases. The global warning potential (GWP) is 1430 at a time horizon of 100 years, CO<sub>2</sub> equivalent = 1.0.

Device 230 V 50 Hz, Refrigerant R-134a	Hei-CHILL 250	Hei-CHILL 350	Hei-CHILL 600	Hei-CHILL 1200
Max. filling weight (kg)	0.085	0.095	0.295	0.575
CO <sub>2</sub> equivalent (t)	0.09	0.16	0.15	0.4

Device 220 V 60 Hz, Refrigerant R-134a	Hei-CHILL 250	Hei-CHILL 350	Hei-CHILL 600	Hei-CHILL 1200
Max. filling weight (kg)	0.08	0.095	0.225	0.585
CO <sub>2</sub> equivalent (t)	0.1	0.1	0.3	0.8

## ➤ Contact / Technical Service

### Questions / Repair work

If any aspect of installation, operation or maintenance remains unanswered in the present manual, please get in touch with the following address.

For repairs please call Heidolph Instruments or your local authorized Heidolph Instruments distributor.

### **Warning!** **Danger of poisoning!**

Contaminated devices can lead to severe injury or death of our employees!

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact substance
- take adequate protective action towards our parts receiving and service personnel
- mark the pack in accordance with Ordinance on Hazardous Substances

A "Confirmation of Condition" form can be found at the end of this operating manual.

→ Prior to shipping the device for repair, complete a copy of this form and submit it in advance.

### Contact details:

#### **Heidolph Instruments Germany**

Heidolph Instruments GmbH & Co. KG  
Technical Service  
Walpersdorfer Str. 12  
D-91126 Schwabach / Germany

Tel.: +49 – 9122 - 9920-74

Fax: +49 – 9122 - 9920-84

E-Mail: [service@heidolph.de](mailto:service@heidolph.de)  
[www.heidolph.com](http://www.heidolph.com)

#### **Heidolph Instruments North America**

Phone: 1-866-650-9604  
E-mail: [service@heidolph.com](mailto:service@heidolph.com)  
[www.heidolphNA.com](http://www.heidolphNA.com)

#### **Heidolph Instruments United Kingdom**

Phone: 01799 - 5133-20  
E-mail: [service@radleys.co.uk](mailto:service@radleys.co.uk)  
[www.heidolph-instruments.co.uk](http://www.heidolph-instruments.co.uk)

#### **All other countries**

You will find contact details of your local Heidolph distributor at  
[www.heidolph.com](http://www.heidolph.com)

## ➤ Warranty

---

Heidolph Instruments provides a 3 year warranty for the products described here (excluding glass and wearing parts) if you register using the warranty card enclosed or online ([www.heidolph.com](http://www.heidolph.com)). The warranty period begins with the date of registration. The serial number is also valid without registration.

This warranty covers material and manufacturing defects.

---

In the event of a material or manufacturing defect, the device shall be repaired or replaced free of charge under the terms of the warranty.

Heidolph Instruments shall not assume liability for any damage incurred as a result of improper handling or transport.

Warranty claim?

- ➔ Please notify Heidolph Instruments or your local Heidolph distributor if you wish to make a warranty claim.

# ➤ Certificate of decontamination

Enclose the Certificate of decontamination, duly completed, with your device return. Submissions without a Certificate of decontamination cannot be processed!

## Certificate of decontamination

### IN CASE OF RETURNS

Please fill in the required fields.

**Note: The sender must package the goods properly and appropriately for transport.**

Heidolph Instruments GmbH & Co. KG  
Walpersdorfer Straße 12  
91126 Schwabach, Germany  
Phone: +49 (0) 9122 9920-380  
**Fax: +49 (0) 9122 9920-19**  
Email: service@heidolph.de

### SENDER

Name _____	First name _____
Company/institution _____	Department _____
_____	Workgroup _____
Address _____	
ZC/City _____	
Country _____	Phone _____
Email _____	

### DEVICE DETAILS

Article number \_\_\_\_\_ Serial no. \_\_\_\_\_

Ticket number \_\_\_\_\_

Reason for sending in \_\_\_\_\_  
\_\_\_\_\_

Has the device been cleaned, decontaminated/disinfected? **Yes** **No** (Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous? **Yes** **No** (Please mark as applicable)

If yes, with which substances did the device come into contact?

### LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product (from serial number on S180003250) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

Hei-CHILL 250	591-00600-00
Hei-CHILL 350	591-00610-00
Hei-CHILL 600	591-00620-00
Hei-CHILL 1200	591-00630-00
Hei-CHILL 3000	591-00700-00
Hei-CHILL 5000	591-00700-01

Maschinenrichtlinie / Machinery Directive 2006/42/EG  
EMV-Richtlinie / Electromagnetic Compatibility Directive 2014/30/EU  
Delegierte (RoHS-) Richtlinie / Delegated (RoHS) Directive 2015/863/EU  
Angewandte (harmonisierte) Normen / (Harmonized) Standards applied:  
EN ISO 12100:2010, EN 61326-1:2013, EN IEC 63000:2018

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## UKCA Declaration of Conformity

# UK CA

### Declaration of Conformity in accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,  
Heidolph Instruments GmbH & Co. KG  
Walpersdörfer Straße 12  
91126 Schwabach / Germany

Product: laboratory chiller

Model:

Hei-CHILL 250	591-00600-xx
Hei-CHILL 350	591-00610-xx
Hei-CHILL 600	591-00620-xx
Hei-CHILL 1200	591-00630-xx
Hei-CHILL 3000	591-00700-xx
Hei-CHILL 5000	591-00710-xx

Description:

Device which cools liquid medias

The object of the declaration described above is in conformity with the relevant  
UK Statutory Instruments (and their amendments):

2008 No. 1597

2008 2016 No. 1091

2012 No. 3032

*The Supply of Machinery (Safety) Regulations*

*The Electromagnetic Compatibility Regulations 2016*

*The Restriction of the Use of Certain Hazardous Substances in  
Electrical and Electronic Equipment Regulations 2012*

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN IEC 63000:2018

UK Authorised Representative (for authorities only):

ProductIP ( UK ) Ltd.  
8. Northumberland Av.  
London WC2N 5BY

Signed for and on behalf of Heidolph Instruments GmbH & Co. KG  
Walpersdorfer Straße 12, 91126 Schwabach / Germany

Schwabach, 01.09.2021

Wolfgang Jaenicke  
Managing Director

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices by Heidolph Instruments GmbH & Co. KG fulfill the requirements of the EU directive 2011 / 65 / EU (RoHS) and the amended of directive 2015 / 863.

All devices are compatible with the requirement of the RoHS for the following materials:  
Max. 0,01% of the weight in homogeneous material for cadmium and max. 0,1 % of the weight in homogeneous material for lead, mercury, hexavalent chromium, polybrominated biphenyl, polybrominated diphenyl ether, Di (2-ethylhexyl) phthalate, butyl benzyl phthalate, dibutyl phthalate, diisobutyl phthalate.

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 22.07.2019

Stefan Peters  
Vice President Marketing, Innovation & Technologie

Marcell Sarré  
Vice President Quality Management & Technical Service

# China RoHS Certification

## China RoHS DECLARATION OF CONFORMITY

Heidolph Instruments GmbH & Co.KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

### Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Instruments GmbH & Co.KG products is 25 years.

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This table is created according to SJ/T 11364-2014

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**Note:** Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
- O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
- X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Instruments GmbH & Co.KG may enter into further devices or can be used together with other appliances .

With these products and appliances in particular, Heidolph Instruments GmbH & Co.KG will not take responsibility for the EFUP of those products and appliances.

Schwabach, 06.08.2021

Wolfgang Jaenicke  
Chief Executive Officer CEO

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Instruments GmbH & Co. KG PRODUCTS
© Heidolph Instruments GmbH & Co. KG