Hei-CHILL

Chiller

Leading Safety Standards

Superior Ease of Use

Reduced Cost of Ownership

Operating Manual

Page 2

Circulation chiller Hei-CHILL 3000, Hei-CHILL 5000
Operating manual must be read before initial start-up.
Please follow the safety instructions provided.
Please keep for future reference.

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About this Document	4
Version and variants	4
About this manual	4
Symbols and keywords.....	5
Safety Instructions	6
General safety instructions.....	6
EU Declaration of Conformity.....	6
Intended use.....	6
Installation	7
Qualifications of personnel	7
Operating company's obligations.....	8
Installation site	8
Modifications to the device.....	8
Safety of personnel.....	8
California Residents	9
Safety during use	9
Disposal	9
Device Description	10
Device overview	10
Device overview Hei CHILL 3000, Hei CHILL 5000	10
Hei CHILL control panel	10
Start-Up	11
Setting up the device	11
Connecting/disconnecting the power cord	12
Switching the device on and off.....	12
Selecting the language.....	13
Filling with thermal fluid	14
Menu navigation.....	16
Presettings	18
Selecting the language	18
Setting the temperature limit	18
Setting the temperature setpoint	19
Setting the date and time	20
Operation	21
Menu structure Overview.....	21
Temperature control operation.....	23
Adjusting the pump pressure.....	24
Special functions	24
Standby operation and autostart.....	24
Setting the operating mode for the cooling unit.....	25
Setting the setpoint offset.....	26
Adjusting the basic settings	27

Control parameters	29
Temperature control via programs	31
Preparing a program	31
Specifying program runs.....	32
Starting and ending a program	33
Device status	34
Interfaces.....	35
Connecting the alarm output	36
Restoring the factory settings.....	37
Cleaning and Maintenance	38
Cleaning.....	38
Cleaning the device.....	38
Maintenance	39
Maintenance	39
Troubleshooting	41
Malfunctions and troubleshooting.....	41
Assembly	46
Electrical connections	46
Peripheral connections	47
Overview of connections	47
Temperature control hoses	48
Connecting the chiller to an external consumer.....	49
Calibrating the temperature sensor	50
Disassembly and Storage	51
Dismantling, transportation and storage	51
Emptying the device	51
Transportation and storage	51
Accessories and Spare Parts	52
Scope of delivery.....	52
Accessories.....	52
Attachments	53
Technical data	53
Performance range.....	54
Thermal fluids	55
Cooling unit	56
China RoHS Declaration of Conformity	57
Service	59
Contact / Technical Service	59
Warranty	60
Confirmation of condition.....	61

About this Document

Copyright

This operating manual is protected by copyright and is exclusively intended for the internal use of the purchaser.

Handing over this operating manual to third parties, reproducing it in any type or form – in part or in full – and utilizing and/or communicating its contents is prohibited without prior written approval from the manufacturer.

Any violations will result in an obligation to provide compensation for damages. Other claims remain unaffected.

➤ Version and variants

Version

This manual describes the function, operation and maintenance of the circulation chillers High-CHILL 3000 and High-CHILL 5000.

Version	Alteration date
1.0	11/2019

Variants

The device is available in different variants. Certain features or functions are only available in specific product variants. These variants are described in this operating manual.

➤ About this manual

This operating manual is an integral part of the device described here.

- Carefully read through this operating manual and observe all of the safety instructions and warnings notices.
- Ensure that every user has carefully read the operating manual before using the device for the first time.
- Store this operating manual so that it is accessible to all users at all times.
- Pass the operating manual on to any subsequent owner.

The current version of this operating manual can also be found at:
<https://heidolph-instruments.com/en/service/downloads/operation-manuals#Chillers> (High-CHILL 3000 - 5000)

In addition, also observe the local regulations at the installation site.

➤ Symbols and keywords

Standardized symbols and signal words are used in this operating manual to warn against any dangers and provide important instructions. These instructions must be strictly observed to avoid accidents and damage.

The following terms and basic symbols are used:

Symbol	Additional signal words / explanation
<p>Warning signs</p> 	<p>The yellow triangle indicates hazardous situations. It is used in combination with the following signal words:</p> <p>DANGER: Indicates a hazardous situation which, if not avoided, will result in serious injury or death.</p> <p>WARNING: Indicates a hazardous situation which, if not avoided, may result in serious injury or death.</p> <p>CAUTION: Indicates a hazardous situation which, if not avoided, may result in property damage and minor or moderate injury.</p>
<p>Prohibitory signs</p> 	<p>Prohibited:</p> <p>The red circle indicates a situation that should be avoided under all circumstances and which, if not avoided, may result in serious injury or death.</p>
<p>Mandatory signs</p> 	<p>Must be observed:</p> <p>The blue circle indicates important information. Please observe this information to avoid any property damage.</p>

Other symbols used:

Symbol	Description
→	Handling instruction, action required
✓	Result of action
▪	List of information
▪	
▪	
a.	List of variants
b.	
c.	

➤ General safety instructions

EU Declaration of Conformity

This device complies with the following EC-Directives*:

- 2006/42/EC Machinery Directive
- 2014/30/EU Electromagnetic Compatibility Directive

*Also see the attachment "EU Declaration of Conformity".

The device has been manufactured according to state-of-the-art technology and in compliance with recognized safety regulations. However, risks may still arise during installation, operation and maintenance.

→ Please ensure that the operating manual is available at all times.

The device may only be used under the following conditions:

- Only operate the device if it is in a technically perfect condition.
- Only operate the device if you are properly aware of the risks and required safety measures.
- Only operate the device in accordance with the instructions given in this operating manual.
- If any information is ambiguous or missing, ask your superior or contact the manufacturer.
- Do not operate anything on the device without authorization.
- Only use the device in accordance with its intended use.

Intended use

The device is intended for use by trained and authorized personnel only.

The device is suitable for the following uses:	<ul style="list-style-type: none"> • Cooling of rotary evaporators: Hei-VAP Industrial • Tempering and delivering non-combustible heat transfer liquids in a closed circuit • Economical cooling as alternative to tap water systems
The device may be operated in research laboratories, other laboratories and production facilities in the following sectors:	<ul style="list-style-type: none"> • Chemistry • Pharmacy • Biology • Environmental analytics • Basic research • Similar research laboratories

Reasonable foreseeable misuse

Any use which deviates from the device's intended use is considered to be improper. The manufacturer is not liable for damage that occurs as a result. The risk is borne by the operator alone.

Amongst other things, the following uses are considered to be foreseeable misuse:

- Operation of the device without heat transfer liquid
- Incorrect connection of tubes
- Placement of the device on a tabletop surface
- Setting an incorrect pump pressure
- Operation of the device with a glass reactor without overpressure protection
- Chilling foodstuffs

Installation

Electrical safety

- The device may only be connected up if the available power supply voltage corresponds to the information stated on the rating plate for the device.
- The power supply connection must be easily accessible at all times.
- An RCD circuit breaker (residual current circuit breaker) must be fitted to the electrical system in the building so that the device is disconnected from the power circuit in the event of a fault.
- Repairs may only be performed by electricians authorised by Heidolph Instruments.
- Never operate the device with a damaged power cord.
- Always turn the device OFF and disconnect the power cord before carrying out any maintenance or repair work.

Thermal fluids

- The device is exclusively designed for nonflammable thermal fluids in Class I according to DIN 12876-1.
- Each thermal fluid covers a specific temperature range. This temperature range must match the temperature range of your application. If these temperatures are exceeded or undercut, serious hazards to health or environment can be caused (also see chapter "Attachment, Technical data", section "Thermal fluids").
- We recommend to use the Heidolph thermal fluid. This has been tested and approved by Heidolph Instruments.
- The safety data sheet must be consulted for the intended use of the device.

Tubes

- Besides material durability please note, that permissible temperature range and maximum pressure must match with your application when selecting tubes for connecting your periphery devices.
- We recommend to use the Heidolph tubes. These have been tested and approved by Heidolph Instruments.

Qualifications of personnel

- The device may only be operated by trained personnel.
- The device may only be operated by persons who have been instructed and supervised in its proper use by trained specialist personnel.
- The device may only be operated by specialist personnel who are above the legal minimum age.
- Other persons may only work on the device under the constant supervision of experienced and trained specialist personnel.
- This operating manual must be read and understood by all persons working with the device.
- Personnel must receive safety training that ensures responsible and safe working practices.

➤ Operating company's obligations

Installation site

- The device must be installed in a suitable location.
- The device must be installed on a firm and stable surface.
- Ensure that the device and all of its components and connections are easily accessible at all times.
- All vents of the device must have a minimum distance of 50 cm (20 inches) to walls or vertical surfaces, so air flow is not restricted.
- It is not permitted to place or store any objects such as accessories, tools or chemicals within this clearance area during operation.
- All screw connections must be securely tightened.
- It is not permitted to operate the device near to highly flammable or explosive substances.
- Only operate the device (depending on the substances being used) in conjunction with an extractor hood (see DIN EN 14175 and DIN 12924).
- The device is designed for indoor use only and under the following ambient conditions:

Ambient temperature	5 - 31 °C at 80 % relative humidity 32 - 40 °C decreasing linearly to a maximum 50 % relative humidity
Installation altitude	0 - 2,000 m above sea level
Contamination level	2
Overvoltage category	II
Permissible supply deviations	± 10 %

- If the device is operated in corrosive atmospheres, the service life of the device will decrease based on the concentration, duration and frequency of the exposure to the corrosive atmosphere e.g. concentrated hydrochloric acid (HCl).

Modifications to the device

- It is not permitted to make any unauthorised modifications or changes to the device.
- Do not attach or install any parts that have not been approved by the manufacturer.
- Unauthorized modifications or changes will void the EC Declaration of Conformity for the device and operation of the device will no longer be permitted.
- The manufacturer is not liable for any damage, dangers or injuries that result from unauthorized modifications and changes or due to the non-observance of the instructions in this manual.

Safety of personnel

- ➔ Ensure that the device is only operated by qualified specialist personnel and trained employees.
- ➔ Observe the following instructions to avoid any personal injuries and property damage:
 - Laboratory regulations
 - Accident prevention regulations
 - Hazardous Substance Act
 - Other generally accepted occupational health and safety regulations
 - Local regulations

California Residents

→ Important information for California residents regarding Prop 65. Please visit www.P65Warnings.ca.gov for more information.

➤ Safety during use

- Wear the appropriate clothing when working on the device (protective clothing, protective glasses and, if necessary, safety gloves).
- Do not use the device in potentially explosive areas. The device is not protected against explosion. There is no explosion or ATEX protection available.
- Do not carry out work with naked flames in the vicinity of the device (risk of explosion).
- Do not operate any devices in the vicinity which may exhibit electromagnetic fields in the frequency range 9×10^3 Hz to 3×10^{11} Hz.
- Do not operate or assemble devices in the vicinity which are emission or radiation sources (electromagnetic radiation) for the frequency range 3×10^{11} Hz to 3×10^{15} Hz (in the optical spectral range wavelengths from 1.000μ to $0,1\mu$).
- Do not operate or assemble appliances in the vicinity of the device which constitute emission or radiation sources for ionizing radiation or in the ultrasonic range.
- Do not operate the device where adiabatic compression or shock waves might occur (shock wave combustion).
- Attach all tube connections securely to ensure safe operation.
- Connect all cables and tubes without kinks and locate them outside the operating and danger zone.
- Do not spill liquids over the device or any parts of it.
- Remove any accidentally spilled liquids immediately.
- Eliminate errors immediately.
- Avoid putting pressure on the display when you are not operating the device.
- Do not use abrasive material to clean the surface. Only wipe with damp cloth.
- Always switch the device OFF after use.

➤ Disposal

Refrigerant

- The device may only be disposed of by a qualified refrigeration technician.
- Dispose of the refrigerant in accordance with EC Directive 215/2067 in combination with 517/2014/EC. (Type and quantity of the refrigerant can be seen on the rating plate; also see chapter „Attachement, Technical data“, section „Refrigerant unit“.)

Device and packaging

- Check the device components for hazardous substances and solvents.
- Clean all components before disposal.
- Dispose of the device in accordance with the relevant national regulations.
- Dispose of the packaging material in accordance with the appropriate national regulations.

Device overview

Device overview Hei CHILL 3000, Hei CHILL 5000

Hei CHILL control panel

➤ Setting up the device

Danger: Risk of poisoning!

Toxic vapours may be emitted depending on the thermal fluid being used and the type of operation.

Ensure that these vapours are extracted adequately.

Warning: Danger crushing!

The chillers are equipped with castors for easier mobility. This allows the device to roll away and possibly tip over.

The devices may only be set up on the floor. Set up the device on a level, non-slip surface with sufficient load-bearing capacity! Secure the device against rolling away by applying the locking brake!

Caution: Possible damage to the cooling system

The liquid from the compressor may leak into the circulation system for the thermal fluid if the device is not set up in an upright position.

The device must not be tilted. Only set up, move and operate the device in an upright position.

Caution: Device damage

The device may overheat in the event of insufficient airflow.

Maintain a minimum distance of 50 cm between the ventilation grille and walls or other vertical surfaces.

- The chiller can be set up lower than the apparatus to be cooled. When refilling the thermal fluid, the pump must be switched on in this case, and all hose connections must be checked for leaks, otherwise the thermal fluid can leak or overflow.
- You can position multiple devices next to one another.

- ➔ Set up the device on a solid, smooth, level and temperature-resistant surface.
- ➔ Release the locking brake on the castors by pressing the lever upward.
- ➔ Lock the castors of the device by pressing the lever downwards.

Locking brake released

Locking brake applied

- ➔ Do not place heavy objects on the device.
- ➔ Clean the installation surface and the castors on the device at regular intervals.

➤ Connecting/disconnecting the power cord

The power cord is permanently connected to the device.

The connection is located on the rear of the device.

- ✓ The main switch is switched OFF ("power off" position).
- ➔ Connect the device connector to your earthed socket.
- ✓ The connector is securely inserted when it is pushed in to the stop.
- ➔ Disconnect the device connector when the device is not operated.

Hei-CHILL 3000

Hei-CHILL 5000

➤ Switching the device on and off

Hei-CHILL 3000


Switching the device on

- ➔ Press the ON/OFF switch upward to "I".
- ✓ The device is switched on.

Switching the device off

- ➔ Press the ON/OFF switch downward to "O".
- ✓ The device is switched off.

Hei-CHILL 3000

Hei-CHILL 5000


Switching the device on

- ➔ Rotate the ON/OFF switch clockwise to "I".
- ✓ The device is switched on.

Switching the device off

- ➔ Rotate the ON/OFF switch anticlockwise to "O".
- ✓ The device is switched off.


Hei-CHILL 5000


Selecting the language


- Switch the device ON.
- ✓ A signal tone sounds.
- ✓ The "Fill mode" submenu is automatically called up when the device is switched on for the first time.
- ✓ If the "Fill mode" submenu is not displayed in a language that you can understand, the correct language must be selected first.
 - Press the "ESC" softkey.
 - ✓ The base menu will be displayed.
 - Press the "Select" button.

- ✓ The main menu will be displayed.
 - Press the down ▼ arrow button to select the "Setup" menu.
 - Confirm the selection using the "Select" button.
- ✓ The "Setup" menu will be displayed.
 - Press the down ▼ arrow button to select the "Basic Setup" menu.

- Confirm the selection using the "Select" button.
- ✓ The "Basic Setup" menu will be displayed.
 - Press the down ▼ arrow button to select the "Language" menu.

- ✓ The Language menu will be displayed.
 - Press the up ▲ arrow button ▼ to select the desired language.
 - Press on the "Select" button to confirm the selected language.
 - Press the "ESC" softkey.
 - ✓ You return to the base menu.

"Language" Menu

➤ Filling with thermal fluid

Danger:
Fire hazard

If you fill the device with a thermal fluid with too low a flash point, it may catch fire.

Only use thermal fluids with a flash point that is at least 20 K above the temperature range for the application. Do not use ethanol as an additive due to its highly flammable properties.

Warning:
Risk of electric shock

If the device is overfilled, thermal fluid may overflow into the device and cause the device to short circuit.

Ensure that the device is not overfilled. Carefully monitor the fill level indicator and the thermal expansion of the thermal fluid.

Warning:
Risk of electric shock

If the device is filled carelessly, thermal fluid may splash into the device and cause the device to short circuit.

Use a funnel for filling. Avoid splashing the thermal fluid.

Warning:
Risk of poisoning! Risk of injury!

When handling chemically prepared thermal fluids, any splashes could cause injuries to areas of bare skin.

Always wear safety goggles, protective gloves and protective clothing when filling the chiller.

The thermal fluid from Heidolph is suitable for all standard applications. It covers a temperature range from minus 30 °C to plus 90 °C. (See chapter "Spare parts and accessories, Accessories")

Start-Up


Filling for the first time

The reservoir for the thermal fluid can hold between 20 and 33 litres (see chapter "Attachments, Technical data").

Filling of the device is supported by a software program. If the fill level is too low:


- the pump and cooling unit will not start
- the Fill mode starts automatically.

→ Close the drain valve by turning the valve Emptying cock closed fully to the left.


- ✓ The device is switched on.
- ✓ The " Fill mode" submenu appears on the display.
- Carefully remove the cover for the filling pipe by pulling it upwards.
- Carefully fill the reservoir with thermal fluid via the filling pipe using a funnel.

Filling pipe for reservoir

Sealing cap for reservoir

- ✓ The display shows the value "0" as the fill level to the right of "started".
- Observe the fill level indicator.
- An acoustic signal sounds at intervals once the value has risen to 4-5.
- ✓ The device is filled sufficiently.

- Connect your external consumer (see chapter "Assembly, Peripheral connections").
- Press the "Standby" softkey to fill the external consumer.
- Continue to refill thermal fluid until the fill level rises to 4-5 again and an acoustic signal will sound.
- End the filling process and close the reservoir with the cover.
- Press the down ▼ arrow button to select "Exit filling".
- Press the "ESC" softkey to exit the "Fill mode" menu.
- ✓ You return to the base menu.

The intervals between the signal tones are shortened until a continuous tone is heard if the filling process does not end. A continuous tone indicates that the device is overfilled and cannot be started. Before starting the pump, some thermal fluid must be drained first


➤ Menu navigation

From the base menu, you reach the main menu, which is divided into 6 submenus. Some of these submenus have up to 4 further submenus.

(see the following chapter "Operation, Menu structure Overview")

Base menu


- ➔ Switch the device ON.
- ✓ A signal tone sounds.
- ✓ The base menu will be displayed after the welcome screen.

- ➔ Press the button "Select" .
- ✓ The main menu will appear.

Main menu

- ➔ Press the up ▲ or down ▼ arrow button.
- ✓ The various menu options are selected and highlighted in colour.
- ✓ An arrow ► after a menu option refers to a further submenu.
- ➔ Press the button "Select" to call up a selected submenu.
- ✓ Either another submenu or the set value for the called up parameter will be displayed.


Start-Up

Submenu for setting of values

- ✓ A cursor flashes under the currently set value.
- Press the right ► or left ◀ arrow button to select a number.
- Press the up ▲ or down ▼ arrow button to change the number.
- Press the right ► or left ◀ arrow button to select the next number.
- Confirm the entire value using the button "Select" .
- ✓ The set value will be stored.

- Press the left ◀ arrow button to return to the previous menu.
- Press the left softkey under "ESC" to return to the base menu.


If you press the "ESC" softkey without first confirming the value with the Enter key, the changes will not be stored.

Softkeys

The softkeys always call up the function displayed in the footer above.

- ✓ You are in a submenu.
- Press the left "ESC" softkey.
- ✓ You will be taken directly to the base menu from any submenu.
- Press the right "Standby" softkey.
- ✓ The device will be switched immediately to the standby or the temperature control mode from any menu.

With a few exceptions, the softkeys in each menu and submenu are assigned the "ESC" and "Standby" functions.

➤ Presettings

Selecting the language

(See chapter "Startup, Selecting the language")

Setting the temperature limit

Before placing the chiller into operation, the desired temperature setting needs to be entered depending on the thermal fluid being used.

Warning:
Risk of poisoning! Risk of injury!

Each thermal fluid covers a certain temperature range. If the temperature exceeds or drops below this temperature range, it could damage the device and cause the thermal fluid to escape. Escaping thermal fluids can lead to slipping and poisoning.

Check the permissible temperature range for the thermal fluid being used and carefully set the temperature limits on the device.


The desired cooling temperature is maintained by defining the following two variables:

- **Temperature limit:**
defines the upper and lower maximum values that must not be exceeded. As soon as the temperature exceeds or falls below these values, a warning signal will sound.
- **Temperature setpoint:**
defines the cooling temperature for the thermal fluid. Until this setpoint is reached, the current temperature will fluctuate above and below this setpoint and thus also between the upper and lower temperature limits.

Both values for the temperature limits are dependent on the thermal fluid being used (see also chapter "Spare parts and accessories, Spare parts").

Setting the upper temperature limit


- ✓ The device is switched on.
- ➔ Select the "Setup" menu option from the main menu.
- ➔ Press the button "Select".
- ➔ Select the menu option "Temp. limits" from the "Setup" submenu.


Start-Up

- ✓ The "Temperature limits" submenu appears.
- Press the up ▲ or down ▼ arrow button to select the lower limit "Til" or the upper limit "Tih".
- Confirm "Til" using the button "Select" .

- ✓ The "Lower limit" submenu appears.
- ✓ A cursor flashes under the set value.
- Correct the currently set value using the arrow buttons right ► or left ◀ and up ▲ or down ▼ .
- Confirm the entire desired value using the button "Select" .
- ✓ The set value will be stored.
- Return to the "Temperature limits" submenu using the left arrow button ◀.
- Press the left softkey under "ESC" to return to the base menu.

Setting the upper temperature limit

- Repeat the procedure for the upper limit value

Setting the temperature setpoint

The adjustable temperature range for the setpoint is automatically limited by the temperature limits. It can be set to a maximum of 2 °C below the upper temperature limit and 2 °C above the lower temperature limit.


At a temperature setpoint or current temperature < 5 °C, the yellow LED on the device will light up. It warns you about the possible incorrect use of the thermal fluid and resulting damage to the device.

Start-Up



- ✓ The device is switched on.
- Select the "Setpoint value" menu option from the main menu.
- ✓ A cursor flashes under the set value.
- Use the arrow buttons to enter the desired temperature setpoint.
- ✓ The value must be between the displayed "Max." and "Min." values.
- Confirm your selection using the button "Select" .

If the entered setpoint is outside of a defined temperature limit, a warning signal will sound. The value will not be saved.

Setting the date and time

- ✓ The device is switched on.
 - Select the submenu "Clock" / "Set Time/Date".
 - Use the arrow buttons to enter the current time and the current date and confirm the entry using the button "Select" .
 - Set the desired display for the date in the "Format/Date" submenu.
-
- ✓ Date and time will be displayed on the top right of the base menu in the header according to your setting.

➤ Menu structure Overview

Main menu

Setpoint Value

Setup (see following page)

Programmer	Program 1 - 5	Status	Start
			Stop
			Hold
			Continue
		Edit	
		Loop	

Interface

Graph	Mode	Online graph	
		Record	Start
		Freeze graph	Start
	Displayed Value	Tset Tint Text	
		Tset Tint	
		Tset Text	
		Tint Text	
		Tint	
		Text	
		Tset	
	Sample Time	2 s (max.2h10min)	
		10 s (max.11h5min)	
		30 s (max.33h20min)	
		1 min (max.66h40min)	
		2 min (max.133h)	
	Time axis	Automatic	
		9 min	
		45 min	
		2 h 15 min	
	Time base	absolute	
		relative	
	Temperature scale	automatic	
		manual	
	Temperature limits	Temp.Scale Min -13 °C	
		Temp.Scale Max 40 °C	

Clock

Set time and date	
Format of date	DD.MM.YYYY
	MM.DD.YYYY

Standby

Operation

Sub menu setup

Control	Control variable	intern Pt 1000	
	Control parameter	intern Pt 1000	Xp
			Tn
			Tv manual/auto
			Tv
			Td
	Setpoint offset	Offset source	off
		Diff.set/actual value ² 0.0	
	Correction limit ¹ 50.0		
Cooling autom.	off		
	on		
	automatic		
Limits	Lo. limit (Til)		
	Up. limit ((Tih)		
Fill mode	Started		
	Stop fill mode		
<u>Basic setup (see hereafter)</u>			
Calibration	Calibration		
	Factory Calibration	no/yes	
Factory setting	All default	no/yes	
	Control	Reset all	no/yes
		ctrl. parameter intern	no/yes
		internal Pt1000	no/yes
		miscellaneous	no/yes
Device status	Fehlerspeicher		
	Gerätedaten		
	SW-Version		
	Typ		
	Seriennummer		

Sub menu Basic setup (submenu of setup menu)

Basic setup	Sounds	Alarm	loud
		Warn.	medium
		Error	low
			off
	Display	Brightness	automatic
			Stage 5 ... 1
			off
	Autostart	off	off
			on
		Curr consumption	14.0A
	Warn level niveau	1	
	Alarm output: Standby	Only alarm	
		Alarm and Standby	
Language	English		
	Deutsch		
	Francais		
	Español		
	Italiano		
		русский	

Temperature control operation


Caution:
Possible damage to the device

If the ventilation openings are not kept free, the motors for the pump and compressor could overheat.

Never cover the ventilation grille. Ensure that there is always a distance of 50 cm between the device and wall or other devices.

- ✓ All peripheral devices are securely connected, and the refrigerant can circulate freely.
- ✓ The device is filled with thermal fluid
- Switch the device ON.
- ✓ The display shows the base menu with the current cooling temperature "Tint".
- ✓ The fill level indicator for the thermal fluid is between 4 and 5.
- Press the "Standby" softkey.
- ✓ The pump for circulating the thermal fluid will start.
- ✓ In the header "Standby" will be replaced by the display of the active cooling.

→ Monitor the fill level indicator and refill the device with thermal fluid if necessary.

Caution:
Possible damage to the device

If you set the chiller to cool down to a minus temperature, it is possible that the thermal fluid will freeze.

If you want to cool the device down to minus temperatures, ensure that you are using an appropriate thermal fluid. (See also chapter "Spare parts and accessories, Spare parts")

If there is a large difference in temperature between the chiller and the external consumer, the temperature drop between the two devices could be so high that it will not be possible to maintain the cooling output.

In this case, increase or decrease the setpoint for the cooling temperature accordingly or ensure that there is a lower volume of items to be cooled in your external consumer.

Adjusting the pump pressure

There is an adjusting wheel on the rear of the unit. It allows you to adjust the pressure for the circulation of the thermal fluid according to your application.

The maximally attainable pressure is 2 bar.

You can read the adjusted pressure on the pressure gauge on the front of the device.

- Turn the adjusting wheel clockwise to adjust the pressure.
- ✓ The pump pressure will increase.
- Turn the adjusting wheel anticlockwise to adjust the pressure.
- ✓ The pump pressure will decrease.
- ✓ The pressure gauge displays the pressure in bar.

Adjusting wheel
Pressure adjustment

Pressure gauge

Special functions

Standby operation and autostart

In standby mode, the device is in an idle state when it is switched on. All displays can be called up in standby mode.

If the standby mode is deactivated or autostart is activated and the unit contains sufficient temperature control fluid, the unit starts to cool immediately when it is switched on and after a power interruption.

Standby mode is activated when the device is delivered.

Do not change this setting when using the device to cool the Hei-VAP Industrial rotary evaporator! The chiller is switched on and off via the rotary evaporator.

Warning:
Risk of explosion! Risk of frostbite! Risk of poisoning! Risk of injury! Risk of slipping!

If the device automatically resumes full operation after a power failure, loose screw connections or defective glass containers could result in the uncontrolled escape of liquids. The liquids that may escape as a result pose a risk of frostbite, poisoning and slippage. Burst glassware also poses a significant risk of injury and can result in exothermic reactions.


If you have activated the autostart function, do not leave the chiller and connected peripheral devices unsupervised under any circumstances. Check that the hose connections are secure at regular intervals. Only work with flawless glass containers.

Activating / deactivating the autostart

- ✓ The device is switched on.
- Select the submenu "Setup" / "Basic settings" / "Autostart".
- Use the arrow buttons to select "on" and confirm your selection using the button "Select".
- ✓ Autostart is activated.
- To switch to standby mode, select the "off" option from the "Autostart" menu.


Setting the operating mode for the cooling unit

The cooling unit can be controlled both automatically and manually.

The setting is "automatic" when the device is delivered.

- "off" deactivates the cooling unit until it is activated again
- "on" activates the cooling unit permanently until another option is selected
- "automatic" switches the cooling unit automatically on or off when needed.
- ✓ The device is switched on.
- Select the submenu "Setup" / "Cooling".
- Use the arrow buttons to select one of the presented options and confirm your selection using the button "Select".
- ✓ The changed operating mode is activated.

The flow line temperature for automatic control of the cooling unit is ± 2 K from the setpoint. Under certain ambient conditions, this may lead to longer downtimes of the cooling unit. This does not impair the cooling output.



Setting the setpoint offset

- ✓ The device is switched on.
- Select the submenu "Setup" / "Control" / "Setpoint offset".
- Use the arrow buttons to select one of the presented options and confirm your selection using the button "Select".

Difference setpoint / actual value

- Select the "Diff. set/actual value" option to determine the upper and lower limit for the setpoint offset.
- Confirm your changes using the button "Select".


➤ Adjusting the basic settings


Adjusting the volume of the signal tones

- ✓ The device is switched on.
- Select the submenu "Setup" / "Basic settings" / "Sounds".
- Use the arrow buttons to select one of the presented options and confirm your selection using the button "Select".

Adjusting the brightness of the display


- ✓ The device is switched on.
- Select the submenu "Setup" / "Basic settings" / "Display" / "Brightness".
- Use the arrow buttons to select one of the presented options and confirm your selection using the button "Select".

Restricting the power consumption

If your mains fuse is less than 16 A, the power consumption can be reduced gradually from 16 A to 8 A.

- ✓ The device is switched on.
- Select the submenu "Setup" / "Basic settings" / "Power consumption".
- Use the arrow buttons to adjust the power consumption and confirm your changes using the button "Select".


Operation

Specifying the minimum fill level warning stage

On delivery, the device is configured such that a low level warning is output at a fill level beginning at level 2. This lower warning level can be adjusted between 0 and 3.

- ✓ The device is switched on.
- Select the submenu "Setup" / "Basic settings" / "Warning level niveau".
- Use the arrow buttons to adjust the warning stage and confirm your changes using the button "Select" .

No warning will be output if you set the warning stage to "0". Cooling mode is switched off as soon as the lowest fill level is reached. A corresponding alarm message will be shown on the display.

➤ Control parameters

The control parameters for cooling the device are set to a thermal fluid with water when the device is delivered. Depending on the application and the thermal fluid used, it may become necessary to adjust the control parameters.

The quality of the control system deteriorates at increasingly lower temperatures due to the change in viscosity of the thermal fluid at low temperatures. This means that if you work with a normal cooling temperature of -20 °C, for example, the control parameters should be adjusted at an operating temperature of -10 °C.

Please keep this in mind when setting the control parameters.

Internal control parameters

The internal control parameters compare the internal setpoint temperature with the flow line temperature. This calculation results in the control variable used for cooling.



Parameter	Designation	Explanation	Unit
Xp	Proportional range	Specifies the temperature range in which the proportional component (P component) of the controller is 0 ... 100 % of the maximum control variable. For example, if the control deviation is 2 K with an Xp of 10 K set, the P component is 20 % of the control value. For a control deviation of 10 K and more, the P component is 100 % of the control variable.	K
Tn	Reset time	Is decisive for the integral component (I component) of the control variable. It specifies the interval during which an existing control deviation is integrated. The larger Tn is, the slower the control deviation is integrated. This makes the control slower. A smaller Tn makes the control more dynamic and finally leads to oscillations.	s
Tv	Rate time	This forms the differential component (D component) of the control variable. It influences the approach speed of the actual value to the setpoint and counteracts the P and I components. The higher the setting for the rate time Tv, the stronger the output signal is dampened. The rule of thumb is: $Tv = Tn \times 0,75$.	s.

Operation

Adjusting the internal control parameters

- ✓ The device is switched on.
- Select the submenu "Setup" / "Control" / "Control parameters" / "Internal Pt1000".
- Use the arrow buttons to select one of the presented options and confirm your selection using the button "Select".
- ✓ By selecting "Tv manual/auto" you can specify the option of automatic or manual control.
- ✓ When selecting any other option, you will be able to adjust the highlighted values.
- ✓ Confirm the changed value using the button "Select".

If the control parameters "Tv manual/auto" is set to automatic control, the subsequent parameters "Tv" and "Td" are derived using factors from Tn. They are blocked from changes and indicated by a lock symbol.

➤ Temperature control via programs


In the "Programmer" submenu, you can define and save the course of a temperature control over a defined period of time.

They can save up to 5 different program sequences. Up to 150 individual steps can be programmed in each program.

Each of the 5 programs can be repeated between 1 and 250 times or you can loop it in continuous operation.

Preparing a program

- ✓ The device is switched on.
- ➔ Select the "Programmer" submenu.
- ➔ Select one of the 5 programs "Program n".
- ➔ Select the "Edit" menu option.

- ➔ To edit, confirm the first line with the button "Select" .
- ➔ Using the arrow buttons, enter for each line the following values:
 - "Tend" the desired end temperature
 - "hh" the desired number of hours
 - "mm" the desired number of minutes
 - "Tolerance" the permissible deviation from "Tend".

No.	Tend	hh	mm	Tolerance
1	30.00	--	--	0.1
2	50.00	0	20	0.0
3	50.00	0	20	0.0
4	70.00	0	20	0.1
5	60.00	0	30	0.0
6	30.00	0	0	0.0

- ➔ To be able to edit it, confirm each new line using the button "Select" .
- ➔ Confirm your changes using the button "Select" .

No time "hh" and "mm" can be entered in the first program line. The temperature entered here will be approached as fast as possible. As soon as the specified temperature is reached, the next program line will be processed.

Changing program parameters

- ➔ Call up the desired program.
- ➔ Use the arrow buttons to select a program parameters to be changed.
- ➔ Confirm the parameter using the button "Select" and change the value.
- ➔ If necessary, change further parameters and confirm all changes using the button "Select" .
- ➔ Use the left arrow button to exit the program.

Operation

Inserting a program line

- Call up the desired program.
- Use the arrow buttons to select a program line to be changed.
- Confirm the line number using the button "Select" .
- ✓ A blank program line is inserted below the line.

Deleting a program line

- Call up the desired program.
- Use the arrow buttons to select a program line to be deleted.
- Confirm the line number with the "delete" softkey.
- ✓ The line is removed from the program.

Change current program

- ✓ You are in the base menu.
- ✓ Press the "Prog.x/y" softkey.
- Use the arrow buttons to select a program parameters to be changed.
- ✓ Entire program lines cannot be inserted or deleted during operation.
- Confirm the desired parameter using the button "Select" and change the value.
- If necessary, change further parameters and confirm all changes using the button "Select" .
- Exit the program using the "ESC" softkey.

The smaller the permissible tolerance is programmed, the more precisely the desired temperature is maintained. However, as a certain amount of time is required before a temperature is reached, unwanted delays may occur if the time specification is tight, since the next program step will only be processed after the desired temperature has been reached.


A greater tolerance should therefore be allowed for large temperature jumps and short time specifications.

The duration of a single program step can be a maximum of 999 hours and 59 minutes. Longer times would have to be distributed to several program steps.

Specifying program runs


- ✓ The device is switched on.
- Select the "Programmer" submenu.
- Select one of the 5 programs "Program n".
- Select the menu option "Loops".



Operation

- ✓ Using the value "0" repeats the program indefinitely.
- Press the left arrow button to enter a 2 or 3 digit number.
- Use the arrow buttons up/down to enter the desired number of cycles and confirm your change using the button "Select" .

Starting and ending a program

- ✓ The device is switched on.
- Select the "Programmer" / "Status" submenu.
- ✓ Only the "Start" option will be displayed if no program is active.
- ✓ If a program is activated, the "Pause" and "Stop" options will be displayed.
- ✓ If a program has been interrupted using "Pause," the "Next" and "Stop" options will be displayed.
- Press the button "Select" to confirm the respective desired option.


➤ Device status

You can call up various data stored in the device from the "Device status" menu. Your device version will be displayed directly under the "Type ..." menu option.

- ✓ The device is switched on.
- ➔ Select the submenu "Setup" / "Device Status".
- ➔ Use the arrow buttons to select one of the presented options and confirm your selection using the button "Select".
- ➔ Use the left arrow button to return to the "Device status" menu.


Error memory

(For an interpretation of the error memory, see chapter: "Troubleshooting, Malfunctions and troubleshooting".)

No.	Source	Code	Type	Date	Time
5	Control	29	Error	30.10.09	10:32
4	Safety	3	Alarm	30.10.09	10:32
3	Control	4	Warn	29.10.09	16:41
2	Safety	29	Error	28.10.09	17:02
1	Control	36	Error	28.10.09	08:04


Device data

Here, the actual and setpoint temperatures are displayed at the various temperature sensors.


SW Version

Following "Control" the software version will be displayed.

Serial number

Following "Control" the serial number for the device will be displayed.

› Interfaces

Warning:
Risk of electric shock

If the interface inputs and outputs are energised too high and are not sufficiently insulated, metallic parts, e.g., the housing, may be energised in the event of a fault.

Isolate low voltage inputs and outputs over 25 V AC or 60 V DC reliably according to DIN EN 61140, or by double or enhanced insulation according to DIN EN 60730-1 or DIN 60950-1.

Use only shielded connection cables. Connect the shielding to the connector.

Caution:
Possible damage to the cooling system

If the connections on the device are not covered, liquids could reach the live components. Ingress of liquids can lead to malfunctions of the interfaces and can lead to a short circuit.

The interface connections on the rear of the device must always remain covered.

In which fault situations a signal is output via the interface can be set via the display.

Connecting the alarm output

The alarm output can be used to output faults to an external device via a potential-free contact.

Warning: Risk of electric shock

If the voltages received by the interface inputs and outputs are too high and not sufficiently insulated, metal components such as e.g., the housing could become live in the event of a fault.

Isolate low voltage inputs and outputs over 25 V AC or 60 V DC reliably according to DIN EN 61140, or by double or enhanced insulation according to DIN EN 60730-1 or DIN 60950-1.


Use only shielded connection cables. Connect the shielding to the connector. Cover unused connectors with protective caps.

Alarm output

The alarm output is located on the front of the unit in the handle recess of the ventilation grille.

Configuring the alarm output

The device is configured in the factory to always output safety-relevant faults to the alarm output. In addition, you can configure the alarm output so that an alarm is also output in standby mode.


A maximum 30 V DC at 0.2A may be connected to the alarm output.

Idle state

If the device is switched off or there is a fault after switching it on, the alarm output is in an idle state:

- Pin 1 and 2 are open.
- Pin 3 and 2 are closed.

Good state


If the device is switched on and operating without any faults, the alarm output is in a good state:

- Pin 1 and 2 are closed.
- Pin 3 and 2 are open.

Operation


- ✓ The device is switched on.
- Select the submenu "Setup" / "Basic settings" / "Alarm output".
- Use the arrow buttons to select one of the presented options and confirm your selection using the button "Select".

Restoring the factory settings


All settings that you have changed in one of the submenus can be restored to the original factory settings under the "Factory settings" menu option.

- ✓ The device is switched on.
- Select the submenu "Setup" / "Factory setting".
- Use the arrow buttons to select one of the presented options and confirm your selection using the button "Select".


Resetting all settings

- ✓ The device is switched on.
- Select the submenu "Setup" / "Factory setting" / "all default".
- Use the arrow buttons to select one of the presented options and confirm your selection using the button "Select".

Resetting settings individually

- ✓ The device is switched on.
- Select the submenu "Setup" / "Factory setting" / "Control".
- Use the arrow buttons to select one of the presented options and confirm your selection using the button "Select".
- ✓ For each setting you will be given the opportunity to confirm with "yes" or "no".

- ✓ You can reset the setpoint and the current consumption under "miscellaneous".

➤ Cleaning

Cleaning the device

The housing and surfaces of the device can be cleaned using a damp cloth and, if required, a mild soap solution.

Warning:
Risk of electric shock

Liquids penetrating into the device could cause the device to become live if it is still connected to the power supply.

Before cleaning the device:

- Switch the device off.
- Disconnect the power cord plug.

Caution:
Device damage

Liquids penetrating into the device could damage the electronics.

Only wipe the surfaces with a damp cloth and not a wet cloth.

Caution:
Damage to surfaces

The surfaces could be damaged by improper cleaning.


Do not use the following under any circumstances:

- Chlorine bleach or cleaning agents containing chlorine
- Substances containing solvents (e.g., acetone)
- Ammonia
- Abrasive cleaning agents such as cleaning wool, scouring agents or cleaning agents with metallic particles

Cleaning the air filter

The ventilation grille in front of the air filter is held by two magnets in addition to the mechanical bracket.

- ✓ The device is switched off and disconnected from the power grid.
- ➔ Use both hands to grab the top grip recess of the ventilation grille.
- ➔ Pull the ventilation grille off the upper magnets.
- ➔ Pull the ventilation grille slightly upwards and pull it off the retaining lugs at the bottom of the device.
- ➔ Brush or vacuum the air filter and, if necessary, the ventilation grille.
- ➔ Slide the ventilation grille at the bottom with the slots over the two retaining lugs.
- ➔ Grab the grip recess and carefully push the ventilation grille back onto the device at the top.

➤ Maintenance

Maintenance

Danger!
Risk of poisoning! Risk for the environment!

Improper repairs to the internal cooling circuit may result in leaks and the refrigerant could escape in an uncontrolled manner. This could result in employees being poisoned and also in serious environmental damage.

Repair and maintenance work on the internal cooling circuit may only be carried out by specialist personnel authorised by Heidolph Instruments.

Warning:
Risk of electric shock Risk of injury!

You could come into contact with live or moving parts when carrying out repairs. This could result in electric shocks, collisions, cutting or crushing

The device must be disconnected from the power grid before all maintenance work.

Repairs may only be performed by specialist personnel authorised by Heidolph Instruments.

Warning!

Risk of burns! Risk of frostbite

Touching hot or cold device components, accessories and thermal fluid could result in burns and frostbite.

Always wear protective gloves and protective clothing when carrying out maintenance work.

Allow the device components, accessories and thermal fluid to reach room temperature before touching them.

If the device has come into contact with hazardous materials, the device needs to be decontaminated before completing any maintenance work.

Maintenance intervals

The maintenance intervals stated in the following table must be observed. The following maintenance tasks are obligatory before all prolonged unsupervised operation.

Interval	Maintenance work
monthly	<ul style="list-style-type: none">▪ Check the drain screw by inspecting it from the outside▪ Check the external hoses for material fatigue▪ Clean the air filter
Six monthly	Check the thermal fluid

Any repairs that may be required may only be performed by a specialist authorized by Heidolph Instruments.

In this case, contact your Heidolph Instruments dealer or a representative of Heidolph Instruments.

Checking the thermal fluid

Contaminated or degenerated thermal fluid must be replaced. Continuing to use the thermal fluid in this case is only permitted after obtaining corresponding test results. The thermal fluid must be checked in accordance with DIN 51529.

➤ Malfunctions and troubleshooting

Alarms are relevant to safety. The device will be switched to standby operation.

- ✓ An acoustic signal will sound rhythmically.
- ➔ After correcting the cause, press the button "Select" to reset the alarms.
- ✓ The alarm has been resolved.

Alarms

Code	Error	Troubleshooting
02	Low level	<ul style="list-style-type: none"> • Check if all connections are securely tightened. • Refill thermal fluid • Press button "Select"

Warnings are not relevant to safety. The device continues to operate.

- ✓ A prolonged acoustic signal will sound that is then repeated periodically.
- ➔ After correcting the cause, press the button "Select" to reset the warnings.
- ✓ The warning has been resolved.

1. Warnings control system

Code	Error	Troubleshooting
001	CAN-receive overflow	Please contact your local Heidolph Instruments distributor
002	Watchdog-Reset	Please contact your local Heidolph Instruments distributor
003	Til limit active	Please contact your local Heidolph Instruments distributor
004	Tih limit active	Please contact your local Heidolph Instruments distributor
005	Corrupt parameter	Please contact your local Heidolph Instruments distributor
006	Corrupt program	Check the parameters of your running program and change the settings
007	Invalid parameter	Please contact your local Heidolph Instruments distributor
008	CAN system	Please contact your local Heidolph Instruments distributor
009	Unknown module	Please contact your local Heidolph Instruments distributor
010	SW Control too old	Please contact your local Heidolph Instruments distributor
011	SW Safety too old	Please contact your local Heidolph Instruments distributor
012	SW Command too old	Please contact your local Heidolph Instruments distributor

Troubleshooting

Code	Error	Troubleshooting
013	SW Command too old	Please contact your local Heidolph Instruments distributor
014	SW Analog too old	Please contact your local Heidolph Instruments distributor
017	SW Valve 0 too old	Please contact your local Heidolph Instruments distributor
018	SW Valve 1 too old	Please contact your local Heidolph Instruments distributor
019	SW Valve 2 too old	Please contact your local Heidolph Instruments distributor
020	SW Valve 3 too old	Please contact your local Heidolph Instruments distributor
021	SW Valve 4 too old	Please contact your local Heidolph Instruments distributor
022	SW Pump 0 too old	Please contact your local Heidolph Instruments distributor
023	SW Pump 1 too old	Please contact your local Heidolph Instruments distributor
024	SW Pump 2 too old	Please contact your local Heidolph Instruments distributor
025	SW Pump 3 too old	Please contact your local Heidolph Instruments distributor
026	SW HTC too old	Please contact your local Heidolph Instruments distributor
028	SW Ethernet too old	Please contact your local Heidolph Instruments distributor
029	SW EtherCAT too old	Please contact your local Heidolph Instruments distributor
033	Clock wrong time	Please contact your local Heidolph Instruments distributor
034	Tset: Prog. is running	Please contact your local Heidolph Instruments distributor
041	Wrong mains voltage	Please contact your local Heidolph Instruments distributor
042	No VC type	Please contact your local Heidolph Instruments distributor
043	No VC voltage	Please contact your local Heidolph Instruments distributor
050	Level very low	Refill thermal fluid
051	Level high	Drain thermal fluid
055	CAN buff. overflow	Please contact your local Heidolph Instruments distributor

Troubleshooting

2. Warnings safety system

Code	Error	Troubleshooting
101	CAN receive overflow	Please contact your local Heidolph Instruments distributor
102	Watchdog-Reset	Please contact your local Heidolph Instruments distributor
103	Heating not correct	Please contact your local Heidolph Instruments distributor
104	Heat 1 failed	Please contact your local Heidolph Instruments distributor
105	Heat 2 failed	Please contact your local Heidolph Instruments distributor
106	Heat 3 failed	Please contact your local Heidolph Instruments distributor
107	Invalid parameter	Please contact your local Heidolph Instruments distributor
108	CAN system	Please contact your local Heidolph Instruments distributor
109	Unknown module	Please contact your local Heidolph Instruments distributor
110	SW Control too old	Please contact your local Heidolph Instruments distributor
111	SW Safety too old	Please contact your local Heidolph Instruments distributor
112	SW Command too old	Please contact your local Heidolph Instruments distributor
113	SW Cool too old	Please contact your local Heidolph Instruments distributor
114	SW Analog too old	Please contact your local Heidolph Instruments distributor
116	SW Contact too old	Please contact your local Heidolph Instruments distributor
117	SW Valve 0 too old	Please contact your local Heidolph Instruments distributor
118	SW Valve 1 too old	Please contact your local Heidolph Instruments distributor
119	SW Valve 2 too old	Please contact your local Heidolph Instruments distributor
120	SW Valve 3 too old	Please contact your local Heidolph Instruments distributor
121	SW Valve 4 too old	Please contact your local Heidolph Instruments distributor
122	SW Pump 0 too old	Please contact your local Heidolph Instruments distributor
123	SW Pump 1 too old	Please contact your local Heidolph Instruments distributor
124	SW Pump 2 too old	Please contact your local Heidolph Instruments distributor

Troubleshooting

Code	Error	Troubleshooting
125	SW Pump 3 too old	Please contact your local Heidolph Instruments distributor
126	SW HTC old	Please contact your local Heidolph Instruments distributor
128	SW Ethernet old	Please contact your local Heidolph Instruments distributor
129	SW EtherCAT old	Please contact your local Heidolph Instruments distributor
155	CAN buff. overflow	Please contact your local Heidolph Instruments distributor

3. Warnings refrigeration unit (when set on "automatic")

Code	Error	Troubleshooting
301	CAN receive overflow	Please contact your local Heidolph Instruments distributor
302	Watchdog-Reset	Please contact your local Heidolph Instruments distributor
303	Missing SM adaptation	Please contact your local Heidolph Instruments distributor
304	Pressure switch activated	Please contact your local Heidolph Instruments distributor
305	Clean condenser	Clean air filter (see chapter "Cleaning & Maintenance", section "Cleaning the air filter".)
306	TO1 out of range (Klixon)	Please contact your local Heidolph Instruments distributor
307	Invalid parameter	Please contact your local Heidolph Instruments distributor
308	CAN system	Please contact your local Heidolph Instruments distributor
309	Unknown module	Please contact your local Heidolph Instruments distributor
310	SW Control too old	Please contact your local Heidolph Instruments distributor
311	SW Safety too old	Please contact your local Heidolph Instruments distributor
312	SW Command too old	Please contact your local Heidolph Instruments distributor
313	SW Cool too old	Please contact your local Heidolph Instruments distributor
314	SW Analogue too old	Please contact your local Heidolph Instruments distributor
316	SW Contact too old	Please contact your local Heidolph Instruments distributor
317	SW Valve 0 too old	Please contact your local Heidolph Instruments distributor

Troubleshooting

Code	Error	Troubleshooting
318	SW Valve 1 too old	Please contact your local Heidolph Instruments distributor
319	SW Valve 2 too old	Please contact your local Heidolph Instruments distributor
320	SW Valve 3 too old	Please contact your local Heidolph Instruments distributor
321	SW Valve 4 too old	Please contact your local Heidolph Instruments distributor
322	SW Pump 0 too old	Please contact your local Heidolph Instruments distributor
323	SW Pump 1 too old	Please contact your local Heidolph Instruments distributor
324	SW Pump 2 too old	Please contact your local Heidolph Instruments distributor
325	SW Pump 3 too old	Please contact your local Heidolph Instruments distributor
326	SW HTC old	Please contact your local Heidolph Instruments distributor
328	SW Ethernet old	Please contact your local Heidolph Instruments distributor
329	SW EtherCAT old	Please contact your local Heidolph Instruments distributor
341	sm0 min too small	Please contact your local Heidolph Instruments distributor
344	Chiller missing	Please contact your local Heidolph Instruments distributor
345	Valve not closed	Please contact your local Heidolph Instruments distributor
347	Configure EEV0	Please contact your local Heidolph Instruments distributor
348	Configure EEV1	Please contact your local Heidolph Instruments distributor
349	Preheat unit	Please contact your local Heidolph Instruments distributor
355	CAN buff. overflow	Please contact your local Heidolph Instruments distributor

Faults

If the fault occurs again after restarting the device or you receive the error message "EXXX", note down the message shown on the display and contact Heidolph Instruments or your local Heidolph Instruments distributor.

Electrical connections

Repairs to the device must **only be carried out by a qualified electrician approved by Heidolph Instruments**. Improper repairs could result in serious dangers. Please contact your specialist dealer if the device needs to be repaired.

Warning: Risk of electric shock

The device must be connected to an earthed power socket.

The rated voltage for the device and the supply voltage must be the same. The rated voltage for the device can be found on the rating plate on the rear of the device.

The device is supplied with a grounded power cord. In case the original plug is replaced, the new plug must have a protective conductor!

Colour code for electrical connections:

1-phase

GREEN/YEL-LOW	PE: Protective conductor (Earth)
----------------------	--

BLUE	N: Neutral conductor
-------------	----------------------

BROWN	P: Phase
--------------	----------

3-phases

GREEN/YEL-LOW	PE: Protective conductor (Earth)
----------------------	--

BLUE	N: Neutral conductor
-------------	----------------------

BROWN	P: Phase
--------------	----------

BLACK	P: Phase
--------------	----------

GREY	P: Phase
-------------	----------

The device is connected to the mains with the power cord supplied. The power connection is located on the rear of the housing.

For countries where a different plug than the standard plug supplied is used:

- The plug supplied may only be changed by a professional electrician.
- If you wish to use an adapter, ensure that it is grounded and approved in accordance with the local regulations.
- Wiring, adapter and plug must at least meet the electrical power stated on the rating plate.
- Do not use extension cables.

Peripheral connections

Overview of connections

Hei-CHILL 3000, Hei-CHILL 5000: Rear of device

Temperature control hoses

**Warning:
Risk of frostbite!**

If you use unsuitable hoses at very low cooling temperatures, the thermal fluid may escape, and you could experience frostbite if you come into contact with it.

Only use hoses with a temperature resistance that corresponds to the operating temperature range of the device.

Warning: Risk of frostbite!

If you come into contact with cold hoses, you could experience frostbite if you are working with very low cooling temperatures.

Only use insulated hoses at operating temperatures $< 0\text{ }^{\circ}\text{C}$.

All hose connections are located on the top rear of the device.

- Pump flow line
- Pump return

All hose connections are fitted with a hose connector suitable for hoses with an internal diameter of $\varnothing 8\text{ mm}$.

Pump connection
return line

Pump connection
flow line

Carefully select hoses and connections for the chiller depending on the fluid, temperature and pressure of your application. We recommend to use hoses with textile reinforcement (also see chapter "Scope of delivery and accessories, Accessories").

The user is responsible for any damage that results from improper connections and the selection of unsuitable hoses.

Connecting the chiller to an external consumer

Warning:

Risk of electric shock Risk of frostbite!

If the hoses between the chiller and the external consumer are not leak-tight and connected without kinks, any escaping thermal fluid could result in frostbite and electric shock in the event of contact with live components.

- Lay all hoses without kinks and protected against kinks.
 - Check the leak-tightness of all hose connections and hoses.
 - Secure all hoses against slippage.
 - Maintain a clearance distance between hoses and power connections or live components.
-

Warning: Risk of explosion! Risk of frostbite! Risk of poisoning! Risk of injury! Risk of slipping!

If the hoses are improperly connected or glassware is exposed to high pressure, the hoses and glassware could burst due to overpressure. The liquids that may then escape as a result pose a risk of frostbite, poisoning and slippage. Burst glassware also poses a significant risk of injury and can result in exothermic reactions.

The maximum pressure at the outlet from the chiller is stated in the Attachments (see chapter "Attachments, Technical data"). If necessary, reduce the pressure in accordance with the connected peripheral devices by fitting a pressure regulator upstream at the outlet from the chiller. Ensure that all hose connections are leak-tight and secured against kinks and slippage to avoid overpressure.

If air penetrates into the external consumer and it is higher than the chiller, it could cause the external circuit to empty into the chiller when the pump is switched off. This could cause the chiller to overflow.

- Shorten the temperature control hoses as much as possible while ensuring that they can be laid without any kinks and without any tension.
- Secure the temperature control hoses against slippage with the aid of hose clamps.

Calibrating the temperature sensor

The internal thermometer in the device is precisely calibrated upon delivery.

If you detect a possible temperature deviation when checking the thermal fluid using a reference thermometer, the measurement values on the temperature sensor can be adjusted.

For this purpose, the reference thermometer needs to be fitted to the flow line on the device.


The factory calibration is overwritten when adjusted. Only use a reference thermometer with a high level of precision for the calibration process.

- ✓ The chiller is switched off and disconnected from the power grid.
- ✓ Allow the device and thermal fluid to warm up or cool down to room temperature.
- Attach a T-piece matching the inside hose diameter (Ø 8 mm) between the connection for the pump flow line and your connection hose.
- Slide your reference thermometer into the still open end of the T-piece and attach the thermometer.
- ✓ The calibrated reference thermometer has been installed in the flow line in a process-capable manner.



Adjusting the temperature value.

- ✓ The device is switched on.
- Start the temperature control.
- Select the submenu "Setup" / "Calibration" / "Calibration".
- Use the arrow buttons to adjust the displayed value to the temperature of your reference thermometer and confirm your changes using the button "Select".

Resetting the temperature value.

- ✓ The device is switched on.
- Select the submenu "Setup" / "Calibration" / "Factory calibration".
- Use the arrow buttons to select an option and confirm your changes using the button "Select".
- ✓ Selecting the "no" option returns you to the previous menu without changes.
- ✓ Selecting the "yes" option resets the calibration to the factory setting.

➤ Dismantling, transportation and storage

Emptying the device

Warning! Risk of frostbite

You could experience frostbite if you come in contact with the thermal fluid while emptying it.

Switch off the device at an early stage before emptying it so that the thermal fluid can warm up to room temperature. Always wear protective gloves and protective clothing when emptying the device.

- Switch the device OFF.
- Allow the device and the thermal fluid to cool down or warm up to room temperature.
- Position a container with sufficient capacity directly underneath the drain screw.

The thermal fluid will immediately flow out of the device when the drain screw is opened.

- Turn the drain screw anticlockwise to open it.
- Observe the regulations for the disposal of used thermal fluid (also see chapter "Safety, Disposal").

Emptying the refrigerant

Danger!: Risk of injury! Risk of poisoning! Environmental hazard!

The container of the cooling unit is under considerable pressure (19-28 bar) and contains environmentally harmful substances. Improper emptying of the refrigerant could result in a serious risk of injury or poisoning as well as environmental hazards.

The refrigerant in the internal cooling circuit may only be emptied by a qualified refrigeration technician.

Transportation and storage

- Only store and transport the device and its components after emptying, cleaning and dismantling it down to the basic device.
- When lifting the device, wear protective gloves and grip underneath the device.
- Store and transport the device and its components in the original packaging or in another suitable container to prevent any damage.
- Carefully seal the packaging against any unauthorized or unintentional opening.
- Store the device in a dry and frost-free location.

Improper handling during storage and transport may result in damage to the device and the mechanical components of the device.

Avoid jolts and vibrations when transporting the device.

➤ Scope of delivery

Component	Quantity	Order number
Hei-CHILL3000*	1	591-00700-00
Hei-CHILL 5000*	1	591-00700-01
Operating manual English /German	1	01-005-005-84
Warranty registration / confirmation of condition	1	01-006-002-78
EU Declaration of Conformity	1	01-001-025-17

* The product numbers are valid for 230 V devices in the European Community. Please contact Heidolph Instruments for further product numbers on country specific varieties.

➤ Accessories

Component	Temperatur range °C	Quantity	Product number
Tube with textile reinforcement, by the meter, not insulated, maximum pump pressure < 16 bar	- 20 to 60	-	591-38000-00
Thermal fluid (10 L)*	- 30 to 90	1	569-00611-00
Thermal fluid (20 L)*	- 30 to 90	1	569-00610-00

* The safety data sheet for the thermal fluid can be found at: <https://heidolph-instruments.com/en/service/downloads/operation-manuals#Rotary-Evaporators/Chillers>

Additional accessories can be found in our general catalogue or at: www.heidolph-instruments.com

➤ Technical data

Overvoltage category	II
Protection class (IEC 61140)	1 Ⓢ
Protection class (IEC 60529)	IP 32
Sound pressure level (dB(A)) (in accordance with IEC 61010)	< 70
Alarm output connection	max. 30 V DC, 0.2 A
Contamination level	2
Distance to surroundings all around (cm)	50
Classification according to DIN 12876-1 (class/designation)	I/NFL
Display	3.5" TFT display (320 x 240 pixels)
Cooling	Air-cooled
Refrigerant	R-449A
Refrigerant fill weight max. (kg)	1,1
Thermal fluid fill volume (max./ min. (L))	33/20
Max. flow rate (L at 1 bar, H ₂ O 20°C)	20
Max. flow pressure (bar, H ₂ O 20°C).	2
Operating temperature range (°C)	-20 to 40
Temperature stability (K)	±0.05
Storage temperature range (°C)	-20 to 50
Dimensions (W x D x H) (mm)	550 x 650 x 970

Model	Hei-CHILL 3000	Hei-CHILL 5000
Standard supply voltage	1 x 230 V (50 Hz) or 1 x 208 - 220 V (60 Hz)	1 x 400 V 3/N/PE (50 Hz) or 1 x 208 - 220 V 3/PE (60 Hz)
Cooling output (W)	3000	5000
GWP _(100a) *	1397	1397
CO ₂ equivalent (t)	1.3	1.5
Rated power	230 V 1800 W	400 V 3300 W
	208 - 220 V 2300 W	208 - 220 V 3600 W
Exhaust air m ³ /h	1300	2500
Weight (kg)	93	98


* Time horizon 100 years - according to IPCC IV

Performance range

Hei-CHILL 3000

Hei-CHILL 5000

Thermal fluids

Each individual thermal fluid covers a recommended temperature range and must be suitable for the temperature range for your application.

→ Never use contaminated or degenerated thermal fluid.

The thermal fluid becomes more viscous at the lower limit of the temperature range and influences the temperature stability, pump output and cooling output. Especially in water, ice may form that could destroy the device. Vapours and odours may develop at the upper limit of the range. Therefore, utilise the full temperature range only when required.

Permissible recommended thermal fluids

Heidolph designation *	Chemical designation	Temperature range in °C	Viscosity (kin) in mm ² /s (at 20 °C)	Viscosity (kin) in mm ² /s at temperature
Thermal fluid	Monoethylene glycol / water	-30 to 90	4	50 at -25 °C

* Safety data sheets for the thermal fluid can be requested at any time if required.

In the case of the Heidolph thermal fluid, the water content will fall after prolonged operation at high temperatures and the mixture will become flammable (flash point 128 °C). Use a hydrometer to check the mixing ratio.

Caution: Damage to the device, damage to the environment

If water is used as the thermal fluid, incorrect additives could damage the device and also the environment if the additives escape from the device.

Do not use the following thermal fluids or additives under any circumstances:

- Distilled, deionised, fully demineralised water or sea water (due to their corrosive properties)
- Water containing iron (due to rusting)
- Untreated river water (due to the formation of algae)
- Inorganic salts (due to flocculation)
- Chlorine or ammonia as additives (e.g., in cleaning agents and disinfectants)

The following water is permissible as a thermal fluid:

- Water free of impurities
- Water with a pH value between 6.0 and 8.5
- Ultra pure water and distilled water after the addition of 0.1 g sodium carbonate (Na₂CO₃) per litre of water

The proportion of alkaline earth ions in the water must lie between 0.71 mmol/L and 1.42 mmol/L (corresponds to 4.0 and 8.0 °dH). Hard water will cause limescale to develop in the device.

Cooling unit

Danger!: Risk of injury! Risk of poisoning! Environmental hazard!

The container of the cooling unit is under considerable pressure (19-28 bar) and contains environmentally harmful substances. Improper emptying of the refrigerant could result in a serious risk of injury or poisoning as well as environmental hazards.

The refrigerant in the internal cooling circuit may only be emptied by a qualified refrigeration technician.

The compressor in the cooling unit is hermetically sealed and equipped with overload protection that responds to the temperature and power consumption of the compressor.

The cooling unit contains refrigerant. Observe the special regulations for its disposal (see also chapter "Safety, Disposal").

Cooling output in kW at ambient temperatures in °C

Cooling output kW	Hei-CHILL 3000	Hei-CHILL 5000
at 20 °C	3,00	5,00
at 10 °C	2,40	3,90
at 0 °C	1,68	2,75
at -10 °C	1,03	1,70
at -20 °C	0,60	1,00

* The cooling output is measured with the thermal fluid maintained at a certain temperature. The ambient temperature for the measurements is 20 °C.

Refrigerant and fill level

The cooling unit contains refrigerant. Observe the special regulations for its disposal (see chapter "Safety, Disposal").

The refrigerant R-449A contains fluorinated greenhouse gases. The global warming potential (GWP) is 1430 at a time horizon of 100 years, CO₂ equivalent = 1.0.

China RoHS DECLARATION OF CONFORMITY

Heidolph Instruments GmbH & Co.KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Instruments GmbH & Co.KG products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
 O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
 X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Instruments GmbH & Co.KG may enter into further devices or can be used together with other appliances.

With these products and appliances in particular, Heidolph Instruments GmbH & Co.KG will not take responsibility for the EFUP of those products and appliances.

Place, date Schwabach, 19.09. 2019

Wolfgang Jaenicke
Chief Executive Officer CEO

Marcell Sarré
Vice President
Quality Management &
Technical Service

➤ Contact / Technical Service

Questions / Repair work

If any aspect of installation, operation or maintenance remains unanswered in the present manual, please get in touch with the following address.

For repairs please call Heidolph Instruments or your local authorized Heidolph Instruments distributor.

Warning! **Danger of poisoning!**

Contaminated devices can lead to severe injury or death of our employees!

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact substance
- take adequate protective action towards our parts receiving and service personnel
- mark the pack in accordance with Ordinance on Hazardous Substances

A "Confirmation of Condition" form can be found at the end of this operating manual.

→ Prior to shipping the device for repair, complete a copy of this form and submit it in advance.

Contact details:

Heidolph Instruments Germany

Heidolph Instruments GmbH & Co. KG
Technical Service
Walpersdorfer Str. 12
D-91126 Schwabach / Germany

Tel.: +49 – 9122 - 9920-74
Fax: +49 – 9122 - 9920-84

E-Mail: service@heidolph.de
www.heidolph.com

Heidolph Instruments North America

Phone: 1-866-650-9604
E-mail: service@heidolph.com
www.heidolphNA.com

Heidolph Instruments United Kingdom

Phone: 01799 - 5133-20
E-mail: service@radleys.co.uk
www.heidolph-instruments.co.uk

All other countries

You will find contact details of your local Heidolph distributor at
www.heidolph.com

➤ Warranty

Heidolph Instruments provides a 3 year warranty for the products described here (excluding glass and wearing parts) if you register using the warranty card enclosed or online (www.heidolph.com). The warranty period begins with the date of purchase. The serial number is also valid without registration.

This warranty covers material and manufacturing defects.

In the event of a material or manufacturing defect, the device shall be repaired or replaced free of charge under the terms of the warranty.

Heidolph Instruments shall not assume liability for any damage incurred as a result of improper handling or transport.

Warranty claim?

- ➔ Please notify Heidolph Instruments or your local Heidolph distributor if you wish to make a warranty claim.

➤ Confirmation of condition

→ In the case of repair, copy and complete the Confirmation of condition form and send it to Heidolph Instruments.

1. Details about the device

Product number _____

Serial number _____

Reason for repair _____

2. Has the device been cleaned, decontaminated/sterilized?

Yes _____ No _____

3. Is the device in a condition which does not represent any health threats for the staff of our service department?

Yes _____ No _____

If not, which substances has the device come into contact with?

4. Legally binding declaration

The customer is aware of being legally liable to Heidolph Instruments for any damages arising from incomplete and incorrect information.

Date

Signature

Company stamp

Please note

The shipper is responsible for the return of the goods in well packed condition, suitable for the mode of transport.

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This is a translation of the original operating manual. We reserve the right to make any technical modifications without prior notification. This documentation is not subject to revision; the latest version of this documentation can be found on our homepage. The copyright for all texts and images is held by Heidolph Instruments GmbH & Co. KG.
© Heidolph Instruments GmbH & Co. KG