RotaChill Large

Chiller

Leading Safety Standards

Superior Ease of Use

Reduced Cost of Ownership

Operating Manual

Page 2

Chiller RotaChill Large

Operating manual must be read before initial start-up.
Please follow the safety instructions provided.
Please keep for future reference.

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About this document	4
Version and variants	4
About this manual	4
Terms and signs used.....	5
Safety Instructions	6
General safety instructions.....	6
EC Declaration of Conformity.....	6
Intended use.....	6
Installation / Electrical safety	6
Personnel qualification.....	7
Operating company's obligations.....	7
Installation site	7
Changes to the unit	8
Safety for the personnel	8
Safety during use	8
Disposal	8
Device Description	9
Device overview	9
RotaChill Large overall view.....	9
RotaChill Large control panel	9
Start Up	10
Set up.....	10
Power cord	10
Main switch.....	11
Process Coolant.....	11
Filling the Reservoir	12
Switch unit on and off	12
Starting process fluid flow.....	13
Operation	14
Normal operation.....	14
Selecting the Temperature Unit (°C or °F)	14
Displaying and Adjusting the Set Point.....	14
Setting Operational Parameters	15
Operational Parameters	15
Enabling/Disabling the Local Lockout (LLO)	17

Cleaning and Maintenance	18
Cleaning.....	18
Maintenance	18
Troubleshooting	20
Possible errors and how to resolve	20
Error codes	20
Other problems	21
Assembly	23
Electrical connections	23
Set up chiller	24
Hoses and connectors	24
Closed System or Cooling Coil Setup	25
Disassembly and Storage	26
Dismantling, transportation and storage	26
Dismantling	26
Transportation and storage	26
Accessories and Spare Parts	27
Scope of delivery	27
Attachments	28
Technical data	28
Performance range.....	29
Contact / Technical Service	30
Warranty	31
Confirmation of condition of unit.....	32

➤ Version and variants

Version

This manual describes the function, operation and maintenance of the compact chiller RotaChill Large.

Version	Alteration date
0.1	11/2017

Variants

The units are available in different configurations. Certain features and functions apply to specific configurations only. The varieties are documented in this manual.

➤ About this manual

This operating manual is a component of the device described here.

- Please read this manual carefully and obey all safety and warning notices.
- Ensure that every operator reads this manual.
- Ensure that this manual is accessible for every operator.
- Pass on the operating manual to the subsequent owner.

For the current version of this operating manual in pdf format, please go to:

<http://www.heidolph-instruments.com/Support/Operation Manuals/Rotary Evaporators/Chillers/RotaChill Large>

In addition please observe the regional regulations.

➤ Terms and signs used

In this manual you will find standardized terms and signs to warn you of possible dangers or give important hints. Please take special notice to these signs to avoid personal injury or damage to the appliance.

The following terms and type of signs (basic drawing elements) are used:

Symbol	Additional term / Description
Warning signs 	The yellow triangle indicates hazardous situations. It is used in combination with the following terms: DANGER: Indicates a hazardous situation which, if not avoided, will result in irreversible injury or death. WARNING: Indicates a hazardous situation which, if not avoided, could result in serious injury or material damage. CAUTION: Indicates a hazardous situation which, if not avoided, may result in a minor injury or material damage.
Prohibitory signs 	Not to do: The red circle indicates a situation that should be avoided under all circumstances to avoid injury or damage.
Mandatory signs 	To do: The blue circle indicates important information. Please obey to avoid property damage.

Other signs used

Symbol	Description
→	Handling instruction, action required
✓	Result of action
▪ ▪ ▪	List of information
a. b. c.	List of variants

➤ General safety instructions

EC Declaration of Conformity

This device complies with the following EC-Directives*:

- 2006/42/EC Machinery Directive
- 2006/95/EC Low-Voltage Directive
- 2004/108/EC Electromagnetic Compatibility Directive

*Also see attachment "EC Declaration of Conformity".

The device has been constructed according to state-of-the-art technology and recognized safety regulations. However, risks may still arise during installation, operation and maintenance.

→ Please ensure the operating manual is available at all times.

The device may only be used under the following circumstances:

- Only operate the device, if it is in full working order.
- Ensure all operators of the device possess the necessary safety and risk awareness.
- Operate the device according to the instructions stipulated in this manual only.
- If there is something you do not understand, or certain information is missing, ask your manager or contact the manufacturer.
- Do not do anything on the device without authorization.
- Only use the device according to its intended use.

Intended use

The device is intended for use by trained and authorized personnel only.

The device is suitable for the following use:	<ul style="list-style-type: none"> ▪ economical cooling as alternative to tap water systems
The operation for the device is suitable in one of the following locations:	<ul style="list-style-type: none"> ▪ in any research application

Improper use

Any use which deviates from the device's intended use is considered to be improper. The manufacturer does not accept liability for any damages resulting from non-permitted uses. The risk is carried by the operator alone.

Installation / Electrical safety

- The device may only be connected when the mains voltage corresponds to the information on the type plate of the unit.
- The mains connection must be accessible at all times.
- Repairs may only be performed by a qualified electrician.
- Never operate the unit with a damaged power cord.
- Always turn the unit OFF and disconnect mains power before performing any maintenance or service.

Personnel qualification

- The device may only be operated by qualified persons.
- The device may only be operated by individuals who have been instructed in its proper use by qualified persons.
- The device may only be operated and maintained by persons who are of legal age.
- Other personnel may only work with the unit under continuous supervision of an experienced qualified operator.
- This manual must be read and understood by all persons working with the device.
- The personnel must have received special safety instructions in order to guarantee responsible and safe work procedures.

➤ Operating company's obligations

Installation site

- The device must be positioned in a suitable location.
- The device must be installed sufficiently stable on a strong and level surface.
- All screw connections must be securely tightened.
- The device should be located as close to possible to the process requiring cooling.
- The device should not be installed closer than 1.4 meters (4 feet) to a heat-generating source, such as heating pipes, boilers, etc.
- If possible, the device should be located near a suitable drain to prevent flooding in the event of leaks.
- Do not place the device where corrosive fumes, excessive moisture, excessive dust, or high room temperatures are present.
- The device may only be operated in enclosed spaces and under the following environmental conditions:

Ambient temperature	5 - 31 °C at 80 % relative humidity 32 - 40 °C decreasing linearly to max. 50 % relative humidity
Installation altitude	0 - 2,000 m (6,500 feet) above sea level
Degree of pollution	2
Overvoltage category	II

- Adequate clearance should be allowed on the front, sides, top, and rear of the device for access to connections and components.
- The front and side vents of the device must be a minimum of 21 cm (8 inches) away from walls or vertical surfaces so air flow is not restricted.
- If the equipment is used in corrosive atmospheres, the service life of the equipment will decrease based on concentration, volumes, and frequency of exposures to these corrosive materials, for example concentrated Hydrochloric Acid (HCl).

Changes to the unit

- No unauthorized changes may be made to the unit.
- No parts may be used which have not been approved by the manufacturer.
- Unauthorized changes result in the EC Declaration of Conformity losing its validity, and the appliance may no longer be operated.
- The manufacturer is not liable for any damage, danger or injuries that result from unauthorized changes or from operating the unit other than described in this manual.

Safety for the personnel

- Ensure that only qualified personnel operate the device.
- Observe the following regulations:
 - Laboratory guidelines
 - Accident prevention regulation
 - Ordinance on Hazardous Substances
 - Other generally accepted rules of safety engineering and occupational health
 - Local regulations

➤ Safety during use

- Wear the appropriate protective clothing when working on the device (clothing, protective glasses and if necessary, safety gloves).
- Do not use the device in potentially explosive areas. The device is not protected against explosion. There is no explosion or ATEX protection available.
- Do not operate or assemble devices in the vicinity which are emission or radiation sources (electromagnetic waves) for the frequency range (3×10^{11} Hz to 3×10^{15} Hz).
- Avoid putting pressure on the display when you do not operate the device.
- Eliminate errors immediately.
- Do not use abrasive material to clean the glass surfaces. Only wipe with damp cloths.
- Always switch the device off after use.


➤ Disposal

- Check the device components for hazardous substances and solvents.
- Clean all components before disposal.
- Dispose of the device according to the valid national legal regulations.
- Dispose of the packaging material in accordance to the valid national legal regulations.

➤ Device overview

RotaChill Large overall view

RotaChill Large control panel

➤ Set up

Power cord

Warning:
Electrical shock!

As long as the unit is not completely set up, electrical voltage is a potential risk to your health.

Do not plug the Chiller into the electrical outlet until the unit is ready for Start up.


Connect power cord

The power cord comes with a three-pole plug and a recess on the bottom side.

The socket is located on the back side of the unit.

- Connect plug to the socket with the recess facing left.
- Push plug up against socket until securely attached.
- Connect plug to power socket.

Disconnect power cord

- Disconnect plug from power socket first.
- Then disconnect plug from socket on unit.

Main switch

The **main switch** is located on the rear side of the unit.

With the main switch the unit is entirely switched on to the electric system and off.

Connect to power supply


→ Press On/Off switch at the upper side.

✓ The unit is now connected.

Disconnect from power supply

→ Press the On/Off switch at the lower side.

✓ The unit is now disconnected.

On/Off switch: On

On/Off switch: Off

Process Coolant

The Chiller is designed to accommodate a variety of coolant fluids (water, glycol mixtures, etc).

- For most applications above 20°C (68°F), distilled water is satisfactory.
- Ethylene glycol (laboratory grade) and distilled water in a 50/50 mixture is satisfactory from +20° to -15°C (68° to 5°F).

Caution: Possible damage to the cooling system.

At low temperatures the fluids might freeze and damage the chiller.

For operation at or below 20°C (68°F), the Chiller must be protected with an antifreeze solution.

Warning: Hazardous to health and environment!

Caustic, corrosive, or flammable fluids might endanger your health.

Do not use the following fluids:

- Automotive antifreeze with additives**
- Hard tap water**
- Deionized water with a specific resistance > 1 meg ohm
- Any flammable fluids
- Concentrations of acids or bases
- Solutions with halides: chlorides, fluorides, bromides, iodides or sulfur
- Bleach (Sodium Hypochlorite)
- Solutions with chromates or chromium salts
- Glycerin
- Syltherm fluids

Only use fluids that will satisfy safety, health, and equipment compatibility requirements.

** At temperatures above 40°C, additives or mineral deposits can adhere to the heater. If deposits are allowed to build up, the heater may overheat and fail. Higher temperatures and higher concentrations of additives can hasten deposit build up.

Start Up

Select a fluid that is compatible with the Chiller's wetted parts (brass, copper, stainless steel, EPDM rubber, Alcryn, nylon, PVC, and polyethylene).

- To prevent algae growth in the system, shield tubing from light. Insulation around the tubing will suffice.
- Use an algacide in the cooling water to minimize algae growth.

Caution:
Possible damage to the unit

Chlorine might damage metall surfces.

Do not use chlorine bleach as an algacide.

Filling the Reservoir

The reservoir holds approximatly 3 liters.

You will find the reservoir inlet at the top of the unit.

- Unscrew the reservoir cap.
- Carefully fill the reservoir with fluid.
- Leave the cap off and add additional fluid as instructed under section "Starting Process Fluid Flow", below.

Reservoir cap

Switch unit on and off

Caution:
Possible damage to the unit!

If proper fluid level is not maintained, the cooling system might get damaged or the unit might not operate.

Always ensure sufficient fluid level, before switching on the device


Switch unit on

- Press the On/Off button at the control panel.
- ✓ The unit is switched on.

Switch unit off

- Press the On/Off button at the control panel again.
- ✓ The unit powers off.

On/Off button

➤ Starting process fluid flow

- ✓ The reservoir is filled with cooling fluid up to the mark "Full".
- ✓ The power cord is connected.
- Place the Circuit Breaker/Power Switch on the rear of the instrument enclosure in the "On" position.
- ✓ Three decimal points will appear on the Temperature display.
- Press the Power Button on the front panel.
- ✓ The fluid level indicator illuminates blue.
- ✓ The set point temperature will appear briefly on the Temperature display.
- ✓ After a few seconds, it will be replaced by the actual fluid temperature.
- Check for leaks.
- ✓ The pump will turn on and fluid will begin circulating through the system.
- ✓ Ninety (90) seconds after power up, the compressor will begin operating.
- ✓ With the pump running, the reservoir's fluid level will drop as the process and/or process cooling lines fill with fluid.
- Slowly add fluid to the reservoir until the liquid level as shown on the front panel fluid level gauge remains stable at the "Full" mark.
- Replace the reservoir cap.

An EFL error (see chapter "Troubleshooting") upon initial start-up or after a routine fluid change may be caused by an air bubble trapped in the pump. Proceed as described below.

EFL Error

- Turn the unit Off and On three times.
- ✓ This will start and stop the pump to remove trapped air.
- Repeat this operation if large amounts of air are still present in the system.
- Monitor the fluid level in the reservoir as fluid level may drop.
- During this process, always verify that there are no fluid leaks and no fluid line restrictions.

The Chiller is equipped with the WhisperCool™ Environmental Control System, which increases/decreases fan speed based on heat load.

You will notice that the fan speed changes during operation. This is particularly beneficial in an environment where low noise level is desirable.

➤ Normal operation

Selecting the Temperature Unit (°C or °F)

The LEDs adjacent to the Temperature Display indicate the unit (°C or °F) in which temperature is displayed. To change from °C to °F or vice versa, proceed as follows:

To change to °F

- Place the Circuit Breaker/Power Switch on the rear of the instrument in the "Off" position.
- Press and hold the Down arrow button while returning the Circuit Breaker/Power Switch to the "On" position.

To change to °C

- Place the Circuit Breaker/Power Switch on the rear of the instrument in the "Off" position.
- Press and hold the Up arrow button while returning the Circuit Breaker/Power Switch to the "On" position.

Caution: Malfunction

User settings (except local lockout, baud rate, and calibration offset) return to the original factory defaults when the unit in which temperature is displayed is changed.

The Chiller's temperature set point and various alarm settings should be reset to the desired values.

Displaying and Adjusting the Set Point

Possible range	-20 °c and +60 °C
Factory setting	+10 °C

- Press the Set button on the front panel.
- ✓ The current set point temperature will be displayed and the decimal point at the bottom right of the display will flash, indicating the temperature can be changed.
- Press the Up and Down arrow button to adjust the set point temperature.
- ✓ The setting is accepted after either pressing the Set button a second time or will be accepted automatically after a few seconds of inactivity

While it is possible to set a temperature set point outside the temperature ranges outlined above, the Chiller will not be able to control at that temperature reliably.

Setting Operational Parameters

The following parameters are user-adjustable.

Display*	Description	Choices / Ranges / Comments	Default Setting
H-71	High Temperature Limit Set Point	1° to 71°C (33° to 159°F)	71°C (159°F)
L-46	Low Temperature Limit Set Point	-46° to -1°C (-50° to 31°F)	-46°C (-50°F)
A 30	Auto-Refrigeration Temperature	20°C – 90°C (always displayed/set in °C)	30°C
P.80	Specific Heat of Fluid	0.10 to 2.00	0.8
PC96	Communication Baud Rate x 100	24, 48, 96, 192	96
C0.0	Temperature Calibration	2.0° to -2.0°C (2.0° to -2.0°F)	0.0°C (0.0°F)
Pd00	Password	For factory use only	00

* Display Values; factory settings

Accessing Operational Parameters Menu Items:

- You can access and scroll through the operational parameter menu items by pressing and releasing the Up and Down arrow buttons simultaneously.
- ✓ The first menu item that will appear on the display is the High Temperature Limit value (shown as an H followed by the present setting; (e.g., H-71).
- To scroll to the next menu item, press and release the Up and Down arrow buttons again.

Changing an Operational Parameter:

- When the operational parameter you wish to change appears, wait until the decimal point on the right of the display flashes and then press the appropriate arrow button to change the displayed value.

Accepting an Operational Parameter Value:

- To accept the displayed value and return to the main operational display, either press the Set button or allow the display to timeout (approximately 15-20 seconds).
- To accept the displayed value and advance to the next menu item, press the Up and Down arrow buttons simultaneously.

Operational Parameters

High Temperature Limit (H-##)

High temperature limit defines:

- how high the temperature set point may be set
- the highest process fluid temperature at which the Chiller will operate

When this value is exceeded, the appropriate alarm or fault message will appear on the display (see chapter "Troubleshooting").

Operation

Low Temperature Limit (L-##)

Low temperature limit defines:

- how low the temperature set point may be set
- the lowest process fluid temperature at which the Chiller will operate

When this value is exceeded, the appropriate alarm or fault message will appear on the display (see chapter "Troubleshooting").

Auto-Refrigeration Temperature (A ##)

This menu item allows you to select the temperature at which refrigeration is activated. When the set point exceeds the auto-refrigeration temperature by more than 1.0°C, the cooling and the fan are turned off.

Caution: Possible wrong setting

This value is always displayed/set in °C.

If your normal operation unit is °F, always ensure with this parameter to convert the value to Celsius

Specific Heat (P.###)

This menu item allows you to adjust chiller performance to the specific heat of the cooling fluid, thus optimizing temperature control. It should be set to 0.80 for a 50% glycol / 50% water solution (the specific heat value for water is 1.00).

Temperature Calibration (C #.#)

The Chiller's temperature probe is accurate to $\pm 0.25^\circ\text{C}$. Therefore, if the displayed temperature reads 15°C, the actual temperature will be between 14.75°C and 15.25°C.

There may be situations where you want the displayed temperature to match a particular value to have standardization between different instruments.

This menu item allows you to adjust the Chiller's temperature reading to match that of a traceable standard or another piece of equipment. It allows you to offset the displayed temperature value by as much as $\pm 2.0^\circ\text{C}$ ($\pm 2.0^\circ\text{F}$).

Positive offsets are denoted with a blank space between C and the offset value (e.g., C 0.3).

Negative offsets are denoted with a minus sign between C and the offset value (C-0.3).

Password (Pd00)

This menu item is reserved for factory use only.

Enabling/Disabling the Local Lockout (LLO)

This feature is used to prevent unauthorized or accidental changes to set point and other operational values. When enabled, the temperature set point and operational parameters can be displayed, but not changed.

- To enable the local lockout, press and hold the Set button until LLO is displayed (approximately 5 seconds).
- ✓ Once enabled, LLO will appear momentarily when the Set button is pressed to display the set point.
- To disable the local lockout, press and hold the Set button until CAn appears, indicating that the local lockout feature has been disabled.

➤ Cleaning

All surfaces of the unit allow for cleaning with a damp cloth and if required a mild soap lotion.

Caution:
Damage to the unit

Liquid detergents might enter the control panel and damage the electronic of software.

Do not allow cleaning liquids or sprays to come in direct contact with the digital display.

Caution:
Damage to the surfaces

All surfaces can be damaged by improper cleaning.

In any case do not use:

- Chlorine bleach or any cleaning agent including chlorine
- Ammonia
- Abrasive cleaning agents such as cleaning rags, scrubbing agents or any other agents which include metal components

➤ Maintenance

The Chiller is designed to require a minimum of periodic maintenance.

Magnetic Drive Centrifugal Pump

When used under continuous operating conditions, this pump should be oiled every six (6) months with SAE 20 oil. The pump incorporates two oil ports for this purpose.

To access the pump:

- Turn both power switches off and unplug the power cord.
- Remove the two screws securing the left side panel (as viewed from the front) of the housing to the frame.

Fluid Level

The fluid level gauge on the front panel of the Chiller should be periodically checked to determine if the fluid level needs to be topped off. Generally, fluid should be added whenever the level in the reservoir is below the "Full" mark.

Condenser, Air Vents, and Reusable Filter

To keep the system operating at optimum cooling capacity, the condenser, the air vents, and reusable filter should be kept free of dust and dirt. They should be checked on a scheduled basis and cleaned as required.


→ To access the air filter, open the filter compartment door by placing your index fingers in the cutouts on the sides of the compartment door (upper corners) and pulling towards you.

→ Remove the filter and remove accumulated dust and dirt.


→ It may be cleaned by either vacuuming, blowing out with a clean air source, and/or washing with a mild detergent and water solution.

→ If using a detergent/water solution, be sure to rinse and dry thoroughly before reinstalling.

Access cutouts
compartment door

Washable
air filter

Any necessary repair must be performed by an authorized Heidolph distributor. Please contact Heidolph Instruments or your local Heidolph Instruments distributor.

➤ Possible errors and how to resolve

Error codes

When certain conditions are detected, message codes flash on the display. Depending on the nature of the condition, power to various systems components is removed. When the condition is rectified, press the front panel Power button or turn power switch/circuit breaker Off then On to clear the fault or error.

Error Message	Possible reason	Troubleshooting
EHL E-HI	<ul style="list-style-type: none"> ▪ Error High Limit: Set point temperature has been set higher than the high limit temperature. ▪ Chiller Operation Chiller continues normal operation; display alternates between EHL, E-HI and the fluid temperature. 	Decrease set point temperature or increase high limit temperature
FHL	<ul style="list-style-type: none"> ▪ Fault High Limit: Fluid temperature has exceeded the high limit temperature for more than 25 seconds. ▪ Chiller Operation: Normal operation stops, with heating and cooling turned off and the pump remaining on. FHL displayed continuously. 	NOTE: Power to the unit must be turned Off before making adjustment. Correct high fluid temperature problem or increase high limit temperature.
ELL E-LO	<ul style="list-style-type: none"> ▪ Error Low Limit: Set point temperature has been set lower than the low limit temperature. ▪ Chiller Operation: Chiller continues normal operation; display alternates between ELL, E-LO and the fluid temperature. 	Increase set point temperature or decrease low limit temperature.
FLL	<ul style="list-style-type: none"> ▪ Fault Low Limit: The fluid temperature has been below the low limit temperature for more than 25 seconds. ▪ Chiller Operation — Normal operation stops, with heating and cooling turned off and the pump remaining on. FLL displayed continuously. 	Correct low fluid temperature problem or decrease low limit temperature. (Power to the unit must be turned Off before making adjustment).
EFL	<ul style="list-style-type: none"> ▪ Error Fluid Low: Fluid flow is low. ▪ Chiller Operation — Chiller continues normal operation; display alternates between EFL and the fluid temperature. 	Add fluid to the reservoir; check fluid lines for restrictions and correct as required. Switch unit On / Off to remove trapped air bubbles (also see chapter "Starting process fluid flow").

Troubleshooting

Error Message	Possible reason	Troubleshooting
FfL	<ul style="list-style-type: none"> Fault Fluid Flow: The fluid flow has been low for more than 20 seconds. Chiller Operation: Normal operation stops, with heating, cooling, and pump turned off. FfL displayed continuously. 	Add fluid to the reservoir; check fluid lines for restrictions and correct as required. (Power to the unit must be turned Off before making adjustment).
Ft3	<ul style="list-style-type: none"> Fault 3: One or all instrument settings are out of range. Chiller Operation: Normal operation stops, with heating, cooling, and pump turned off. Ft3 displayed continuously. 	Contact your local Heidolph Instruments distributor
FOT	<ul style="list-style-type: none"> Fault OTP: Heater over-temperature fault. Chiller Operation: Normal operation stops, with heating, cooling, and pump turned off. FOT displayed continuously. 	Contact your local Heidolph Instruments distributor
Ft5	<ul style="list-style-type: none"> Fault 5: Heater Triac failure. Chiller Operation: Normal operation stops, with heating, cooling, and pump turned off. Ft5 displayed continuously. 	Contact your local Heidolph Instruments distributor
Ft7	<ul style="list-style-type: none"> Fault Temperature Probe Fault: The main temperature probe has failed. Chiller Operation: Normal operation stops, with heating, cooling, and pump turned off. Ft7 displayed continuously. 	Contact your local Heidolph Instruments distributor

Other problems

Many problems can be resolved by restoring the factory defaults. If this solves the problem, be careful when restoring your operational settings in order not to repeat the problem.

To restore the factory default settings:

- Place the Power Switch/Circuit Breaker on the rear of the unit in the OFF position.
- Press and hold either the Up (for °C) or Down (for °F) arrow button while returning the Circuit Breaker/Power Switch to the On position.

Warning: **Electrical shock**

When electrical power is On, dangerous voltages exist within chassis components.

Refer servicing to qualified service personnel.

Troubleshooting

Problem	Possible reason	Troubleshooting
Unit does not run (digital displays blank)	No power to unit	<ul style="list-style-type: none"> Check that the electrical cord is secure and connected to an operating electrical outlet. Check that Power Switch / Circuit Breaker on rear of unit is ON.
Unit does not run (three decimal points appear on temperature display, two decimal points on pressure/ flow rate display)	Unit in Standby mode	Press Power Button on front panel.
No fluid circulation	Insufficient fluid in reservoir Blockage in circulating system	<ul style="list-style-type: none"> Add fluid to reservoir Remove blockage
	Pump is not operating	Replace pump
Pump cycles On and Off at startup; EFL flashes on display	No flow condition detected due to insufficient fluid in reservoir or blockage in circulating system	<ul style="list-style-type: none"> Add fluid to reservoir as required Remove blockage as required
Insufficient circulation	Fluid viscosity too high External tubing diameter too small	<ul style="list-style-type: none"> Replace with lower viscosity fluid Replace with larger diameter tubing
	Restrictions in fluid lines	Check and correct as required.
	Low line voltage	Check and correct as required
Unit does not cool or cooling is insufficient	Dust build up on air filter or condenser	Clean air filter and/or condenser
	Blocked air ventilation screens	Remove blockages
	<ul style="list-style-type: none"> Excessive heat load Ambient air temperature too high 	<ul style="list-style-type: none"> Check that heat load does not exceed capacity of chiller; correct as required. Decrease ambient air temperature.
	Low or high line voltage	Check and correct as required.

If you experience a malfunction which can not be resolved, please contact your authorized Heidolph Instruments distributor immediately.

➤ Electrical connections

- This unit must **only** be connected **by a professional electrician**.
- **Any repair** on this unit must **only** be carried out **by a qualified specialist**.

**Warning:
Electric shock!**

This unit must only be connected to a grounded electrical socket.

Repairs are allowed by qualified professionals only. Any improper repair can result in a dangerous situation. Contact your local Heidolph Instruments distributor for any repair you may have.

The unit and voltage must match. The type label on the back side of the unit provides all voltage specifics.

When shipped, the unit is grounded. In case the original plug is changed, the new plug must have a protective conductor!

Color code for electrical connections:

Green/Yellow	PE: Protective conductor (Earth)
Blue	N: Neutral conductor
Brown	P: Phase
Black	P: Phase
Grey	P: Phase

The unit is connected with the power cord supplied to the electrical socket. It plugs into the housing on the rear side of the unit.

For countries utilizing a plug other than the standard one supplied:

- The plug supplied may only be changed by a professional electrician.
- If utilizing an adaptor ensure that it is approved by the local regulations.
- To help prevent voltage drops, position the Chiller as close as possible to the power distribution panel. Avoid voltage drops by using a properly grounded power outlet wired with 14 gauge or larger diameter wire.
- The use of an extension cord is not recommended.

➤ Set up chiller

Caution:
Possible damage to the cooling system

Fluid from compressor might leak into cooling circulation when unit is not in upright position.

Keep unit upright when moving.

The Chiller may be located at a level below that of the equipment being cooled. As long as the process remains closed, overflow will not occur when adding cooling fluid to the Chiller reservoir.

Hoses and connectors

To maintain a safe workplace and avoid leaks, special care should be taken when choosing hoses and connectors for the Chiller. It is the user's responsibility to ensure that the tubing and fittings connected to the Chiller are compatible with the fluid, temperature, and pressure being used.

Pressure Ratings	Hoses should be able to withstand the largest pressure that they will encounter.
Flexible Tubing	Avoid tubing that will expand and take up fluid volume when operating at the desired pressure.
Hose Diameter	Process piping/hosing with a diameter smaller than 0.5 inch ID can be used if desired. However, keep in mind that using smaller diameter hosing increases pressure and reduces flow in the circulating system and may result in a flow alarm. The minimum recommended tubing size is 3/8 inch ID.
Couplings and Clamps	The use of screw-tightened hose clamps is necessary on all joints to insure good, tight connections. Quick connectors are not recommended as they have the potential for restricting flow rate.
Fluid Filter	If particulate matter may be present in the fluid, the installation of a 50 micron filter on the chiller inlet is recommended (see section "Accessories").

Reservoir Drain

A 0.5 inch NPT connection is provided for the reservoir's gravity drain.

- It should be piped to a drain or receptacle positioned below the bottom of the drain.
- If a receptacle is used, be sure it is of sufficient volume to hold all the water in the reservoir, process, and process lines.

Closed System or Cooling Coil Setup

- Connect the Chiller's inlet and outlet to the external apparatus with hoses or pipes.
- The direction of the flow through the system can be controlled by the way the connections are made.
Fluid is drawn into the Chiller through the "Inlet" connection; fluid is pumped out of the Chiller through the "Outlet" connection.

➤ Dismantling, transportation and storage

Dismantling

Warning:
Risk of injury!

While the unit is plugged in you might be at harm if the device is accidentally switched to ON mode.

Prior to dismantling always disconnect the plug from the mains.

- ✓ The unit is switched off.
- Disconnect plug from power socket.
- Disconnect plug from socket on the unit.
- Empty complete coolant fluid.

Transportation and storage

- Store and transport the unit and its components only if they were emptied and cleaned in the original packing material.
- Alternatively store and transport the unit in a container which prevents damages during transportation.
- Seal the packing carefully against unauthorized or mistaken opening.
- Store the unit in a dry place.

Improper transportation may result in damages to the system and the mechanical components.

Avoid any kind of shocks during transportation.

➤ Scope of delivery

Component	Variant	Quantity	Product number
RotaChill Large		1	591-00240-00
Operating Manual English / German		1	01-005-005-66
Warranty registration / Confirmation of condition		1	01-006-002-78
EC Declaration of Conformity		1	01-001-025-16
Inlet/Outlet Adaptor	1/2" NPT to 3/8" (10 mm) hose connector	1	23-30-01-07-64
Power cord		1	

Additional accessories are shown in the general catalog or on our website at

- www.heidolph-instruments.com

➤ Technical data


Standard supply voltage	1 x 230 V (50/60 Hz)
Operational voltage	216-264 V (50/60 Hz)
Supply power (W)	1,610
Pump	centrifugal drive pump
Refrigerant type	R-404a
Heat exchanger	Stainless steel copper-brazed plate
Noise rating	53 dBA
Temperature stability (°C)	±0.1
Overheat protection	yes
Auto-restart on power failure	yes
Low flow alarm and power cut off	yes
Reservoir volume (L)	2.65
Dimensions (mm) (WxHxD)	254 x 483 x 607
Weight (kg)	46.3
Permissible ambient conditions (°C)	80 % relative humidity, no condensation 5 - 35 °C
Maximum altitude (m)	2,000
Fluid inlet / outlet	0.5 inch female brass NPT
Ø Recommended tubing	3/8" (10 mm) minimum


Performance range

Cooling capacity	Ethylene Glycol & Water (50/50 mix)
Working temperature range 1	-20° to +40°C / -4° to +104°F
Operating temperature range 2	-20° to +60°C / -4° to +140°F

Pump	M2 Centrifugal		MY Centrifugal		TY Turbine		
	@	Watts	BTU/Hr	Watts	BTU/Hr	@	Watts
-20 °C	230	785	240	820	-20 °C	140	478
-10 °C	435	1485	460	1571	-10 °C	330	1127
0 °C	680	2322	700	2391	0 °C	500	1708
+10°C	1030	3518	1070	3996	+10°C	750	2561
+20 °C	1160	3962	1190	4064	+20 °C	830	2835
+30 °C	1380	4713	1420	4850	+30 °C	990	3381
+40 °C	1545	5276	1580	5396	+40 °C	1105	3774

Pump option	M2 Centrifugal	MY Centrifugal	TY Turbine
Maximum flow	12.9 lpm / 3.4 gpm	11.7 lpm / 3.1 gpm	8.3 lpm / 2.2 gpm
Maximum Pressure	5.5 psi / 0.38 bar	10.5 psi / 0.72 mbar	32 psi / 2.2 bar
Maximum Head	4 m / 12.8 ft H2O	7 m / 23 ft H2O	17 m / 56.5 ft H2O

* Performance data based on 240V, 50hz input power, 20°C ambient temperature and a 50/50 mix of Ethylene Glycol and distilled water as coolant.

➤ Contact / Technical Service

Questions / Repair work

If any aspect of installation, operation or maintenance remains unanswered in the present manual, please get in touch with the following address.

For repairs please call Heidolph Instruments or your local, authorized Heidolph Instruments distributor.

**Warning:
Danger of poisoning!**

Contaminated units can lead to severe injury or death of our employees!

When shipping items for repair that may have been contaminated by hazardous substances, please

- advice exact substance
 - take adequate protective action towards our parts receiving and service personnel
 - mark the pack in accordance with Ordinance on Hazardous Substances
-

Included at the end of this manual you will find a „Confirmation of condition of unit“.

→ Kindly copy and fill in this form and submit it prior to shipping the unit for repair.

Contact details:

Heidolph Instruments Germany

Heidolph Instruments GmbH & Co. KG
Technical Service
Walpersdorfer Str. 12
D-91126 Schwabach / Germany

Tel.: +49 - 9122 - 9920-74
Fax: +49 - 9122 - 9920-84

E-Mail: service@heidolph.de
www.heidolph.de

Heidolph Instruments North America

Phone: 1-866-650-9604
E-mail: service@heidolph.com
www.heidolphNA.com

Heidolph Instruments United Kingdom

Phone: 01799 - 5133-20
E-mail: service@radleys.co.uk
www.heidolph-instruments.co.uk

All other countries

You will find contact details of your local Heidolph distributor at
www.heidolph.com

➤ Warranty

Heidolph Instruments provides a 3-year warranty for the products described here (excluding glass and wearing parts) if you register using the warranty card enclosed or online (www.heidolph.com). The warranty is valid from the point of registration. The serial number is also valid without registering.

The warranty covers part and manufacturing defects.

In the case of a part or manufacturing defect, the device shall be repaired or replaced free of charge under the terms of the warranty.

Heidolph Instruments shall not assume liability for any damages incurred as a result of improper handling or transport.

Warranty claim?

- ➔ Please inform Heidolph Instruments or your local distributor should you wish to make a warranty claim.
-

➤ Confirmation of condition of unit

→ In the case of repair, copy and complete the Confirmation of condition of unit and send it to Heidolph Instruments.

1. Details about the unit

Product number _____

Serial number _____

Reason for repair _____

2. Has the device been cleaned, decontaminated/sterilized?

Yes _____ No _____

3. Is the unit in a condition which does not represent any health threats for the staff of our service department?

Yes _____ No _____

If not, which substances has the unit come into contact with?

4. Legally binding declaration

The customer is aware of being legally liable to Heidolph Instruments for any damages arising from incomplete and incorrect information.

Date _____ Signature _____

Company stamp

Please note

The shipper is responsible for the return of the goods in well packed condition, suitable for the mode of transport.

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We reserve the right to make any technical modifications without prior notification. This documentation is not subject to revision; the latest version of this documentation can be found on our homepage. The copyright for all texts and images is held by Heidolph Instruments GmbH & Co. KG.
© Heidolph Instruments GmbH & Co. KG