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Assembly instructions
Page 50 – 94

Contents

About this document	50
Typographic conventions	50
Copyright protection	50
Work safety	51
Personal protective equipment (PPE)	51
Scope of delivery	52
General notes	53
Part A: Hei-VOLUME Distimatic Pro with Hei-VAP Industrial	54
Vacuum unit	55
Glassware	57
Buffer vessel	58
Collector vessel	59
Educt sensor vessel	60
Vapor tube distributor unit	61
Condenser connection tube	62
Distimatic Pro support	63
Tube connections	64
Tube connection vacuum unit – compressor Distimatic Pro – buffer vessel inlet	65
Tube connection buffer vessel outlet – vacuum connection socket	66
Tube connection terminal block mounting plate 2 – collector vessel	67
Tube connection collector vessel – condenser	68
Tube connection evaporation flask – condenser – collector unit – residue container – distillate container	69
Tube connection educt container – educt sensor vessel – evaporation flask	70
Cable duct	71
Retaining clips for the base	72
Conversion of existing devices	73
Vacuum sensor block	73
Replacement of the vacuum connection socket	75

Part B: Hei-VOLUME Distimatic Pro with Hei-VAP Expert/Ultimate..... 76

Assembly description77

- Condenser 78
- Collector..... 80

Valve matrix82

Tube connections.....84

- Vacuum unit (valve V5) – buffer vessel inlet 85
- Outlet, buffer vessel – condenser 86
- Collector – valve block vacuum 87
- Residue draining, valve V3 – condenser 88
- Educt inlet, valve V1 – condenser..... 89
- Educt supply & residue draining 90
- Distillate draining..... 91
- Condenser coolant supply 92

Part C: Cable connections Hei-VOLUME Distimatic Pro – Hei-VAP Industrial/Expert/Ultimate 93

Cable connections.....93


- Connections at rear of Distimatic Pro 94

About this document

These instructions describe the assembly, wiring and tubing of the components of an automatic module, type Hei-VOLUME Distimatic Pro on large-scale rotary evaporators of the type Hei-VAP Industrial (Part A) and on rotary evaporators of the type Hei-VAP Expert and Hei-VAP Ultimate (Part B).

Typographic conventions

Standardized symbols, highlighting elements, and signal words are used in this document to identify warnings, cautions, important information, and special text contents.

Symbol	Signal word / explanatory note
	<p>Warning symbols in combination with a signal word indicate dangers:</p> <p>DANGER Indicates an immediate dangerous situation which, if not avoided, will result in death or serious injury.</p> <p>WARNING Indicates a potential danger which, if not avoided, may result in death or serious injury.</p> <p>CAUTION Indicates a potential hazard which, if not avoided, damage to property and minor to moderate injuries can occur.</p>
	<p>Mandatory signs are used to indicate important and useful information on handling a product.</p> <p>This information is used to ensure operational safety and to maintain the value of the product.</p>
	<p>The arrow symbol indicates instructions to be followed in order to ensure the operational safety when handling the product.</p>

Copyright protection

This publication is protected by copyright and intended for internal use by the purchaser of the product only.

No part of this publication may be transmitted or reproduced in any form, by any means, without the prior written consent of the copyright owner Heidolph Instruments GmbH & Co. KG. Any violation is subject to compensation for damages.

Work safety

CAUTION: Material damage, loss of production

In the event of improperly carried out assembly work, improper system tubing and wiring and unauthorized changes to the system setup, there is a risk of direct and indirect damage to property!

- Pay attention to all notes and follow all relevant information and instructions in this document.
- Do not make any unauthorized changes to the device, system, tubing, nor the existing cable connections.
- If necessary, contact our technical service (see next page).

- All installation and assembly work may only be carried out by qualified and specially trained personnel.
- Pay attention to all notes and follow all instructions and regulations in this document and all general and statutory safety and health regulations applicable in the place of use.
- Only use suitable tools and suitable aids for the assembly/disassembly of components.
- Always carry out the described assembly operations with at least two people.
- Before carrying out assembly work, make sure that all units concerned and all peripheral devices are completely disconnected from the power supply and are secured against reconnection.

The professional mounting and installation of the device including the supplied accessories is within the sole responsibility of the operator!

Noncompliance will invalidate any warranty against Heidolph Instruments.

The operator is solely liable for all damage resulting from unauthorized changes or modifications to the unit, from the use of unauthorized or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Personal protective equipment (PPE)

The operator must determine and provide the necessary PPE, depending on the respective application and the tools and chemicals used.

The corresponding instruction of the personnel is solely within the operator's responsibility.

Scope of delivery

- The contents of the automatic module, Hei-VOLUME Distimatic Pro for large-scale rotary evaporators of the type Hei-VAP Industrial, include the following components:
 - Distimatic Pro
 - Distimatic Pro support
 - Vacuum unit, consisting of mounting plate 1 and mounting plate 2
 - Glassware, consisting of buffer vessel, collector vessel, educt sensor vessel, distributor, condenser coupling incl. fastening clamps
 - Tubing set

- The contents of the automatic module, Hei-VOLUME Distimatic Pro for rotary evaporators of the type Hei-VAP Expert and Hei-VAP Ultimate, include the following components:
 - Distimatic Pro
 - Complete valve matrix
 - Glassware, consisting of condenser, Woulff bottle, collector vessel, educt sensor vessel
 - Tubing set

Check the delivery immediately for visible damage and completeness.

In case of need, please contact an authorized dealer or our technical service:

Heidolph Instruments GmbH & Co. KG

Walpersdorfer Str. 12
D-91126 Schwabach/Germany

Phone: +49 9122 9920--0

Fax: +49-9122-9920- 84

Email: service@heidolph.de

Representations

To find your local distributor please visit www.heidolph.com

General notes

If necessary, Heidolph offers a professional assembly service for the mechanical setup and installation of the device including all components supplied with it and including commissioning the Hei-VAP Industrial or Hei-VAP Expert/Ultimate and Hei-VOLUME Distimatic Pro system.

If necessary, contact our technical service:

Heidolph Instruments GmbH & Co. KG

Technical Service
Walpersdorfer Str. 12
D-91126 Schwabach/Germany
Phone: +49 - 9122 - 9920-380
Fax: +49-9122-9920- 84
Email: service@heidolph.de

Representations

To find your local distributor please visit www.heidolph.com

If you assemble it yourself, follow the specific safety instructions below:

WARNING

Risk of electric shock.

Risk of direct and indirect personal injuries and/or damage to property due to contact with chemical residues in the device components and tubing!

- Completely disconnect the rotary evaporator and – in case of a conversion – the Distimatic module from the power supply, disconnect the devices from the power supply system and secure them against reconnection.
 - Before dismantling, empty all system containers properly (heating bath, glassware).
 - Clean and dry dismantled components and store them properly in the respective original packaging.
-

Part A of these mounting instructions describes the procedure for setting up a system consisting of a Hei-VOLUME Distimatic Pro automatic module and a large-scale rotary evaporator of the type Hei-VAP Industrial.


Vacuum unit

The vacuum unit of the Hei-VOLUME Distimatic Pro automatic module is supplied pre-assembled and wired ready for operation on mounting plates:

Fasten the mounting plates of the vacuum unit with screws on the glassware protection housing

- Screw the pre-assembled mounting plate 1 of the vacuum unit on the rear panel of the glassware protection housing: Insert the supplied mounting screws from the inside through the existing mounting holes in the rear panel of the glassware protection housing of the rotary evaporator.

- Screw the pre-assembled mounting plate 2 of the vacuum unit underneath the mounting plate 1 on the rear panel of the glassware protection housing.
- Connect the terminal blocks of the two mounting plates as shown in the figure above, using a suitably cut tube piece (use the supplied tube material!).


Use the supplied tubing material for all tube connections. Be sure to cut the tube pieces at both ends at right angles and flat, using a suitable tool.

Each cut surface of the tube must be in contact with the counterpart over the entire diameter.

This is to ensure that every connection closes sufficiently tightly!


Glassware

The following figure shows the mounting position of all glass vessels on the rotary evaporator (maximum configuration):

Buffer vessel

The buffer vessel (Woulff bottle) is supplied with a pre-assembled mounting bracket and two matching mounting screws:

Mounting bracket
(pre-assembled)

Mounting screws
(scope of delivery)


Mounting position on HVI:
left panel, below the
glassware protection housing

- Using a suitable screwdriver, mount the buffer vessel in the intended position on the left side panel, below the glassware protection housing.
- To do this, use the existing threaded holes.

Collector vessel


The collector is supplied as a complete unit with a pre-assembled mounting bracket and two matching mounting screws:

- Using a suitable screwdriver, mount the collector in the intended position on the left side panel, below the glassware protection housing (next to the buffer vessel).
- To do this, use the existing threaded holes.

Educt sensor vessel

The educt sensor vessel is supplied as a complete unit with a pre-assembled mounting bracket and two matching mounting screws:

Mounting position on HVI:
left panel, below the
glassware protection housing
(next to the collector)


- Using a suitable screwdriver, mount the educt sensor vessel in the intended position on the left-hand side panel, below the glassware protection housing (next to the collector).
- To do this, use the existing threaded holes.

Vapor tube distributor unit

The distributor unit for the vapor tube is supplied as a pre-assembled unit. Observe the instructions contained in the operating instructions of the rotary evaporator when assembling the distributor unit!

Condenser connection tube

Observe the instructions contained in the operating instructions of the rotary evaporator when replacing the condenser connection tube!


Use the supplied clamps to mount the connection tube.

Distimatic Pro support

The Distimatic Pro support is supplied as a pre-assembled unit. Insert the support bars into the openings for transport tubes of the rotary evaporator and mount the fastening lugs using the supplied screws on the side panel of the rotary evaporator.

The base plate of the Distimatic Pro support is clamped with two hand screws on each support bar. Thus, the position of the Distimatic Pro to the rotary evaporator may be adjusted at any time

Tube connections

Make the required tube connections between the vacuum unit, the buffer vessel, the collector and the educt sensor vessel as shown in the following sections.

Only use the pieces of tube supplied as raw material for all tube connections.

Always select the tube diameter to match the respective connections.

Be sure to cut the tube pieces at both ends at right angles and flat, using a suitable tool.


Each cut surface of the tube must be in contact with the counterpart over the entire diameter.

This is to ensure that every connection closes sufficiently tightly!

Tube set (included)

Tube connection vacuum unit – compressor Distimatic Pro – buffer vessel inlet

- 1 Tube connection terminal block mounting plate 1 – compressor OUT/control box (pneumatic tube, blue)


- 2 Tube connection terminal block mounting plate 1 – buffer vessel inlet (PTFE, transparent)

Buffer vessel inlet

Inner tube
buffer vessel
(pre-mounted)

Tube connection buffer vessel outlet – vacuum connection socket

- 3 Tube connection vacuum connection socket – buffer vessel outlet (PTFE, transparent)

Buffer vessel outlet

Tube connection terminal block mounting plate 2 – collector vessel


- 4 Tube connection terminal block mounting plate 2 – collector vessel inlet (PTFE, transparent)

Collector vessel inlet


Tube connection collector vessel – condenser

Connect the two outer (lateral) connection sockets of the collector to the return and inlet sockets from/to the condenser of the rotary evaporator. The flow direction is arbitrary!


Tube connection evaporation flask – condenser – collector unit – residue container – distillate container

- Connect the residue drain tube, coming from the evaporation flask to the inlet of valve V3.
- Attach the tube to the residue tank to the outlet of valve V3.
- Establish the tube connection between the condenser and the distillate inlet of the collector.
- Attach the tube to the distillate container to the outlet of valve V2.

Tube connection educt container – educt sensor vessel – evaporation flask


- Connect the inlet of the distributor (connection educt valve V1) to the outlet of valve V1 on the educt sensor vessel.
- Attach the educt suction tube to the inlet of the educt sensor vessel.

Educt valve V1 connection at distributor unit, connected to educt sensor vessel

Cut-off sensor S3, to be connected to the signal amplifier of the Distimatic

Educt valve V1 at educt sensor vessel, outlet to educt valve V1 connection at distributor unit



Educt sensor vessel

Educt sensor S1, to be connected to any given coded Distimatic terminal

Educt inlet, from educt container




Cable duct

- Guide all cables through the cable holders provided on the mounting plates of the assembled vessels, see the following figures (arrows)!
- Group the routed and connected cables with cable ties.
- Be careful not to bend or pinch the cables! Use the existing cut-outs on the rotary evaporator housing and the Distimatic Pro support.

Retaining clips for the base

To stabilize the overall system, the four rotary evaporator feet are fixed to the base with clip-on retaining clips [1]:

- Remove the rotary evaporator front panel to install or remove the front two retaining clips [2] (the rear clips are directly accessible).
- [3a]: Retaining clip attached; [3b]: Retaining clip in plugged position, front left.
- [4]: Retaining clip in plugged position, front right.

Before commissioning, make sure that all four clamps are properly inserted to ensure the stability of the entire system!


Conversion of existing devices

Vacuum sensor block

When retrofitting existing devices, the existing vacuum sensor block must be converted accordingly first. This step is omitted for new devices!

Remove the existing vacuum sensor unit from the rear panel of the glassware protection housing on the rotary evaporator.


- Disconnect the connecting cable from the vacuum sensor.
- Loosen the hose fitting on the sensor block and remove the vacuum tube.
- Loosen the mounting screws of the sensor block (accessible from the inside of the glassware protection housing).
- Remove the sensor block from the rear panel of the glassware protection housing.
- Remove the vacuum sensor from the existing sensor block.

→ Mount the existing vacuum sensor on the sensor block of the new vacuum unit (see following section "Vacuum unit" on page 55).


1. Unscrew the vacuum sensor from the existing sensor block (turn anti-clockwise)

2. Screw the vacuum sensor into the new vacuum unit (turn clockwise)

3. Connect the vacuum sensor →



Replacement of the vacuum connection socket


Only for conversion. This step is omitted for new devices!

If converting existing devices, the existing vacuum connection socket on the connection block at the rear of the rotary evaporator must be replaced with the new vacuum connection socket supplied. The latter has a visibly larger cross-section!

Replace the existing with the **new** vacuum connection socket

Mounting position vacuum connection socket

Part B of these mounting instructions describes the procedure for setting up a system of a Hei-VOLUME Distimatic Pro automatic module and a rotary evaporator of the type Hei-VAP Expert and Hei-VAP Ultimate.

Assembly description

A Hei-VOLUME Distimatic Pro automatic module & Hei-VAP Expert or Hei-VAP Ultimate rotary evaporator system consists of different assemblies, which are delivered pre-assembled and must be put together and connected in the place of use as described in the following sections.


During the assembly, pay attention to cleanliness and good accessibility of the assembly units and comply with the following specific instructions:

- Make sure that the assembly/work area is clean and dry.
- Where possible, maintain a distance of approx. 30 cm between the back of the rotary evaporator and other devices or elements of the building, to ensure easy access to the plug-in connections at the rear.
- We strongly recommend that all assembly/installation work is always carried out by at least two persons.
- In particular, note that the valve matrix is not tightly connected to the rotary evaporator. This means that the whole system cannot be moved as a compact unit after assembly/installation. Instead, the valve matrix and the rotary evaporator must be moved individually and at the same time so that the existing tube and cable connections are not exposed to any mechanical loads during transport.
- Only use the components supplied and the recommended tools!

Condenser


The following figure shows the preassembled G9 XL condenser in the as-delivered condition:

Component / Function	
1	Condenser G9 XL
2	Level sensor S4
3	Emergency stop sensor S3

Condenser assembly


- Feed the lance of the level sensor [2] through the vapor tube [4] of the rotary evaporator.
- Position the condenser on the drive unit [1] of the rotary evaporator and fix the condenser with the union nut [3].

Component / Function	
1	Rotary evaporator drive unit
2	Level sensor S4, inserted in the vapor tube
3	Union nut
4	Vapor tube

Collector

The following figure shows the preassembled collector in the as-delivered condition:


Component / Function

1	Collector vessel
2	Retaining and stabilization element
3	Strain relief, condenser fitting
4	Intermediate valve V4
5	Collector sensor S2
6	Distillate valve V2

Collector assembly


- Preload the strain relief by turning the two wing screws upwards to the right-hand stop as shown [3a].
- Undo the four knurled screws in the adjusting plates of the retaining and stabilization element [2].
- Position the collector on the bottom screw thread of the condenser [1] and fix the collector with the union nut of the strain relief [4].
 - At the same time, make sure that the retaining and stabilization element [2] is free from mechanical stress between the union nut, which fixes the condenser on the drive unit (see condenser assembly), and the flange of the drive unit grips around the vapor tube.
- Align both adjusting plates of the retaining and stabilization element [2] so that the contact closure in the area of the drive unit is as uniform as possible, and tighten the four knurled screws finger tight.
- Slacken the strain relief, by turning the two wing screws downwards to the left-hand as shown [3b].


Valve matrix

With the exception of the buffer vessel (Woulff bottle), the valve matrix is delivered completely preassembled. To assemble the buffer vessel (delivered in separate packaging with two suitable knurled screws [1]), proceed as follows:

- Unscrew the buffer vessel from the bracket.
- Use the two knurled screws to screw the bracket of the buffer vessel into the threaded holes on the left side of the valve matrix [2]. Photo [3] shows the correct positioning.
- Screw the buffer vessel back into the bracket, photo [4]:


Thanks to its special design, the valve matrix has the inherent stability required for safe operation. For this reason, after assembly of the buffer vessel, the valve matrix is only pushed towards the rotary evaporator so that the recess in the bottom [1] grips around the foot of the rotary evaporator [2], see following figure. A direct mechanical connection between the valve matrix and rotary evaporator is not planned or required.

Valve matrix components

The following figure shows the valve matrix in the standard configuration:

Component / Function

	Component / Function
1	Valve V5, connection between vacuum unit or compressor and condenser (via switching valve)
2	Inlet valve V1, from the educt container
3	Valve V6 to the collector (collector draining)
4	Switching valve UV1, switching between vacuum pump (negative pressure) and compressor (overpressure)
5	Venting valve BV1 (Distimatic)
6	Educt sensor vessel with educt sensor S1
7	Valve V3 to the residue container (residue draining, variant with automatic residue draining)
8	Vacuum vent valve, solenoid valve DN3 (evaporation flask)
9	Mounting plate for optional vacuum valve (INTEGRATED variant)


Additional vacuum valve option (INTEGRATED variant)

Variants with valve-controlled vacuum pump must be equipped with an additional vacuum valve, to which the optional diaphragm vacuum pump is connected directly.

In this variant the corresponding connection to the Distimatic Pro is also omitted!

- Assemble the optional vacuum valve [1] supplied on the mounting plate [9] of the valve matrix provided (see figure in section “Valve matrix components” on page 83):

Tube connections

All tube connections are push-fit connections. The assembly only requires a suitable tube cutter (see previous section for designations of the components “Valve matrix components” on page 83).

CAUTION

Especially in the area of glass threads/glass connections, there is a risk of fracture due to the application of excessive force when inserting the tube ends.

- When inserting the tube ends into the connectors, always exercise the required care.
- To insert the tube ends, where possible, dismantle the tube connectors from the glass threads/glass components.
- Always insert the tube ends at right-angles to the connector, avoid canting.


Vacuum unit (valve V5) – buffer vessel inlet

Preferably carry out this step BEFORE positioning the valve matrix on the rotary evaporator:

A preassembled piece of tube including a suitable angle clip is supplied for the tube connection between the vacuum unit [1] with valve V5 and the inlet of the buffer vessel [2] (distributor unit with inner tube [3] in the vessel, see figure).

- Firstly, fit the piece of tube in the angle clip supplied (clamped connection, no tools required).
- Connect the distributor unit [2] (with inner tube [3]) of the buffer vessel to the tube connection of the vacuum unit [1]:


Outlet, buffer vessel – condenser


Preferably carry out this step BEFORE positioning the valve matrix on the rotary evaporator:

- Cut a piece of tube approx. 65 cm long (PTFE, transparent, sold by the meter, included).
- Dismantle the tube connector [2] from the condenser, insert the tube end into the tube connector and screw the tube connector back onto the condenser [3].
- Connect the tube end to the unused outlet of the buffer vessel [1].

Collector – valve block vacuum

- Required tube: PTFE, transparent (included, sold by the meter).
- Connect the distributor unit of the valve block [1] to the distributor unit of the collector [2].


Residue draining, valve V3 – condenser


Only variant with residue draining. Variant without residue draining does not have valve V3.

- Required tube: PTFE, transparent (included, sold by the meter).
- Connect the distributor unit of the residue valve V3 [1] to the bottom left-hand tube connection of the condenser [2].

Educt inlet, valve V1 – condenser

- Required tube: PTFE, transparent (included, sold by the meter).
- Connect the distributor unit of the educt valve V1 [1] to the bottom right-hand tube connection of the condenser [2].


Educt supply & residue draining

- Required tube: PTFE, transparent (included, sold by the meter).
- Position the educt container and the residue container in the place of use and cut the two pieces of tube to the required lengths between the container and valve matrix.
- Fit the tube ends, which are plugged into the valve matrix, in the angle clips supplied for this purpose (see figure, Item [1] and [2]).


Make sure that the free end of the tube protrudes from the angle clip far enough to ensure a secure connection to the tube connector!

- Educt supply: Tube connection [1] is used to supply the system with educt. Route the tube [1] from the educt sensor vessel into the educt container.
- Residue draining: Tube connection [2] is used to drain the residue from the system. Route the tube [2] from the residue valve into the residue container.

Distillate draining

- Required tube: PTFE, transparent (included, sold by the meter).
- Position the distillate container and cut the piece of tube to the required length.
- Insert one end of the tube into the distributor unit provided on the distillate valve V2 and route the tube into the distillate container.

Condenser coolant supply

- Required tube: PTFE, transparent (included, sold by the meter).
- Position the external condenser and cut the pieces of tube required to the necessary length (inlet, outlet).
- Connect the inlet and return tube to the condenser with connections [1] and [2] and make the connections to the inlet and outlet of the external condenser. The flow connection of the coolant can be set as required, i.e. the connections are interchangeable.

This part of the mounting instructions describes the cable connections between any Hei-VAP Industrial/Expert/Ultimate rotary evaporator and the control module of the Distimatic Pro automatic module.

Cable connections

All outgoing cables and connectors from the valve matrix and all connection sockets of the Distimatic Pro are indelibly labeled and can be easily assigned with no further tools.


All cable connections to the Distimatic Pro automatic module are fitted with coded connectors. Thanks to these, the connection cables can be connected to the plug-in cards as required, see following section "Connections at rear of Distimatic Pro" on page 94.

Connections at rear of Distimatic Pro

The following figure shows the rear of the Distimatic Pro with standard plug-in cards and optional expansion cards. The actual assignment depends on the scope of the order!

The coded connectors of the system's sensors and valves are plugged into the connection sockets of the sensor cards, are identified automatically and are assigned in the device software.

In addition, observe the marking of the existing connection sockets and all information on replacing existing or installing optional expansion cards in the supplied assembly documentation!

The **NOT-HALT EXTERN/EMERGENCY STOP** connection socket on the automatic module must be assigned, i.e., either the integrated connector (as shown in the following figure) or an external emergency stop switch must be plugged in here. Otherwise, the automatic module is in the emergency stop state and no processes can be started!

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© Heidolph Instruments GmbH & Co. KG