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Translation of the original instructions
Page 80 – 149

Certifications

Content

Introduction

| | |
|------------------------------|----|
| About this document..... | 80 |
| Typographic conventions..... | 80 |
| Copyright protection..... | 80 |

General information

| | |
|--|----|
| Notes on the product | 81 |
| Directives applied, product certification..... | 81 |
| California Residents..... | 81 |
| Residual risk..... | 81 |
| Intended use..... | 81 |
| Reasonably foreseeable misuse | 81 |
| Transportation..... | 82 |
| Storage..... | 82 |
| Acclimatization | 82 |
| Permissible ambient conditions | 82 |

Safety

| | |
|--|----|
| General safety instructions..... | 83 |
| Electrical safety..... | 83 |
| Data security | 83 |
| Operational safety | 83 |
| Overpressure operation | 84 |
| Biohazard..... | 84 |
| Special hygiene measures for the use of laboratory equipment in food, cosmetics and pharmaceutical production | 84 |
| General measures..... | 84 |
| Device-specific measures | 85 |
| Other regulations..... | 85 |

Assembly

| | |
|----------------|----|
| Assembly | 86 |
|----------------|----|

Device description

| | |
|---|----|
| Principle of operation | 87 |
| Mechanical design | 88 |
| Control elements | 88 |
| Legend | 88 |
| Connections rear side Distimatic Pro Basis | 89 |
| Hei-VOLUME Distimatic Pro & Hei-VAP Industrial | 90 |
| Hei-VOLUME Distimatic Pro & Hei-VAP Expert/Ultimate | 91 |
| Heating tape (optional) | 92 |
| Signaling of the operating states (optional) | 93 |
| Structure of the user interface | 94 |
| Control and display elements | 95 |
| System time, timer, stopwatch | 97 |
| Color coding | 98 |
| Symbols | 98 |
| Input keys | 99 |
| User input | 99 |

Commissioning

| | |
|--|-----|
| Commissioning | 100 |
| Switch the device on/off | 100 |
| Calibrate the flask level sensor | 101 |
| Initial filling of the evaporation flask | 104 |
| Device configuration | 105 |
| Main menu | 105 |
| Applications menu | 107 |
| Settings menu | 108 |
| Screen lock | 108 |
| Display brightness | 109 |
| Calibrate filling level | 110 |
| Maximum pressure | 110 |
| Sensitivity cut-out sensor S3 | 110 |
| Factory settings | 110 |
| Selecting the user language | 110 |
| Select formats | 111 |
| Automatic drainage of residue | 112 |
| Select evaporator type/glassware | 112 |
| Data recording | 113 |
| Connectivity | 115 |

| | |
|------------------------------------|-----|
| System check | 116 |
| Testing valves and signaling | 117 |
| Perform overpressure test | 118 |
| Testing the safety switch | 119 |

Operation

| | |
|--|-----|
| Selecting the operating mode | 120 |
| Sensor mode | 120 |
| Parameter sensor mode | 124 |
| Parameters Hei-VAP | 127 |
| Data logging | 127 |
| Recipes | 128 |
| Recipe list | 129 |
| Time mode | 130 |
| Time mode parameters | 132 |
| Rinse mode | 135 |
| Rinse mode parameters | 137 |
| Manual mode | 139 |
| Opening the error list | 141 |
| Resetting the device to factory settings | 141 |
| Opening system information | 142 |
| Updating the firmware and software | 142 |

Troubleshooting

| | |
|-----------------------|-----|
| Troubleshooting | 143 |
|-----------------------|-----|

Appendix

| | |
|--|-----|
| Technical specifications | 145 |
| Scope of delivery | 146 |
| Device service | 147 |
| General cleaning instructions | 147 |
| Repairs | 147 |
| Maintenance | 147 |
| Dismantling | 147 |
| Disposal | 148 |
| Heidolph International contact information | 148 |
| Warranty Statement | 148 |
| Certificate of decontamination | 149 |

About this document

These operating instructions describe the functions and operation of the automatic Hei-VOLUME Distimatic Pro module. The operating instructions are an integral part of the described device!

Typographic conventions

Standardized symbols, signal words and highlighting are used in this document to warn of hazards, and to identify important information and special text contents.

| Symbol | Signal word / explanation |
|---|--|
|  | <p>Warning symbols in combination with a signal word indicate dangers:</p> <p>DANGER Indicates an immediate dangerous situation which, if not avoided, will result in death or serious injury.</p> <p>WARNING Indicates a potential danger which, if not avoided, may result in serious injury.</p> <p>CAUTION Indicates a potential hazard which, if not avoided, can result in damage to property and minor to moderate injuries.</p> |
|  | <p>Mandatory signs indicate important and useful information on handling a product.</p> <p>This information is used to ensure operational safety and to maintain the value of the product.</p> |
| [GUI] | Parameter designations, display texts , and device labels are typographically highlighted in the text and in tables to facilitate their assignment on the device. |
| → | The arrow indicates specific instructions to be followed to ensure operational safety when handling the product. |

Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, and utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Instruments GmbH & Co. KG. Any violation is subject to compensation for damages.

Notes on the product

Directives applied, product certification

| | |
|---|---|
|  | <p>CE marking</p> <p>The device meets all requirements of the following directives:</p> <ul style="list-style-type: none"> • Low voltage Directive 2014/35/EU • EMC Directive 2014/30/EU |
|  | <p>NRTL certification</p> <p>The device has been tested in accordance with the following standards:</p> <ul style="list-style-type: none"> • UL 61010-1:2012/R:2018-11 CAN/CSA-C22.2 No. 61010-1:2012/A1:2018-11 • UL 61010-2-010:2019 CSA-C22.2 No. 61010-2-010:2019 • Supplemented by UL 61010-2-081:2019 • Supplemented by CSA-C22.2 No. 61010-2-081:2019 |

California Residents

Important information for California residents regarding Prop 65. Please visit www.P65Warnings.ca.gov for more information.

Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

Intended use

The automatic Hei-VOLUME Distimatic Pro module has been specially developed for the automatic filling and draining of Heidolph rotary evaporators. It is partly completed machinery as defined by the Machinery Directive 2006/42/EC.

The following processes can be controlled with the automatic module:

- Dosed addition of a starting material
- Drainage of distillate
- Drainage of residue (optional)

Reasonably foreseeable misuse

Additional measures may be necessary, and/or specific directives and safety regulations may have to be observed for use under conditions or for purposes deviating from the intended use (see, for example, section „Special hygiene measures for the use of laboratory equipment in food, cosmetics and pharmaceutical production“ on page 84). Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

Transportation

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device.

Keep the original packaging in a dry and protected place for later use.

Storage

Always store the device in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

Acclimatization

After each transportation and after storage under critical climatic conditions (e.g. high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

Permissible ambient conditions

The device is designed for indoor use only. Permissible ambient conditions for operation:

- 5 °C – 31 °C at up to 80 % rel. humidity
- 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
- Installation altitude: 2,000 m asl

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

The device **IS NOT** suitable for outdoor use!

The device **IS NOT** suitable for use in potentially explosive areas!

General safety instructions

- Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.
- Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself and, where applicable, on connected devices or the supply connections.
- If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.
- Only use the device in accordance with the regulations on intended use („Intended use“ on page 81).

Electrical safety

- Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.
- Ensure that the power socket is protected by a residual-current device (RCD).
- Always use the 3-pole power cord supplied with the device (phase, neutral, ground).
- Prior to each start-up, check that the device and the power cord are free from visible damage.
- Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Instruments.
- Always switch the device OFF and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

Data security

- The user is responsible for ensuring data security when transferring data between the described device and other devices.
- Only use secure networks for the data transfer and avoid use of critical infrastructure.
- Only use high-quality shielded data cables for the data transfer.
- For data transfer via a USB B connection, an industrial standard USB hub should be preferably used to ensure the most stable connection possible.

Operational safety

- Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (see EN 14175 and DIN 12924).
- Do not make any unauthorized changes or modifications to the device!
- Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!
- Rectify malfunctions or faults on the device immediately.
- Shut down the device and disconnect it from the power supply if it is not possible to eliminate the malfunction or rectify the fault immediately.
- Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

Overpressure operation

The distillate and residue are drained applying overpressure, i.e. the used distillate and residue containers are pressurized accordingly!

When using unsuitable distillate or residue containers, there is a risk that the containers will rupture or that liquid will escape uncontrolled and contaminate the environment.

There is also a high risk of injury, especially when processing toxic or corrosive substances!

Only use the distillate and residue containers recommended by Heidolph (see section „Scope of delivery“ on page 146) or others that are sufficiently pressure resistant and equipped with an adequately dimensioned pressure relief valve.

Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with the biohazard warning symbol.

The evaluation of corresponding measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Special hygiene measures for the use of laboratory equipment in food, cosmetics and pharmaceutical production

When laboratory equipment is used in the production processes of the food, cosmetics or pharmaceutical industry, special hygiene measures must be taken by the user to avoid sample contamination and to minimize any risk to humans and the environment as far as possible.

Please observe the following manufacturer's recommendations:

General measures

- Ensure a clean working and storage environment when handling substances and materials.
- Train all employees in the field of occupational hygiene, document all training measures and check the implementation of all required hygiene measures regularly during operation.
- Use a hygiene control concept such as HACCP (Hazard Analysis and Critical Control Points). The HACCP comprises the following criteria:
 - Hazard analysis
 - Identification of critical control points
 - Definition of critical limit values
 - Establishment of a system for monitoring and controlling critical control points (CCP)
 - Corrective actions for uncontrollable CCP
 - Establishment of a system to verify the implementation of all HACCP measures
 - Establishment of a system for documenting all associated procedures and records

The evaluation of the applicability of the mentioned rules and regulations is within the sole responsibility of the operator!

Device-specific measures

- Regularly clean components that come into contact with the product, such as flasks, seals, tubes, etc. in the autoclave (if available or possible) or chemically (e.g. with ethanol) to sterilize all surfaces.
- Make sure that even products that are intended for single use only are of sufficient purity.
- Do not use open containers.
- Avoid contamination caused by careless handling of contaminated vessels, apparatus or aids.

Contact information

For further information, please contact our after sales service at any time.
Email: sales@heidolph.de

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Instruments.

The operator is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unauthorized or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Assembly

If required, Heidolph offers a professional assembly service for the mechanical setup and installation of the device, including all accessories supplied and commissioning .

The majority of the sensors and valves of the automatic module are plugged into the coded connectors on the back of the Distimatic Pro (sensor cards).

Thanks to the coding of the connectors, all units are automatically detected and assigned in the device software!

All other units must be plugged into the correspondingly marked connection sockets. For a detailed description of all connections, please refer to the assembly instructions supplied.

Faulty or disconnected units are detected and displayed when performing a system test (integrated software function).

The assembly instructions supplied describe all steps for the correct connection of all tubes and cables between the Distimatic Pro, the components supplied (containers, sensors, valves, glassware) and the evaporation system in detail!

CAUTION

Due to improper system tubing and wiring as well as unauthorized changes to the system, there is a risk of direct and indirect damage to property!

For system tubing and wiring, follow all notes and instructions in the assembly instructions provided (part number 01-001-009-12).

Do not make any unauthorized changes to the device, system, tubing, nor the existing cable connections.

If necessary, contact our technical service (see „Heidolph International contact information“ on page 148).

The professional mounting and installation of the device including the supplied accessories is within the sole responsibility of the operator!

Noncompliance will invalidate any warranty against Heidolph Instruments.

The operator is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unauthorized or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Principle of operation

The automatic Hei-VOLUME Distimatic Pro module has been specially developed for use in combination with rotary evaporators. The automatic module enables the automatic filling and draining of the system in continuous operation.

The following processes are controlled with the automatic Hei-VOLUME Distimatic Pro module:

- Filling the evaporation flask
- Drainage of distillate
- Drainage of residue

The following operating modes are available to the user:

- **SENSOR MODE**
Standard mode for automatic filling and draining in continuous operation (see section „Sensor mode“ on page 120).
- **TIME MODE**
The time mode is used specifically for processing media that form sticky residues or that tend to foam (see section „Time mode“ on page 130).
- **RINSE MODE**
The rinse mode is used to flush the supply lines and drains with a solvent for cleaning purposes (see section „Rinse mode“ on page 135).
- **MANUAL MODE**
In manual mode, the connected rotary evaporator can be operated for processing small quantities (see section „Manual mode“ on page 139).

In automatic operation (sensor and time mode), the distillate and – depending on the variant (with or without automatic drainage of residue) – the residue is drained from the collector or the evaporation flask by applying overpressure. Distillate and residue can also be drained manually by the push of a button. The required overpressure is generated by a compressor integrated in the Distimatic Pro.

The evaporation flask is filled from the educt container using the vacuum in the rotary evaporator system. The required vacuum is generated by an external vacuum pump.

A pressure mode valve is used to switch between the compressor (overpressure operation) and the vacuum pump (underpressure operation).

Mechanical design

This section provides an overview of the mechanical design of the automatic Hei-VOLUME Distimatic Pro module. The Distimatic Pro Basis is housed in a sheet steel housing with protection class IP 42 with removable side covers and a removable control panel. The Distimatic Pro Basis contains the main circuit board, all function modules and the compressor (without illustration) for generating the overpressure to drain the residue and the distillate.

Control elements

CAUTION

Tip or sharp objects can damage the surface of the touchscreen.
To operate the device, touch the touchscreen only with your fingertips or a special stylus with a soft-touch tip.

Legend

- | | |
|---|---------------------------|
| 1 | On/off switch |
| 2 | Control panel, detachable |
| 3 | Touchscreen |
| 4 | Push-and-rotary controls |
| 5 | Standby button |

Connections rear side Distimatic Pro Basis

The following figure shows the rear side of the Distimatic Pro with standard plug-in cards and optional expansion cards. The actual assignment depends on the scope of the order!

The sensors and valves of the system are plugged into the coded connection sockets of the sensor cards. Thanks to the coding of the sockets, all units are automatically detected and assigned in the device software. In addition, observe the labels of the existing connection sockets, as well as all information on installing optional expansion cards in the supplied assembly instructions for each module!

Legend

| | |
|----|--|
| 1 | Valve and sensor terminals |
| 2 | Expansion cards |
| 3 | Sub-D connection chiller |
| 4 | Hei-VAP Control terminal |
| 5 | External cut-off |
| 6 | Overpressure connection (HVI sensor unit) |
| 7 | Aspiration connection: filter or inert gas |
| 8 | Measuring amplifier terminal level sensor S4 |
| 9 | IEC appliance inlet |
| 10 | Switchable built-in power socket |

Hei-VOLUME Distimatic Pro & Hei-VAP Industrial

Use of the automatic Hei-VOLUME Distimatic Pro module requires a special set of glassware consisting of a Woulff bottle [1], a collector vessel [2] and an educt sensor vessel [3].

The following figure shows the maximum configuration with Woulff bottle, collector vessel and educt sensor vessel, including all valves and sensors, completely mounted on the housing of a rotary evaporator:

Observe the instructions contained in the assembly instructions delivered (article no. 01-001-009-12), as well as the information given in chapter „Assembly“ on page 86.

For the exact scope of delivery, please refer to the order/delivery documentation.

Hei-VOLUME Distimatic Pro & Hei-VAP Expert/Ultimate

A Hei-VOLUME Distimatic Pro automatic module & Hei-VAP Expert or Hei-VAP Ultimate rotary evaporator system consists of different assemblies, which are delivered preassembled and must be put together and connected in the place of use as described in the following sections.

The following figures show the pre-assembled G9 XL condenser with level sensor S4 and emergency stop sensor S3 (left) and the valve matrix in standard configuration (right):

Observe the instructions contained in the assembly instructions delivered (article no. 01-001-009-12), as well as the information given in chapter „Assembly“ on page 86.

For the exact scope of delivery, please refer to the order/delivery documentation.

Heating tape (optional)

If necessary, the automatic Hei-VOLUME Distimatic Pro module can be operated with a heating tape (not included in the standard scope of delivery, see section „Scope of delivery“ on page 146, Heat control expansion card).


Principle of operation

The heating tape is wound around the tubing and preheats it during operation so that viscous media are passed reliably out of the system.

The heating tape is connected to the Distimatic Pro and is integrated in the system software via the corresponding heat control expansion card.

- As soon as the heating tape has been detected correctly, a corresponding icon appears on the software interface. Touch the associated icon to switch on/off the heating tape:

The operating behavior is defined by the process parameters of the selected mode (automatic switching on, on-period/pause time in percent).

PRACTICAL EXAMPLE

Period 100 %, defined heating tape on-period 30 %, pause time 70 %, for a 100-second cycle, equates to an on-period of 30 seconds per cycle.

Signaling of the operating states (optional)

The different operating states of the automatic module can also be made visible over large distances by an optional signal tower and the associated signal light expansion card (signal tower accessories set, article number 591-66000-00).

The installation of the additional card is described in detail in the enclosed sheet supplied.


The following operating states are indicated:

| Color code/state | Meaning |
|-------------------------|--|
| flashing green | Process is finished |
| continuously lit green | Process is in progress |
| flashing yellow | Educt is empty or collection container is full |
| continuously lit yellow | Note the warning on the display |
| flashing red | Emergency stop has been triggered |
| continuously lit red | Device in fault condition, note the error message on the display |

Structure of the user interface

All device functions are controlled via the integrated control panel and the graphical user interface.

After switching on the device and initializing the device control, the Start window of the last active operating mode appears on the display (in the example: Sensor mode):

From the start window, you can access the main menu, select an operating mode directly, and select the time display.

The function keys in the footer can be used to drain the distillate, to drain the residue, and to start the process directly.

Control and display elements

| Control element | Function |
|---|---|
| Main menu function button | <p>Touch this function button to open the main menu with the following options:</p> <ul style="list-style-type: none"> ▪ Applications → operating mode ▪ Settings → system parameters ▪ Languages → menu language ▪ Formats → system configuration and display formats ▪ Data logging → data recording ▪ System check → complete system check ▪ Error list → open error list ▪ Info → open system info window |
| Function button Mode | <p>Touch this function button to open the Applications sub-menu directly. Operating modes:</p> <ul style="list-style-type: none"> ▪ Sensor mode ▪ Time mode ▪ Rinse mode ▪ Manual mode |
| Function button System time/Timer | <p>Touch this function button to toggle between the following functions:</p> <ul style="list-style-type: none"> ▪ System time display ▪ Timer ▪ Stopwatch <p>See the following section „System time, timer, stopwatch“ on page 97 for detailed information.</p> |
| Drain distillate button | <p>Touch this button to drain the distillate from the collector → only when no process is active!</p> <p>By touching this button, valves V2 and V4 are opened, see section „Principle of operation“ on page 87.</p> |
| Drain residue button | <p>Touch this button to drain the residue from the evaporation flask → only when no process is active!</p> <p>By touching this button, valve V3 is opened, see section „Principle of operation“ on page 87.</p> |
| Start process button | Touch this button to start the selected process. |
| Process visualization main display | The current process step is visualized in the main display. |


| | |
|-------------------------------|---|
| Left push-and-rotary control | <p>In open lists, you can switch between entries by turning the left push-and-rotary control. A selected value (highlighted in orange) can be accepted by pressing the push-and-rotary control (similar to touching the entry on the display).</p> <p>The color of the LED ring indicates process status:</p> <ul style="list-style-type: none">▪ green → overpressure▪ blue → underpressure (vacuum)▪ white → ventilate & atmospheric pressure |
| Right push-and-rotary control | <p>Turn the right push-and-rotary control to switch between the display screens (in the case of multi-page displays, similar to the wiping movement on the display).</p> <p>When the on-screen keyboard is open, the display value can be increased or decreased by turning the right push-and-rotary control. The step width (1, 25 or 50) varies depending on the rotation speed of the push-and-rotary control.</p> <p>Color code LED ring:</p> <ul style="list-style-type: none">▪ LED ring lights up white when the on-screen keyboard is open▪ LED ring flashes red when the system is in a safe operating condition <p>Press and hold both push-and-rotary controls simultaneously for approximately three seconds to lock the display.</p> |
| Standby button | <p>Press the standby button to put the device in standby mode. Press the button again to put the device into the previously active operating mode.</p> <p>Pressing the standby button will stop the current process!</p> |

System time, timer, stopwatch

System time/timer/stopwatch selection

Briefly touch the **System time/Timer** function button in the header of the start page to open the selection menu (see following figure) with the following options:

- Display of current time (12 or 24-hour display, depending on the selection in the **Settings** menu).
- Stopwatch: In this mode, the actual time elapsed since the start of the process is displayed.
- Timer: In this mode, the time remaining until the scheduled end of the process is displayed.

This setting is applied in all operating modes and can be changed at any time!

Set time

- Switch to **Time** mode.
- Touch the **System time/Timer** function button for approx. two seconds to open the on-screen keyboard (see section „User input“ on page 99) and correct the system time (date, time).

Program timer

- Switch to **Timer** mode.
- Touch the **System time/Timer** function button for approx. two seconds to open the on-screen keyboard (see section „User input“ on page 99) and to program the timer.


Color coding

In the process visualization area (see section „Structure of the user interface“ on page 94), all process lines and function elements are color-coded to illustrate the operating states:

| Color | Meaning |
|------------|---|
| Gray | PASSIVE or off state (e.g. non-active/closed line, non-switched valve, sensor low) |
| Orange | ACTIVE or on state (e.g. active/open line, switched valve, sensor high) |
| White | VENTILATED condition, atmospheric pressure present |
| Light blue | VACUUM condition, vacuum pump active, evacuation |
| Green | OVERPRESSURE condition, compressor on, draining |



Symbols

In the process visualization (see section „Structure of the user interface“ on page 94), the following symbols are used to illustrate the process sequences and the integrated function elements (for color coding, see section above):

| Symbol | Meaning |
|---|---------------------|
|  | Compressor |
|  | Vacuum pump |
|  | Pressure mode valve |
|  | 2-way valve |
|  | Sensor |
|  | Heating tape |


Input keys

The following input keys are available in any menu:

| Symbol | Meaning |
|---|---|
|  | Delete button. Press this button to correct a numeric or alphanumeric value entered. |
|  | ok / confirm button. Press this button to confirm and accept a selection or input value. |
|  | Discard button. Press this button to discard a selection or input value and return to the previous view. |

User input

User input for configuration, parameterization and programming is carried out via the numerical and the alphanumeric on-screen keyboard. The on-screen keyboard appears automatically as soon as an input field is touched:

Commissioning

Observe the following steps when commissioning the automatic Hei-VOLUME Distimatic Pro module.

Switch the device on/off

To switch on and off the automatic module, use the main switch on the front of the Distimatic Pro (see section „Control elements“ on page 88).

CAUTION

If the Distimatic Pro is switched on and off incorrectly, there is a risk that peripheral devices will not be detected correctly or will switch in an uncontrolled manner!

Observe the order given when switching the Distimatic Pro and the connected peripheral devices on and off.

In particular, before turning off the Distimatic Pro, make sure that all processes have been completed properly.

Observe the following sequence when turning the automatic module and attached peripherals on and off:

- **TURNING ON THE SYSTEM:** First, turn on all peripheral devices so that they are recognized when initializing the automatic module software.
- **TURNING OFF THE SYSTEM:** First, turn off all peripherals to put the entire system in a safe operating condition.

Calibrate the flask level sensor

The flask level sensor (S4) must be calibrated during the initial commissioning and every time the solvent is changed.

WARNING

In case of inadequate calibration or failure of a level sensor (S4), overflow/leaking of the flask contents is prevented by an emergency stop sensor (S3)!

As soon as this sensor responds in the cooler, the Distimatic Pro and the connected peripheral devices are switched to a safe operating state and no further medium is fed in.

In extreme cases, there is a danger of the emergency stop sensor (S3) not clearly detecting the overflow of very low polar media such as hexane or heptane. If you change to such a medium, contact our application specialists if necessary (see „Heidolph International contact information“ on page 148).

CAUTION

The function of the level sensor in the evaporation flask (S4) may be affected by electromagnetic fields. This can lead to a loss of the operational quality of the entire system.

During calibration, commissioning and operation, ensure that there are no sources of interference in the immediate vicinity of the evaporation flask to ensure proper operation of the level sensor at all times.


To calibrate the flask level sensor, follow these steps:

- Make sure that the sensor is connected to the Distimatic Pro correctly.
- Switch on the Distimatic Pro and wait approximately 10 minutes for the electronics to warm up to operating temperature. In this way, measurement inaccuracies are avoided during calibration!

Calibration to maximum sensitivity, low-polar solvents:

- Remove the sensor from the evaporation flask.
- Hold the sensor in the open space.
- Make sure there are no objects or bodies in a hemispherical space one (1) meter in front of the sensor head. These would be detected as a signal and falsify the calibration result!

- Start the calibration process in the device software:
 - Open the main menu using the **Main menu** function button in the header of the start window:


- Open the **Settings** menu and select **Calibrate filling level** (see also „Main menu“ on page 105).
- Follow the instructions on the display.

OR

Calibration on high-polar solvents and mixtures:

- Make sure the sensor is in the flask.
- Fill the evaporation flask with at least 500 ml of the product.
- Slide the sensor into the stopped flask until the sensor head is just above the fluid level or – for example, in the case of mixtures with water – until the sensor head touches the fluid level slightly.
- Start the calibration process in the device software:
 - Open the main menu using the **Main menu** function button in the header of the start window:

- Open the **Settings** menu and select **Calibrate filling level** (see also „Main menu“ on page 105).
- Follow the instructions on the display.


PRACTICAL TIP

When processing solvents with a tendency toward boiling retardation (e.g. ethanol), it is recommended to calibrate the sensor slightly less sensitively.

This prevents the boiling retardation from being detected as a fluid and the fluid level in the rotary flask from decreasing excessively before adding new medium.

The sensor should only respond to liquids, never to boiling retardation or foam!


- Start the calibration process with **ok**.
- The calibration process can be canceled at any time with **discard**.
- The progress is visualized by the traffic light display. When the calibration is complete, the green dot on the right lights up.
- Confirm the completion of the calibration process with **ok**.

Initial filling of the evaporation flask

Complete the following steps for filling the evaporation flask to ensure continuous supply of the educt during operation:


- Turn on the rotary evaporator and wait for the software to initialize (start screen appears on the control panel).
- Switch on the heating bath heater and the recirculation chiller.
- Turn on Distimatic Pro and wait for the software to initialize (start screen appears on the control panel).
- Start any given process by touching the **Start process** button on the control panel of the Distimatic Pro.
- Open the inlet valve on the educt sensor vessel.


- Close the inlet valve when the displayed level is reached.

Device configuration


The basic settings for operation are made in the main menu. After switching on the Distimatic Pro and initializing the software, the start page of the last active mode appears. To open the main menu, touch the **Main menu** function button in the header of the display:

Main menu


The following options are available from the main menu:

| Button | Function/Meaning |
|---|--|
|  | <p>Home button. Touch this function button to return to the previously opened view (start window).</p> |
| Operating modes | <p>Touch this function button to open the Applications menu and select an operating mode, see section „Applications menu“ on page 107.</p> |
| Settings | <p>Touch this function button to open the Settings menu. Various basic settings can be adjusted here and the automatic module can be reset to the factory settings, see section „Settings menu“ on page 108.</p> |
| Languages | <p>Touch this function button to open the Languages menu and load a desired user language, see section „Selecting the user language“ on page 110.</p> |
| Formats | <p>Touch this function button to open the Formats menu. Here you can select the desired format for the display of numerical values, time, date and temperature. At the same time, the connected rotary evaporator / glassware is assigned and the operation with or without residue drainage is selected, see section „Select formats“ on page 111.</p> |
| Data logging | <p>Touch this function button to start the Data logging function, see section „Data recording“ on page 113.</p> |
| Hei-CONTROL Pro | <p>Touch this function button to open the screen form for definition of the network parameters, see section „Connectivity“ on page 115.</p> |
| System check | <p>Touch this function button to perform a full system check, see section „System check“ on page 116.</p> |
| Error list | <p>Touch this function button to open the list of all error messages, see section „Opening the error list“ on page 141.</p> |
| Info | <p>Touch this function button to open the information window of the system, including firmware and hardware data, see section „Opening system information“ on page 142.</p> |

Applications menu

The desired operating mode can be selected in the **Applications** menu.

To select and activate an operating mode, proceed as follows:


- Complete all running processes before changing the operating mode.
- Touch the entry of the desired operating mode. The background color of the selected entry changes to blue, as shown.
- Press **ok** to confirm the selection and activate the selected mode.
 - The start window of the activated mode appears on the display.
- Press **discard** to discard the selection and to return to the previous view.

The **Applications** menu can also be opened directly by touching the **Mode** function button in the header of the start window.

Settings menu

In the **Settings** menu, various device settings can be adjusted:

The various options are described in detail in the following sections.

Screen lock

The **Automatic screen lock** menu item can be used to activate the automatic screen lock for the display (**active**) or to deactivate it (**inactive**).

When the screen lock is active, the display locks after a certain period of inactivity (no operator actions on the control panel).

The process progress window is displayed in the locked screen. As soon as the display is touched, the unlock slider appears.

The screen lock prevents accidental operation during operation, but does not provide protection against unauthorized access.

To protect the device from unauthorized access, the **PIN for Screen lock** option must be enabled!

Unlocking the screen

To unlock the screen, touch any point on the screen and swipe the slider from left to right:

Time until automatic screen lock is activated

Using the option **Automatical lock after ... min** you can specify after which period of time the automatic screen lock should be activated.

- Touch the function button to open the numeric on-screen keyboard and set the desired time interval in minutes.
- Confirm the input with **ok**.
- Touch the **discard** button to discard the entry and return to the previous menu without making any changes.

From a continuous on-time of approx. 60 minutes, there is a risk of the display becoming burnt in and static display elements remaining visible as a visible pattern temporarily, even after switching off. It is therefore advisable to set a maximum time of 60 minutes before the screen lock is activated.

Screen lock PIN

The **PIN for Screen lock** option allows you to set a security code that must be entered to unlock the screen. This feature can be used to protect the device from unauthorized access.

- Touch the **PIN for Screen lock** function button to open the on-screen numeric keypad.
- Enter a four-digit PIN and press **ok** to confirm the entry.
- Re-enter the PIN to confirm (confirmation prompt!) and end with **ok** to activate the security code → **PIN for screen lock active**.
- Touch the **PIN for Screen lock** function button again to disable the security code → **PIN for screen lock inactive**.

The security code is automatically disabled when the screen lock is disabled.

Display brightness

The brightness of the display can be adjusted via the menu item **Brightness (Main menu → Settings)**.

- Touch the **Brightness** function button to open the on-screen numeric keyboard.
- Enter a percentage value between 10 and 100 (in increments of ten) and confirm the input with **ok**.
- Touch the **discard** button to discard the selection and return to the **Settings** menu without making any changes.

Calibrate filling level

The flask level sensor S4 can be calibrated via the menu item **Calibrate filling level (Main menu → Settings)**, see section „Calibrate the flask level sensor“ on page 101.

Maximum pressure

Via the **Max. pressure** menu item, the maximum permissible absolute target pressure (corresponding to atmospheric pressure + compressor pressure) can be set within a range of 1,100 –1,400 mbar.

The higher the target pressure, the better the performance of the system. An adjustment may also be necessary if larger distances or heights (e.g. tanks) have to be overcome.

Please note that exceeding a target pressure of 1,400 mbar during operation activates the overpressure shutdown of the evaporator and the system is ventilated.

Sensitivity cut-out sensor S3

Via the **Sensitivity cut-out sensor S3** menu item, the sensitivity of the emergency stop sensor can be adjusted in a range of 0 –100 % in steps of [1 %].

The **Sensitivity cut-out sensor S3** parameter may only be changed by authorized technical personnel during installation or after consultation with the Heidolph application engineering department.

The standard setting covers most applications, an adjustment is only necessary for extremely low- or high-polarity media.

Factory settings

The device can be reset to the factory settings via the menu item **Factory settings (Main menu → Settings)**, see section „Resetting the device to factory settings“ on page 141.

Selecting the user language

In the **Languages** menu (**Main menu → Languages**), the user language can be changed at any time. Available languages are English, German, Mandarin, French, Spanish, Portuguese, Italian, Russian and Polish.

- Touch the **Languages** function button.
- The list of all available menu languages opens.
- Touch the entry for the desired language.
- The background color of the selected entry changes to orange.
- Press **ok** to confirm the selection and to return to the Start window of the active mode.
- Touch the **discard** button to discard the entry and to return to the start window of the active mode without making any changes.

Select formats

The following system settings can be made in the **Formats** menu (**Main menu** → **Formats**):

- Set date format (US or EU)
- Set time format (12/24-hour)
- Set number separator (US or EU)
- Set temperature unit (degrees Celsius, Fahrenheit, Kelvin)

- The active selection is highlighted in orange. Touch the desired function button to switch between options. Several changes can be made at the same time!
- Confirm the new settings with **ok**. All changes are immediately accepted without a confirmation prompt!
- To discard any changes and to return to the previous menu, touch the **discard** button.

Automatic drainage of residue

In the **Formats** menu, you can define whether the residue is automatically drained from the evaporation flask (blowing out in overpressure operation), or whether the flask is drained by hand.

A system **with automatic drainage of residue**, in sensor or time mode can be operated automatically in continuous operation (condition: continuous supply with educt!).

Operation with automatic drainage of residue

In operation **with automatic drainage of residue**, the residue is automatically drained under the following conditions:

- Stagnation reached in sensor mode.
- Distillate drain counter expired in sensor mode.
- Refill counter expired in sensor mode.
- Additional evaporation time for distillate drain or refill counter has expired in sensor mode.
- Number of defined steps (filling/evaporation) reached in time mode.

Operation without automatic drainage of residue

In operation **without automatic drainage of residue**, the evaporation flask must be drained by hand.

Drainage of residue selection

- The active selection is highlighted in orange when the **Formats** menu is opened. To change, touch the desired function button **with / without automatic Residue Extraction**.
- Confirm the new setting with **ok**. The change is accepted immediately without a confirmation prompt!
- Touch the **discard** button to discard the change.

Select evaporator type/glassware

In the **Formats** menu, you can define the connected evaporator type including the glassware used. The selection stored in the system appears in the form of function buttons:

- The active selection (in the example: Hei-VAP Industrial with set of glassware R) is highlighted in orange. Touch the desired function button to switch between options.
- Confirm the new setting with **ok**. The change is accepted immediately without a confirmation prompt!
- Touch the **discard** button to discard the change.


The selection evaporator type/glassware is necessary for display purposes only!

Data recording

The **Data logging** function can be used to record the individual steps of a process. The function records each individual step with start and end time and stores this information in the Data logging memory. The function must be activated before the process starts:


- To activate data recording, switch from the start/overview page of the active mode to the parameters page.
- Touch the **Data logging** button in the footer:

- When the data logging function is activated, the red recording dot appears in the function button and in the process progress display.

The recorded values can be read out, saved and deleted using the **Data logging** option


- Open the **Main menu** and touch the function button **Data logging**.
- The Data logging list opens:

- The recorded values show when a process was started and stopped again.

In the event of a power interruption or when the device is switched off during the current data recording, the recording process is aborted and the corresponding entry marked with "(!)". The displayed data set contains the recorded process data that was recorded up to the power interruption:

In the data log export file, the corresponding file names are also given the addition "(!)" before the file extension: "(!).csv", or "(!).json".

Deleting entries

To delete an entry from the list, touch the associated Delete button.

Saving entries on a USB stick

Recorded entries can be saved on an USB-stick. To do this, connect a VFAT-formatted USB stick to the USB port of the display and open the data list as described.

As soon as a USB stick is connected, a **Save** button (diskette symbol) is assigned to each list entry.

Touch the associated **Save** button to save an entry to the attached USB stick.

PRACTICAL TIP

Large amounts of data are generated during the long-term recording of process data (several days).

It is therefore recommended to save the current data recording memory on a USB stick and to remove the existing entries from the internal memory before starting the process.

Connectivity

All necessary network parameters and basic settings for time synchronization can be defined in the **Hei-CONTROL Pro** menu:

- The connection status of the Hei-CONTROL Pro application is displayed in the [Status] field.
- All associated parameters can be defined manually in [static IP configuration] mode (checkbox deselected, figure on the left). Select the checkbox (figure on the right) to obtain the IP address of the device automatically.
- Select the [NTP server] checkbox to synchronize the system time of the device with a time server. This function requires a valid IP configuration!
 - Place the cursor in the corresponding input field and use the on-screen keyboard to enter the required server address.

If the time is synchronized via a time server, the **System time/Timer** function button in the display header can only be used to select a time zone.


If an ethernet connection exists, the system time is updated automatically via the defined time server. To define the system time manually, the [NTP server] option must be deselected!

The full functionality of the Hei-CONTROL Pro application is not available in the current device version.

System check

All connected sensors and actuators or peripherals are displayed in the **System check** menu. Various tests can be carried out via the **System check** sub-menu to ensure the functional and operational safety of the system:

- Open the **Main menu** and touch the **System check** function button in the footer. The **System check** menu appears:

- **Error-free** components are marked by a green checkmark
- **Faulty** components are marked by a red X Please note the corresponding error messages.
- Swipe up or down on the touch panel to scroll through the list.
- Touch the **TEST** function button in the header of the **System check** menu to open the **System check** sub-menu.

CAUTION

Faulty connections or signaling can lead to malfunctions or uncontrolled switching operations!

Leaks reduce the performance of the entire system.

Correct any indicated errors and/or detected malfunctions immediately.


In the event of a continued fault indication/malfunction, contact your local sales representative or our technical service department, see „Heidolph International contact information“ on page 148.

Testing valves and signaling

The **Test all actuators** function can be used to test the function of all connected valves and signal lights.

→ Open the **System check** menu:


Test all valves and signal lights

- In the footer of the **System check** overview, touch the **Test all actuators** function button to activate the test function (the button icon changes from gray to orange).
- Start the test with **Start test**:
 - During the test, all valves and the signal lights are switched twice in succession (note the color change of the corresponding symbol!). At the end of the test, all switching elements are back in their initial state.
- The test can be stopped at any time by **Stop test**.

Perform overpressure test

With this test function, the system can be tested for leaks.

- Open the **System check** menu.
- Touch the **Overpressure** function button in the footer to activate the test function (button changes to orange).

WARNING

Due to faulty components, faulty assembly of components or faulty connections, there is a risk that pressurized components will suddenly loosen or rupture during an overpressure test!

Before performing an overpressure test, ensure that all components are properly installed. In particular, make sure that all screw connections are properly tightened.

Before starting the overpressure test, close all safety covers, hoods and doors of the entire system.

Ensure that third parties are clear of the danger zone.

Use the required personal protective equipment (eye protection, hand protection, etc.).

- Touch the **Start test** function button.
- Confirm the subsequent confirmation prompt with **ok** to carry out the overpressure test:
 - During the test, the system pressure is built up to the defined maximum pressure (**Settings** menu). This pressure level must be reached within a maximum of 30 seconds!
 - When the maximum pressure is reached, a measurement is made for a period of 240 seconds to determine whether the system pressure drops below a defined limit value.
 - After this time has elapsed, the system message **Overpressure test successfully completed** is shown on the display, provided that all required values have been achieved.
- The **Overpressure test failed** error message appears if the defined minimum pressure is not reached during the overpressure test or the system pressure falls below a defined limit value. In this case, follow the instructions on the display.
- The test can be stopped at any time by **Stop test**.
- Close the menu with **ok**

Testing the safety switch

The function of the integrated safety switch can be tested with the test function **Safety switch**.

The safety switch is activated when a defined system overpressure is reached during compressor operation. When the safety switch is activated, the compressor is shut down so that the system pressure can be released.

- Open the **System check** menu.
- In the footer, touch the function button **Safety switch** to activate the test function (button changes to orange).
- Touch the **Start test** function button.
- Confirm the subsequent confirmation prompt with **ok** to perform the safety switch test:
 - During the test, the compressor is used to build up pressure in the system. After a certain time has elapsed, the pressure mode valve is automatically switched to vacuum mode by software and the compressor line is blocked abruptly.
- A corresponding error message is displayed if the safety switch does not trigger.
- If the safety switch triggers correctly, **The safety switch works properly** appears.

If the safety switch triggers correctly, the standard error message for switching off the system after the safety switch has been triggered also appears!

Follow the instructions on the display.

- The test can be stopped at any time by **Stop test**.
- Close the menu with **ok**.


Selecting the operating mode

After switching on the device and initializing the device control, the start window of the last active operating mode appears on the display. Touch the **Mode** function button to open the **Applications** menu and to activate the desired mode (see section „Applications menu“ on page 107).

Sensor mode

Activating sensor mode

- Open the menu **Applications** and touch the function button **Sensor mode**. The entry is highlighted in orange.
- Confirm the selection with **ok** to activate the mode and to open the sensor mode start window:

In sensor mode, the filling quantity in the evaporation flask is regulated by the flask level sensor S4. The following sensor statuses are possible:

- The filling process is started and runs as long as the sensor does not have contact with the medium → S4 in the status **low**.
- The filling process is stopped when the sensor comes into contact with the medium → S4 in the state **high**.

Preconditions

The following conditions must be met for operation in sensor mode:

- Heating bath heater of the rotary evaporator is switched on
- Recirculation chiller is switched on
- Vacuum pump is switched on
- Rotary evaporator is under filling vacuum
- Level sensor S4 calibrated
- Educt sensor vessel is filled


Filling the evaporation flask

After the initial filling of the evaporation flask, the filled educt is brought to process temperature.

For the initial filling, it is essential to observe the instructions in section „Initial filling of the evaporation flask“ on page 104.

- Touch the **Start process** button in the footer of the start window. The following message appears on the display:


- The filled educt is saved within a preset time (factory setting: 300 s) to the heating bath temperature.
- Observe the progress bar in the message box.
 - Press **ok** to skip the heating process.

Note that skipping the heating process will automatically start the evaporation process. This activates the refill, the emergency stop and the educt sensor!

- After the pre-heating time has elapsed or after the step has been skipped, the filling pressure is automatically regulated to the set process pressure. The following message appears on the display:

- Start automatic mode with **ok** as soon as the evaporation process begins (distillate precipitates in the condenser).

Automatic mode

The process time starts when automatic mode is started.

- As the evaporation starts, the filling level in the evaporation flask decreases:
 - The flask level sensor S4 loses contact with the medium and changes from **high** to **low**. With this status change, the **filling time** is started (→ **low** = start filling time).
- As soon as the defined setpoint for the filling time has been reached (**filling time** expired), the evaporation flask is automatically filled with new educt.
 - The **stagnation time** starts when the level sensor S4 changes from low to high during filling (→ **high**= start stagnation time).
- The filling level in the evaporation flask rises again. As soon as the level sensor S4 detects media contact, the status changes again from **low** to **high**(→ **high** = start stagnation time).
- The continuous evaporation process reduces the filling level in the evaporation flask again:
 - The level sensor S4 loses contact with the medium and changes from **high** to **low** (→ **low** = start filling time).
 - When the status changes again, the **filling time** is started and the **stagnation time** is set to zero.
 - As soon as the defined **filling time** has been reached again, the evaporation flask is filled with new educt.
- The filling level in the evaporation flask rises again and the described process steps are repeated until the end of the process.

The **filling time** defines the time until the next filling of the evaporation flask and thus the amount of educt that is evaporated until the next refilling.

The **stagnation time** runs as soon as the level sensor S4 changes to the high status after filling . If the stagnation time expires completely, the automatic residue emptying is initiated (depending on the set of glassware and the parameterization) or the flask must be emptied manually.

End of process in sensor mode

In sensor mode, the process sequence is stopped under one of the following conditions:

- Educt container runs empty
 - After the **filling time** has elapsed, the evaporation flask is no longer filled → Status change of the feed sensor S1 (**high to low**).
 - Unplanned end is initiated: Additional evaporation time is started (setting via process parameters, see section „Parameter sensor mode“ on page 124).
 - After the additional evaporation time has elapsed, the system switches to a safe operating state: Peripheral shutdown (recirculation chiller and vacuum pump off, evaporator in standby).
- Maximum filling rate reached
 - Using the process parameter **max. filling: Additional Evaporation time** a maximum filling rate can be set in order to achieve a defined concentration of the residue (setting via process parameters, see section „Parameter sensor mode“ on page 124).
- Maximum number of distillate drainage steps reached
 - Use the process parameter **to empty max. distillate: Additional Evaporation time**, a maximum number of distillate drainages can be specified in order to achieve a defined concentration of the residue (setting via process parameters, see section „Parameter sensor mode“ on page 124).
- Residue level in the evaporation flask reached, flask level sensor S4 remains permanently in the status **high** (medium contact), nominal **stagnation time** is reached.
 - In this case, the sensor control is stopped and the residue is drained automatically or manually. The draining process depends on the glassware (with/without drainage of residue) and the **Drain residue on/off** default in the parameters screen form (setting via process parameters, see section „Parameter sensor mode“ on page 124).
- Maximum filling level in distillate/residue container → Level sensor S8 (distillate container) or S9 (residue container) indicate by a status change from **low to high** that the maximum level has been reached (option).


Both of the level sensors S8 (distillate container) and S9 (residue container) are available as optional accessories and are not included in the standard scope of delivery!

The level sensors prevent overflow of the monitored vessel. When a level sensor is triggered, the system is switched to a safe operating state.


Parameter sensor mode

General parameters



- To adjust the **general parameters**, switch from the start page of sensor mode to the corresponding parameters screen form (tap the parameters overview or swipe from right to left):

- Use the two navigation keys to return to the start window or to continue to the **Recipes** view:

The following parameters can be adjusted:

| Parameters | Function |
|-------------------------|---|
| Cycles | <p>Maximum number of evaporation cycles</p> <ul style="list-style-type: none"> Selection switch in position  → no limit. Selection switch in  position → Define a certain number of cycles in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard). <p>Touch the  /  icon to switch between the two options.</p> |
| Fill | <p>Maximum filling rate</p> <ul style="list-style-type: none"> Selection switch in position  → no limit. Selection switch in  position → Define a maximum number of filling operations in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard). <p>Touch the  /  icon to switch between the two options.</p> |
| Drain residue | <p>Residue draining time</p> <ul style="list-style-type: none"> Selection switch in position off → no defined draining time. Selection switch in on position → Define a residue draining time in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard). <p>Touch the on / off icon to switch between the two options.</p> |
| Drain distillate | <p>Distillate draining time</p> <ul style="list-style-type: none"> In the selection field, define a distillate draining time (adjustment via plus/minus buttons or tap/on-screen keyboard). |
| Distillate | <p>Maximum number of distillate drainage steps</p> <ul style="list-style-type: none"> Selection switch in  position → no limit. Selection switch in  position → Define a maximum number of distillate drainage steps in the selection window (adjustment via plus/minus buttons or tap/on-screen keyboard). <p>Touch the  /  icon to switch between the two options.</p> |
| End of process | <p>Specify here whether the distillate and/or residue are drained at the end of the process. Touch the corresponding checkbox to activate (checkmark) or deactivate the selection.</p> |

The **Fill** and **Distillate** parameters are mutually exclusive!
 For an evaporation cycle, you can specify either a maximum filling rate **OR** a maximum number of distillate drainage steps.


→ Confirm or discard all changes with **Accept** or **Cancel**.

Process parameters

- Touch the **Process parameters** button in the footer to open the menu.
- Tap the button of the desired parameter in the menu and adjust the value using the on-screen keyboard, see also „User input“ on page 99.

| Parameters | Function |
|---|---|
| Heating tape | On-period of the heating tape in operation. Range: 10 – 100% in increments of ten. (see „Heating tape (optional)“ on page 92) |
| Heating tape automatic | Automatically switch on heating tape yes/no. Range: on / off. |
| Evacuation after draining the collector | Evacuation period of the bypass line after draining the collector vessel. Range: 1 – 9999 s. |
| Warm-up time, filling in sensor mode | Warm-up time when filling the evaporation flask in sensor mode. Range: 1 – 9999 s. |
| Educt empty: additional evaporation time | When the educt container runs empty: Additional evaporation time for evaporation of the remaining educt in the evaporation flask. Range: 1 – 9999 s. |
| | Maximum filling rate: Selectable when a specific number of filling operations is specified in the start window: |

Max. filling: additional evaporation time


In the example: five fillings, then optionally an additional evaporation time runs. This allows a desired concentration to be achieved!

The parameters **max. filling: Add. evaporation time** and **max. drain distillate: additional evaporation time** are mutually exclusive and cannot be activated at the same time!

Selectable if a specific number of distillation drainage steps is specified in the start window:

Max. drain distillate: Additional evaporation time

In the example: drain distillate six times, then the optional additional evaporation time runs. This allows a desired concentration to be achieved!

The parameters **max. filling: Add. evaporation time** and **max. drain distillate: additional evaporation time** are mutually exclusive and cannot be activated at the same time!

Forced draining of the collector

Forced draining of the collector active: Collector will be drained after 1 – 9999 seconds, regardless of the filling level.

Range: 1 – 9999 s.

Activate/deactivate forced draining of the collector (determines previous parameter!).

Range: active / inactive

Standby timer after educt container runs empty

Time until the standby mode is activated after the educt container runs empty.

Range: 1 – 9999 s.

Parameters Hei-VAP

In the sub-menu **Parameters Hei-VAP** various parameters of the connected Hei-VAP can be stored (vacuum, hysteresis, rotational speed, heating bath temperature, recirculation chiller temperature).

These values are for informational purposes only and may be included in a new recipe.

Touch the **Hei-VAP parameters** button in the footer of the Parameter screen form to open the submenu and to enter the required values for the connected rotary evaporator by hand.


Data logging

See section „Data recording“ on page 113.

Recipes

By touching the **Save** button, the actual settings and values can be stored in the system as a recipe.

→ Touch the **Save** button to open the on-screen keyboard (see also section „User input“ on page 99):

→ Use the on-screen keyboard to enter any recipe name.

It is not possible to assign recipe names twice! If a desired recipe name is already stored in the system, a corresponding warning message appears. In this case, you can enter a different name or update the existing recipe with the new values.

→ Confirm your entry with **ok** to store the new recipe or the new settings in the recipe list.

Recipe list

To open the recipe list, switch from the start window of sensor mode to the corresponding **Recipes** view (two swipes from right to left, or use the navigation button from the parameters window).

Here, new recipes can be created (**New** button) and existing recipes can be loaded (**Load** button) or deleted (**Delete** button).


The entire recipe list can be saved on a USB stick. Recipe lists can also be loaded from a USB stick to the device memory.

To do this, connect a USB stick to the USB port on the control panel and follow the instructions on the display.

Time mode

Activating time mode

- Open the menu **Applications** and touch the function button **Time mode**. The entry is highlighted in orange.
- Confirm the selection with **ok** to activate the mode and open the time mode start window:

In the time mode, the filling of the evaporation flask, the distillation and the draining of the distillate and residue are carried out in a purely time-controlled manner:

- Within a cycle, a desired amount of concentrate is built up as a residue in the evaporation flask in the time mode.
- A cycle can consist of 1 – 9999 steps.
- One step consists of the **filling** and **evaporation** processes.
- In each step, the desired end concentration is reached.

Preconditions

The following conditions must be met for operation in time mode:

- Heating bath heater of the rotary evaporator is switched on
- Recirculation chiller is switched on
- Vacuum pump is switched on
- Rotary evaporator is under filling vacuum
- Process parameters are set correctly on the rotary evaporator
- Level sensor S4 calibrated

Before the initial filling of the evaporation flask or before the first process step, it must be ensured in time mode that all process parameters (vacuum, temperature, speed, cooling) are set correctly on the rotary evaporator.

Otherwise, there may be inaccuracies in the recipe processing!

Procedure

- Initial filling in time mode: This step is only necessary to prevent drying or crystallization of the educt. Activate with the **Initial filling** in the parameters screen form, see following section „Time mode parameters“ on page 132. For the initial filling, it is essential to observe the instructions in section „Initial filling of the evaporation flask“ on page 104.
- In time mode, a defined quantity is filled into the evaporation flask and immediately vaporized within a specified time:
 - The fill quantity is defined by the **filling time**.
 - The **evaporation time** is defined by the user.
- Due to the **pre-evacuation** of the evaporation flask, the evaporation process starts immediately while filling. After processing a step, the evaporation flask is filled again and the following step is processed:
 - The number of steps in a cycle is parameterized by the user. This value depends on the desired target quantity, which is to form in layers in the evaporation flask.
- After processing the specified number of steps, a cycle is run through and the residue is pressed out of the evaporation flask by means of overpressure.
 - The number of cycles is parameterized by the user.
- At the end of the last cycle, all connected devices are switched off (operating safe state, end of process).


In time mode, the system containers are monitored by sensors! The end process can be initiated and/or an ongoing process can be terminated and the evaporator can be put into a safe operating state when a corresponding signal is emitted by a sensor.


Time mode parameters

General parameters

- To adjust the **general parameters**, switch from the time mode start page to the corresponding parameters screen form (tap parameters overview or swipe from right to left):

- Use the two navigation keys to return to the start window or to continue to the **Recipes** view:

The following parameters can be adjusted:

| Parameters | Function |
|--------------------------|--|
| Pre-evacuation | <p>Adjustable time for the pre-evacuation phase</p> <p>Before starting the process loops evaporation – filling – draining, the system is pre-evacuated for the time set here.</p> |
| Initial filling | <p>Adjustable time for the initial filling</p> <p>Before starting the process loops consisting of evaporation – filling – draining, the system can be filled initially.</p> <p>An initial fill ensures that the process is started with minimum quantity, thus avoiding possible recrystallization.</p> <ul style="list-style-type: none"> ▪ Selection switch in position off → no initial filling. ▪ Selection switch in on position → Define a time period for the initial filling in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard). <p>Touch the on / off icon to switch between the two options.</p> |
| Fill | <p>Filling time setting</p> <p>Adjust the filling time as required using the plus/minus buttons.</p> |
| Evaporate | <p>Evaporation time setting</p> <p>Adjust the evaporation time as required using the plus/minus buttons.</p> |
| End concentration | <p>Additional evaporation step to reach a specific end concentration</p> <p>Activate this option to pass an additional evaporation step after passing through the process loops to achieve a desired end concentration or to remove any residual moisture from the residue in solid-liquid mixtures.</p> <ul style="list-style-type: none"> ▪ Selection switch in position off → no additional evaporation step. ▪ Selection switch in on position → Define a time period for one additional evaporation step (adjustment via plus/minus buttons or tap/on-screen keyboard). <p>Touch the on / off icon to switch between the two options.</p> |
| Drain residue | <p>Residue draining time</p> <ul style="list-style-type: none"> ▪ Selection switch in position off → no defined draining time. ▪ Selection switch in on position → Define a residue draining time in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard). <p>Touch the on / off icon to switch between the two options.</p> |

| | |
|-------------------------|--|
| Drain distillate | Distillate draining time <ul style="list-style-type: none"> In the selection field, define a distillate draining time (adjustment via plus/minus buttons or tap/on-screen keyboard). |
| End of process | Specify here whether the distillate and/or residue are drained at the end of the process. Touch the corresponding checkbox to activate (checkmark) or deactivate the selection. |

→ Confirm or discard all changes with **Accept** or **Cancel**.

Process parameters

- Tap the **Process parameters** button in the footer.
 - The process parameters menu opens.
- Touch the function button of the desired parameter and adjust the value using the screen keyboard, see also „User input“ on page 99.

| Parameters | Function |
|---|---|
| Heating tape | Performance of the heating tape during operation. Range: 10 – 100% in increments of ten. (see „Heating tape (optional)“ on page 92) |
| Heating tape automatic | Automatically switch on heating tape yes/no. Range: on / off. |
| Evacuation after draining the collector | Evacuation period of the bypass line after draining the collector vessel. Range: 1 – 9999 s. |
| Educt empty: additional evaporation time | When the educt container runs empty: Additional evaporation time for evaporation of the remaining educt in the evaporation flask. Range: 1 – 9999 s. |
| Forced draining of the collector | Forced draining of the collector active: Collector will be drained after 1 – 9999 seconds, regardless of the filling level. Range: 1 – 9999 s. Activate/deactivate forced draining of the collector (determines previous parameter!). Range: active / inactive |
| Standby timer after educt container runs empty | Time until the standby mode is activated after the educt container runs empty. Range: 1 – 9999 s. |


- **Hei-VAP parameters**, see section „Parameters Hei-VAP“ on page 127.
- **Data logging**, see section „Data recording“ on page 113.
- **Save**, see section „Recipes“ on page 128.

Rinse mode

Activating rinse mode

- Open the **Applications** menu and touch the **Rinse mode** function button. The entry is highlighted in orange.
- Confirm the selection with **ok** to activate the mode and to open the rinse mode start window:

The rinse mode allows the fill and concentrate line to be flushed with a suitable solvent. The flushing process removes residues from previous processes and impurities. When rotation is activated, the evaporation flask is also cleaned.

Multiple cycles should be carried out, depending on the degree of contamination!

To clean the complete system including the distillation area and the condenser, an evaporation process with a suitable solvent is necessary. Adjusting the parameters appropriately, the distillate produced this way will moisten all glass parts, thereby cleaning them.

PRACTICAL TIP

- Select a vacuum value adequate for the solvent and set the hysteresis to 50 mbar.
- A large pressure range is passed through by the large hysteresis.
- All glass parts are first wetted, after which drops are formed, which then run off.
- As a result, the glass parts in the distillate area are cleaned within a short time.

Preconditions

The following conditions must be met for rinsing the distillate tube:

- Heating bath heater of the rotary evaporator is switched on
- Recirculation chiller is switched on
- Vacuum pump is switched on

Procedure

- Depending on the selection, the various lines of the system are flushed.
- To fill the evaporation flask, see also section „Commissioning“ on page 100.


Rinse mode parameters

General parameters

- To adjust the **general parameters**, switch from the rinse mode start page to the corresponding parameters screen form (tap the parameters overview or swipe from right to left):

- Use the two navigation keys to return to the start window or to continue to the **Recipes** view:

The following parameters can be adjusted:

| Parameters | Function |
|-------------------------|---|
| Pre-evacuation | Adjustable time for the pre-evacuation phase Before starting the process loops evaporation – filling – draining, the system is pre-evacuated for the time set here. |
| Fill | Filling time setting Adjust the filling time as required using the plus/minus buttons. |
| Evaporate | Evaporation time setting Adjust the evaporation time as required using the plus/minus buttons. |
| Drain residue | Residue draining time <ul style="list-style-type: none"> ▪ Selection switch in position off → no defined draining time. ▪ Selection switch in on position → define a residue container draining time in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard). Touch the on / off icon to switch between the two options. |
| Drain distillate | Distillate draining time <ul style="list-style-type: none"> ▪ In the selection field, define a distillate draining time (adjustment via plus/minus buttons or tap/on-screen keyboard). |
| End of process | Specify here whether the distillate and/or residue are drained at the end of the process. Touch the corresponding checkbox to activate (checkmark) or deactivate the selection. |

Process parameters

Tap the **Process parameters** button in the footer.

– The process parameters menu opens.

→ Touch the function button of the desired parameter and adjust the value using the screen keyboard, see also „User input“ on page 99.

| Parameters | Function |
|---|---|
| Heating tape | Performance of the heating tape during operation. Range: 10 – 100% in increments of ten. (see „Heating tape (optional)“ on page 92) |
| Heating tape automatic | Automatically switch on heating tape yes/no. Range: on / off. |
| Standby timer after educt container runs empty | Time until the standby mode is activated after the educt container runs empty. Range: 1 – 9999 s. |


→ **Hei-VAP parameters**, see section „Parameters Hei-VAP“ on page 127.

→ **Save**, see section „Recipes“ on page 128.

Manual mode

Activating manual mode

- Open the **Applications** menu and touch the **Manual mode** function button. The entry is highlighted in orange.
- Confirm the selection with **ok** to activate the mode and to open the manual mode start window:


In manual mode, single batches can be processed. For automatic filling of the evaporation flask and draining of the distillate, the **Drain distillate** and **Start process** buttons are available on the operating panel.

Operating parameters must be individually defined on the connected rotary evaporator.

Procedure

- Touch the **Start process** button.
 - The vacuum pump is switched on, the pressure in the evaporation flask is lowered to process pressure.
 - In the footer, the **Fill educt** and **Stop process** buttons are displayed

- Touch the **Fill educt** button to start the automatic filling of the evaporation flask.

By briefly tapping the **Fill educt** button, a small amount can be injected into the evaporation flask.
 Touch the **Fill educt** button for approximately 2 seconds (observe button progress) to fill the evaporation flask continuously. The filling has to be stopped manually!

- Touch the **Fill educt** button again to stop filling.
- Touch the **Stop process** button to stop the process (end process).

When the process is ended, the **Drain distillate** button appears on the control panel.

- Touch the **Drain distillate** button → Distillate is drained (default: approx. 30 s).
- Touch the **Drain distillate** button → Drainage of the distillate is stopped.

The evaporation flask must be drained by hand.
 Please observe the relevant information and instructions in the operating instructions for the rotary evaporator!


Opening the error list

- To display the error list, open the **Main menu** and touch the **Error List** function button in the footer.
- To view an individual error message, tap the corresponding entry:
 - The error message appears in plain text on the display.
 - The entries **cannot** be deleted for service reasons!
- Close the error list with **ok**. The display returns to the **Main menu**.

Resetting the device to factory settings

Use the following procedure to reset the device to the factory settings:

- From the **Main menu**, open the **Settings** menu.
- Touch the **Factory settings** function button.
- The following prompt appears on the display.

- Confirm the confirmation prompt with **ok** to reset the device to the factory settings. The device must be restarted!
- Touch the **discard** button to discard the reset and return to the previous menu without making any changes.

Opening system information


The **Information** view displays all version and license information for your automatic module.

- To view system information, open the **Main menu** and touch the **Info** function button in the footer.
- Use the **Versions** and **Licenses** buttons to access the desired information about the software and firmware for your product.
- Scan the stored QR code with a suitable device to access additional product information and relevant documents via our homepage.
- Close the view with **ok**. The display returns to the **Main menu**.

Updating the firmware and software

The system is updated using an update file stored on a USB stick. The upgrade includes the system's software and firmware by default.

- Open the safety cover of the USB port on the control panel and insert the USB stick with the update file.
- Open the **Main menu**:

- Tap the **System update** button at the bottom of the main menu (this button appears only if a USB stick with a valid update file is detected on the USB port!):

- Confirm the confirmation prompt with **ok** to update the device software.
- Follow the instructions on the display
- Touch the **discard** button to discard the operation and to return to the previous menu without making any changes.
- Restart the device after the software update.
- After the restart, the firmware of the device can be updated. To do this, also confirm the corresponding confirmation prompt with **ok**,
- After the update, unplug the USB stick and put the security cover back into the USB port.

Troubleshooting

| Error message | Possible remedy |
|--|--|
| Cable breakage or sensor defect | Check the cable of the displayed component. |
| Short circuit | Check the cable connections of the displayed component. |
| Values out of range | Correct the values of the displayed parameter. |
| Compressor has been stopped | Excessive pressure from compressor, reduce pressure. |
| Pressure sensor failed! | Check pressure sensor for function, system check. |
| Error due to barometer! | Check the function of the barometer. |
| Temperature sensor failed! | Check temperature sensor for function. |
| Emergency stop has been triggered! | Eliminate the cause of the emergency stop from the external emergency stop switch. |
| Educt container S1 is empty! | Refill educt. |
| Overflow sensor S8 to the distillate container has triggered! | Drain the distillate container. |
| Overflow sensor S9 to the residue container has triggered! | Drain the residue container. |
| S3 cut out sensor has triggered! | Eliminate the cause of the emergency stop on the rotary evaporator. |
| Failure message from recirculation chiller | Observe the indicator on the recirculation chiller. |
| Invalid sensor plugged in! | Check the coding of the sensor, check that the sensor is in perfect condition. |
| Invalid valve plugged in! | Check the coding of the valve, check that the valve is in good condition. |
| Communication error! | Check the cable connection and connectors of the displayed components. |
| Operation mode cannot be started! | Open the error list, remove the indicated error cause(s). |

| | |
|---------------------------------------|---|
| ... has been plugged in twice. | Identical sensors or actuators plugged in, remove one unit. |
| Set switchpoint failed! | Switching error, contact technical service! |
| Calibration failed! | Repeat calibration procedure. |
| Update failed! | Check USB stick for correct connection, restart update. |

All error messages and warning messages are shown in plain text on the display.

Follow the instructions on the display.

In case of recurring errors, please contact the responsible sales department or our technical service. Contact address see „Heidolph International contact information“ on page 148.

Technical specifications

| General device data | |
|---------------------------------------|--|
| Model | Automatic Hei-VOLUME Distimatic Pro module |
| Dimensions (W × H × D) | 186 × 429 × 521 mm, including control panel |
| Weight | approx. 17 kg |
| Housing degree of protection | IP42 |
| Control panel degree of protection | IP42 |
| Control panel | detachable, 7 inch touchpanel, two control knobs |
| Electrical data | |
| Rated voltage | 100 – 240 V AC; 50/60 Hz |
| Connection | L+N+PE |
| Protection class | I |
| Overvoltage category | II |
| Degree of pollution | 2 |
| Power input | 1,500 W max. |
| Fuse | 15 A |
| Permissible ambient conditions | |
| Operating temperature | 5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly) |
| Installation altitude | up to 2,000 m asl |

Scope of delivery

| Component | Quantity | Product no. |
|---|-----------------|--------------------|
| Hei-VOLUME Distimatic Pro Basis | 1 | 591-60001-00 |
| Operating instructions | 1 | 01-005-006-51 |
| Assembly instructions | 1 | 01-001-009-12 |
| Safety guide | 1 | 01-005-006-53 |
| Guarantee registration / Certificate of decontamination | 1 | 01-006-002-78 |

Further information, in particular about the available expansion modules and the available accessories, can be found on our website at www.heidolph.com!

Device service

When carrying out service work on the device (cleaning, maintenance, repair), observe the general instructions and safety information described in this section.

WARNING

Live components are installed inside the device.

When opening the device, there is a risk of touching live components.

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

Penetrating liquid poses the danger of an electric shock.

When cleaning, avoid the penetration of liquids.

General cleaning instructions

Wipe all surfaces and the control panel with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

CAUTION

Improper cleaning can damage the surfaces of the device.

Penetrating liquid can damage the electronic components inside the device.

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

Repairs

Repairs to the device may only be carried out by authorized experts!

Unauthorized repairs during the warranty period will result in the loss of the warranty claim.

The owner is solely liable for damage caused by unauthorized repairs.

In case of a need for repair, contact an authorized dealer or our technical service, see „Heidolph International contact information“ on page 148.

Include the completed certificate of decontamination with every device return, see „Certificate of decontamination“ on page 149.

Maintenance

The device contains no user-serviceable components. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example), contact our technical service, see „Heidolph International contact information“ on page 148.

Dismantling

Observe the instructions contained in the assembly instructions supplied when dismantling.

Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
- Properly remove and dispose of residues of substances that are hazardous to health, the environment and are biohazardous!

Heidolph International contact information

Heidolph Instruments North America

Phone: 1-866-650-9604
 Email: service@heidolph.com
www.heidolphNA.com

Heidolph Instruments United Kingdom

Phone: +44 (0) 1799- 5133-20
 Email: service@radleys.co.uk
www.heidolph-instruments.co.uk

Local distributors

To find your local distributor please visit www.heidolph.com

Warranty Statement

Heidolph Instruments provides a three-year warranty for material and manufacturing defects.

The warranty does not cover glass and wear parts, transportation damage and damage resulting from improper handling or non-intended use of the product.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our website - www.heidolph.com.

The warranty period for unregistered products begins on the date of serial production (determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

Certificate of decontamination

Enclose the declaration of no objection, duly completed, with your device return. Submissions without a declaration of no objection cannot be processed!

CERTIFICATE OF DECONTAMINATION

IN CASE OF RETURNS

Please fill in the required fields.

Note: The sender must package the goods properly and appropriately for transport.

Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12
91126 Schwabach, Germany
Phone: +49 (0) 9122 9920-380
Fax: +49 (0) 9122 9920-19
Email: service@heidolph.de

SENDER

Name _____ First name _____
Company/institution _____ Department _____
_____ Workgroup _____
Address _____
ZC/City _____
Country _____ Phone _____
Email _____

DEVICE DETAILS

Article number _____ Serial no. _____

Ticket number _____

Reason for sending in _____

Has the device been cleaned, decontaminated/disinfected? Yes No (Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous? Yes No (Please mark as applicable)

If yes, with which substances did the device come into contact?

LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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The image displays two side-by-side screenshots of the Hei-CONTROL Pro configuration interface. Both screens show the 'Statut' section at the top with a green checkmark and the text 'Aucun clé d'appareil disponible'. Below this, there are fields for 'ID d'appareil :', 'URL du serveur :', and 'Configuration IP :'. In the left screenshot, 'Configuration IP' is set to 'Statique' and the 'Serveur NTP' checkbox is checked. In the right screenshot, 'Configuration IP' is set to 'Automatique (DHCP)' and the 'Serveur NTP' checkbox is also checked. Both screenshots show the same network parameters: Adresse IP (127.0.0.1), Masque réseau (255.255.255.0), Passerelle (192.168.2.1), and DNS (192.168.2.1).

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Heidolph Instruments GmbH & Co. KG

Walpersdorfer Straße 12

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hereby declare, that the product designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

Following essential requirements of of the Machinery Directive 2006/42/EC, Annex I have been fulfilled:

| | |
|-------------------------|----------------|
| DISTIMATIC HVI M.ABS. A | 591-61110-xx-x |
| DISTIMATIC HVI M.ABS. R | 591-61120-xx-x |
| DISTIMATIC HVI O.ABS. A | 591-61210-xx-x |
| DISTIMATIC HVI O.ABS. R | 591-61220-xx-x |

Machinery Directive 2006/42/EG
Electromagnetic Compatibility Directive 2014/30/EU
Delegated RoHS Directive 2015/863/EU

(Harmonized) Standards applied:
EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2020, EN IEC 63000:2018

The special technical documents were prepared in accordance to Annex VII, Part B. In response to a reasoned request by the national authorities, relevant information on the partly completed machinery can be transmitted via mail or email

The partly completed machinery must not be put into service until the final machinery into which it is to be incorporated has been declared in conformity with the provisions of this Directive, where appropriate.

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UKCA-Einbauerklärung - UKCA Declaration of Incorporation

UK DECLARATION OF INCORPORATION

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,
Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12
91126 Schwabach / Germany

Product: Laboratory automatic device

Model:

DISTIMATIC PRO HVI W. DRAINAGE A. ADAPT. A
DISTIMATIC PRO HVI W. DRAINAGE A. ADAPT. R
DISTIMATIC PRO HVI W/O DRAINAGE A. ADAPT. A
DISTIMATIC PRO HVI W/O DRAINAGE A. ADAPT. R

591-6111x-xx-x
591-6112x-xx-x
591-6121x-xx-x
591-6122x-xx-x

Description:

Device for filling and emptying of evaporation equipment.

The object of the declaration described above is classified as Partly Completed Machinery in accordance with UK Statutory Instrument 2008 No. 1597 Supply of Machinery (Safety) Regulations 2008.

The object of the declaration described above also is in conformity with the relevant UK Statutory Instruments (and their amendments):

2008 2016 No. 1091 The Electromagnetic Compatibility Regulations 2016
2012 No. 3032 The Restriction of the Use of Certain Hazardous Substances
 in Electrical and Electronic Equipment Regulations 2012

This product is supplied as an optional accessory for incorporation to a Machine Tool. The device has to be installed in accordance with the instructions contained in the Installation Guide.

The partly completed machinery must not be put into service until it has been established that the machinery into which the partly completed machinery is to be incorporated complies with the provisions of a relevant product directive as a whole.

And complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):

ProductIP (UK) Ltd.
8. Northumberland Av.
London WC2N 5BY

Signed for and on behalf of Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12, 91126 Schwabach / Germany

Schwabach, 17.08.2021

Wolfgang Jaenicke
Managing Director

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices by Heidolph Instruments GmbH & Co. KG fulfill the requirements of the EU directive 2011 / 65 / EU (RoHS) and the amended of directive 2015 / 863.

All devices are compatible with the requirement of the RoHS for the following materials:

Max. 0,01% of the weight in homogeneous material for cadmium and max. 0,1 % of the weight in homogeneous material for lead, mercury, hexavalent chromium, polybrominated biphenyl, polybrominated diphenyl ether, Di (2-ethylhexyl l) phthalate, butyl benzyl phthalate, dibutyl phthalate, diisobutyl phthalate.

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 22.07.2019

Stefan Peters
Vice President Marketing, Innovation & Technologie

Marcell Sarré
Vice President Quality Management & Technical Service

China RoHS Certification

China RoHS DECLARATION OF CONFORMITY

Heidolph Instruments GmbH & Co.KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Instruments GmbH & Co.KG products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
- O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.

- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
- X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Instruments GmbH & Co.KG may enter into further devices or can be used together with other appliances .

With these products and appliances in particular, Heidolph Instruments GmbH & Co.KG will not take responsibility for the EFUP of those products and appliances.

Schwabach, 06.08.2021

Wolfgang Jaenicke
Chief Executive Officer CEO

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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© Heidolph Instruments GmbH & Co. KG