Industrial
evaporation

Betriebsanleitung
Operating instructions
Notice d'instructions

Hei-VAP Industrial

 **heidolph**
research made easy

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Contents

Introduction

About this document.....	70
Typographic conventions.....	70
Copyright protection.....	70
Safety instructions in the official languages of the European Union	70

General product information

Directives applied, product certification.....	71
Customs declaration	71
California Residents.....	71
Residual risk.....	71
Intended use.....	71
Compliant use.....	71
Reasonably foreseeable misuse	72
Transportation.....	72
Storage.....	72
Acclimatization	72
Permissible ambient conditions	72

Safety

General safety instructions.....	73
Electrical safety.....	73
Data security	73
Operational safety	73
Occupational safety	74
Personal protective equipment (PPE).....	74
Environmental protection.....	74
Biohazard.....	75
Other regulations.....	75

Device description

Mechanical design.....	76
Main Switch	77
Emergency Stop.....	77
Control panel	78
Interfaces.....	79
Sensor/peripherals connections	79

Cooling water supply	80
Flask clamping system EASY LOCK	81
Mechanical design	81
Work positions	82
Handling	82
Adjust pressure	83
Flask support	84
Mechanical design	84
Handling	85
Heating Bath	85
Moving the heating bath manually	86
Automatic moving	86
Heating bath liquids	87
Automatic water bath refill	87
Safety features	88
Guard hood (optional)	88
Overpressure cut-out	89
Heating bath overheat protection	89
Level monitoring (heating bath liquid)	89

Commissioning

Setting up the device	90
Power supply	90
Switching the device on/off	91
Insert/remove the evaporation flask	91
Filling/emptying the heating bath	92
Bath liquid selector	92
Change the heating bath liquid	94

Operation

Activating/deactivating the device control	95
Structure of the user interface	95
Buttons and symbols	96
Numeric on-screen keyboard	97
Alphanumeric on-screen keyboard	97
[Settings] menu	98
Device parameters – preferences	99
Vacuum control mode	105
Operating mode SETpressure	105
Operating mode AUTOaccurate	106

Distillation under atmospheric pressure	108
Set rotation speed	110
Recommended rotation speed	110
Set the heating bath temperature	111
Timer function	111
Favorites and Profiles	112
Favorites mode	112
Create/edit favorites	113
Load favorite.....	113
Live view.....	114
Create/edit process profiles	115
Load process profile	117
Data logging	118
Activating/deactivating data logging	118
Setting the logging rate	119
Logging capacity	119
Data list	120
Deleting measurement series	120
Saving measurement series on a USB flash drive	120
Data output	121

Troubleshooting

Troubleshooting.....	122
----------------------	-----

Appendix

Technical specifications	124
Scope of delivery	125
Accessories	125
Device service.....	125
General cleaning instructions.....	126
Repairs	126
Maintenance.....	126
Disposal.....	127
Warranty statement.....	127
Contact information	127
Certificate of decontamination.....	128

About this document

This operating instructions manual describes the features and operation of Hei-VAP Industrial large-scale evaporators. The operating instructions are an integral part of the described device!

Typographic conventions

The following symbols, signal words and highlighting elements are used in this document:

Symbol	Signal word / explanation
	<p>Warning symbols in combination with a signal word indicate dangers:</p> <p>DANGER Indicates an immediate dangerous situation. Failure to respect the indications will result in death or serious injury.</p> <p>WARNING Indicates a potential danger. Failure to respect the indications will result in serious injuries.</p> <p>CAUTION Indicates a potential hazard which, if not avoided, damage to property and minor to moderate injuries can occur.</p>
	<p>Mandatory signs are used to indicate important and useful information on handling a product. This information is used to ensure operational safety and to maintain the value of the product.</p>
[GUI]	<p>Parameter designations, display texts, and device labels are typographically highlighted in the text and in tables to facilitate their assignment on the device.</p>
→	<p>The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.</p>

Copyright protection

This document is protected by copyright and is generally intended for use by the purchaser of the product.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Instruments GmbH & Co. KG. Any violation is subject to compensation for damage.

Safety instructions in the official languages of the European Union

A summary of all safety instructions in the official languages of the European Union can be found in our Safety Guide for the product group rotary evaporators (Ref. 01-005-006-67). This document is available for download on our homepage in the most up-to-date version.

Directives applied, product certification

CE marking

The device meets all requirements of the following directives:

- European Machinery Directive, 2006/42/EC
- EMC Directive, 2004/108/EC

NRTL Certification

The device has been tested in accordance with the following standards:

- UL 61010-1:2012/R:2016-04
CAN/CSA-C22.2 No. 61010-1:2012/U2:2016-04
- UL 61010-2-010:2015
CAN/CSA-C22.2 No. 61010-2-010:2015
- UL 61010-2-051:2015
CAN/CSA-C22.2 No. 61010-2-051:2015

Customs declaration

Rotary evaporators, as combustion and cleaning equipment, may be subject to notification to a competent customs authority in the country of destination.

The evaluation of a customer obligation to notify and, if applicable, the notification to a competent customs authority in the country of destination is generally the responsibility of the user!

California Residents

Important information for California residents regarding Prop 65. Please visit www.P65Warnings.ca.gov for more information.

Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

Intended use

The described device was designed by the manufacturer for the separation of substances, powder drying, concentration, crystallization of substances, and recycling of solvents under vacuum.

Due to its design, the device in its delivery condition may generally be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

Any other use of this device is not considered as intended!

Compliant use

The user is generally responsible for evaluating the conformity of his application and, if necessary, for taking additional measures.

Reasonably foreseeable misuse

Additional measures may be necessary for use under conditions or for purposes deviating from the intended use, and/or specific guidelines and safety regulations must be observed. Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is generally the responsibility of the operator.

All risks resulting from improper use are borne generally by the operator.

The device may generally be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the general responsibility of the operator!

Transportation

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device. Keep the original packaging in a dry and protected place for later use.

Storage

Always store the device in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

Acclimatization

After each transportation and after storage under critical climatic conditions (e.g. high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

Permissible ambient conditions

The device is designed for indoor use only. The device is **NOT** suitable for outdoor use! The device is **NOT** suitable for use in potentially explosive areas!

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

General safety instructions

- Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.
- Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself and, where applicable, on connected devices or the supply connections.
- If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.
- Only use the device in accordance with the regulations on intended use („Intended use“ on page 71).

Electrical safety

- Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.
- Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).
- Always use the power supply cord supplied with the device. All installation work for the power supply of the device may only be carried out by an authorized electrician or by the technical service of Heidolph instruments!
- Prior to use, check that the device and the power supply cord are free of visible damage.
- Always have repairs and/or maintenance work on the device carried out by an authorized electrician or by the technical service department of Heidolph Instruments.
- Always switch off and disconnect the device from the power supply, preventing reconnection, before carrying out maintenance work, cleaning, or repairs.

Data security

- The user is responsible for ensuring data security when transferring data between the described device and other devices.
- Always use secure networks for the data transfer and avoid use of critical infrastructure.
- Always use high-quality shielded data cables for the data transfer.
- For data transfer via a USB B connection, an industrial standard USB hub should be preferably used to ensure the most stable connection possible.

Operational safety

- Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (see EN 14175 and DIN 12924).
- Do not make any unauthorized changes or modifications to the device!
- Always use genuine spare parts and accessories, or those expressly approved by the manufacturer!
- Rectify malfunctions or faults on the device immediately.
- Switch off and disconnect the device from the power supply, preventing reconnection, if it is not possible to eliminate the malfunction or rectify the fault immediately.
- Observe all relevant general and safety instructions for the connected peripheral devices (observe the supplied documentation!).
- Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

Occupational safety

- Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.
- Do not operate any other devices in the immediate vicinity of the device ...
 - which can generate electromagnetic fields in the frequency range between 9×10^3 Hz to 3×10^{11} Hz,
 - which generate emission or radiation sources in the frequency range 3×10^{11} Hz to 3×10^{15} Hz (in the optical spectral range wavelengths from 1,000 μm to 0,1 μm),
 - which generate ultrasonic or ionizing waves.
- Do not operate the unit when adiabatic compression or shock waves may occur (pressure wave ignition).
- Do not use substances that could release energy in an uncontrolled way and cause a pressure increase (exothermic reaction, spontaneous ignition of dusts).
- Do not process hard, brittle materials such as stones, soil samples, etc., that could destroy the evaporation flask.
- Only use heating bath media that guarantee sufficient heat transfer.
- Do not operate the device with overpressure.
- Do not expose the glass components to a pressure difference of more than 2 bar.
- Ensure that the coolant overpressure does not exceed a level of 2 bar.
- Ensure that the flow velocity does not exceed 1 m/s when aspirating liquids with flammable components (electrostatic charge, danger of ignition!).
- Avoid the formation of explosion group IIC gases and potentially explosive distillation residues.

Personal protective equipment (PPE)

The operating company must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The corresponding instruction of the personnel is generally within the operating company's responsibility.

Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid hazards to the environment.

The evaluation of corresponding measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the general responsibility of the operator!

Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with the biohazard warning symbol.

The evaluation of corresponding measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the general responsibility of the operator!

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Instruments.

The operator is generally liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Mechanical design

The following figure shows an example of a device structure with glassware RC and optionally available base cart.

Legend

1	Base (option)
2	Control panel
3	Evaporation flask
4	Guard hood
5	Emergency stop
6	Rotational axis interlock
7	Glassware RC incl. condenser and receiving flask
8	Glassware protection housing with PMMA panel

Main Switch

The main switch for switching the device on and off is located on the rear of the device.

Turn the main switch clockwise to the ON position to switch on the device.

Note that the mains voltage is also present at the device when it is switched off (main switch in the OFF position). In order to de-energize the device for maintenance or repair work, the device must be disconnected from the power socket!

Emergency Stop

The emergency stop button of the device is located on the front, above the control panel.

By pressing the emergency stop button, all device functions are switched off and the running process is stopped. The device switches to the safe operating state (rotation and heating off, heating bath moves to the lower end position), but the mains voltage is still present!

In the event of an emergency stop, a warning signal sounds for approximately 5 seconds.

The emergency stop button engages when pressed. Correct the cause of the error and proceed as follows to unlock the emergency stop switch:

- Turn off the device
- Unlock the emergency stop button, turning it slightly to the left.
- Switch the device back on after a wait of at least 10 seconds

As soon as the cause of the error has been corrected and the emergency stop button has been released, the device is ready for operation again.

Before switching the device back on after an emergency stop, it is essential to correct the cause(s) of the fault or disturbance.

Use the emergency stop button only to switch off the device in emergency situations or in the event of a fault.

Always use the main switch to switch off the device during normal operation!

Control panel

CAUTION

Pointed or sharp objects can damage the surface of the display on the control panel!

To operate the device functions, always touch the display with your fingertips or a special stylus with a soft-touch tip.

The control elements of the central display are touch-sensitive (touch panel) and can be operated by finger or with a suitable stylus.

Legend

1	Menu button
2	Mode button (operating mode)
3	Button [Datalogging Start/Stop]
4	Button [Timer Start/Stop]
5	Button Up/Down heating bath lift
6	[Process parameters] area with mode button, Set/Actual display and Start/Stop button for single processes
7	Button [General process stop]: stop all processes simultaneously
8	Button [Warm-up]
9	Button [General process start]: start all processes simultaneously
10	Control [Vacuum override]
11	On/Off switch

Interfaces

The device interfaces are located on the bottom of the control panel:

- 1 × USB 2.0, max. 500 mA
- 1 × LAN RJ-45 (not working in the current device version)

WARNING

If the interfaces are assigned too high voltages and are not sufficiently insulated, metallic parts such as the housing can be live in the event of a fault.

Safely separate low-voltage inputs and outputs via 25 V AC or 60 V DC according to DIN EN 61140, or by double or reinforced insulation according to DIN EN 60730-1 or DIN 60950-1.

Always use shielded connection cables. Connect the shield to the connector housing.


CAUTION

Make sure that the device interfaces on the bottom of the display are always properly covered!

Penetrating liquid can cause interface malfunctions and short-circuiting.

Sensor/peripherals connections

The connections for the sensors and the peripheral devices are located on the rear panel of the device (further information, see mounting instructions O1-001-009-12):

Connections

1	Vacuum sensor
2	Vacuum valve
3	Venting valve
4	Heating bath water valve
5	Vacuum pump
6	Domestic water for automatic water bath refill of the heating bath
7	Return of the cooling liquid in the recirculation chiller
8	Supply of the cooling liquid in the recirculation chiller
9	Emergency stop sensor for Distimatic Pro (optional)
10	Temperature sensor PT (cooling liquid of the recirculation chiller)

Cooling water supply


CAUTION

Keep the [Supply] and [Return condenser] connections closed if no external condenser is connected!

In order to not exceed the maximum permissible pressure in the glass condenser (2 bar) during operation, open the [Supply] and [Return condenser] connections when an external condenser is connected!

The connections for the cooling water supply are located on the left-hand side panel of the device, above the receiving flask(s). The feed of the cooled condenser is made via the connections [2] and [4] (main cooling, cooling of the medium). The connections [1] and [3] are equipped with shut-off valves and can be used for optional aggregates:

Connections

1	Supply of cooling liquid (optional aggregates)
2	Return of cooling liquid from condenser (main cooling)
3	Return of cooling liquid (optional aggregates)
4	Flow of cooling liquid to condenser (main cooling)

Flask clamping system EASY LOCK

The device is equipped with the patented flask clamping system EASY LOCK for Hei-VAP evaporation flasks.

CAUTION

Incorrect or careless handling can damage the flask clamping system and/or the flask inserted!

Before inserting a flask for the first time, familiarize yourself with the handling of the EASY LOCK flask clamping system (see „Handling“ on page 82).

Before inserting a flask for the first time, familiarize yourself with the handling of the flask support (see „Handling“ on page 85).

Unsuitable flasks can damage the flask clamping system and/or the flask inserted!

The EASY LOCK flask clamping system is only suitable for flasks with an opening diameter of 150 mm. The use of flasks with a smaller opening diameter (e.g. LR 20) is not permitted, as these cannot be fixed in a positive manner due to the geometric properties.

Mechanical design

The following figure shows the structure of the EASY LOCK flask clamping system:


Components

1	Clamping lever
2	Clamping lever lock
3	Outer ring
4	Sliding segments to fix the flask

Work positions

The following figure shows the different working positions of the flask clamping system:

- "closed-locked" position [1]: In this position, the evaporation flask is clamped in the EASY LOCK, the evaporator can be operated normally.
- "unlocked" position [2]: The internal sliding segments of the EASY LOCK are open enough that the evaporation flask inserted can be rotated and positioned for a change.
- "open" position [3]: The evaporation flask can be inserted or removed. Furthermore, the sliding segments and the outer ring can be removed.

Handling


Before opening the EASY LOCK, the rotation axis of the device must be locked using the latching rotary knob!

To lock the rotational axis, follow these steps:

- Pull out the latching rotary knob and turn it up to the [Lock] position to lock the rotation axis.
- Depending on the rotation position, it may not be possible to engage the locking knob directly in the locking position! In this case, turn the outer ring of the EASY LOCK and the latching rotary knob at the same time until it can be turned to the locking position and locks into place.
- Turn the latching rotary knob into the [Unlock] position to unlock the rotation axis.

After locking the rotational axis, the following message appears on the device display:


A process start is only possible when the rotational axis is unlocked and the warning message on the display is confirmed with [OK] !

- To open the EASY LOCK, first press in the clamping lever lock to unlock the clamping lever.
- Pull the unlocked clamping lever up to the stop until it is held in the open position (2) by the pin on the threaded bolt of the pretensioning mechanism. The sliding segments are released, the flask can be moved but not removed!
- Grip the clamping lever with your hand and turn the outer ring of the EASY LOCK downwards up to the stop (3). When doing so, the clamping lever must not be pressed against the outer ring. The EASY LOCK is now fully open and the flask and the PTFE seal can be inserted or removed.

- Insert a suitable evaporation flask (see also „Handling“ on page 85)
- Grip the clamping lever with your hand and use it to pull the outer ring of the EASY LOCK upwards up to the stop. The sliding segments are pre-tensioned, the flask can be moved, but is securely held.
- Then press the clamping lever against the outer ring with normal manual force until the clamping lever lock audibly engages.

Adjust pressure


Due to industrial manufacturing tolerances, the contact pressure of the EASY LOCK must be adjusted individually for each evaporation flask.

To do this, proceed as follows:

- Lock the rotational axis.
- Unlock and open the EASY LOCK clamping lever to the stop position (see „Handling“ on page 82).
- Insert the set screw completely, turning clockwise (marked in the following figures): In this position, minimum contact pressure is applied to the flask.
- Insert the flask and lock the EASY LOCK (see „Handling“ on page 85).
- Check if the clamped flask can be moved in the locked EASY LOCK.

If the flask still rotates:

- Unlock and open the EASY LOCK clamping lever to the stop position.
- Pull back the inserted flask.
- Turn the set screw one to two turns counter-clockwise to increase the contact pressure step by step.
- Insert the evaporation flask and lock the EASY LOCK.
- Check if the clamped flask can be moved in the locked EASY LOCK.
- Repeat the procedure until the clamping lever closes and locks against moderate resistance and the flask can no longer be moved by hand in the closed EASY LOCK

CAUTION

Excessive contact/closing pressure can damage the EASY LOCK and/or the flask inserted!

Always follow the instructions for the correct adjustment of the contact pressure in this section!

Flask support

The flask support facilitates the insertion of the evaporation flask and prevents the flask from falling into the heating bath during removal.

Mechanical design

The flask support consists of the following mechanical components (see also the following figure):

- Two guide bars, fixed to the chassis
- Two support bars with support element, movable on the two guide bars
- Supported locking lever with handle and defined locking positions

Components

1	Guide bars
2	Support bars with support element
3	Locking lever

Handling

The two support bars of the flask support are connected with the locking lever so that both can be moved parallel along the guide bars to and from the EASY LOCK.

The distance between the support bars is defined by the locking positions for three flask sizes.

To set the appropriate distance, do the following:

- Uncouple the support bars by pulling up on the locking lever on the handle.
- Specify the appropriate locking position (refer to the volume mark on the evaporation flask!).
- Move the two support bars so that the locking lever engages in the appropriate locking position when moving it down.
- Using the locking lever, re-couple the support bars.
- Ensure that the support bars are securely coupled. To do this, move the unit along the guide bars.

The evaporation flask is held by the support bars in such a way that it can be rotated in all axles.

Due to the mechanical design of the flask supports, the evaporation flask cannot fall into the heating bath or roll from the flask support when the locking lever is correctly set and engaged.

For normal operation, the two support bars must be decoupled and moved outwards to their respective end position. The end position of the support bars is monitored electronically.

As long as the support bars are not in their respective end position, the following message appears on the display:

A process start is only possible when the support bars are in their respective end position in the flask support and the warning message on the display has been confirmed with [OK]!

Heating Bath

The heating bath is moved up or down with a motor lift. The control is carried out via the device software and the control panel of the device.


In normal operation, the heating bath is moved to the upper end position in order to immerse the flask into the heating bath liquid.

To replace the evaporation flask, to fill the heating bath and to carry out cleaning, maintenance or repair work, etc., the heating bath is moved to the lower end position.

The heating bath lift only works if the level sensor (float switch) is correctly mounted in the heating bath.

Moving the heating bath manually

The [Lift UP/DOWN] arrow keys can be used to manually move the heating bath lift. The buttons are available on the start window after the device is switched on and initialized:

- Touch the [Lift DOWN] button to move the heating bath down manually.
- Touch the [Lift UP] button to move the heating bath up manually.
 - The movement stops as soon as you release the button.
 - The movement stops automatically at the upper or lower end position.

CAUTION

When the evaporation flask is immersed, the already filled heating bath liquid is displaced and the liquid level in the heating bath rises!

Especially when using oil, the device and the surrounding floor surface can be heavily contaminated by overflowing heating bath liquid.

When filling the heating bath, make sure that the maximum filling quantity is not exceeded.

Automatic moving

The heating bath is automatically moved up or down when one of the following device functions is performed:

- Warm Up: when the heating process is started, the heating bath is automatically moved to the lower end position - the flask is lifted out of the heating bath.
- Start All: when starting all processes in automatic mode, the heating bath is automatically moved to the upper end position - the flask is immersed in the heating bath fluid.
- Stop All: when stopping all processes, the heating bath is automatically moved to the lower end position - the flask is lifted out of the heating bath.
- Emergency Stop: when the emergency stop button is pressed, the heating bath is automatically moved to the lower end position - the flask is lifted out of the heating bath.
- Timer:Lift: after a set time period in timer mode, the heating bath is automatically moved to the lower end position - the flask is lifted out of the heating bath.

Heating bath liquids

DANGER

When using oils with a flash point < 285 °C as a heating bath liquid, uncontrolled thermal reactions can occur at high heating bath temperatures and an explosive atmosphere can arise: Always use oils with a flash point > 285 °C as heating bath liquid.

From a heating bath temperature of approx. 75 °C – 80 °C, the use of polyethylene glycol or oil is recommended as heating bath liquid.

Reference values for heating bath liquids

- Tap water; temperature range 20 °C – 90 °C
- Polyethylene glycol; temperature range 20 °C – 170 °C
- Silicone oil, low viscosity (max. 40 cP), flash point > 285 °C; temperature range 20 °C – 170 °C

Automatic water bath refill

This function is switched to active when the bath liquid selector for the heating bath liquid is set to [Water] and the relevant setting has been adjusted: see „Menu item Bath liquid“ on page 100.

Functional principle

Due to the heat development in the heating bath, the fluid level gradually decreases due to evaporation.

In order to ensure a constant immersion depth of the flask and thus a continuous evaporation process, the heating bath is automatically lifted up step by step toward the upper end position. This lifting movement is controlled by the level sensor.

Operating behavior

As soon as the heating bath reaches the upper end position (minimum level), water is automatically refilled until the maximum level in the heating bath is reached again.

Monitoring

When refilling the heating bath via the automatic water bath refill, the filling duration is limited to 30 seconds in the as-delivered condition due to safety reasons. This function can be deactivated or the refilling duration can be adjusted in the menu [Settings], [Water Timeout] parameter (see section „Menu item Water Timeout“ on page 100)!


CAUTION

Possible overflow during the refilling of the heating bath! The monitoring of the refill duration during the automatic water bath refill should only be deactivated in exceptional cases.

If the maximum level during an active automatic water bath refill is not reached within the defined time period, the automatic refill is shutdown and the device switches to the safe operating state.

The following prompt appears on the display:

- Confirm this prompt and restart the device. Check the water supply from the house connection to the supply above the heating bath for correct connection and continuity as well as the level sensor for correct function.

If the maximum level is reached within the defined time period, the automatic refill is switched off normally and the device is ready for use.

Safety features

Guard hood (optional)

Depending on the equipment variant, the device is equipped with a guard hood that completely covers the evaporation flask and the heating bath, see section „Mechanical design“ on page 76.

WARNING

Even when the guard hood is closed, there is a risk that the heating bath liquid will spill out of the heating bath at high rotational speeds. There is a risk of scalding and/or contamination.

Also refer to the recommendations in section „Set rotation speed“ on page 110.

CAUTION

Always use the handle provided to open and close the guard hood.

When closing the guard hood, keep sufficient distance to the closing edges with your free hand.

The guard hood is equipped with an interior lighting which allows a better view of the ongoing distillation process. The guard hood fulfills the following safety functions:

- Basic protection against touching the rotary/heated evaporation flask
- Basic protection against direct contact with the (heated) bath liquid
- Basic protection against flying glass parts and spriting liquid in the case of implosion or thermal reaction

An error message appears on the display when the guard hood is opened during operation. At the same time, the flask rotation is stopped and all running processes are stopped. As soon as the guard hood is fully closed again, the flask rotation starts again and the stopped processes continue.

Due to design and function, the guard hood is not designed to tightly seal the heating bath area!

Overpressure cut-out

The device is equipped with an electronic overpressure monitoring system. As soon as the integrated vacuum sensor detects an impermissible pressure level, an error message appears on the display; at the same time, the device is switched off and put into an operationally safe state.

The vacuum limit is set to 1,200 mbar at the factory. Depending on the prevailing ambient conditions, this limit value can be adjusted in a range of 900 – 1.400 mbar.

The overpressure limit is set in the menu [Settings] – Safety Settings, Menu item Safety Settings* on page 99.

PRACTICAL TIP

The higher the application location is above sea level, the lower the vacuum limit value:


- Measure the ambient pressure on the ventilated system and at high air pressure.
- Add about 50 mbar to the reading.

Heating bath overheat protection

The device is equipped with an autonomous overheat protection function that prevents the heating bath from overheating. When reaching the maximum heating bath temperature, the heating function is switched off by a mechanical switch

The overheat protection has to be reset manually:

- Make sure that the heating bath is in the lowest position and allow the heating bath liquid to cool down.
- Remove the safety assembly cover on the right side of the rotary evaporator housing, see figure.
- Operate the reset button with a suitable tool (blunt tip, insulated or non-conductive material).
- Mount the safety assembly cover on the right side of the rotary evaporator housing.

Level monitoring (heating bath)

The level in the heating bath is monitored during operation by means of a permanently installed level sensor (float switch in the heating bath). The movement of the heating bath lift is stopped when the level sensor is activated. Limit values of the level sensor:

- minimum liquid level in the heating bath: 50 mm from bottom of the heating bath
- maximum liquid level in the heating bath: 145 mm from bottom of the heating bath

Operating behavior when refilling tap water automatically

During the initial filling or the refilling of the heating bath via the automatic water supply of the device, the refill duration is limited to five (5) minutes due to safety reasons. This value cannot be changed.

If the maximum level is not reached within the defined time period, the automatic refill is shutdown and the device switches to the safe operating state.

→ In this case, check the water supply from the house connection to the supply above the heating bath for correct connection and continuity as well as the level sensor for correct function.

If the maximum level is reached within the defined time period, the automatic refill is switched off normally and the device is ready for use.

Setting up the device

Ensuring proper system and/or device mounting, including all accessory parts, is the general responsibility of the operating company.

If required, Heidolph offers a professional system and/or device mounting service, including all supplied components and including commissioning. The associated mounting instructions (Ref. 01-001-009-12) describe all steps for proper system and/or device mounting in detail, including wiring and tubing.

CAUTION

Due to improper mounting and/or installation of the system/device, as well as unauthorized changes to the mounting, there is a risk of direct and indirect damage to property!

If necessary (system mounting and installation carried out by the operator or third parties, necessary changes to the existing mounting), contact the manufacturer's technical service (see section „Contact information“ on page 127).

To use the device, place it on a clean, stable, level and horizontal surface. Before switching on the device, make sure that it is sufficiently stable and use all available locking means on the base.

The operator is generally liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions! At the same time, in this case, any warranty claims against Heidolph Instruments will invalidate.

Power supply

DANGER

Please consider all notes in section „General safety instructions“ on page 73, especially section „Electrical safety“ on page 73.

Switching the device on/off

CAUTION

Before each operation of the device, check the caps and connecting elements of all glass components for correct and tight fit!

At loose/leaky connections, ambient air is sucked in during the vacuum build-up and the required vacuum pressure cannot be reached in a stable manner!

In the event of excess pressure, there is a risk that liquid will escape from the system! The device and surrounding areas can be severely contaminated by leaking media!

Tighten loose caps/fasteners by hand or correct the seating of the retaining clips.

To switch on and off, use the main switch on the rear of the device, see section „Main Switch“ on page 77.

Insert/remove the evaporation flask

CAUTION

In the event of improper handling, the evaporation flask that is to be used may fall to the ground.

- Glass breakage!
- Risk of injury/contamination by leaking medium!

Use only the evaporation flasks approved by the manufacturer.

When inserting and removing an evaporation flask, please consider all instructions for handling the „Flask clamping system EASY LOCK“ on page 81 as well as all instructions for the „Handling“ on page 85.

WARNING

The smallest damage to the evaporation flask can lead to glass breakage when pressurized.

- Glass breakage!
- Risk of injury/contamination by leaking medium!

Before each operation of the device, check the evaporation flask for visible damage.

Insert the evaporative flask

To insert the evaporation flask, follow these steps:

- Stop all running processes (if applicable).
- Move the heating bath to the lower end position („Heating Bath“ on page 85).
- Stop the rotation movement of the flask (if applicable).
- Open the glass guard hood.
- Lock the rotation axis and fully open the flask clamping system.
- Position the flask support over the heating bath according to the flask size.

- Align the evaporation flask on the flask support so that the opening can be inserted into the open EASY LOCK.
- For new evaporation flasks, adjust the contact pressure of the EASY LOCK.
- Close the EASY LOCK.
- Make sure that the EASY LOCK clamping lever is properly engaged.

Remove the evaporation flask

To remove the evaporation flask, follow these steps:

- Stop all running processes (if applicable).
- Move the heating bath to the lower end position („Heating Bath“ on page 85).
- Stop the rotation movement of the flask (if applicable).
- Lock the rotational axis.
- Open the glass guard hood.

CAUTION

The heating bath is designed for a maximum operating temperature of 210°C (using silicone oil). From a liquid temperature of 50 °C, there is a risk of injury in the event of contact.

Allow heated liquid to cool below 50 °C before removing the flask.

Use appropriate personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to remove the flask.

- Position the flask support according to the flask size.
- Fully open the flask clamping system.
- Using the flask support, pull the evaporation flask out of the open EASY LOCK and rotate the it so that the opening is vertical up.
- Remove the evaporation flask.

Filling/emptying the heating bath

Before the initial filling or refilling of the heating bath, the preselection water or oil has to be made via the [Bath liquid] selector.

Bath liquid selector


Before filling the heating bath for the first time and when changing the bath fluid, the bath fluid selector on the back of the device must be moved to the required position.

The bath liquid selector setting affects the following device functions:

- Process control (automatic mode)
- Heating bath level sensor
- Automatic water bath refill

The selection can only be changed if the device is switched off and the heating bath liquid has cooled down to at least 50 °C!

- Slide the bath liquid selector to the [Oil] or [Water] position.
- [Water] position
 - Maximum temperature setting: 90 °C
 - Automatic water bath refill active
 - Heating bath level sensor active
- [Oil] position
 - Maximum temperature setting: 180 °C
 - Automatic water bath refill inactive
 - Heating bath level sensor deactivated
- The setting of the bath liquid selector must be confirmed in the [Settings] menu, see „Menu item Bath liquid“ on page 100.

Fill the heating bath

The heating bath can always be filled manually and when using tap water it can also be filled automatically. Before the automatic fill, make sure that the device is correctly connected to the house connection. Observe the instructions in section „Operating behavior when refilling tap water automatically“ on page 90.

CAUTION

If the level in the bath is too high, there is a risk that the heating bath liquid will overflow when the flask is immersed and will severely contaminate surrounding areas.

Consider the mark for maximum filling level in the heating bath when choosing the manual refill!

- Set the bath liquid selector on the rear of the device to the correct position, see section „Bath liquid selector“ on page 92.
- Close the drain valve on the back of the device (handle in 90 ° – position to the drain valve).
- Ensure that the required evaporation flask is correctly mounted („Insert/remove the evaporation flask“ on page 91).
- Move the heating bath to the upper end position („Moving the heating bath manually“ on page 86).
- With the flask immersed, fill the heating bath with the desired medium until the liquid level reaches a height of 40 mm below the upper edge of the heating bath (maximum level) or open the valves of the house connection to fill the heating bath via the filling system of the rotary evaporator.

Empty the heating bath

CAUTION

The heating bath is designed for a maximum operating temperature of 210°C (using silicone oil). From a liquid temperature of 50 °C, there is a risk of injury in the event of contact.

Allow heated liquid to cool below 50 °C before removing the flask.

Use appropriate personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to remove the flask.

- Ensure that a suitable drain tube is attached to the drain valve.
- Move the heating bath to the lower end position („Moving the heating bath manually“ on page 86).
- Open the drain valve on the back of the device (handle parallel to the drain valve).
- Clean the heating bath and the drain tube between the heating bath and the drain cock.
- Close the drain valve on the back of the device!

Change the heating bath liquid

To change the heating bath liquid, proceed as described in the preceding section „Remove the evaporation flask.“ on page 92:

- Emptying the heating bath, see „Empty the heating bath“ on page 94.
- Filling the heating bath, see „Fill the heating bath“ on page 93.
- Observe the instructions in section „Bath liquid selector“ on page 92!

When changing the bath liquid, the heating bath and the drain tube between the heating bath and the drain valve must be thoroughly cleaned and dried.


In particular, when changing from water to oil, make sure that the heating bath temperature is only slowly and gradually augmented during restart (recommendation: 10 °C), so that remaining water in the system can evaporate until the process temperature is reached.

Activating/deactivating the device control

Make sure that the device is switched on (main switch in position ON).

To activate/deactivate the device control, touch the on/off button in the lower right corner of the control panel (see also „Control panel“ on page 78).


First, the following dialog box appears on the device display:

- Confirm the prompt with [Yes] to deactivate the device control.
- Confirm the prompt with [No] to activate the device control.

Structure of the user interface

After the device control has been initialized, the user interface start page appears on the display.

The following options are available:


- [Settings] menu: touch the [Menu] button on the top left to open the menu.
- Selection of the operating mode
 - Touch the [Set Pressure] button to operate the device at a specific vacuum pressure, see „Operating mode SETpressure“ on page 105.
 - Touch the Auto Accurate button to activate the AUTOaccurate function and operate the system with automatic vacuum control, see „Operating mode AUTOaccurate“ on page 106.

The software keys/buttons of the active function are highlighted in orange in all menu windows.


Buttons and symbols

The following table provides an overview of all buttons and icons on the control panel and in the menu navigation:

Symbol	Function
	ON/OFF button on the control panel
	[Menu] button
	Start/stop function
	Heating bath lift up / down
	Symbol [Function deactivated]
	Accept value (OK), activate function / symbol [Function activated]
	Open info dialog (checklist)
	Jump to the end (checklist)
	Arrow keys step-by-step back/forward
	[Delete] button
	Toggle keys numeric/alphanumeric characters
	Shift button upper case
	Insert space
	Add new favorite
	Add new row (Profile)
	Remove individual entry
	Remove entire list or operation
	Zoom in/out
	Save Favorite
	Save and apply Favorite

Numeric on-screen keyboard


The numeric on-screen keyboard appears automatically when you touch a numeric entry field within a mask.

- Use the keyboard to enter a desired value within the range of values allowed in each case.
- Touch the [Cancel] button to return to the parent dialog without making any changes.
- Touch the [Delete] button to cancel the last entry.
- Confirm your entry with [OK] to accept the new value. The on-screen keyboard closes.

Alphanumeric on-screen keyboard


The alphanumeric on-screen keyboard appears automatically when you touch an input field for alphanumeric values within a mask.

- Use the keyboard to enter a desired designation.
- Touch the [Cancel] button to return to the parent dialog without making any changes.
- Touch the [123] button to open the highlighted keyboard with numeric values and special characters.
- Touch the [Delete] button to cancel the last entry.
- Confirm your entry with [OK] to accept the new value. The on-screen keyboard closes.

[Settings] menu

From the start window, touch the [Menu] button to open the main menu. Open the [Settings] menu by touching [Settings]. Here you can set various default settings for the operation of the device:

The menu includes several screens. Touch the display and swipe up or down to scroll through the menu.

The operating functions are identical in all screens:

- Touch the arrow button (left) at the bottom of the screen to return to the [Settings] menu without making any changes.
- Touch the [OK] button at the bottom right of the screen to confirm an entry.
- Valid changes to the parameters are accepted without further confirmation prompt!
- Entries outside the respective valid value range are not accepted.

Device parameters – preferences

Menu item Safety Settings

CAUTION

Due to incorrect/insufficient safety settings there is a risk that processes will not be shut down or will be shut down with a time delay when an unacceptable pressure level is reached.

Only qualified personnel may adjust the safety settings of the device.

When required, please contact our technical service, see „Contact information“ on page 127.

To open the [Safety Settings] menu, touch the arrow button [right] in the list entry. Here, you can define the conditions under which a running process is automatically stopped. To adjust the values, follow these steps:

- Touch the input field in the entry [Overpressure].

- The numeric on-screen keyboard appears: Use the keyboard to enter the desired value:
 - Default [Overpressure]: 1200 mbar
 - Setting range [Overpressure]: 900 – 1400 mbar
- Confirm the change with [OK].

The stop criterion [Leakage Pressure] has to be activated manually if required or necessary:


- Touch the icon [Function deactivated] of the entry [Leakage Pressure]. The symbol changes to [Function activated], the stop criterion is active and can be adjusted.

- Touch the input field in the entry [Leakage Pressure].
- The numeric on-screen keyboard appears: Use the keyboard to enter the desired value:
 - Default [Leakage Pressure]: 900 mbar
 - Setting range [Leakage Pressure]: 1 – 1399 mbar
- Confirm the change with [OK].

Menu item Water Timeout

Via this menu item, you can adjust and, if necessary, deactivate the refilling duration for the function [Automatic water bath refill].

→ Adjust refilling duration, symbol [Function activated]: Touch the input field for seconds [s]:

- Use the keyboard to enter the desired wait time in a range of [1 s – 9999 s].
- Deactivate the function: touch the icon [Function activated]. The icon changes to [Function deactivated]:

Menu item Temperature

You can use this menu item to select the unit for temperature displays.

- Touch the button of the desired display unit (°C, °F, or K). The active selection is highlighted in orange.
- All temperature values are displayed in the selected unit.

Menu item Pressure

You can use this menu item to select the unit for pressure displays.

- Touch the button of the desired display unit (mbar, hPa or Torr). The active selection is highlighted in orange.
- All pressure values are displayed in the selected unit.

Menu item Bath liquid

In this menu, the current Bath liquid setting is displayed.

- First, check the position of the bath liquid selector on the back of the device (see section „Bath liquid selector“ on page 92) and adjust it, if necessary.
- Touch the [Change] button to switch between the settings [Water] and [Oil] (according to the position of the selector).
- Note and confirm the following system info with [OK]:

Menu item Lock screen after

You can use this menu item to specify the time after which the display is automatically locked from the last operation.

→ Touch the input field for [Minutes]:

→ Use the keyboard to enter the desired wait time in a range of at least one (1) to a maximum of five (5) minutes.


→ To unlock the screen during operation, drag the onscreen slider to the right.

Menu item Manual vacuum override active while screen locked

As soon as the screen lock is active, the system stops responding to the vacuum control knob of the control panel in its normal state.

This menu item allows you to deactivate the vacuum control knob lock.

→ Touch the icon [Function deactivated] of the entry. The icon changes to [Function activated]. In this state, the vacuum pressure can be adjusted using the vacuum control knob despite the screen lock being active.


→ Proceed in the same way to reactivate the vacuum control knob lock: Touch the icon [Function activated], the icon switches to [Function deactivated].

Menu item Date/Time

Via this menu item, you can adjust the format for the display of the system time.

→ Touch the arrow button [right] in the entry. The [Date/Time] input screen appears on the display:


- Set the date and time format by touching the corresponding input field (active format highlighted in orange).
- Touch the input fields [Date] and [Time] one after the other and entry the desired values via the keyboard.

Menu item Light inside Safety Hood

This menu item allows you to activate or deactivate the automatic activation of the interior lighting of the guard hood.

- Touch the icon [Function deactivated] of the entry. The icon changes to [Function activated]. In this state, the interior lighting of the guard hood is automatically switched on when the device is switched on.


- Proceed in the same way to switch off the automatic activation of the interior lighting again: Touch the icon [Function activated], the icon switches to [Function deactivated].

Menu item Glassware assembly

You can use this menu item to select the mounted glassware.

- Touch the arrow button [right] in the entry. The selection screen [Glassware] appears on the display:

- Touch the entry of the mounted glassware (A, A2 or R) to load the corresponding system settings. The settings are immediately applied in the main menu and in the checklist.

Menu item Show Checklist

You can use this menu item to call up a system checklist. You can use this checklist to check all the relevant components on the device structure: Basic safety instructions, installation of the peripherals, installation of the glassware, installation of the sensors, installation of the tubing, process aid


- Touch the arrow button [right] in the entry and confirm the system prompt [Show Checklist] in the following window with OK to load the checklist stored for the selected glassware:

Menu item Show Checklist when switching on the device?

You can use this menu item to specify whether the checklist is loaded automatically when the device control is switched on (factory setting!), or manually as described above.

- Touch the icon [Function activated] of the entry. The icon changes to [Function deactivated]. In this state, the checklist is not loaded when the device control is switched on.

- Follow the same procedure to enable automatic loading of the checklist. Touch the icon [Function activated], the icon switches to [Function deactivated].

Menu item Language

You can use this menu item to specify the desired menu language.


- Touch the arrow button [right] in the entry and then the icon [Function activated] of the desired language in the list of all available languages. The icon switches to [Function activated], the setting is immediately applied to the system:

Menu item Show System Information

Use this menu item to open the [System information] window and load available updates. The following information is displayed:

- **HMI Version**
 - **HMI Update Version:** This entry changes to **Insert USB stick** as soon as you have inserted a USB stick with an HMI update into the USB slot on the underside of the control panel.
 - **Controller Version**
 - **Controller Update Version:** This entry changes to **Insert USB stick** as soon as you have inserted a USB stick with a system update into the USB slot on the bottom of the control panel.
- Touch the arrow button [right] in the entry to open the [System information] window.

When updating the HMI or system, do not remove the USB stick from the USB slot until the installation is complete. Note the progress indicator on the display.

When the installation is complete, a message appears at the bottom of the display.

After each update, the device control restarts automatically.

Menu item Show Datalog Series

You can use this menu item to display recorded Datalog measurement series. Measurement series are recorded using the Datalogging function. For more information on Datalogging, see „Data logging“ on page 118.


- Touch the arrow button [right] in the entry. The available measurement series are displayed.

Menu item Error History

You can use this menu item to call up all recorded system messages. All messages are displayed with a time stamp.

- Touch the arrow button [right] in the entry to open the error history (figure by way of example):

- Touch the button [Remove individually] of an entry to remove it from the list.
- Touch the button [Remove all] at the bottom edge of the display to remove the whole list.

Vacuum control mode

The device can be operated with manual (SETpressure) and automatic (AUTOaccurate) vacuum control.


For correct vacuum control, the vacuum valve, the vacuum sensor and the ventilation valve must be connected to the rear of the device (see the associated assembly instructions).

Regardless of the vacuum control mode, the vacuum pressure can be adjusted at any time during operation with the [Vacuum override] rotary control (see „Control panel“ on page 78 and „Menu item Manual vacuum override active while screen locked“ on page 101)!

Operating mode SETpressure

In this operating mode, all process parameters are specified manually by the operator. To open the start window of this mode, touch the [Set Pressure] button on the start window:

Legend

1	SETpressure activated
2	Actual value of the coolant temperature of the PT-100 sensor
3	Button [Save]
4	Process parameter [Vacuum], display of actual and set value, Start/Stop button
5	Process parameter [Rotation], display of actual and set value, Start/Stop button
6	Process parameter [Bath temperature], display of actual and set value, Start/Stop button
7	Process parameter [Steam temperature], display of the actual value
8	Button [Stop all]
9	Button [Warm-up]
10	Button [Start all]

- Touch a parameter's field to open the on-screen keyboard and set a desired value. The input values are accepted without a security prompt:
 - Vacuum setpoint (Adjustment Range: 1 – 1,400 mbar)
 - Hysteresis setpoint (1 – 50 mbar)
 - Rotation speed (6 – 160 rpm) (see also section „Set rotation speed“ on page 110)
 - Heating bath temperature water (20 – 100 °C) or oil (20 – 180 °C) (see also section „Set the heating bath temperature“ on page 111)

WARNING

Due to incorrect/insufficient process settings, there is a risk that impermissible operating values are reached: Risk of personal injuries/damage to property caused by overheating and/or by leaking bath liquid/distillation material!

Only qualified personnel may adjust the process settings of the device.

When determining the pressure values, observe the chemical-physical properties of the respective material to be distilled.

Pay particular attention to the specific safety instructions for the processing of toxic substances!

The parameters [Vacuum], [Rotation] and [Bath Temperature] are assigned a Start/Stop pair of buttons. Touch these buttons to individually activate or deactivate the associated process (build vacuum/maintain level, start/stop rotation, bath heating on/off).

Touch the [Warm Up] button to activate the system warm up function.

Touch the [Start All] button to start all processes with the displayed set values at the same time.

Touch the [Stop All] button to stop all processes at the same time.

The current settings can be stored and loaded at any time as Favorite in the system memory using the [Save] button, see section „The timer function can be stopped at any time with [Stop All]. In this case, the running process is canceled!“ on page 112.

Operating mode AUTOaccurate

In this operating mode, the vacuum pressure is regulated automatically and depending on the temperature. The rotation speed and heating bath temperature are manually set by the operator as in the SETpressure mode.

AUTOaccurate function

The temperature-dependent vacuum control feature AUTOaccurate is based on two values:

- Inlet temperature cooling liquid: this value is measured with a temperature sensor (PT 100) at the chiller inlet on the rear of the device (see section „Mechanical design“ on page 76)
- AUTOaccurate temperature: This temperature is measured at the AUTOaccurate sensor in the condenser

In order to use the function, a value ΔT must be defined in the device control (parameter [Auto Accurate], see the following figure). This value is added to the inlet temperature of the cooling liquid.

The sum ΔT + inlet temperature cooling liquid serves as a reference value for the shutdown of the vacuum pump. This means that as soon as this temperature value is measured at the AUTOaccurate sensor in the condenser after the process has started, the evacuation stops.

As soon as the temperature at the AUTOaccurate sensor has subsequently decreased by a defined value again, the vacuum pump is switched on and the system will be evacuated again until the reference value at the sensor is reached again (with relatively increased vacuum pressure) and the vacuum pump switches off again.

In this way, a consistent performance is achieved in the processing of solvents or mixtures throughout the distillation process.

The described process steps are repeated until a defined maximum vacuum pressure is reached (input value [Vacuum & Steam temperature] → [end p]). The distillation process is then switched off.

PRACTICAL TIP

Before starting the process in AUTOaccurate mode, ensure that all set process temperatures are reached in a stable manner (observe display values!). Only then can the best possible and constant distillation results be achieved!

A warm-up time of 15 minutes is recommended as a guide. Use the device's warm-up function (function button [Warm Up] in the footer of the start page).

To open the start window of this mode, touch the [Auto Accurate] button on the start window:

Legend

1	AUTOaccurate activated
2	Actual value of the coolant temperature of the PT-100 sensor
3	Button [Save]
4	Process parameter [Max. vacuum pressure], actual value of the steam temperature
5	Process parameter [Rotation], display of actual and set value, Start/Stop button
6	Process parameter [Bath temperature], display of actual and set value, Start/Stop button
7	Process parameter [ΔT], actual value of the AUTOaccurate
8	Button [Stop all]
9	Button [Warm-up]
10	Button [Start all]

Touch a parameter's field to open the on-screen keyboard and set a desired value. The input values are accepted without a security prompt:

- Maximum vacuum pressure [end p] (setting range: 1 – 1,400 mbar)
- Rotation speed (6 – 160 rpm) (see also section „Set rotation speed“ on page 110)
- Heating bath temperature water (20 – 100 °C) or oil (20 – 180 °C) (see also section „Set the heating bath temperature“ on page 111)
- Input value ΔT (0.1 – 30 °C)

WARNING

Due to incorrect/insufficient process settings, there is a risk that impermissible operating values are reached: Risk of personal injuries/damage to property caused by overheating and/or by leaking bath liquid/distillation material!

Only qualified personnel may adjust the process settings of the device. When determining the pressure values, observe the chemical-physical properties of the respective material to be distilled. Pay particular attention to the specific safety instructions for the processing of toxic substances!

The hysteresis value for the AUTOaccurate operation cannot be adjusted. Factory setting: 0.5 °C.

The parameters [Rotation], [Bath Temperature] and [Auto Accurate] are assigned a Start/Stop pair of buttons. Touch these buttons to individually activate or deactivate the associated process (start/stop rotation, bath heating on/off, start/stop AUTOaccurate).

Touch the [Warm Up] button to activate the system warm up function.

Touch the [Start All] button to start all processes with the displayed set values at the same time.

Touch the [Stop All] button to stop all processes at the same time.

The current settings can be stored and loaded at any time as Favorite in the system memory using the [Save] button (see section „The timer function can be stopped at any time with [Stop All]. In this case, the running process is canceled!“ on page 112).

Distillation under atmospheric pressure

In certain cases, distillation under atmospheric pressure (low boiling) may be necessary. Refer to the specific safety and operating instructions in this section.

During atmospheric pressure distillation, only the flask rotation and the heating bath are started manually in the [SETpressure] operating mode after filling the evaporation flask.

The vacuum function must remain switched off, otherwise the vacuum pump will be switched on and vacuum will be built up in the system.

CAUTION

During distillation under atmospheric pressure, excess pressure may occur in the system. There is a risk that, in a high system pressure situation, gas or liquid will escape into the ambient air in the area of the seals and joints. Risk of personal injuries/damage to property

There is a risk of bursting under extreme pressure conditions!

Carefully compare the process parameters with the prevailing permissible ambient conditions and the chemical-physical properties of the respective material to be distilled (refer to data sheets)!

Pay particular attention to the specific safety instructions for the processing of toxic substances!


Only operate the device with the guard hood closed or under a closed laboratory fume hood.

Use appropriate personal protective equipment (PPE).

Procedure

- In the system settings, deactivate the [Leakage Pressure] stop criterion, see section „Menu item Safety Settings“ on page 99. Touch the icon [Function activated] of the entry [Leakage Pressure]. The icon changes to [Function deactivated].

- Open the start page again and switch to the operating mode SETpressure (see section „Operating mode SETpressure“ on page 105).
- Adjust the parameters [Rotation] and [Bath Temperature] taking into account the characteristics of the educt and the prevailing ambient conditions.
- Start the flask rotation and the heating bath using the associated start/stop buttons one after the other.

Be careful not to start the vacuum function manually!

Do not start the distillation process with the [Start All] function, as this automatically activates the vacuum function.

Check the pressure in the evaporation flask at regular intervals (display value process parameter [Vacuum]).

Set rotation speed

The procedure for adjusting the rotation speed is described in detail in sections „Operating mode SETpressure“ on page 105 and „Operating mode AUTOaccurate“ on page 106. Please also observe the information in this section!

WARNING

When operating the device without a guard hood, there is a risk of contact with rotating parts. Loose clothing, jewelry and open hair can be pulled in!

- If possible, operate the device with the guard hood closed or under a closed laboratory fume hood.
- Never touch rotating parts with your hands. Do not bend over rotating parts.
- Wear tight-fitting work clothes, remove jewelry before starting work, tie long hair under a hair net or under a suitable head covering.

At high rotation speeds, the heating bath liquid can spray away from the evaporation flask by forming film and/or slop out of the heating bath. There is a risk of scalding and/or contamination.

- If possible, operate the device with the guard hood closed or under a closed laboratory fume hood.
- Adjust the rotation speed of the evaporation flask according to the recommendations in this section.
- Use appropriate personal protective equipment (PPE).
- Clean dirty surfaces immediately!

Even when the guard hood is closed, there is a risk that the heating bath liquid will spill out of the heating bath at high rotational speeds. There is a risk of scalding and/or contamination.

- Adjust the rotation speed of the evaporation flask according to the recommendations in this section.
- Use appropriate personal protective equipment (PPE).
- Clean dirty surfaces immediately!

Recommended rotation speed

The speed of the evaporation flask affects the rate of distillation that can be achieved during the distillation process: The higher the rotation speed, the faster and more gentle the distillation material is processed.

Observe the following recommendations when using water as heating bath liquid:

Volume evaporation flask	Maximum recommended speed
20 liters	120 rpm
10 liters	140 rpm
6 liters	160 rpm

Due to the possible viscosity and volume range (depending on the oil type and the set heating bath temperature), no generally valid recommendation for the rotation speed of the evaporation flask can be given when using oil as heating bath liquid.

Set the heating bath temperature

The procedure for adjusting the heating bath temperature is described in detail in sections „Operating mode SETpressure“ on page 105 and „Operating mode AUTOaccurate“ on page 106 . Please also observe the information in this section!

CAUTION

The heating bath is designed for a maximum operating temperature of 210°C (using silicone oil). From a liquid temperature of 50 °C, there is a risk of injury in the event of contact.

At high rotational speeds, highly heated heating bath liquid can spray from the evaporation flask and/or swap out of the heating bath. There is a risk of scalding.


If possible, operate the device with the guard hood closed or under a closed laboratory fume hood.

Use appropriate personal protective equipment (PPE).


Timer function

The timer function can be used to define a specific time period for a process sequence. At the same time, it is possible to specify which individual functions are stopped after a defined time has elapsed.

- Touch the [Timer...] button on the start window to open the [Set Timer] window:

- Touch the [Days], [Hours], [Minutes] and [Seconds] input fields to open the on-screen keyboard and enter the desired value (max. 99 days, 23 hours, 59 minutes, 59 seconds).
- In the selection [When Timer Expires:-], specify which functions should be stopped after the timer expires (in the example above, all).
- Touch the [Close] button to dismiss all entries and return to the start window.
- Touch the [Start] button to start the timer and thus the process flow. The elapsed process time is visualized by means of a progress display on the start page:

When the timer expires, the message **Timer expired** appears. At the same time, the [Set Timer] window is opened again. Here you can read when the last process sequence and which individual functions were terminated.

The timer function can be stopped at any time with [Stop All]. In this case, the running process is canceled!

Favorites and Profiles

The set process parameters can be stored in the system memory as favorite (process parameter data set) in any operating mode. For the SETpressure operating mode, it is also possible to store certain process sequences as profile.


The procedure for storing and managing favorites and profiles is described in detail in the following sections.

Favorites mode

- Switch to the required operating mode and define the process parameters, see sections „Operating mode SETpressure“ on page 105 or „Operating mode AUTOaccurate“ on page 106.
- Touch the [Save] button:

- The [Save favorite] dialog window appears. Touch the [Please enter name] input field and use the on-screen keyboard to specify a name for the favorite.

- Confirm the entry with [Save] to save the new process parameter data record as a favorite.
 - The list of available favorites appears again.
 - If necessary, create additional favorites as described above ([New] button).


or

- To save the new process parameter data record as a favorite and load it immediately, finalize the entry with [Confirm].
- Touch the [Delete] button to eliminate the opened favorite from the system memory.
- Touch the [Return] button to return to the list of available favorites without applying any changes.

Create/edit favorites

The stored process parameter data records can be managed in the [Favorites] menu:

- Switch to the desired operating mode.
- From the start window, touch the menu button to open the main menu.
- Use [Favorites] to open the list of available favorites:

- Touch the [New] button to create a new record.

or

- Touch the arrow button of an entry to open the corresponding record.
- Touch the [Edit] button in the footer to edit the new/existing favorite (see sections „Operating mode SETpressure“ on page 105 or „Operating mode AUTOaccurate“ on page 106).

- Load the record with [Apply].
- Touch the [Return] button to return to the list of available favorites without applying any changes.

Load favorite


- Switch to the desired operating mode.
- From the start window, touch the menu button to open the main menu.
- Use [Favorites] to open the list of available favorites (see previous section).
- Select the desired data record and confirm with [Apply].


- Then start the process with the loaded values with [Start All].

Live view

→ From the start window, touch the [Show diagram] button to open the process diagram:

In the process live view, the parameters of the current process are visualized in real time:

The orange lines indicate the overall course.

The white colored area shows the previous progress (corresponding progress display).


The current process can be manually canceled at any time. To do this, touch the field with the name of the current favorite (at the top left of the start page) and confirm the security prompt.

Create/edit process profiles

User-specific process sequences can be stored as profiles in the SETpressure operating mode:

- Switch to the Set Pressure operating mode.
- From the start window, touch the menu button to open the main menu.
- Open the [Favorites] menu and switch to [Profiles] to open the list of available profiles:

Edit an existing profile

- Touch the arrow button of an entry to open the corresponding record.
- Touch the [Edit] button in the footer to edit the profile (see also the following section „Create a new profile“ on page 115).

- Confirm any changes with [Apply].
- Touch the [Return] button to return to the list of available profiles without applying any changes.

Use existing profile as template

- Touch the arrow button of an entry to open the corresponding record.
- Touch the [Edit] button in the footer, and then touch [Copy] on the input screen to duplicate the profile.
- Use the on-screen keyboard to specify a label for the new profile and adjust the values as desired (see also the following section „Create a new profile“ on page 115).
- Confirm any changes with [Apply].
- Touch the [Return] button to return to the list of available profiles without applying any changes.

Create a new profile

- Touch the [New] button to create a new record.

- In the [Settings] input screen, first define a hysteresis value and the shutdown behavior for the individual functions (see section „Timer function“ on page 111):

Set Pressure Profile
name

Table Settings Graphs

Hysteresis 10 mbar

When Profile ends:

Vacuum Will stop

Rotation Will stop

Bath Will stop

Lift Will stop


→ Switch to the [Table] input screen:

Set Pressure Profile
abc

Table Settings Graphs

Time [hh:mm]	Vacuum [mbar]	Rotation [rpm]	Bath Temperature [°C]
0:00	800	0	20
0:01	<input type="checkbox"/> 800	<input type="checkbox"/> 60	<input type="checkbox"/> 60
0:02	<input type="checkbox"/> 800	<input type="checkbox"/> 80	<input type="checkbox"/> 80
0:03	<input type="checkbox"/> 1200	<input type="checkbox"/> 80	<input type="checkbox"/> 80
0:04	<input type="checkbox"/> 1200	<input type="checkbox"/> 120	<input type="checkbox"/> 90
0:05	<input type="checkbox"/> 1200	<input type="checkbox"/> 120	<input type="checkbox"/> 90


- Define one process step per line here: Touch an input field ([Time], [Vacuum], [Rotation], [Bath Temperature]) to open the on-screen keyboard.
- Set the desired value and confirm each entry with [OK].
- Insert new line: Highlight the line above which you want to add a new line and touch the [New Line] button.
- Delete Row: Highlight the row you want to delete and touch the [Delete Row] button.
- After you have completed the parameterization, switch to the [Diagram] screen:

- Use the [Zoom+/-] buttons to adjust the time axes stretch for the programmed ramp as desired.
- Confirm with [Save] to save the new profile in the system memory without loading it directly.
- Touch the [Apply] button to finalize the entry, to save the new profile in the system memory and to load it directly.
- Touch the [Return] button to return to the list of available profiles without applying any changes.

Load process profile

- Switch to [Set Pressure] operating mode.
- From the start window, touch the menu button to open the main menu.
- Open the [Favorites] menu and switch to [Profiles] to open the list of available profiles:

- Touch the arrow button of the desired profile and confirm with [Apply] to load the profile:

- Touch the [Return] button to return to the list of available profiles without applying any changes.

When the profile is completely processed, the message **Profile processed** appears. Here you can read when the profile was processed and which individual functions were terminated.

The processing of profiles can be canceled manually at any time with [Stop All].

Data logging

Using the data logging function, various process data of the Hei-VAP Industrial can be recorded in the form of measurement series. The internal database offers a total of 20 slots with a capacity of 72,000 measuring points per measurement series for this.

Data is recorded at individually defined intervals. With the shortest recording interval of one second, a maximum total recording time of approximately 27 hours can be achieved. If the Hei-VAP Industrial is operated in combination with an automatic module Hei-VOLUME Distimatic Pro, a recording interval of ≥ 30 seconds is recommended to ensure a recording period of at least three weeks.

When combined with an automatic module Hei-VOLUME Distimatic Pro, the Datalogging function of the Hei-VAP Industrial records the values of the integrated pressure sensor, including the overpressure values reached when draining the distillate out of the evaporator.


While draining the distillate by overpressure, which is controlled by the automatic module Hei-VOLUME Distimatic Pro, the process area (evaporation flask and condenser) is completely isolated by two closing valves and no medium is fed in or discharged.

During this period of up to 30 seconds, the vacuum values of the process area are not recorded!


After draining the distillate, the associated area is brought back to the process vacuum and the bypass is switched on again.

Activating/deactivating data logging

The data logging function is manually activated/deactivated. The function status active/inactive is visualized using the appearance of the [Datalog...] button:

- Touch the [Datalog...] button to open the main window of the data logging function.
 - Depending on the operating state, one of the following status messages is displayed:
Data logging in progress
Data logging stopped

- Status message **Data logging stopped**: activate the data logging function with [Start].
- Status message **Data logging in progress**: Deactivate the data logging function with [Stop].

The actual remaining storage time is displayed in days, hours, minutes, and seconds.

As soon as all slots are occupied, the [Slots] button changes its color. In this case, existing measurement series must be swapped out or deleted, see section „Deleting measurement series“ on page 120.

- Touch the button [Close] to close the window.


Setting the logging rate

The recording rate can be defined in the main window.

- Touch the input fields [Data logging every XX minutes/seconds] one after the other and enter the desired values via the keyboard.

Logging capacity

A total of 20 slots are available for recording a maximum of 100,000 measuring points. The measuring points can be recorded on the available measuring series (slots) as desired. When the maximum recording capacity is reached, the following message appears:


- Confirm the message and proceed as described in the following section "Delete measurement series" to delete measurement series.

Data list

The data list can be opened from the main window of the data logging function/button [Slots] or via the [Settings] menu/**Show Datalog Series**. The data list contains all available measurement series. Each entry consists of the time stamp and the number of measuring points recorded.

Deleting measurement series

- To delete an individual entry, touch the corresponding delete key (symbol ✕) and confirm the security prompt.
- To delete all entries, touch the [Delete All] button in the footer and confirm the security prompt.

To delete entries, the data logging function must be deactivated!

During the deletion process of all entries no individual deletion or saving processes are possible and the data list cannot be closed.

Saving measurement series on a USB flash drive

Recorded measurement series can be saved on a USB stick. To do so, connect a VFAT formatted USB flash drive to the USB interface at the underside of the display and open the data list as described.

As soon as a suitable USB stick is connected, a [Save] button is assigned to each list entry. At the same time, the [Save all] button appears in the footer of the main window.

Select a single entry in the list and touch the associated [Save] button to save that entry individually to the connected USB flash drive.

Alternatively, all list entries can be saved simultaneously. To do this, touch the [Save all] button in the footer:

During the saving process of all entries no individual saving processes are possible and the data list cannot be closed.

Progress display when saving measurement series

When the saving process begins, the [Save] button is highlighted in orange. Once the process is completed successfully, the USB symbol will be replaced by a check mark.

The progress indicator is reset when you exit the mask

Data output

All data are saved in CSV and JSON format on the connected USB flash drive.

When opening the CSV file in Microsoft Excel®, the temperature values may be displayed in date format!

In this case, open the menu **File / Options / Extended** in the opened CSV file and deselect the options **Apply separator from operating system**.

Troubleshooting

System message	Possible remedy
Heating bath heating cannot be started	Water level too low, bath heating cannot heat up: Refill water bath.
Emergency stop activated	Press the emergency stop button again, check the parameters and adjust if necessary, restart the process.
Rotation lock active	Rotary knob of the rotation lock is set to lock. Unlock.
Flask support not in start position	Support bars not in end position. Place support bars in end position.
Maximum pressure reached	In the vacuum control mode AUTOaccurate the value "end p" was reached, process was stopped. Distillation process finished!
Guard hood open	Guard hood was opened during a running process. Close the guard hood.
Communication between panel and controller interrupted	Internal COM error. Contact technical service!
Loading	Device operation while charging is in progress. No operator action until the message disappears.
Short circuit AUTOaccurate sensor	Contact technical service!
Overheat protection triggered	Check the level of the heating bath and the setting of the bath liquid selector on the rear of the device. If the heating bath heating unit is defective, contact authorized dealer or technical service.
Bath lift motor failure	Contact technical service!
Bath lift sensor failure	Contact technical service!
Bath mode (water/oil) cannot be changed during operation!	Bath liquid selector was actuated during a process. Move selector to original position.
Set bath temperature exceeded	Reduce the bath temperature and then increase gradually, clean the heating bath. If the heating bath heating unit is defective, contact authorized dealer or technical service.
Short circuit bath temperature sensor	Contact technical service!
Short circuit cooling water temperature sensor	Contact technical service!
Maximum permissible pressure exceeded	Decrease temperature and/or vacuum and increase gradually again, check setting in the Safety Settings menu.
Vacuum sensor broken	Replace vacuum sensor.
Short circuit vacuum sensor	Contact technical service!
Boiling temperature sensor not connected	Connect the boiling temperature sensor. Replace the boiling temperature sensor if defective.
Short circuit boiling temperature sensor	Contact technical service!
AUTOaccurate sensor not connected	Connect the AUTOaccurate sensor. Replace the AUTOaccurate sensor if it is defective.
Bath heating temporarily suspended. Chiller temperature too high	Reduce bath temperature, decrease hysteresis value, connect more powerful chiller.
Bath temperature sensor broken	Contact technical service!
Coolant temperature sensor not connected	Connect the coolant temperature sensor. Replace coolant temperature sensor if defective.
Rotation lock activated during operation	Unlock the rotation lock.

System message	Possible remedy
System leaking	Check the seals and connections, check the correct fit of the PTFE sealing, check the cleanliness and fit of the joints, grease if necessary, check the vacuum line and pump, check the setting in the Safety Settings menu.
Flask drive motor defective. Error Number:...	Contact technical service!
Actual vacuum value in the ventilated state differs from the air pressure.	Installation location is above normal zero, calibrate vacuum sensor.
Heating bath lift does not move up or down	Flask support engaged. Slide outward.
	Guard hood open. Close.
	Lift blocked. Check for blocking objects.
	Gap between housing and heating bath not free. Check for blocking objects.
No flask rotation	If the lift drive is defective, contact technical service.
	Flask support engaged. Slide outward.
	Rotation lock is in lock position. Unlock.
	Guard hood open. Close.
Ventilation valve does not switch	If the flask drive is defective, contact technical service.
Distillation prematurely stopped	Connect the ventilation valve. If the valve is defective, contact authorized dealer or technical service.
Vacuum is not built up	Check end p value, check Set Time parameter.
	Check PTFE seal, replace if necessary. Check sandwich seal, replace if necessary.

All error messages and warning messages are shown in plain text on the display. Follow the instructions on the display.

In case of recurring errors, please contact the responsible sales department or our technical service. Contact address see „Contact information“ on page 127.

Technical specifications

General device data	
Model	Rotary evaporator Hei-VAP Industrial
Dimensions (W × H × D)	998 × 2,225 × 694 mm
Weight	approx. 120 kg, without glassware
Acoustic pressure	< 85 (dB(A)) (in accordance with IEC 61010)
Flask drive	EC motor
Speed range	6 – 160 rpm
Lift travel	180 mm
Lift speed	45 mm/s
Evaporation rates (L/h) ΔT^* 40 °C (in continuous operation)	Toluene 26.0 Acetone 25.9 Ethanol 12.9 Water 4.7
Permissible leakage rate	1 mbar/min
Maximum cooling surface	1.2 m ²
Control panel	detachable, 7 inch touchpanel, control knobs
Overtemperature protection	integrated, mechanical
Degree of protection, device	IP20
Control panel degree of protection	IP45
Electrical data	
Rated voltage	3 × 400 V (50/60 Hz) or 1 × 230 V (50/60 Hz)
Connection	L1+L2+L3+N+PE or L+N+PE
Protection class	I  (IEC 61140)
Overvoltage category	II
Degree of pollution	2
Power input	4,600 W max.
Heating bath, cooling	
Heating bath	Inner Ø: 397 mm, outer Ø: 400 mm Material: V4A steel (1.4404) X2CrNiMo 17-12-2 Volume: 18 L H ₂ O / 14 L oil
Heating power	4,000 W
Temperature range heating bath	20 – 100 °C H ₂ O / 20 – 180 °C oil
Accuracy heating bath control	± 1 °C
Heating bath temperature control	Microprocessor
Required pump capacity	1,5 – 3 (m ³ /h)
Consumption cooling liquid	150 – 250 (L/h)

Permissible ambient conditions	
Storage conditions (recommendation)	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 m asl

* ΔT = difference between heating bath temperature and boiling temperature

Scope of delivery

Item	Quantity	Product no.
Rotary evaporator Hei-VAP Industrial, complete*	1	518-13100-00
Operating instructions	1	01-005-005-58
Guarantee registration / Certificate of decontamination	1	01-006-002-78

*The actual delivery depends on the respective device variant.

Accessories

Detailed information on the available accessories for your device variant can be found on our website at www.heidolph-instruments.com.

Device service

DANGER

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

When cleaning, avoid the penetration of liquids.

Before replacing the fuses, switch off the device and disconnect the power supply cord.

Always replace the two device fuses in pairs with original manufacturer fuses. Further information on available accessories can be found on our website at www.heidolph.com!

After fuse replacement, check the device for a safe condition according to IEC 61010-1.

General cleaning instructions

CAUTION

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

Repairs

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim.

The owner is generally liable for damage caused by unauthorized repairs.

If repairs are required, contact an authorized dealer or our technical service, see „Contact information“ on page 127.

Include the completed certificate of decontamination with every device return, see „Certificate of decontamination“ on page 128.

Maintenance

There are no user-serviceable components in the unit housing. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example), contact our technical service, see „Contact information“ on page 127.

- Check the caps and fasteners of all glass components for correct and tight fit at regular intervals and especially before each start-up of the device.
- Check the level sensor in the heating bath regularly for proper operation, contamination, and limescale build-up to prevent malfunctions.
- Check the PTFE seal of the EASY LOCK flask clamping system at regular intervals for heavy contamination and visible damage.
- Clean the PTFE seal of the EASY LOCK flask clamping system at regular intervals to prevent excessive deposits or encrustations.

CAUTION

At loose/leaky connections, ambient air is sucked in during the vacuum build-up and the required vacuum pressure cannot be reached in a stable manner!

In the event of excess pressure, there is a risk that liquid will escape from the system! The device and surrounding areas can be severely contaminated by leaking media!

Tighten loose caps/fasteners by hand or correct the seating of the retaining clips.

Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
- Properly remove and dispose of residues of substances that are hazardous to health, the environment and biohazardous!

Warranty statement

Heidolph Instruments GmbH & Co. KG provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage www.heidolph.com.

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

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Heidolph Instruments North America

E-mail: service@heidolph.com
www.heidolphNA.com

Heidolph Instruments United Kingdom

E-mail: service@radleys.co.uk
www.heidolph-instruments.co.uk

Local distributors

To find your local distributor please visit www.heidolph.com

Certificate of decontamination

Enclose the certificate of decontamination, duly completed, with your device return.
Submissions without a certificate of decontamination cannot be processed!

CERTIFICATE OF DECONTAMINATION

IN CASE OF RETURNS

Please fill in the required fields.

Note: The sender must package the goods properly and appropriately for transport.

Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12
91126 Schwabach, Germany
Phone: +49 (0) 9122 9920-380
Fax: +49 (0) 9122 9920-19
Email: service@heidolph.de

SENDER

Name _____ First name _____
Company/institution _____ Department _____
_____ Workgroup _____
Address _____
ZC/City _____
Country _____ Phone _____
Email _____

DEVICE DETAILS

Article number _____ Serial no. _____

Ticket number _____

Reason for sending in _____

Has the device been cleaned, decontaminated/disinfected? Yes No (Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous? Yes No (Please mark as applicable)

If yes, with which substances did the device come into contact?

LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product (from serial number on 200067811) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

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Declaration of Conformity
 In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,
 Heidolph Instruments GmbH & Co. KG
 Walpersdorfer Straße 12
 91126 Schwabach / Germany

Product: Laboratory evaporator
 Model:

Hei-VAP Industrial xxx 518-xxxx-xx-x

Description:
 Rotary Evaporators with heating bath and human machine interface

The object of the declaration described above is in conformity with the relevant
 UK Statutory Instruments (and their amendments):

- | | |
|--------------------|---|
| 2008 No. 1597 | <i>The Supply of Machinery (Safety) Regulations</i> |
| 2008 2016 No. 1091 | <i>The Electromagnetic Compatibility Regulations 2016</i> |
| 2012 No. 3032 | <i>The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012</i> |

and complies with the following technical standards :
 EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):
 ProductIP (UK) Ltd.
 8, Northumberland Av.
 London WC2N 5BY

Signed for and on behalf of Heidolph Instruments GmbH & Co. KG
 Walpersdorfer Straße 12, 91126 Schwabach / Germany

Schwabach, 01.09.2021

Wolfgang Jaenicke
 Managing Director

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices by Heidolph Instruments GmbH + Co. KG fulfill the requirements of the EU directive 2011 / 65 / EU (RoHS) and the amended of directive 2015 / 863.

All devices are compatible with the requirement of the RoHS for the following materials:
Max. 0,01% of the weight in homogeneous material for cadmium and max. 0,1 % of the weight in homogeneous material for lead, mercury, hexavalent chromium, polybrominated biphenyl, polybrominated diphenyl ether, Di (2-ethylhexyl) phthalate, butyl benzyl phthalate, dibutyl phthalate, diisobutyl phthalate.

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 22.07.2019

Stefan Peters
Vice President Marketing, Innovation & Technology

Marcell Sarré
Vice President Quality Management & Technical Service

China RoHS DECLARATION OF CONFORMITY

Heidolph Instruments GmbH & Co.KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Instruments GmbH & Co.KG products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
- O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
- X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Instruments GmbH & Co.KG may enter into further devices or can be used together with other appliances.

With these products and appliances in particular, Heidolph Instruments GmbH & Co.KG will not take responsibility for the EFUP of those products and appliances.

Schwabach, 06.08.2021

Wolfgang Jaenicke
Chief Executive Officer CEO

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Instruments GmbH & Co. KG PRODUCTS
© Heidolph Instruments GmbH & Co. KG