Safe
Heating and Mixing

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Operating instructions
Notice d'instructions

EKT Hei-Con

 **heidolph**
research made easy

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Seite 4 – 22

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Translation of the original instructions
Page Page 26 – 44

Contents

Introduction

About this document	26
Typographic conventions	26
Copyright protection	26

General notes

Basic product information	27
Guidelines applied, product certification	27
California Residents	27
Residual risk	27
Intended use	27
Reasonably foreseeable misuse	27
Transportation	27
Storage	28
Acclimatization	28
Permissible ambient conditions	28

Safety

General safety information	29
Other regulations	29

Device description

Device components	30
Display elements and controls	31
Principle of operation	31
Safety functions	31
Power supply, switching the device on and off	31

Operation

Measurement setup, device connection	32
Operation	34
System test	34
Settings	34
Starting the heating process	36
Control mode	36
Display of the controller status	37
Safety switch-off	37
Measurement range monitoring	38
Sensor monitoring	38

Troubleshooting

Troubleshooting.....	39
----------------------	----

Appendix

Technical Specifications	40
PIN assignment of five pin connection plug	40
List of the factory settings	41
Scope of delivery	41
Accessories.....	41
Device service.....	42
General cleaning instructions	42
Repairs	42
Maintenance.....	42
Disposal.....	43
Contact information Heidolph international.....	43
Warranty Statement	43
Certificate of decontamination	44

About this document

These operating instructions describe all the functions and the operation of EKT Hei-Con type electronic temperature controllers. The operating instructions are an integral part of the described device.

Typographic conventions

Standardized symbols, highlighting elements, and signal words are used in this document to identify warnings, cautions, important information, and special text contents.

Symbol	Signal word / explanatory note
	<p>Warning symbols in combination with a signal word indicate dangers:</p> <p>DANGER Indicates an immediate dangerous situation which, if not avoided, will result in death or serious injury.</p> <p>WARNING Indicates a potential danger which, if not avoided, will result in serious injury.</p> <p>CAUTION Indicates a potential hazard which, if not avoided, damage to property and minor to moderate injuries can occur.</p>
	<p>Mandatory signs are used to indicate important and useful information on handling a product.</p> <p>This information is used to ensure operational safety and to maintain the value of the product.</p>
[GUI]	<p>Parameter designations, display texts, and device labels are highlighted in text and tables in a typographic manner to facilitate the assignment on the device.</p>
→	<p>The arrow symbol indicates instructions to be followed in order to ensure the operational safety when handling the product.</p>

Copyright protection

This publication is protected by copyright and intended for internal use by the purchaser of the product only.

No part of this publication may be transmitted or reproduced in any form, by any means, without the prior written consent of the copyright owner Heidolph Instruments GmbH & Co. KG. Any violation is subject to compensation for damages.

Basic product information

Guidelines applied, product certification

CE Marking

The device complies with the following standards:

- 2014/35/EU (Low Voltage Directive)
- 2014/30/EU (EMC Directive)

California Residents

Important information for California residents regarding Prop 65. Please visit www.P65Warnings.ca.gov for more information.

Residual risk

The device was designed and manufactured in accordance with the latest technical standards at the time of development and the recognized safety regulations. During installation and use, as well as during maintenance work, repairs and cleaning, there are nevertheless certain residual risks associated with the device described.

These are identified and described at the appropriate points in this document.

Intended use

The device has been developed by the manufacturer as an electronic temperature controller for liquids. The device is solely intended for use in combination with magnetic stirrers of the type Hei-Standard and Hei-PLATE Mix'n'Heat Core made by Heidolph Instruments.

Any other use of this device is not considered as intended!

Reasonably foreseeable misuse

Additional measures may be necessary for use under conditions or for purposes deviating from the intended use, and/or specific guidelines and safety regulations must be observed. Corresponding requirements must be evaluated and observed by the owner/operating company in each individual case.

Compliance with and implementation of all relevant guidelines and safety measures for the respective field of application is within the sole responsibility of the owner/operating company.

All risks resulting from improper use are solely borne by the operator.

The device may exclusively be operated by authorized and instructed personnel.

Training and qualification of the operating personnel as well as ensuring that the device is operated with responsibility are the sole responsibility of the operator!

Transportation

During transport, avoid severe shocks and mechanical stresses that can cause damage to the device.

Keep the original packaging in a dry and protected place for later use.

Storage

Always store the device in its original packaging. To protect against damage and unreasonable material aging, store the device in a dry environment that should be as temperature-stable and dust-free as possible.

Acclimatization

After each transport and after storage under critical climatic conditions (e.g. high temperature difference between inside and outside), allow the device to acclimatize at room temperature for a minimum of two hours to prevent possible damage from condensation before putting it into operation at the place of use.

If necessary, extend the acclimatization phase if the temperature differences are very high.

Permissible ambient conditions

The device is designed for indoor use only. When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

General safety information

- Only use the device in accordance with the regulations for intended use ("Intended use" on page 27).
- Before commissioning and using the device, familiarize yourself with all the safety regulations and guidelines for occupational safety applicable at the place of use and observe them at all times.
- Only operate the device if it is in perfect technical condition and there is no visible damage.
- Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Instruments.
- If there is missing or misleading information on the device or on occupational safety, contact the responsible safety specialist or our technical service.
- Do not make any unauthorized changes or modifications to the device!
- Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!
- Rectify malfunctions or faults on the device immediately.
- Shut down the device if it is not possible to eliminate the malfunction or rectify the fault immediately.
- Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.
- Wipe off any fluid that may have spilled on the device immediately.

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty against Heidolph Instruments.


The operator is solely liable for all damage resulting from unauthorized changes or modifications to the unit, from the use of unauthorized or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Device components

The electronic temperature controller is made up of the following components:

- Control unit [1]
- Connection cable to the magnetic stirrer (spiral cable, five-pin, hard-wired, with connection plug in acc. with DIN 45322) [2]
- Bracket [3] for the control unit with two retaining clips \varnothing 10 mm for fixing on a wall grid or similar and two retaining clips \varnothing 13 mm for fixing on the stay bar of a magnetic stirrer:

- Sensor (hard-wired) [4]

Display elements and controls

The control unit has the following display elements and controls:

- LCD multi-functional display [1]
- **max** button [2]: press the button briefly to display the set maximum temperature.
- **Increase value** arrow key [3]: press the button briefly to increase the value in increments; keep the button pressed to increase the displayed value quickly.
- **Reduce value** arrow key [4]: press the button briefly to reduce the value in increments; keep the button pressed to reduce the displayed value quickly.
- **set** button [5]: press the button briefly to confirm and accept a displayed (context dependent) value.

Principle of operation

The electronic EKT Hei-Con type temperature controller was specially designed for applications with liquids and replaces conventional mercury contact thermometers in the laboratory.

The microprocessor-controlled electronics provide the control functions and ensure a high degree of linearity and measurement accuracy. The integrated FUZZY logic ensures an optimum heating rate without overshooting the temperature.

Thanks to the remote sensor, the control unit can be positioned outside the area of high heat development and steam formation.

Safety functions

The electronic temperature controller is equipped with the following safety functions:

- Sensor short-circuit protection
- Sensor break protection
- Safety and control circuit to DIN 12878/Class 1 and 2
- Safety switch-off if the set temperature (MAX temperature) is exceeded by 25 °C
- Safety switch-off if the temperature sensor is not immersed in the liquid within a certain time

During operation, the sensor and the measuring line are tested continuously for interruptions and short circuit. In case of an error, the heating switches off and the error is shown on the display.

Power supply, switching the device on and off

The device is supplied with the necessary 8 – 15 V DC voltage via the five-pole connection cable of the control unit from the connected magnetic stirrer.

The device does not have its own on/off switch but instead is switched on and off at the same time as the connected magnetic stirrer.

Measurement setup, device connection

The control unit can be connected to a wall grid or similar or can be fixed onto the stay bar of the magnetic stirrer using the bracket supplied.

- Fit the retaining clips on the rear rail of the bracket for the control unit.
 - Use the two retaining clips \varnothing 10 mm to fix the bracket onto a wall grid or similar.
 - Use the two retaining clips \varnothing 13 mm to fix the bracket on the stay bar of a magnetic stirrer (supplied with the magnetic stirrer, for assembly of the stay bar see operating instructions of the magnetic stirrer). The following figure shows the back of the bracket for the control unit with fitted retaining clips:

- Insert the control unit into the preassembled bracket.
- Connect the connection plug of the control unit's connection cable with the five-pin connection socket [1] on the back of the magnetic stirrer:


- Insert the sensor in the preassembled sensor clamping system (optional accessory) and route the connection cable along the setup as close-fitting as possible. Make sure that the connection cable is not exposed to any mechanical loads (bending, tension).

- Align all components according to the external circumstances and according to the vessels/attachments used.
- The figure on the right shows the example of a set up consisting of
 - Hei-PLATE Mix'n'Heat Core magnetic stirrer [1]
 - Sensor clamping system [2] (optional accessory)
 - EKT Hei-Con [3], fixed on the stay bar of the magnetic stirrer

When using the optional sensor clamping system:

- Secure the sensor by means of the supplied nozzle in the sensor clamping system:

CAUTION

Damage of the device cables by the effect of heat

Fix and secure the cable of the electronic EKT in such way on the setup that a contact with the heating plate will be avoided during operation.

Operation

WARNING

Risk of overheating and measurement errors

Note that if the EKT is connected, the temperature of the magnetic stirrer's heating plate is adjusted to the sensor automatically. As long as the sensor is not immersed in the sample, the ambient air is used as the reference temperature. The heating plate of the magnetic stirrer can therefore be heated up to 300 °C (maximum heating plate temperature of devices of the type Hei-PLATE Mix'n'Heat Core) without being noticed!

- Make sure that during operation, the sensor is always immersed by at least 20 mm in the medium to be measured.
- Always use the necessary personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to process samples from a temperature of 50 °C.
- When processing flammable samples, make sure that the set maximum heating plate temperature is at least 25 °C below the flashpoint of the sample.

WARNING

Risk of burns at the sensor

In operation, the sensor can heat up to the medium temperature, depending on the intensity and the duration of exposure.

- Avoid direct skin contact with the metallic lance of the sensor when removing it from the medium.
- Always use the necessary personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to process samples from a temperature of 50 °C.
- Leave the sensor to cool down sufficiently before further use or cleaning.

System test

Directly after switching on, all segments of the control unit's digital display are shown for two seconds. The device is then ready and the temperature currently measured at the sensor is shown on the display.

Settings

The unit settings are accessed via the device's menu navigation. Open the main menu as follows:

- Make sure that the device is connected to a magnetic stirrer.
- Make sure that the magnetic stirrer and thus the device is switched off.
- Press the **set** button on the control unit and keep it pressed.
- Switch on the magnetic stirrer:
 - The **Hei/Con/Vx.xx/User/Menu/Unit** value is shown on the control unit's display.
- Release the **set** button.
 - The first menu item appears on the display, see following table.
- Use the two **Increase value** and **Reduce value** arrow keys to switch between the menu items.
- Change the respective setting by pressing repeatedly the **set** button.

Menu item	Setting options
Unit	<p>This menu is used to select the temperature display: Press set to switch between the options °C (temperature display in degrees Celsius, factory settings) and °F (temperature display in degrees Fahrenheit).</p>
SSOP	<p>In this menu the sensitivity of the safety switch-off can be adjusted (see also section "Safety switch-off" on page 37): Open the menu and press repeatedly the set-button to adjust the minimum required slope of the temperature curve per time interval during the heating phase in a range of 2-20 [equals 0,2 K – 2,0 K]. The time interval can be adjusted in the menu SINT. The value is increased or reduced by 1 [equals 0.1 K] with each press of the button.</p>
FAST / FINE	<p>This menu is used to select the control mode: In the control mode FINE the set temperature can be regulated gradually and with minimal overshoot. This mode is particularly suitable for temperature-sensitive samples and for processing samples in small containers (beakers and heating baths smaller than 1 L). In the control mode FAST the set temperature can be regulated fastly which can result in a stronger overshoot of the temperature curve. This mode is particularly suitable for temperature-insensitive samples and for processing samples in larger containers (beakers and heating baths larger than 1 L).</p>
SINT	<p>In this menu the time interval for the safety switch-off can be adjusted, cf. parameter SSOP. Open the menu and press repeatedly the set-button to adjust the time interval in a range of 10 – 60 seconds. The value increases or reduces by 10 s with each press of the button.</p>
SSOT	<p>In this menu, the safety switch-off of the magnetic stirrer can be adjusted (see section "Control mode" on page 36) Open the menu and press repeatedly the set button to adjust the safety temperature in a range of 10 – 25 °C. The value increases or reduces by 1 °C with each press of the button.</p>


→ Press the **max** button two times to exit the device menu and to switch to the normal mode.

Starting the heating process

The control process can only be started if the temperature of the magnetic stirrer's heating plate is below 50 °C. Allow the magnetic stirrer to cool if necessary!

After switching on the magnetic stirrer, in normal mode the actual [1] and the last set temperature [2] to be set are shown on the display:

Procedure

- Make sure that the sensor is immersed in the sample by at least 20 mm.
- Switch on the heating function of the magnetic stirrer.
- If necessary, use the two arrow keys to adjust the set temperature. This value can be changed at any time.
- Confirm the displayed/adjusted set temperature with **set** to start the control process with the selected parameter settings.

After switching on the device, after a voltage interruption (power outage, error switch-off) and after each adjustment, the displayed set value must be re-confirmed and the heating of the magnetic stirrer must be switched on again.

Otherwise the control process is not started!

Control mode

To avoid overshooting of the temperature when heating the sample as far as possible, the EKT Hei-Con uses FUZZY logic as control optimization: the heating of the magnetic stirrer is switched off before reaching the set temperature so that the heating is then only provided by the heating plate's heat supply.

The control range of the EKT Hei-Con is limited by the maximum heating plate temperature set at the magnetic stirrer: if the heating plate temperature is limited, e.g. to 90 °C at the magnetic stirrer, a higher temperature cannot be set at the EKT Hei-Con (error message: **NO SENSR**).

Set the temperature to around twice the value of the set temperature to heat small quantities (< 500 mL). The maximum temperature should only be used for large quantities.

At set values above 150°C, a slight control deviation can occur. This can be compensated for by setting a higher set value.

Display of the controller status

The heating function is switched on if the set value is above the actual value. Heating mode is indicated by a small rectangle, which wanders from left to right at the top of the control unit's display.

The symbol goes out as soon as the heating is switched off prematurely by the FUZZY logic.

Safety switch-off

In case of a safety switch-off of the magnetic stirrer, **NO SENSR** or **MAX Temp.** appears as information on the display. A safety switch-off occurs in the following cases:

Switch-off criterion	Conditions
Sensor outside the sample when switched on	<ul style="list-style-type: none"> → $T_{act} < T_{set}$ and $\Delta T > -20$ K → Temperature curve gradient/time interval during the heating process is below the SSOP requirements (0.2 K – 2.0 K) in positive or negative direction in three consecutive time intervals. <p>The monitoring is not activated until 300 seconds after starting due to the control section's own measurement. If the sensor is not immersed, the safety switch-off of the magnetic stirrer occurs three minutes after monitoring is activated. The NO SENSR message is shown on the display of the control unit.</p>
Sensor removed from heated sample	<ul style="list-style-type: none"> → $T_{act} < T_{set}$ and $\Delta T > -5$ K → After three consecutive time intervals of three seconds each, the temperature drop is > 0.7 K in the negative direction. <p>If the sensor is removed from the sample, the safety switch-off of the magnetic stirrer occurs after three minutes have expired. The NO SENSR message is shown on the display of the control unit.</p>
MAX Temp. safety temperature reached	<ul style="list-style-type: none"> → $T_{act} < T_{set}$ and $\Delta T > +25$ K → After reducing the set temperature, this criterion does not become active again until $T_{act} < T_{set}$ and $\Delta T > +5$ K. <p>On reaching the set safety temperature (query with max button), the safety switch-off of the magnetic stirrer occurs. The MAX Temp. message is shown on the display of the control unit.</p>

The error messages on the electronic contact thermometer (EKT) can only be acknowledged if the heating plate temperature of the connected magnetic stirrer is below 50 °C (residual heat indicator OFF, if necessary, see operating instruction of the magnetic stirrer!).

To acknowledge an error message on the EKT independently, the cable connection between the EKT and the magnetic stirrer has to be disconnected.

Measurement range monitoring

If the measurement range (+350 °C) is exceeded, the **HI** message appears on the control unit's display.

If the temperature falls below the measurement range (-80 °C), the **Lo** message appears on the control unit's display.

Sensor monitoring

In the event of a sensor interruption, the **HI** message appears on the control unit's display.

In case of a short-circuit in the sensor, the **Lo** message appears on the control unit's display.

Troubleshooting

The following table includes possible failures and corresponding corrective measures:


Malfunction	Possible cause/ remedy
No heating function	→ EKT Hei-Con not or incorrectly fitted: Check connection plug for correct fit.
	→ Safety switch-off after sensor defect (HI or Lo display), check the sensor and cable for visible damage/short circuit, contact technical service.
	→ Max ERROR display after safety switch-off on reaching the set safety temperature: Check the device for damage, contact the technical service. → Maximum allowable heating plate temperature exceeded at the magnetic stirrer: Leave the magnetic stirrer to cool, check settings for plausibility, restart.
NO SENSR display	→ Sensor not immersed, immerse sensor in the sample (at least 20 mm).
	→ The set value for the magnetic stirrer set temperature is lower than the set EKT Hei-Con set temperature: Check settings for plausibility and adjust.

If a fault cannot be rectified with the described suggestions, please contact an authorized sales representative or our technical service (see section "Contact information Heidolph international" on page 43).

Technical Specifications

Model	EKT Hei-Con
Control unit dimensions (W × H × D)	126 × 54 × 22 mm
Weight	250 g
Protection class	IP54
Temperature measurement range	-80 – +350 °C
Resolution	1 °C
Measurement accuracy	± digit
Measuring interval	1.0 s
Safety switch-off temperature	25 K above set value
Sensor	Pt1000, V2A tube, Ø 3 mm
Sensor length	210 mm
Sensor cable length	700 mm
Sensor minimum immersion depth	20 mm
Switching power at the output	12 V / 3 mA
Supply voltage	8 – 15 V
Connection cable to magnetic stirrer	Spiral cable, max. length approx. 1.5 m, with five pin DIN connector
Ambient temperature	-5 – +60 °C
Storage temperature	-30 – +70 °C
Housing material	Polyamide

PIN assignment of five pin connection plug

List of the factory settings

Parameter	Factory settings	Range
Temperature control nom.	15 °C	0 – 300 °C 32 – 572 °F
Temperature [UNIT]	°C	°C / °F
Control parameter SSOP	6	2 – 20
Control mode	FINE	FINE / FAST
Control parameter SINT	60	10 – 60
Safety switch-off temperature SSOT	Solltemperatur +25 K	10 – 25

The factory settings of the device should only be adjusted if, in some cases, it is unavoidable from a process-specific view!

Scope of delivery

Item	Quantity	Product no.
EKT Hei-Con complete with bracket	1	509-88000-00
Retaining clip Ø 10 mm for bracket	2	11-008-010-19
Retaining clip Ø 13 mm for bracket	2	11-008-010-18
Nozzle	1	23-09-03-01-47
Operating instructions	1	01-005-004-61
Warranty registration	1	01-006-002-78

Accessories

Detailed information on the available accessories can be found on our website at www.heidolph-instruments.com.

In case of need, contact an authorized dealer or our technical service, see "Contact information Heidolph international" on page 43.

Device service

When carrying out service work on the device (cleaning, maintenance, repair), observe the general instructions and safety information described in this section.

General cleaning instructions

Wipe all surfaces of the device with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

CAUTION: Damage to the device

Improper cleaning can damage the surfaces of the device.

Penetrating liquid can damage the electronic components inside the device.

- Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.
 - Never use any aggressive or abrasive cleaning agents or aids.
-

Repairs

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim.

The owner is solely liable for damage caused by unauthorized repairs.

In case of repair contact an authorized dealer or our technical service, see "Contact information Heidolph international" on page 43.

Include the completed declaration of no objection with every device return, see "Certificate of decontamination" on page 44.

Maintenance

There are no user-serviceable components in the unit housing. If necessary, in the event of abnormal operating behavior such as excessive noise or heat generation, for example, contact our technical service, see "Contact information Heidolph international" on page 43.

Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment, and biohazardous before disposing.
- Properly remove and dispose residues of substances that are hazardous to health, the environment and biohazardous!

Contact information Heidolph international

Heidolph Instruments North America

Phone: 1-866-650-9604
 E-mail: service@heidolph.com
www.heidolphNA.com

Heidolph Instruments United Kingdom

Phone: 01799 - 5133-20
 E-mail: service@radleys.co.uk
www.heidolph-instruments.co.uk

Local distributors

To find your local distributor please visit www.heidolph.com

Warranty Statement

Heidolph Instruments provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage www.heidolph.com.

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

Certificate of decontamination

Enclose the Certificate of decontamination, duly completed, with your device return. Submissions without a Certificate of decontamination cannot be processed!

CERTIFICATE OF DECONTAMINATION

IN CASE OF RETURNS

Please fill in the required fields.

Note: The sender must package the goods properly and appropriately for transport.

Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12
91126 Schwabach, Germany
Phone: +49 (0) 9122 9920-380
Fax: +49 (0) 9122 9920-19
Email: service@heidolph.de

SENDER

Name _____ First name _____
Company/institution _____ Department _____
_____ Workgroup _____
Address _____
ZC/City _____
Country _____ Phone _____
Email _____

DEVICE DETAILS

Article number _____ Serial no. _____
Ticket number _____
Reason for sending in _____

Has the device been cleaned, decontaminated/disinfected? Yes No (Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous? Yes No (Please mark as applicable)

If yes, with which substances did the device come into contact?

LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

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Authorized to compile the technical file:

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**UK
CA**

Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,
Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12
91126 Schwabach / Germany

Product: Laboratory Temperature Controller
Model:
EKT-HeiCON-2 509-88000-00-x

Description:
Electronic Contact Thermometer for Magnetic Stirrer

The object of the declaration described above is in conformity with the relevant
UK Statutory Instruments (and their amendments):

2008 2016 No. 1091	<i>The Electromagnetic Compatibility Regulations 2016</i>
2012 No. 3032	<i>The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012</i>

and complies with the following technical standards :
EN 61326-1:2013, EN IEC 63000:2018

UK Authorised Representative (for authorities only)
ProductIP (UK) Ltd.
8. Northumberland Av.
London WC2N 5BY

Signed for and on behalf of Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12, 91126 Schwabach / Germany

Schwabach, 21.12.2022

Wolfgang Jaenicke
Managing Director

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices and components by Heidolph Instruments GmbH & Co. KG fulfill the requirements of the EU directive 2011/65/EU (RoHS) and the extension 2015/863. All devices are compatible with the requirements of the RoHS Directive for the following materials:

Max. 0,01% of weight in homogenous material for cadmium and max. 0,1% of the weight in homogenous material for Lead, Mercury, Hexavalent chromium, Polybrominated biphenyls (PBB), Polybrominated diphenyl ethers (PBDE), Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP).

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 06.Okt.2022

W. Jauho

China RoHS Certification

China RoHS DECLARATION OF CONFORMITY

Heidolph Instruments GmbH & Co.KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Instruments GmbH & Co.KG products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
- O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
- X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Instruments GmbH & Co.KG may enter into further devices or can be used together with other appliances.

With these products and appliances in particular, Heidolph Instruments GmbH & Co.KG will not take responsibility for the EFUP of those products and appliances.

Schwabach, 06.08.2021

Wolfgang Jaenicke
Chief Executive Officer CEO

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Instruments GmbH & Co. KG PRODUCTS
© Heidolph Instruments GmbH & Co. KG